

# MB-10 Tablets®

A robust disinfectant, virucidal, tuberculocidal, sterilant in a convenient, easy to use oxidizing tablet.



## Ingredients

### Active ingredients:

Sodium chlorite.....	20.8%
Sodium dichloroisocyanurate dihydrate.....	7.0%
Other ingredients:.....	72.2%
Total.....	100.0%

EPA Reg. No. 70060-19-46269

Always read and follow the label directions and recommendations on the package before using.

## Directions for Use

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Apply a chlorine dioxide use-solution to hard, non-porous surfaces and/or instruments thoroughly wetting surfaces with a cloth, mop, sponge or sprayer or by immersion. Treated surfaces must remain wet for 10 minutes (30 minutes for Foot & Mouth Disease Virus). Wipe dry with a cloth, sponge or mop or allow to air dry. For sprayer applications, use a coarse spray device. Spray 6-8 inches from the surface and spray until totally wet. Do not breathe spray. Allow to air dry.

See label for specific use directions.

## Presentation:

Item#	Description
190452	6g tablets 35-2pks per carton 4 cartons/case

## \*Description:

Use MB-10 Tablets to disinfect pre-cleaned surfaces and to decontaminate instruments and veterinary exam rooms, kennels, waiting rooms, livestock trailers, surgery & necropsy rooms, and any hard non-porous surfaces needing sanitation.

A 100 ppm chlorine dioxide solution with a contact time of 10 minutes has been shown to be effective against:

- Vaccinia virus
- Influenza A virus (Hong Kong)
- Porcine Respiratory & Reproductive Syndrome virus (PRRSV)
- Infectious bursal disease virus (IBDV)
- Marek's disease virus (MDV)
- Avian Influenza A (H3N2) virus
- Porcine circovirus type 2
- Canine Parvovirus
- Hantavirus
- Foot and Mouth Disease Virus after 30 minutes of contact

A 200 ppm chlorine dioxide solution with a contact time of 10 minutes has been shown to be effective against:

### Viruses

- Hepatitis B Virus
- Hepatitis C Virus
- Newcastle Disease virus after 5 minutes of contact
- Norovirus
- Human Coronavirus
- HIV-1 (the virus that causes AIDS) after 1 minute of contact

### Bacteria

- Escherichia coli O157:H7
- Pseudomonas aeruginosa
- Staphylococcus aureus
- Salmonella enterica
- Bordetella bronchiseptica
- Corynebacterium bovis at 100 ppm use-solution
- Helicobacter pylori at 100 ppm use-solution
- Clostridium difficile (vegetative)
- Staphylococcus aureus (MRSA)
- Enterococcus faecalis (VRE) after 5 minutes of contact
- Mycobacterium bovis after 5 minutes of contact at 20 degrees C

\*This technical sheet is intended only as a guide in providing general information associated with the use of this product. Always read and follow the label instructions on the package before using.

# MB-10 Tablets®

Disinfecting & Sterilizing Tablets

*A robust disinfectant in a convenient, easy to use tablet.*

Designed for use in veterinary offices and exam rooms, kennels, animal rearing and confinement facilities.

## Features:

- Effectively removes bio-film
- Convenient, easy to use tablet form
- Rapidly dissolves in water
- Eliminates the hazards associated with mixing
- Broad spectrum efficacy
- Two tablets per package (one tablet per gallon for 100 ppm use-solution)
- Reduces shipping cost and storage space needed
- Delivers consistent solutions time and time again
- Environmentally safe



## Benefits of MB-10 Tablets®:

- Non-corrosive
- Penetrates, kills and strips bio-film
- Powerful and selective deodorizer
- Tolerates organic material
- Fast acting - 100 & 200 ppm solutions proven to kill most viruses and bacteria after 10 minutes



## Also Available...

### MB DeoTrex®-30 Sachets

- Controls odor causing bacteria, mold, mildew and chemical odors in confined spaces where moisture is present.
- Recommended for use in such places as refrigerators, freezers, trashcans, trashbags and basements.



Item#	Description
190453	MB DeoTrex-30 Sachets 50g pk 18pks/case

