




Veratox[®] for Total Milk

Note: Please read kit instructions completely before performing test.

Questions? Call 800/234-5333 or 517/372-9200



1. Add 150 μ L controls and extracted samples to transfer wells.



2. Transfer 100 μ L to the antibody wells. Incubate for 10 minutes. Mix for 20 seconds by sliding back and forth on a flat surface.



3. Dump liquid from antibody wells.



4. Wash wells thoroughly. Repeat wash step wash 10 times times with washing solution.



5. Tap out water on absorbent paper towel.



6. Transfer 100 μ L conjugate from reagent boat to antibody wells using 12-channel pipettor. Incubate for 10 minutes. Mix for 20 seconds by sliding back and forth on a flat surface.



7. Repeat steps 3–5 by dumping out the liquid, thoroughly washing the wells and tapping dry.



8. Transfer 100 μ L substrate from reagent boat to antibody wells using 12-channel pipettor. Incubate for 10 minutes. Mix for 20 seconds by sliding back and forth on a flat surface.



9. Transfer 100 μ L Red Stop from reagent boat to antibody wells. Mix.



10. Read results in a microwell reader.

Product Specifications

Format: Sandwich ELISA

Range of quantitation: 2.5–25 ppm NFDM

Controls provided: 0, 2.5, 5, 10 and 25 ppm NFDM

Testing time: 30 minutes

Ordering Information

8470 Veratox for Total Milk Allergen

Call **800/234-5333** to order or visit us online at www.neogen.com



800-234-5333 (USA/Canada) or 517-372-9200
foodsafety@neogen.com • foodsafety.neogen.com

