*Read instructions carefully before starting test*

# Veratox<sup>®</sup>

## **for Aflatoxin M<sub>1</sub>**

### **THE TOXIN**

Aflatoxin M<sub>1</sub> is a toxin metabolite of aflatoxin B<sub>1</sub>, which is produced by certain strains of the molds *Aspergillus flavus* and *A. parasiticus*. There are four principle types of aflatoxin: B<sub>1</sub>, B<sub>2</sub>, G<sub>1</sub> and G<sub>2</sub>. Aflatoxin B<sub>1</sub> is the most frequently encountered of the group and the most toxic.

Aflatoxin is carcinogenic and may be present in grains and animal feeds. When animals are fed grain or feed contaminated with the toxin, aflatoxin B<sub>1</sub> is converted by hydroxylation to aflatoxin M<sub>1</sub>, which is subsequently secreted in the milk of lactating cows.

Many countries have set maximum allowable levels of aflatoxin M<sub>1</sub> in dairy products. Therefore, accurate determination of the toxin's presence is of major importance to those monitoring the quality of products in which aflatoxin may occur. Testing these commodities for the toxin requires careful sampling, chemical extraction, sanitation, and quantitative analysis.

The U.S. Food and Drug Administration has issued regulatory levels for aflatoxin at 500 parts per trillion (ppt) in milk. In the European Union, the limit for aflatoxin M<sub>1</sub> is 50 ppt in raw milk, heat treated milk, and milk for the manufacture of milk-based products.

### **INTENDED USE**

Veratox<sup>®</sup> for Aflatoxin M<sub>1</sub> is intended for the quantitative analysis of aflatoxin M<sub>1</sub> in commodities such as milk, milk powder, butter and cheese.

### **INTENDED USER**

The test kit is designed for use by quality control personnel and others familiar with commodities possibly contaminated with aflatoxin M<sub>1</sub>. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed training.

### **ASSAY PRINCIPLES**

Veratox for Aflatoxin M<sub>1</sub> is a competitive enzyme-linked immunosorbent assay (ELISA) that allows the user to obtain exact concentrations of aflatoxin M<sub>1</sub> in dairy products in parts per trillion (ppt). Free aflatoxin M<sub>1</sub> in the samples and controls is added to microwells that have been coated with anti-aflatoxin M<sub>1</sub> monoclonal antibody. A wash step will remove any unbound aflatoxin.

During a second incubation, an aflatoxin-HRP conjugate will bind to the remaining free binding site on the antibody. After a wash step, substrate is added and allowed to react with the bound conjugate to produce blue color. More blue color means less aflatoxin M<sub>1</sub>. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form the standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of aflatoxin M<sub>1</sub>.

## STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F).

## MATERIALS PROVIDED

1. Microwell holder
2. 96 antibody-coated microwells
3. 96 red-marked mixing wells
4. 6 yellow-labeled bottles of 0, 5, 15, 30, 60, and 100 ppt aflatoxin controls
5. 1 bottle of sample diluent
6. 1 bottle of 25X wash buffer. Each bottle is enough to prepare 1 L in distilled or deionized water (pH 7.4)
7. 2 blue-labeled bottles of aflatoxin-HRP conjugate solution
8. 1 green-labeled bottle of K-Blue® Substrate solution
9. 1 red-labeled bottle of Red Stop solution

## MATERIALS REQUIRED BUT NOT PROVIDED

1. Scale capable of weighing 1–11 grams
2. Microwell reader with a 650 nm filter (Neogen item 9303)
3. 12-channel pipettor (Neogen item 9273)
4. 100 µL pipettor (Neogen item 9272)
5. Tips for 100 µL and 12-channel pipettors (Neogen item 9410)
6. Paper towels or equivalent absorbent material
7. Plastic bucket for use as waste receptacle
8. Plate rotary shaker
9. Timer (Neogen item 9426)
10. Vortex
11. Waterproof marker
12. Wash bottle (Neogen item 9400)
13. 3 reagent boats for 12-channel pipettor (blue, green and red-labeled) (Neogen item 9435)
14. Centrifuge with rotor for 50 cc tubes
15. 50 cc tubes with caps
16. Distilled or deionized water
17. Dichloromethane (cheese testing only)
18. Hexanes (cheese testing only)
19. 100% methanol (cheese testing only)
20. Chemical fume hood (cheese testing only)
21. High-speed blender (cheese testing only)
22. Grater (cheese testing only)
23. Water bath evaporator with nitrogen gas stream (cheese testing only)
24. Incubator capable of maintaining 50°C (cheese testing only)
25. Glass tube for evaporation (cheese testing only)

## PRECAUTIONS

1. Methanol, dichloromethane and hexanes are highly flammable. Keep container tightly closed and keep away from heat, sparks, open flame and those smoking. They are toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
2. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze.
3. Do not use kit components beyond expiration date.
4. Do not mix reagents from one kit lot with reagents from a different lot.
5. Do not run more than 24 wells per test.
6. Run controls and samples in duplicate.
7. Follow proper pipetting technique (e.g., prime tips and use clean tips).
8. Use of incubation times other than those specified may give inaccurate results.
9. Kits and components should be brought to room temperature 18–30°C (64–86°F) prior to use.
10. Avoid prolonged storage of kits at ambient temperatures.
11. Treat all used liquids, including sample extract, and labware as if contaminated with aflatoxin. Gloves and other protective apparel should be worn at all times.
12. To avoid cross-contamination, use clean pipette tips and glassware for each sample, and thoroughly detoxify and wash all glassware between samples.

## PROCEDURAL NOTES

1. **Substrate.** K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until it is needed.
2. **Conjugate.** The conjugate supplied with this kit is ready to use. The two bottles provided are enough for 96 wells. Cover the reagent boat to keep conjugate protected from direct light and contaminants.
3. **Wash Buffer.** Prepare the wash buffer solution by pouring all the 25X wash buffer concentrate into an empty 1 L container. Rinse the wash buffer concentrate bottle with distilled or deionized water and pour into the 1 L container to ensure all the concentrate is used. Fill the 1 L container with additional distilled or deionized water, and swirl to ensure thorough mixing. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).  
**NOTE:** Discard unused portions of extraction solution and wash buffer when the test kit has been used completely.
4. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.

## SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques (see Neogen's Food Allergen Handbook). Cheese samples should be ground and thoroughly mixed prior to proceeding with the extraction. Store samples at 2–8°C (35–46°F) until analyzed.

### Liquid milk

1. Centrifuge milk samples to separate fat for **10 minutes** at 3500 rpm (2740 x g) at 10°C (50°F).  
**NOTE:** If a refrigerated centrifuge is not available, pre-cool the sample in an ice bath for **10–15 minutes**.
2. After centrifugation, collect the defatted supernatant (bottom layer).  
**NOTE:** Be sure to avoid collecting the upper fatty layer as this will add an additional matrix to the sample, which could adversely impact the outcome of the assay.
3. The supernatant (skimmed milk) will be used directly in the test (100 µL per well).  
**NOTE:** No dilution factor is used.

## Dried milk powders

1. Weigh 10 g of milk powder in a flask. Add 100 mL of distilled water.
2. Stir to dissolve and homogenize.
3. Extract in a rotary shaker for **30 minutes**.  
**NOTE:** If whole milk powder is used, continue with steps 4 and 5. If not testing whole milk powder, proceed to step 6.
4. Centrifuge milk samples to separate fat for **10 minutes** at 3500 rpm (2740 x g) at 10°C (50°F).  
**NOTE:** If a refrigerated centrifuge is not available, pre-cool the sample in an ice bath for **10–15 minutes**.
5. After centrifugation, collect the defatted supernatant (bottom layer).  
**NOTE:** Be sure to avoid collecting the upper fatty layer as this will add an additional matrix to the sample, which could adversely impact the outcome of the assay.
6. The supernatant (skimmed milk) will be used directly in the test (100 µL per well).  
**NOTE:** Because the milk powders were diluted 1:10 during extraction, predicted values need to be multiplied by 10 in order to obtain ppt of aflatoxin M<sub>1</sub> per gram of the powdered milk.

For example, if the predicted reading is 10 ppt, this refers to the level in reconstituted milk. To calculate the amount in milk powder, multiply by 10:

$$10 \text{ ppt} \times 10 = 100 \text{ ppt/g powdered milk}$$

## Cheese

**NOTE:** Carefully grate a representative sample. In cheeses with high amounts of external mold, shave off the outer layer prior to sampling, as the mold will hinder the performance of the assay.

1. Weigh out 2 g of finely grated cheese in a 50 cc centrifuge tube.
2. Add 10 mL of 100% dichloromethane.
3. Vortex for **1 minute**.
4. Extract in a rotary shaker for **15 minutes**.
5. Incubate at 50°C (122°F) for **30 minutes**.
6. Centrifuge at 3500 rpm (2740 x g) for **15 minutes** at 10°C (50°F).
7. Withdraw 5 mL of the dichloromethane supernatant and place in a glass tube.
8. Evaporate dichloromethane to dryness using a water bath evaporator set at 60°C (140°F) with the use of a nitrogen gas stream to help speed evaporation.
9. Re-dissolve the oily residue in 0.5 mL of 100% methanol and 0.5 mL of distilled water, then add 4 mL of 100% hexanes.
10. Vortex for **1 minute**.
11. Centrifuge at 3500 rpm (2740 x g) for **10 minutes** at 10°C (50°F).
12. Remove hexanes (upper layer) completely using a Pasteur pipette.
13. Take a 0.5 mL aliquot of the remaining methanol-water layer and dissolve in 2 mL of the sample diluent.
14. Vortex for **1 minute**.
15. Use this reconstituted material directly in the test (100 µL per well).

**NOTE:** Because the cheeses were diluted 1:5 during extraction, predicted values need to be multiplied by 5 in order to obtain ppt of aflatoxin M<sub>1</sub> per gram of the cheese.

For example, if predicted reading is 10 ppt, this refers to the level in cheese. To calculate the amount in cheese, multiply by 5:

$$10 \text{ ppt} \times 5 = 50 \text{ ppt/g cheese}$$

## TEST PROCEDURE

Allow all reagents to warm to room temperature 8–30°C (64–86°F) before use. All standards and samples will be run in duplicate wells. This will be achieved by transferring standards and samples from 1 red-marked mixing well into 2 antibody-coated microwells.

1. Remove 1 red-marked mixing well for each sample to be tested plus 6 red-marked wells for controls, and place in the well holder.
2. Remove 2 antibody-coated wells per 1 red-marked mixing well to be used. Return antibody wells that will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of strip with a “1” and place strip in the well holder with the marked end on the left. Do not mark the inside or bottom of the wells.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Using a new pipette tip for each, transfer 250 µL of controls and samples to the red-marked mixing wells as described below.

0	5	15	30	60	100	S1	S2	S3	S4	S5	S6	Strip 1
---	---	----	----	----	-----	----	----	----	----	----	----	---------

5. Using a 12-channel pipettor, transfer 100 µL from the above mixing wells to the 2 antibody-coated wells as described below.

0	5	15	30	60	100	S1	S2	S3	S4	S5	S6	Strip 1
---	---	----	----	----	-----	----	----	----	----	----	----	---------

0	5	15	30	60	100	S1	S2	S3	S4	S5	S6	Strip 2
---	---	----	----	----	-----	----	----	----	----	----	----	---------

6. Place the microwell holder on an automatic plate shaker (~600 rpm) and allow continuous shaking for **20 minutes** at room temperature (18–30°C, 64–86°F). Discard the red-marked mixing wells.
7. Shake out the contents of the antibody wells. Fill the wells with diluted wash buffer and dump them out. Repeat this step a total of 5 times, then turn the wells upside-down and tap out on a paper towel until the remaining wash buffer has been removed.
8. Pour the needed volume of conjugate from the blue-labeled bottle into the blue-labeled reagent boat.
9. With new tips on the 12-channel pipettor, prime and pipette 100 µL of conjugate into the wells.
10. Place the microwell holder on an automatic plate shaker (~600 rpm) and allow continuous shaking for **10 minutes** at room temperature (18–30°C, 64–86°F).
11. Shake out the contents of the antibody wells. Fill the wells with diluted wash buffer and dump them out. Repeat this step a total of 5 times, then turn the wells upside-down and tap out on a paper towel until the remaining wash buffer has been removed.
12. Pour the needed volume of substrate from the green-labeled bottle into the green-labeled reagent boat.
13. With new tips on the 12-channel pipettor, prime and pipette 100 µL of substrate into the wells.
14. Place the microwell holder on an automatic plate shaker (~600 rpm) and allow continuous shaking for **15 minutes** at room temperature 18–30°C (64–86°F).
15. Discard remaining substrate and rinse the reagent boat with water.
16. Pour Red Stop Solution from the red-labeled bottle into the red-labeled reagent boat.
17. Eject the excess substrate from the 12-channel pipettor, prime the tips, and pipette 100 µL of Red Stop to each well. Mix by sliding back and forth on a flat surface. Discard the tips.
18. Wipe the bottom of the microwells with a dry cloth or towel and read in a microwell reader using a 650 nm filter. Air bubbles should be eliminated, as they could affect analytical results. Results should be read within **20 minutes** after the addition of Red Stop.
19. Read and calculate results using Neogen’s microwell reader, or an equivalent strip reader. If using another strip/plate reader, calculate results using Neogen’s Veratox for Windows® software.

## RETESTING

If positives occur in commodities not previously tested, confirm with an additional approved method prior to taking action.

## PERFORMANCE CHARACTERISTICS

**Cross-reactivity:** Aflatoxin M<sub>1</sub> 100%. Aflatoxin B<sub>1</sub>, B<sub>2</sub>, G<sub>1</sub>, G<sub>2</sub>, and M<sub>2</sub> < 1%.

**Limit of detection:** 4.3 ppt (Determined by the mean average of 10 aflatoxin M<sub>1</sub>-free samples plus 2 standard deviations.)

**Limit of quantitation:** 5 ppt (Described as the lowest concentration point on the calibration curve that this test can reliably detect aflatoxin M<sub>1</sub>.)

**Range of quantitation:** 5–100 ppt (For quantitating samples above 100 ppt, contact Neogen Technical Services for dilution instructions.)

**Validated matrixes:** Liquid raw milk, nonfat dried milk, dried whole milk powders, cheese and butter.

## CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

## SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at [www.neogen.com](http://www.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

## TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit [www.neogen.com/Corporate/termsconditions.html](http://www.neogen.com/Corporate/termsconditions.html).

## WARRANTY

Neogen makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## TEST KITS AVAILABLE FROM NEOGEN

### Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

### Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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