Read instructions carefully before starting test

Veratox[®]

for Gliadin R5

Refrigerate at 2–8°C (35–46°F) • Do not freeze

GLIADIN/GLUTEN

Gliadin is an alcohol-soluble protein found in wheat that belongs to a group of proteins called prolamins. Other prolamins include secalin, found in rye, and hordein, found in barley. Gluten consists of two groups of proteins (prolamins and glutelins) that are found in differing amounts in wheat, barley, rye and oats.

Gliadin and other prolamins have been identified as major causal agents in a number of disorders, including wheat allergy and gluten intolerance (celiac disease). Wheat allergy is a specific immune response to a number of wheat proteins, including gliadin, albumin, globulin, and glutenin. Celiac disease is a chronic reaction to gluten proteins that results in the poor absorption of nutrients in the small intestine.

Those who must avoid gluten rely upon the correct labeling of food to make appropriate, safe food choices. Testing for the presence of gluten components ensures food manufacturers that an unlabeled—and potentially dangerous—ingredient did not make its way into a food product.

INTENDED USE

Veratox[®] for Gliadin R5 is intended for the quantitative analysis of ingredients, clean-in-place solutions, and finished food products intended to be gluten-free for the presence of gliadin and prolamins found in wheat, barley and rye. Rice flour, cereal, bread and cooked hamburger were validated by AOAC Research Institute. Clean-in-place solutions and other foods were validated internally by Neogen.

INTENDED USER

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by gliadin. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed the Neogen training.

ASSAY PRINCIPLES

Veratox for Gliadin R5 is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Gliadin is extracted from samples with a 60% ethanol solution by shaking in a shaker or rotator. The extract is diluted in Phosphate Buffered Saline (PBS) and diluted samples are added to R5 antibody-coated wells (capture antibody) where gliadin will bind to the antibody during an incubation period. Any unbound gliadin is washed away and a second antibody, which is R5 enzyme-labeled (detector antibody) is added. The detector antibody binds to the gliadin during another incubation period. Unbound enzyme-labeled antibody is washed away and a one step substrate is added. Color develops as a result of the presence of bound-labeled antibody. A stopping reagent is added and the color of the solution is observed. Blue indicates samples containing high levels of gliadin while purple or red samples contain little or no gliadin. The optical densities of the controls form a standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of gliadin in parts per million (ppm).

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F).

MATERIALS PROVIDED

1. 48 antibody-coated microwells
2. 48 red-marked transfer wells
3. 6 yellow-labeled bottles of 0, 2.5, 5, 10, 20, and 40 ppm gliadin controls (5, 10, 20, 40, and 80 ng/mL gliadin)
4. 2 blue-labeled bottles of enzyme-labeled antibody conjugate
5. 1 green-labeled bottle of K-Blue[®] Substrate
6. 1 red-labeled bottle of Red Stop solution
7. 40 mL of 10 mM PBS-Tween washing reagent in a widemouth bottle. Each bottle can prepare 1 L in distilled or deionized water (pH 7.4)
8. 1 foil pouch sample diluent concentrate of 10 mM PBS dry powder containing enough powder to prepare 1 L of dilution buffer
9. 50 g of extraction additive in 2 specimen cups
10. Plastic scoop to measure extraction additive

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Orbital rotator or shaker to hold 50 cc centrifuge tubes for a 1 g sample, or shaker water bath with clamps adjusted to hold 125 mL (4 oz) extraction bottles for a 2 g sample
2. Oven or water bath adjustable to 50°C if analyzing heat-processed samples
3. Scale capable of weighing 0.25 ± 0.01 g if using gliadin renaturing cocktail solution to analyze heat-processed samples
4. Gliadin renaturing cocktail solution for heat processed samples (Neogen item 8515, 8515S, 8515B)
5. Laboratory grade ethanol (190 proof)
6. Microwell reader with a 650 nm filter (Neogen item 9303)
7. Pipettor, 50–200 µL adjustable (Neogen item 9276)
8. Pipettor, 12-channel (Neogen item 9273)
9. Pipette tips (Neogen item 9410, 9417, 9407)
10. Two 1 L bottles to prepare washing solution and sample extract dilution solution (Neogen item 9472)
11. Test tubes to perform sample extract dilution
12. Timer (Neogen item 8426, 9452)
13. 3 reagent boats for 12-channel pipettor (Neogen item 9435)
14. Microwell strip holder (Neogen item 9402)
15. Wash bottle (Neogen item 9400)
16. Paper towels or equivalent absorbent material
17. Waterproof marker
18. Distilled or deionized water
19. Centrifuge (optional)
20. Vortex (Neogen item 9494)

PRECAUTIONS

1. Ethanol solution is highly flammable. Keep container tightly closed, and keep away from heat, sparks, open flame and those smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
2. Components of Veratox for Gliadin R5, such as controls and extraction additive, may contain one or more of the following potentially allergic materials: gluten, casein and soy protein. If allergic to any of these compounds, use caution when using this product.
3. Concentrated food additives, colors and flavors may cause interferences on ELISA test methods. Contact Neogen's technical services for validation information.
4. Sponges should not be used for sample collection and gliadin testing. Sample collection swabs other than Neogen swabs should be validated prior to use. General sponges and swabs may contain solutions or materials that may interfere with the test kits.
5. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze test kits, and avoid prolonged storage at ambient temperatures.
6. Bring kits to room temperature (18–30°C, 64–86°F) prior to use.
7. Do not use kit components beyond expiration date.
8. Do not mix reagents from one kit lot with reagents from a different lot.

9. Do not run more than 24 wells per test.
10. Follow proper pipetting techniques (e.g., prime tips and use clean tips).
11. Use only incubation times specified; others may give inaccurate results.
12. Use clean pipette tips and glassware for each sample to avoid cross-contamination. Thoroughly wash all glassware between samples.

PROCEDURAL NOTES

1. **Sample extract dilution solution (PBS).** Prepare extract dilution solution by adding a foil pouch of sample dilution concentrate, 10 mM PBS, to 1 L distilled or deionized water. Swirl to mix thoroughly.
2. **Controls.** Six controls are provided with this kit. Neogen recommends using all six or a combination of at least five controls with each assay. This combination may vary. One possible combination is to eliminate the 2.5 ppm control to yield results in the range of 5–40 ppm. If needed, contact Neogen Technical Services for more information on the use of these controls.
3. **Wash buffer.** Prepare the wash buffer solution by pouring the wash buffer concentrate into a 1 L container. Add 960 mL of distilled or deionized water. Swirl to ensure thorough mixing.
NOTE: Discard unused portions of extract dilution solution and wash buffer when the test kit has been used completely.
4. **Extraction additive.** To be used with the provided scoop with all the samples utilizing extraction procedure A or B (non-heat processed samples). Heat processed samples (procedure C) should **only** use the extraction additive if samples contain buckwheat, chestnut flour or tannins/phenolic compounds such as chocolate, coffee, cocoa, wine, herbs or fruits.
5. **Substrate.** K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until needed.
6. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.

SAMPLE PREPARATION AND EXTRACTION

For analyzing heat-processed samples or samples of an unknown origin, follow extraction procedure C. For all commodities that were not heat-processed, follow either extraction procedure A or B.

A. Extraction of non-heat processed samples with orbital shaker or rotator

1. Prepare 60% ethanol extraction solution by combining 6 parts ethanol with 4 parts distilled water. Prepare sample extract dilution solution (PBS) as detailed in procedural note 1.
2. Add 1 g ground sample, or 1 mL liquid sample, to a clean 50 cc tube.
3. Add 1 level scoop of extraction additive to the tube.
4. Add 10 mL (9 mL for liquid samples) of 60% ethanol to the tube, cap tightly, then shake the tube vigorously by hand for about **1 minute**, or vortex for **30 seconds**, to ensure complete mixing.
5. Extract by shaking (150 rpm) in an orbital shaker or rotator by laying the tube down on its side over the flat pad of the instrument, and holding it tightly using a rubber band or tape. Rotate or shake for **10 minutes** at room temperature.
6. Centrifuge sample (if necessary) for **10 minutes** at ≥ 2500 g at room temperature.
7. Dilute each sample 1:50 by withdrawing 100 μ L of the upper layer of the extract and transferring it to a small tube or vial containing 4.9 mL of sample extract dilution solution (PBS).
8. To mix, vortex the tube for **5 seconds**.
9. Test diluted samples within **2–3 hours** of extraction.

B. Extraction of non-heat processed samples with shaker or shaker water bath

1. Prepare 60% ethanol extraction solution by combining 6 parts ethanol with 4 parts distilled water. Prepare sample extract dilution solution (PBS) as detailed in procedural note 1.
2. Add 2 g ground sample, or 2 mL liquid sample, to a 125 mL clean extraction bottle.
3. Add 1 level scoop of extraction additive to the bottle.
4. Add 20 mL (18 mL for liquid samples) of 60% ethanol, cap the bottle tightly, then shake vigorously by hand for about **20 seconds** to ensure complete mixing.
5. Extract by shaking (150 rpm) in a shaker for **10 minutes** at room temperature (a shaker water bath can work, but do not turn the heat on). Remove the bottle from shaker or bath.
6. Centrifuge sample (if necessary) for **10 minutes** at ≥ 2500 g at room temperature.
7. Dilute each sample 1:50 by withdrawing 100 μ L of the upper layer of the extract and transferring it to a small tube or vial containing 4.9 mL of sample extract dilution solution (PBS).
8. To mix, vortex the tube for **5 seconds**.
9. Test diluted samples within **2–3 hours** of extraction.

C. Extraction of heat-processed or unknown origin commodities

Heat-processed commodities require the gliadin renaturing cocktail solution (Neogen item 8515) that renatures the heated samples and allows the accurate detection of any possible gliadin in a sample. To extract gliadin from heat-processed samples:

1. Prepare 80% ethanol extraction solution by combining 8 parts ethanol with 2 parts distilled water.
2. Prepare sample extract dilution solution (PBS) as detailed in procedural note 1.
3. Weigh out 0.25 g sample into a 50 cc screw cap centrifuge tube.
4. Add 2.5 mL of renaturing cocktail solution.
5. If samples contain buckwheat, chestnut flour or tannins/phenolic compounds such as chocolate, coffee, cocoa, wine, herbs or fruits, add 1 level scoop of extraction additive. For all other commodity types, do not add extraction additive.
6. Cap and vortex **30 seconds** to homogenize cocktail and sample.
7. Incubate **40 minutes** at 50°C (water bath or oven).
8. Remove samples and let cool for **5–10 minutes**.
9. Add 7.5 mL of 80% ethanol and vortex again for **10–20 seconds**.
10. Shake (150–200 rpm) for **1 hour** at room temperature on a rotator (tube on its side).
11. Centrifuge sample (if necessary) for **10 minutes** at ≥ 2500 g at room temperature.
12. Dilute the sample 1:12.5 into PBS (200 μ L sample into 2.3 mL PBS).
13. To mix, vortex the tube for **5 seconds**.
14. Test diluted samples within **2–3 hours** of extraction.

TEST PROCEDURE

Allow the test kit and all reagents to warm to room temperature (18–30°C, 64–86°F) before using.

1. Remove 1 red-marked mixing well for each sample to be tested plus 1 red-marked well for each control, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells that will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of the strip with a "1," and place strip in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. The controls (see procedural note 2) are supplied ready to use—do not dilute. Using a new pipette tip for each, transfer 150 µL of controls and diluted samples to the red-market mixing wells as shown in one of the templates below. Only run up to two 12-well strips at a time.

If using all 6 controls for a range of 2.5–40 ppm (AOAC method):

| | | | | | | | | | | | |
|----|-----|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 0 | 2.5 | 5 | 10 | 20 | 40 | S1 | S2 | S3 | S4 | S5 | S6 |
| S7 | S8 | S9 | S10 | S11 | S12 | S13 | S14 | S15 | S16 | S17 | S18 |

If using 5 controls for a range of 5–40 ppm:

| | | | | | | | | | | | |
|----|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 0 | 5 | 10 | 20 | 40 | S1 | S2 | S3 | S4 | S5 | S6 | S7 |
| S8 | S9 | S10 | S11 | S12 | S13 | S14 | S15 | S16 | S17 | S18 | S19 |

5. Place tips on the 12-channel pipettor and transfer 100 µL of the controls and sample extracts to the antibody-coated wells. Mix for 20 seconds by sliding the well holder back and forth on a flat surface.
6. Incubate microwells 10 minutes at room temperature (18–30°C, 64–86°F). Discard the red-marked transfer wells.
7. Empty the contents of the wells into a sink. With a wash bottle, fill each antibody well with the wash buffer solution and dump out. Repeat the washing 5 times, then turn the wells upside down and tap out on a paper towel until all washing solution is removed.
8. Pour the needed volume of conjugate from the blue-labeled bottle into a clean reagent boat.
9. Using the 12-channel pipettor and new tips, transfer 100 µL of the conjugate into all the wells and mix for 20 seconds by sliding the well holder back and forth on a flat surface.
10. Incubate for 10 minutes at room temperature (18–30°C, 64–86°F).
11. Wash all wells with the wash buffer solution as described in step 7.
12. Pour the needed volume of substrate solution from the green-labeled bottle into a clean reagent boat.
13. Place new tips on the 12-channel pipettor and transfer 100 µL of substrate into each well and mix for 20 seconds. Do not eject tips.
14. Incubate for 10 minutes at room temperature (18–30°C, 64–86°F).
15. Pour the needed volume of Red Stop solution from the red-labeled bottle into a clean reagent boat.
16. With the same tips used to dispense the substrate, transfer 100 µL of Red Stop into each well and mix for 20 seconds.
17. Wipe the bottom of the microwells and read in a microwell reader with a 650 nm filter.
18. Interpret the test's results using the Neogen 4700 microwell reader, or an equivalent strip reader. If using a strip reader, calculate the results using Neogen's Veratox Software for Windows.

INTERPRETATION OF RESULTS

Standard controls were made from wheat gliadin and calculated as gliadin. Approximately 50% of the gluten is available as gliadin. Therefore, to calculate the gluten value of the samples, multiply the ppm gliadin results by 2.

PERFORMANCE CHARACTERISTICS

Limit of quantitation: From AOAC validation: Cereal–5.71; cooked hamburger–3.15; bread–1.84; and rice flour–3.31. Neogen internally validated commodity LOQs were all under 2.5 ppm (described as the lowest concentration point on the calibration curve that this test can reliably detect gliadin).

Limit of detection: From AOAC validation: Cereal–2.16; cooked hamburger–1.21; bread–0.63; and rice flour–1.23. Neogen internally validated commodity LODs were all under 2.5 ppm (defined as the mean value of 10 negative samples + 3.3 standard deviations).

Range of quantitation: 2.5–40 ppm (For quantitating samples above 40 ppm, contact a Neogen representative for dilution instructions).

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/en/terms-and-conditions.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustacea, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnut, multi-treenut

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