Clean-Trace™ Test

-  **(EN)** Water Positive Control (30ng)
-  **(FR)** Témoin positif pour l'eau (30 ng)
-  **(DE)** Wasser-Positivkontrolle (30 ng)
-  **(IT)** Controllo positivo dell'acqua (30 ng)
-  **(ES)** Control positivo para agua (30 ng)
-  **(NL)** Waterpositieve controle (30 ng)
-  **(SV)** Vatten positiv-test (30 ng)
-  **(DA)** Vand Positiv Kontrol (30 ng)
-  **(NO)** Positiv kontroll til vann (30 ng)
-  **(FI)** Water positiivinen kontrolli (30 ng)
-  **(PT)** Controle Positivo de Água (30 ng)
-  **(EL)** Θετικής Επιβεβαίωσης (θετικός μάρτυρας) για Νερό (30 ng)
-  **(PL)** Kontrola dodatnia wody (30 ng)
-  **(HU)** Víz pozitív kontroll (30 ng)
-  **(CS)** Pozitivní kontrola vody (30 ng)
-  **(RO)** Control pozitiv al apei (30 ng)
-  **(RU)** Положительный контроль воды (30 нг)
-  **(TR)** Su Pozitif Kontrolü (30 ng)
-  **(JA)** 水中ポジティブコントロール (30ng)
-  **(ZH)** 水质采样棒阳性对照 (30ng)

Clean-Trace™ Test

Product Instructions

Water Positive Control (30ng)

Product Description and Intended Use

The measurement of the standard concentration of Adenosine Triphosphate (ATP) by 3M™ Clean-Trace™ Water ATP Test Kits is used as an assurance of reagent and instrument function. In addition, the degree of quenching caused by interfering substances in the water sample can be determined using both the 3M Clean-Trace Water ATP Test Kits and 3M™ Clean-Trace™ Biomass Detection Kit. The test kit contains 10 ready to use reagent vials. Each vial contains 30 ng or 5×10^{-11} moles of ATP. Also included are 10 vials of ATP free water to re-hydrate the ATP.

3M Clean-Trace Water ATP Test Kits – A 3M™ Clean-Trace™ Water ATP Sampling Stick is used to collect a sample of re-hydrated ATP positive control. Upon activation of the test, reagent in the cuvette of the device reacts with ATP collected on the sample collection rings to produce light. The intensity of light is proportional to the amount of ATP. Measurements are performed in the 3M™ Clean-Trace™ Luminometer and results displayed on the digital read out in Relative Light Units (RLU).

3M Clean-Trace Biomass Detection Kit – The vial of ATP positive control reagent is re-hydrated with 500 μ L of ATP free water. A reagent and instrument assurance test is then undertaken by mixing ATP free water, extractant, enzyme and ATP standard. The intensity of light is proportional to the amount of ATP. Measurements are performed in the 3M Clean-Trace Luminometer and results displayed on the digital read out in relative light units (RLU).

The signal obtained should be above a specified level, (see below) which indicates that both the reagents and instrument are operating correctly. The unique format of this product provides the user with a convenient method of assuring that all components used during testing are operating correctly, and that the operator is using the correct method. In some cases, the water samples under test can contain components that reduce or “quench” the activity of the reagents, thus reducing the bioluminescent signal. The degree of quenching can be determined by repeating the assurance test detailed in the descriptions above by replacing the ATP free water with a sample of the water under test. The test indicates quenching if the result obtained is lower than the result achieved during the assurance test. Alternatively, some substances can enhance light output.

This kit has been designed for industrial and environmental applications and contains a higher level of ATP to give an appropriate signal for this function.

⚠ Warnings

Consult the Safety Data Sheet for additional information and local regulations for disposal. ATP is a common substance. It is present on skin, hair and on many surfaces. Do not touch the swab or the interior surfaces of the positive control vials. Do not use devices past expiration date.

Notice

To reduce the risk of inaccurate results, 3M™ Clean-Trace™ Water Positive Control Kit should be stored at 2°-8°C (36°-46°F). To avoid inaccurate results, keep reconstituted 3M Clean-Trace Water Positive Control Kit at refrigerated temperatures 2°-8°C (36°-46°F). To reduce the risks associated with cross contamination while preparing the assay the use of gloves is recommended. To reduce the risk of inaccurate results do not shake the reconstituted reagents. Follow instructions for reconstituting reagents.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.3M.com/foodsafety, or contact your local 3M representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user’s responsibility in selecting any test method to evaluate a sufficient number of samples with the appropriate matrices and challenges to satisfy the user that the chosen test method meets the user’s criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of this product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, 3M DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any 3M Food Safety Product is defective, 3M or its authorized distributor will, at its option, replace, or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify 3M within sixty days of discovery of any suspected defects in a product and return it to 3M. Please call Customer Service (1-800-328-1617 in the U.S.) or your official 3M Food Safety representative for a Returned Goods Authorization.

Limitation of 3M Liability

3M WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall 3M's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage and Disposal

For maximum shelf life, store between 2°-8°C (36°-46°F). After use, may contain microorganisms that may be a potential biohazard. Follow current industry standards for disposal. Do not freeze. Protect from direct sunlight.

Instructions for Use

Before the 3M Clean-Trace Water Positive Control Kit is used, it is important to ensure your 3M Clean-Trace Luminometer is switched on and initialized. Refer to the manual provided with the instrument for full details.

3M™ Clean-Trace™ Water Plus - Total ATP and 3M™ Clean-Trace™ Water - Free ATP Test Applications

1. Place the 3M Clean-Trace tests devices and positive control vials at room temperature for at least 10 minutes before use.
2. Gently tap the positive control vial to ensure that the reagent pellet is at the bottom.
3. Carefully remove the stopper from the vial, avoid contaminating the opening to the vial.
4. Use a disposable Pasteur pipette or laboratory micropipette to transfer 2 mL of the ATP free water into the ATP vial.
5. Swirl the vial gently for 5 seconds.
6. Use a disposable Pasteur pipette or laboratory micropipette to transfer ALL the solution from the ATP vial back into the bottle containing the ATP free water.

Note: Most 1 mL laboratory micropipettes and some disposable 3 mL Pasteur pipettes will not be able to transfer all the solution because they cannot reach the bottom of the ATP vial due to their diameter. Therefore, a smaller volume micropipette or Pasteur pipette would also be required.

7. Swirl to mix the sample. Remove a 3M Clean-Trace Water Plus - Total ATP test device from the foil pouch. Remove the sample stick from the test device and immerse the sample rings into the liquid under test, tapping the handle gently if bubbles form. Immediately remove the sample stick from the solution to be tested and carefully return the sample stick to the device such that the handle is at its starting position.
8. Process the sampling stick as per the 3M Clean-Trace Water Plus - Total ATP or 3M Clean-Trace Water - Free ATP product instructions. A maximum of five devices may be tested using a single reconstituted vial of ATP positive control. The reconstituted ATP positive control should be used within 6 hours if stored at ambient temperature (21°C, 70°F) or 24 hours if stored at 2-8°C (36° - 46°F).

Test Procedure: Quench check

- Take a bottle of ATP free water and discard its contents.
- Fill this bottle with 10 mL of the water sample under test.
- Use this sample water to reconstitute an ATP vial, as above.
- Test with a 3M Clean-Trace Water ATP test device.

3M Clean-Trace Biomass Detection Kit Applications

1. Prepare the 3M Clean-Trace Biomass Detection Kit reagents as per the product instructions.
2. Gently tap the positive control vial to ensure that the reagent pellet is at the bottom.
3. Carefully remove the stopper from the vial, avoid contaminating the opening to the vial.
4. Use a disposable Pasteur pipette or laboratory micropipette transfer 500 μL of the ATP free water into the ATP vial.
5. Swirl the vial gently for 5 seconds.
6. Mix together in a cuvette:
 - 100 μL of ATP Free water
 - 100 μL of extractant
 - 100 μL of enzyme
 - 20 μL of ATP standard
7. Read the level of bioluminescence in your instrument.

Test Procedure: Quench check

1. Repeat the Assay Assurance test then mix together:
 - 100 μL of the water sample under test
 - 100 μL of extractant
 - 100 μL of enzyme
2. Read the level of bioluminescence in your instrument.
3. Add of 20 μL ATP standard.
4. Re-read the level of bioluminescence in your instrument.

Interpretation of Results

Assay Assurance Test

Testing using the 3M Clean-Trace Water ATP Test Kits should give readings greater than 10,000 RLU.

Testing using the 3M Clean-Trace Biomass Detection Test Kit should give readings greater than 20,000 RLU. If the tests fail to reach these levels repeat the test procedure using a fresh vial of ATP positive control reagent. If low results are again achieved it may indicate a possible problem with either the reagent and/or instrument.

Quench check

Comparing the results of the 3M Clean-Trace Water ATP Test Kits and 3M Clean-Trace Biomass Detection Kit with the result for the respective quench check can derive the degree of quenching. The degree of quenching should be taken into account when interpreting your readings for that particular water sample.

Trouble Shooting Guide

The following section is a list of potential problems that may cause a low signal:

- Incomplete mixing - shake the device and read again
- Reagents have passed their expiration date
- Reagents stored incorrectly
- Reconstituted ATP is contaminated
- Reagents not allowed to reach room temperature before testing
- Contaminated luminometer
- Damaged luminometer

For further information, contact 3M Food Safety Technical Services.

Explanation of Symbols

www.3M.com/foodsafety/symbols

3M Food Safety

3M United States

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-328-6553

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-563-2921

3M Latin America

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-954-340-8263

3M Europe and MEA

3M Deutschland GmbH
Carl-Schurz-Strasse 1
D41453 Neuss/Germany
+49-2131-14-3000

3M United Kingdom PLC

Morley Street, Loughborough
Leicestershire
LE11 1EP
United Kingdom
+(44) 1509 611 611

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Euro Plaza
Gebaude J, A-1120 Wien
Kranichberggasse 4
Austria
+(43) 1 86 686-0

3M Asia Pacific

No 1, Yishun Avenue 7
Singapore, 768923
65-64508869

3M Japan

3M Health Care Limited
6-7-29, Kita-Shinagawa
Shinagawa-ku, Tokyo
141-8684 Japan
81-570-011-321

3M Australia

Bldg A, 1 Rivett Road
North Ryde, NSW 2113
Australia
61 1300 363 878

3M Health Care

2510 Conway Ave
St. Paul, MN 55144 USA
www.3M.com/foodsafety

© 2018, 3M. All rights reserved.
3M and Clean-Trace are trademarks of 3M. Used under license in Canada.
34-8723-5573-9

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3M Österreich GmbH