Read instructions carefully before starting test

Veratox[®]

for Aflatoxin

Quantitative Test

FGIS-2018-113

Refrigerate at 2–8°C • Do not freeze

THE TOXIN

Aflatoxin is a toxic and carcinogenic substance produced by certain strains of the molds *Aspergillus flavus* and *A. parasiticus*. There are four principal types of aflatoxin: B₁, B₂, G₁ and G₂. Aflatoxin B₁ is the most frequently encountered of the group and the most toxic. The commodities most affected by aflatoxin are corn, peanuts, cottonseed, milo and the majority of tree nuts.

The effects in animals of ingesting excessive amounts of the toxin range from chronic health and performance problems to death. Aflatoxin has been shown to cause liver damage or cancer, decreased milk and egg production, immune suppression and interference with reproductive efficiency.

The U.S. Food and Drug Administration (FDA) has set maximum allowable levels of aflatoxin in food and feed. Therefore, accurate determination of the presence of the toxin is of major importance to those monitoring the quality of food and feed in which aflatoxin may occur. Testing these commodities for the toxin requires careful sampling, chemical extraction, sanitation and quantitative analysis.

The FDA has issued regulatory levels for aflatoxin as follows:

For	Level	Commodities
Humans	20 ppb	All food except milk
All animal species	20 ppb	All feed (exceptions below)
Exceptions:		
Breeding cattle, breeding swine, mature poultry	100 ppb	Corn
Finishing swine (>100 lbs)	200 ppb	Corn
Finishing beef cattle	300 ppb	Corn
All animal species	300 ppb	Cottonseed meal used in feed

INTENDED USE

Veratox[®] for Aflatoxin is intended for the quantitative analysis of aflatoxin in commodities such as barley, corn, corn meal, corn gluten meal, corn/soy blend, cottonseed, cottonseed meal, DDGS, hominy, milo, peanuts, pet food, popcorn, rice, soy meal, walnuts and wheat.

INTENDED USER

This test kit is designed for use by quality control personnel and others familiar with food and feed possibly contaminated by aflatoxin. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed the Neogen training.

ASSAY PRINCIPLES

Veratox for Aflatoxin is a direct competitive ELISA in a microwell format, which allows the user to obtain exact concentrations in parts per billion (ppb). Free aflatoxin in the samples and controls is allowed to compete with enzyme-labeled aflatoxin (conjugate) for the antibody binding sites. After a wash step, substrate is added, which reacts with the bound conjugate to produce blue color. More blue color means less aflatoxin. The test is read in a microwell reader to yield optical densities.

The optical densities of the controls form the standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of aflatoxin.

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F).

MATERIALS PROVIDED

1. 48 antibody-coated microwells
2. 48 red-marked mixing wells
3. 4 yellow-labeled bottles of 0, 5, 15 and 50 ppb aflatoxin controls (see precautions for handling of methanol solution)
4. 1 blue-labeled bottle of aflatoxin-HRP conjugate solution
5. 1 green-labeled bottle of K-Blue[®] Substrate solution
6. 1 red-labeled bottle of Red Stop solution

MATERIALS REQUIRED BUT NOT PROVIDED

1. Extraction materials (items c through e available in kit form from Neogen, item 8052):
 - a. 70% ACS Grade methanol (Neogen item 8055, 8056)
 - b. 250 mL graduated cylinder (Neogen item 9368)
 - c. Container with 125 mL capacity (Neogen item 9428)
 - d. Neogen filter syringes, Whatman #1 filter paper, or equivalent (Neogen item 9420, 9430)
 - e. Sample collection tubes (Neogen item 9421)
2. High-speed blender (Neogen item 9493, 9477)
3. Agri-Grind grinder or equivalent (Neogen item 9401, 9453)
4. Scale capable of weighing 5–50 grams (Neogen item 9427)
5. Microwell reader with a 650 nm filter (Neogen item 9303)
6. Pipettor, 12-channel (Neogen item 9273)
7. Pipettor, 100 µL (Neogen item 9272, 9290)
8. Pipette tips for 100 µL and 12-channel pipettors (Neogen item 9410, 9407, 9417)
9. Paper towels or equivalent absorbent material
10. Plastic bucket for use as waste receptacle
11. Microwell holder (Neogen item 9402)
12. Timer (Neogen item 9426)
13. Waterproof marker
14. Wash bottle (Neogen item 9400)
15. 2 reagent boats for 12-channel pipettor (Neogen item 9450)
16. Distilled or deionized water

PRECAUTIONS

1. Methanol solution is highly flammable. Keep container tightly closed, and keep away from heat, sparks, open flame and those smoking. It is toxic if swallowed or if vapor is inhaled. Avoid contact with skin.
2. Store test kit refrigerated, between 2–8°C (35–46°F) when not in use. Do not freeze.
3. Kits should be brought to room temperature (18–30°C, 64–86°F) prior to use.
4. Avoid prolonged storage of kits at ambient temperatures.
5. Do not use kit components beyond expiration date.
6. Do not mix reagents from one kit serial with reagents from a different kit serial.
7. Do not run more than 24 wells per test.
8. Follow proper pipetting techniques, including priming of tips.
9. Use of incubation times other than those specified may give inaccurate results.
10. Treat all used liquids, including sample extract, and labware as if contaminated with aflatoxin. Gloves and other protective apparel should be worn at all times.
11. To avoid cross-contamination, use clean pipette tips and glassware for each sample, and thoroughly detoxify and wash all glassware between samples.
12. Commodities tested should have a pH of 6–8. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH, contact your Neogen representative or Technical Services.

PROCEDURAL NOTE

1. K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until it is needed.
2. **Antibody wells:** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test is set to begin.

SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques. The sample should be ground and thoroughly mixed prior to proceeding with the extraction. Store samples at 2–8°C (35–46°F) until analyzed. **NOTE:** If you are using Neogen's Mycotoxin Extraction Kit, follow the instructions in that kit for the extraction procedure. If you are preparing your own extraction solution, continue with the instructions that follow.

1. If not using Neogen's prepared solution, prepare a 70% methanol solution by mixing 7 parts ACS-grade methanol with 3 parts distilled or deionized water for each sample to be tested.
2. Obtain a representative sample. Grind the entire sample so that at least 75% of the ground material passes through a 20 mesh sieve, the particle size of a fine instant coffee.
3. Vigorously shake, using hand or mechanical means, 5 grams of ground sample in 25 mL of 70% methanol for 3 minutes.

AOAC Method: Vigorously shake, using hand or mechanical means, 50 grams ground sample in 250 mL of 70% methanol for 3 minutes.

USDA/GIPSA Method: Blend 50 grams of ground sample with 250 mL of 70% methanol for 1 minute in a high-speed blender.

4. Filter the extract by pouring at least 5 mL through a Whatman #1 filter (or Neogen filter syringe) and collecting the filtrate as a sample.
5. The sample is now ready for testing.

TEST PROCEDURE

Allow all reagents to warm to room temperature (18–30°C, 64–86°F) before use.

1. Remove 1 red-marked mixing well for each sample to be tested plus 4 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells that will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of strip with a “1,” and place strip in the well holder with the marked end on the left. Do not mark the inside or bottom of the wells.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Place 100 µL of conjugate from the blue-labeled bottle in each red-marked mixing well.
5. Using a new pipette tip for each, transfer 100 µL of controls and samples to the red-marked mixing wells as described below.

0	5	15	50	S1	S2	S3	S4	S5	S6	S7	S8	Strip 1
---	---	----	----	----	----	----	----	----	----	----	----	---------

S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19	S20	Strip 2
----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	---------

6. Using a 12-channel pipettor, mix the liquid in the wells by pipetting it up and down 3 times. Transfer 100 µL to the antibody-coated wells. Discard the red-marked mixing wells.
7. Set timer for **2 minutes**, mixing the wells for the first **10–20 seconds** of the room temperature incubations by sliding the microwell holder back and forth on a flat surface without splashing reagents from wells.
8. Shake out the contents of the antibody wells. Fill the wells with distilled or deionized water and dump them out. Repeat this step 5 times, then turn the wells upside-down and tap out on a paper towel until the remaining water has been removed.
9. Pour the needed volume of substrate from the green-labeled bottle into the green-labeled reagent boat.
10. With new tips on the 12-channel pipettor, prime and pipette 100 µL of substrate into the wells.
11. Set timer for **3 minutes**, mixing the wells for the first **10–20 seconds** by sliding back and forth on a flat surface. Discard remaining substrate and rinse the reagent boat with water.
12. Pour Red Stop solution from the red-labeled bottle into the red-labeled reagent boat.
13. Eject the excess substrate from the 12-channel pipettor, prime the tips, and pipette 100 µL of Red Stop to each well. Mix by sliding back and forth on a flat surface. Discard the tips.
14. Wipe the bottom of the microwells with a dry cloth or towel and read in a microwell reader using a 650 nm filter. Air bubbles should be eliminated, as they could affect analytical results. Results should be read within **20 minutes** after the addition of Red Stop.
15. Read and calculate results using Neogen’s Stat Fax microwell reader, or equivalent. If using a strip/plate reader, calculate results using Neogen’s Veratox software.

RETESTING

If positives occur in commodities not previously tested, confirm with an additional approved method prior to taking action.

PERFORMANCE CHARACTERISTICS

Limit of detection: 1.4 ppb (Determined by the mean average of 10 aflatoxin free samples plus 2 standard deviations.)

Limit of quantitation: 5 ppb (Described as the lowest concentration point on the calibration curve that this test can reliably detect aflatoxin.)

Range of quantitation: 5–50 ppb (For quantitating samples above 50 ppb contact Neogen Technical Services for dilution instructions.)

AOAC validated matrices: corn

Neogen validated matrices: Alfalfa meal, ammoniated corn, ammoniated cottonseed meal*, barley, barley flour, barley lightly pearled, beet pulp*, canola meal, coconut, copra, corn, corn bran*, corn germ meal, corn gluten meal, corn grits, corn meal, corn oil, corn silage**, corn/soy blend, corn starch, corn steep*, cottonseed, cottonseed meal, DDGs*, DDGs wet cake and syrup*, figs, flaxseed meal, hazelnut, hominy, kamut, lentils, malted barley flour, milled rice, millet, milo, oats, oat fiber, oat hulls*, oats (naked), pea fiber, peanut meal*, peanuts (raw or roasted), peanut hulls, pet food*, pistachio, popcorn, potato (white), potato with skins powder*, pumpkin seeds, quinoa, raw flour, rice, rice bran, rice gluten, rice hulls, rye (raw), rye flour, sorghum (grain), soy flour, soy germ meal, soy hydrolysate, soybean meal, soybeans, soy pallet, sunflower meal, tapioca, wheat, wheat bran*, wheat bran aleurone, wheat flour 2nd clear, wheat middlings and wheat midds

**May require a pH adjustment*

***Contact Neogen for special extraction procedure.*

NOTE: Neogen continues to validate new commodities. Please contact a representative for the latest validated commodity list.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/en/terms-and-conditions.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TESTING KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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