Read instructions carefully before starting test

Quantitative Aflatoxin High Sensitivity Test

REFRIGERATE AT 2-8°C • DO NOT FREEZE

THE TOXIN

Aflatoxin is a toxic and carcinogenic substance produced by certain strains of the molds *Aspergillus flavus* and *A. parasiticus*. There are four principle types of aflatoxin: B₁, B₂, G₁, and G₂. Aflatoxin B₁ is the most frequently encountered of the group and the most toxic. The commodities most affected by aflatoxin are corn, peanuts, cottonseed, milo, and the majority of tree nuts.

The effects in animals of ingesting excessive amounts of the toxin range from chronic health and performance problems to death. Aflatoxin has been shown to cause liver damage or cancer, decreased milk and egg production, immune suppression and interference with reproductive efficiency.

The FDA has set maximum allowable levels of aflatoxin in food and feed. Therefore, accurate determination of the presence of the toxin is of major importance to those monitoring the quality of food and feed in which aflatoxin may occur. Testing these commodities for the toxin requires careful sampling, chemical extraction, sanitation and quantitative analysis.

The FDA has issued regulatory levels for aflatoxin as follows:

For	Level	Commodities
Humans	20 ppb	All food except milk
All animal species	20 ppb	All feed (exceptions below)
Exceptions:		
Breeding cattle, breeding swine, mature poultry	100 ppb	Corn
Finishing swine (>100 lbs.)	200 ppb	Corn
Finishing beef cattle	300 ppb	Corn
Finishing beef cattle, swine, poultry	300 ppb	Cottonseed meal

INTENDED USE

Veratox for Aflatoxin HS is intended for the quantitative analysis of aflatoxin in such commodities as corn, cornmeal, corn germ meal, corn gluten meal, milo, peanuts, peanut butter, popcorn, rice, tree nuts and wheat.

INTENDED USER

The test kit is designed for use by quality control personnel and others familiar with food and feed possibly contaminated by aflatoxin. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed the training.

ASSAY PRINCIPLES

Veratox for Aflatoxin HS is a competitive direct enzyme-linked immunosorbent assay (CD-ELISA) which allows the user to obtain exact concentrations in parts per billion (ppb). Free aflatoxin in the samples and controls is allowed to compete with enzyme-labeled aflatoxin (conjugate) for the antibody binding sites. After a wash step, substrate is added, which reacts with the bound conjugate to produce blue color. More blue color means less aflatoxin. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form the standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of aflatoxin.

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2-8°C (35-46°F).

MATERIALS PROVIDED

- 48 antibody-coated microwells
- 48 red-marked mixing wells
- 5 yellow-labeled bottles of 0, 1, 2, 4 and 8 ppb aflatoxin controls (see precautions for handling of methanol solution)
- 1 blue-labeled bottle of aflatoxin-HRP conjugate solution
- 1 green-labeled bottle of K-Blue® Substrate solution
- 1 red-labeled bottle of Red Stop solution

MATERIALS REQUIRED BUT NOT PROVIDED

- Extraction materials (items c through e available in kit form from Neogen, item #8052):
 - 70% ACS Grade methanol (Neogen item #8055)
 - 250 mL graduated cylinder (Neogen item #9368)
 - Container with 125 mL capacity
 - Neogen filter syringes, Whatman #1 filter paper, or equivalent
 - Sample collection tubes
- High-speed blender
- Agri-Grind grinder or equivalent (Neogen item #9401)
- Scale capable of weighing 5-25 grams (Neogen item #9427)
- Microwell reader with a 650 nm filter (Neogen item #9301/#9302)
- 12-channel pipettor (Neogen item #9273)
- 100 µL pipettor (Neogen item #9272)
- Tips for 100 µL and 12-channel pipettors (Neogen item #9410)
- Paper towels or equivalent absorbent material
- Plastic bucket for use as waste receptacle
- Microwell holder (Neogen item #9402)
- Timer (Neogen item #9426)
- Waterproof marker
- Wash bottle (Neogen item #9400)
- 2 reagent boats for 12-channel pipettor (Neogen item #9435)
- Distilled or deionized water

PRECAUTIONS

- Methanol solution is highly flammable. Keep container tightly closed, and keep away from heat, sparks, open flame and those smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
- Store test kit between 2-8°C (35-46°F) when not in use, do not freeze.
- Do not use kit components beyond expiration date.
- Do not mix reagents from one kit serial with reagents from a different kit serial.
- Do not run more than 24 wells per test.
- Follow proper pipetting techniques, including priming of tips.
- Use of incubation times other than those specified may give inaccurate results.
- Kits should be brought to room temperature (18-30°C, 64-86°F) prior to use.
- Avoid prolonged storage of kits at ambient temperatures.
- Treat all used liquids, including sample extract, and labware as if contaminated with aflatoxin. Gloves and other protective apparel should be worn at all times.
- To avoid cross-contamination, use clean pipette tips and glassware for each sample, and thoroughly detoxify and wash all glassware between samples.
- Commodities tested should have a pH of 6-8. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH contact your Neogen representative or Technical Services.

PROCEDURAL NOTE

K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until it is needed. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.

SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques. The sample should be ground and thoroughly mixed prior to proceeding with the extraction. Store samples at 2-8°C (35-46°F) until analyzed. **Note:** If you are using Neogen's Mycotoxin Extraction Kit, follow the instructions in that kit for the extraction procedure. If you are preparing your own extraction solution, continue with the instructions that follow.

- If not using Neogen's prepared solution, prepare a 70% methanol solution by mixing 7 parts ACS Grade methanol with 3 parts distilled or deionized water for each sample to be tested.
- Obtain a representative sample. Grind the entire sample so that at least 75% of the ground material passes through a 20 mesh sieve, the particle size of a fine instant coffee.
- Blend 25 grams of ground sample with 125 mL of 70% methanol for 2 minutes in a high-speed blender. Alternative method: Add 5 grams of ground sample to 25 mL of 70% methanol and shake vigorously for 3 minutes.
- Filter the extract by pouring at least 5 mL through a Whatman #1 filter (or Neogen filter syringe) and collecting the filtrate as a sample.
- The sample is now ready for testing.

TEST PROCEDURE

Allow all reagents to warm to room temperature (18-30°C, 64-86°F) before use.

- Remove 1 red-marked mixing well for each sample to be tested plus 5 red-marked wells for controls, and place in the well holder.
- Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of strip with a "1", and place strip in the well holder with the marked end on the left. Do not mark the inside or bottom of the wells.

- Mix each reagent by swirling the reagent bottle prior to use.
- Place 100 µL of conjugate from the blue-labeled bottle in each red-marked mixing well.
- Using a new pipette tip for each, transfer 100 µL of controls and samples to the red-marked mixing wells as described below.

0	1	2	4	8	S1	S2	S3	S4	S5	S6	S7	Strip 1
S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19	Strip 2

- Using a 12-channel pipettor, mix the liquid in the wells by pipetting it up and down 3 times. Transfer 100 µL to the antibody-coated wells. Mix by sliding the microwell holder back and forth on a flat surface for 10-20 seconds without splashing reagents from the wells. Incubate for **10 minutes** at room temperature (18-30°C, 64-86°F). Discard the red-marked mixing wells.
- Shake out the contents of the antibody wells. Fill the wells with distilled or deionized water and dump them out. Repeat this step 5 times, then turn the wells upside-down and tap out on a paper towel until the remaining water has been removed.
- Pour the needed volume of substrate from the green-labeled bottle into the green-labeled reagent boat.
- With new tips on the 12-channel pipettor, prime and pipette 100 µL of substrate into the wells. Mix by sliding back and forth on a flat surface for 10-20 seconds.
- Incubate **10 minutes**. Discard remaining substrate and rinse the reagent boat with water.
- Pour Red Stop solution from the red-labeled bottle into the red-labeled reagent boat.
- Eject the excess substrate from the 12-channel pipettor, prime the tips, and pipette 100 µL of Red Stop to each well. Mix by sliding back and forth on a flat surface. Discard the tips.
- Wipe the bottom of the microwells with a dry cloth or towel and read in a microwell reader using a 650 nm filter. Air bubbles should be eliminated, as they could affect analytical results. Results should be read within 20 minutes after the addition of Red Stop.
- Read and calculate results using Neogen's Awareness Stat Fax microwell reader, or equivalent. If using an EL301 reader or other strip/plate reader, calculate results using Neogen's Veratox for Windows software.

RETESTING

If positives occur in commodities not previously tested, confirm with an additional approved method prior to taking action.

PERFORMANCE CHARACTERISTICS

Limit of detection: 0.5 ppb (Determined by the mean average of 10 aflatoxin free samples plus 2 standard deviations.)

Limit of quantitation: 1 ppb (Described as the lowest concentration point on the calibration curve that this test can reliably detect aflatoxin.)

Range of quantitation: 1 - 8 ppb (For quantitating samples above 8 ppb contact Neogen Technical Services for dilution instructions.)

Validated matrixes: Ammoniated corn, barley, beat pulp*, coconut, copra, corn, corn bran*, corn grits, corn/soy blend, corn starch, cornmeal, corn germ meal, corn gluten meal*, cotton seed, cotton seed meal, DDGS*, figs, flaxseed meal, hominy, kamut, lentils, oat fiber, milo, raw peanuts**, roasted peanuts**, peanut butter, peanut meal**, peanut hulls, pet food, popcorn, quinoa, rice, rice hulls, rice bran, soy flour, soy germ meal, soy meal, tree nuts, wheat, and wheat midds.

*May require pH adjustment.

**Contact Neogen for special procedure.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Service can be reached between 8 a.m. and 6 p.m. Eastern time by calling 800/234-5333 or 517/372-9200 and asking for a Neogen sales representative or Technical Services. Assistance is available on a 24-hour basis by calling 800/867-0308. Training on this product, and all Neogen test kits is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's Food Safety test kits, on Neogen's Web site at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TESTING KITS AVAILABLE FROM NEOGEN**Natural Toxins**

- Aflatoxin, DON, Ochratoxin, Zearalenone, T-2 Toxin, Fumonisin, Histamine

Foodborne Bacteria

- E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, Yeast and Mold, Total Plate Count, Generic *E. coli* and Total Coliforms, Protein Residues

Food Allergens

- Peanuts, Milk, Eggs, Almonds, Gluten, Soy Flour, Hazelnut

Genetic Modification

- CP4 (Roundup Ready®)

Ruminant By-products

- Meat and Bone Meal, Feed

620 Leshar Place • Lansing, MI 48912
800/234-5333 (USA/Canada) or 517/372-9200 • fax: 517/372-2006
e-mail: foodsafety@neogen.com • www.neogen.com

©Neogen Corporation, 2008. Neogen, Veratox and K-Blue are registered trademarks of Neogen Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.
1572B V-AflaHS-ENSP_1208

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