Read instructions carefully before starting test

ALERT for Histamine

REFRIGERATE AT 2–8°C. DO NOT FREEZE.

HISTAMINE

High levels of histamine may develop in a variety of fish species as they decompose. These species include tuna, mahi-mahi, marlin, bluefish, sardines, anchovy, bonito, herring and mackerel. Ingestion of histamine may cause scombroid poisoning in humans, which may lead to a variety of symptoms, including rash, nausea, vomiting, diarrhea, hypotension, palpitations and muscle weakness. Paralysis and death have also been reported in cases of scombroid poisoning.

Because of its potential for human illness, the U.S. Food and Drug Administration (FDA) has ruled that extensive refrigeration records and/or histamine testing must be included in a hazard analysis and critical control point (HACCP) program for relevant fish species. The FDA has set an action level of 50 parts per million (ppm) for histamine in domestic and imported fish.

INTENDED USE

Alert® for Histamine is intended for the qualitative analysis of histamine in scombroid species of fish, such as tuna, bluefish, mahi-mahi and fish meal.

INTENDED USER

The test kit is designed for use by quality control personnel and others familiar with histamine analysis in fish. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed Neogen training.

ASSAY PRINCIPLES

Alert for Histamine is a competitive direct enzyme-linked immunosorbent assay (CD-ELISA). Histamine is extracted from a sample in a quick water extraction process. This extract is filtered and then diluted into a buffer solution supplied with the test. Free histamine, in the buffered sample and control, is allowed to compete with enzyme-labeled histamine (conjugate) for the antibody binding sites. After a wash step, substrate is added, which reacts with the bound enzyme conjugate to produce blue color. The color of the sample is visually compared to the color of the control. If the sample has more blue color than the control, it contains less histamine than the control. If the sample contains less blue color (or more red) than the control, it contains more histamine than the control.

ALERT®

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F).

MATERIALS PROVIDED

1. 24 antibody-coated microwells
2. 24 red-marked mixing wells
3. 1 yellow-labeled bottle of 50 ppm histamine control
4. 1 blue-labeled bottle of histamine-HRP conjugate solution
5. 20 bottles of sample extract diluent buffer
6. 1 bottle of 80 mL wash buffer
7. 1 green-labeled bottle of K-Blue® Substrate solution
8. 1 red-labeled bottle of Red Stop solution
9. 12 disposable beakers

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Alert for Histamine extraction kit (Neogen item 9520) — contains 20 bottles, 20 filter syringes and 20 tubes
2. Filter syringes (Neogen item 9420H), Whatman no. 1 filter paper (Neogen item 9430), centrifuge or equivalent
3. Tubes with caps, 5 mL (Neogen items 9421, 9421B)
4. Graduated cylinder
5. Blender
6. Scale capable of weighing 10–50 grams (Neogen item 9427)
7. Pipettor, 100 µL (Neogen items 9272, 9290)
8. Pipette tips (Neogen items 9410, 9407, 9417)
9. Paper towels or equivalent absorbent material
10. Microwell holder (Neogen item 9402)
11. Timer (Neogen item 9426)
12. Waterproof marker
13. Wash bottle (Neogen item 9400)
14. Distilled or deionized water

PRECAUTIONS

1. Store test kit between 2–8°C (35–46°F) when not in use.
2. Do not use kit components beyond expiration date.
3. Glassware should not be used for extraction purposes. Histamine may adhere to glass, therefore using glassware may affect test results.
4. Do not mix reagents from a kit with one serial number with reagents from a kit with a different serial number.
5. Do not run more than six wells per test.
6. Follow proper pipetting techniques, including priming tips by filling and dispensing solution once before use.
7. Use of incubation times other than those specified may give inaccurate results.
8. Kits should be brought to 18–30°C (68–82°F) prior to use.
9. Avoid prolonged storage of kits at ambient temperatures.
10. Do not freeze test kits.
11. To avoid cross-contamination, use clean pipette tips for each sample.
12. Commodity extracts should have a pH of 6–8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH, contact Neogen.

PROCEDURAL NOTES

1. **Substrate.** K-Blue Substrate is pre-activated and ready for use. To use, fill a clean, disposable beaker provided in the kit about 1/4 full. This will provide an excess of substrate. Keep the substrate protected from light until needed. **Do not return unused substrate to the bottle.** If leftover substrate becomes contaminated, and turns dark blue in the bottle, it should be disposed of.
2. **Control.** The control provided with the kit is ready to use. Do not dilute the control in the same manner as sample extracts. This control is equivalent to a concentration of 50 ppm in a properly extracted fish sample.
3. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after the samples are extracted, and the test procedure is set to begin.

SAMPLE PREPARATION

A. Canned tuna: AOAC 937.07b

Place entire contents of can/pouch, meat and liquid, into a blender. Blend until homogenous. Store samples at 2–8°C (35–46°F) until analyzed.

B. Fresh or thawed frozen raw fish: AOAC 937.07a

Clean and eviscerate three fish. Cut three cross-sectional pieces 2.5 cm (1 inch) thick, from back of the pectoral fin, halfway to vent and one posterior to the vent. Debone slices and blend or grind combined samples until homogenous. Store samples at 2–8°C (35–46°F) until analyzed.

SAMPLE EXTRACTION

NOTE: If using Neogen's Alert for Histamine Extraction Kit, follow the instructions in that kit. If preparing your own extraction solution, continue with the following instructions.

1. Add 10 g of the homogenous mixture to a clean disposable extraction bottle containing 190 mL of distilled or deionized water.
2. Tightly cap and vigorously shake the bottle for **15–20 seconds** to suspend the fish tissue in the water.
3. Wait approximately **5 minutes**, then shake the bottle for **15–20 seconds** to resuspend the fish tissue.
4. Wait an additional **5 minutes**, and again shake the bottle for **15–20 seconds** to re-suspend the fish tissue. Allow the tissue to settle to the bottom of the bottle for about **30 seconds**.
5. If necessary, filter the contents through folded filter paper or a filter syringe into a clean container. The sample is now ready for extract dilution.

ALTERNATE: Centrifuge the sample and use the clear supernatant for testing.

SAMPLE EXTRACT DILUTION

Because of the sensitivity of the ELISA format, the sample extract must be diluted before performing the test. Using a clean pipette tip, add 100 µL of the extract to a bottle of sample diluent buffer provided with the kit. Gently swirl to mix. The sample is now ready to test. Repeat for all samples.

TEST PROCEDURE

Allow all reagents to warm to 18–30°C (68–82°F) prior to use.

1. Remove 1 red-marked mixing well for each sample to be tested plus 1 red-marked well for the control, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells that will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of strip with a “1”, and place strip in the well holder with the marked end on the left. Do not mark on the inside or bottom of the wells.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Pipette 100 µL of conjugate from the blue-labeled bottle into each red-marked mixing well.
5. Using a new pipette tip, add 100 µL of the control from the yellow-labeled bottle to the first red-marked mixing well. Mix the control with the conjugate in the well by pipetting it up and down 5 times. Discard the used tip.
6. Using a new pipette tip, add 100 µL of diluted sample to the next red-marked mixing well. Mix as in the previous step. Record the location of each sample for identification. Repeat the process for each sample, using a new pipette tip for each.
7. Using a new pipette tip for each, transfer 100 µL from each red-marked mixing well to an antibody-coated well. Timing is very important with the transfers. All transfers must be completed in less than a minute, or the accuracy of the subsequent results may be jeopardized.
8. Mix by sliding the microwell holder back and forth on a flat surface for **10–20 seconds** without splashing reagents from the wells. Incubate for **10 minutes** at 18–30°C (68–82°F). Discard red-marked mixing wells.
9. Shake out the contents of the antibody wells. Wash the antibody wells by completely filling the wells with wash buffer contained in the dispensing bottle, and dumping them out. Repeat this wash step 3 more times. Following the last wash, turn the wells upside down and tap out the remaining liquid on an absorbent towel.
10. Add 100 µL of substrate from the green-labeled bottle to each antibody well from the beaker of substrate prepared earlier (see **Procedural Notes**).
11. Mix by sliding back and forth on a flat surface for **10–20 seconds**. Incubate for **10 minutes** at 18–30° C (68–82° F).
12. Add 100 µL of Red Stop from the red-labeled bottle to each well in the same order as the substrate.

INTERPRETATION OF RESULTS

If a sample well is **darker blue** than the control well, the sample contains **less than 50 ppm** of histamine. If a sample well shows **less blue** color or **more red** color than the control, the sample contains **more than 50 ppm** of histamine. For optimum observation of color differences, place the wells on a white surface and read looking down through the solution.

ALTERNATIVE: Microwells may be read in a microwell reader. Wipe bottom of wells, blank reader on air using a 650 nm filter, and compare sample readings to 50 ppm control reading.

RETESTING AND LIMITATIONS

Sample wells with equivalent color to the control contain histamine levels close to the action level set by the FDA of 50 ppm. Repeating the test for these samples may be advisable. Confirmation of the results may also be desired. The results obtained when using a 10 gram sample may not be representative of the entire sample. The results are qualitative and should be used for screening only.

VALIDATED MATRIXES

Fresh, canned and pouched tuna in oil or water, mahi-mahi, bluefish and fish meal.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts

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