Read instructions carefully before starting test

Veratox[®] for Fumonisin 5/10 Quantitative Test

REFRIGERATE AT 2-8°C • DO NOT FREEZE

THE TOXIN

Discovered in 1989, fumonisins are a family of mycotoxins produced by different species of the mold *Fusarium*. These molds commonly infect corn (in fact, they are considered ubiquitous in corn) and rice, hence the potential for fumonisins to be found in feed and foodstuffs is high. Fumonisin affects various animals differently and have been linked to esophageal cancer in humans. The Environmental Protection Agency classifies fumonisins as Category II-B carcinogens.

Horses are extremely sensitive to low amounts of fumonisin, which can cause leukoencephalomalacia (liquefaction of the brain). In swine, research has shown fumonisin attacks the cardiopulmonary system causing pulmonary edema, as well as liver and pancreatic lesions.

The FDA has issued a final guidance for total fumonisins (FB₁+FB₂+FB₃) in food and animal feeds:

HUMAN FOODS (total fumonisins)

| | |
|---|-------|
| Degermed dry milled corn products..... | 2 ppm |
| Whole/partially dry milled corn products, dry milled corn bran, cleaned corn for masa production..... | 4 ppm |
| Cleaned corn for popcorn..... | 3 ppm |

ANIMAL FEEDS (corn/corn by-products, total fumonisins)

| | |
|---|------------------------|
| Equids and rabbits | 5 ppm, ≤ 20% of diet |
| Swine and catfish | 20 ppm, ≤ 50% of diet |
| Breeding ruminants, breeding poultry and breeding mink..... | 30 ppm, ≤ 50% of diet |
| Ruminants ≥ 3 months old being raised for slaughter and mink being raised for pelt production..... | 60 ppm, ≤ 50% of diet |
| Poultry being raised for slaughter | 100 ppm, ≤ 50% of diet |
| All other species or classes of livestock and pet animals..... | 10 ppm, ≤ 50% of diet |

INTENDED USE

The Veratox for Fumonisin 5/10 test kit is a competitive direct enzyme-linked immunosorbent assay (CD-ELISA) in a microwell format for the quantitative analysis of fumonisin in such commodities as corn, cornmeal and rice.

INTENDED USER

The test kit is designed for use by quality control personnel and others familiar with food and feed possibly contaminated by fumonisin. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed the Neogen training.

ASSAY PRINCIPLES

Veratox for Fumonisin 5/10 allows the user to obtain exact concentrations in parts per million (ppm). Free fumonisin in the samples and controls is allowed to compete with enzyme-labeled fumonisin (conjugate) for the antibody binding sites. After a wash step substrate is added, which reacts with the bound conjugate to produce blue color. More blue color means less fumonisin. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form the standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of fumonisin.

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2-8°C (35-46°F).

MATERIALS PROVIDED

1. 48 antibody-coated microwells
2. 48 red-marked mixing wells
3. 5 yellow-labeled bottles of 0, 0.5, 1, 3 and 6 ppm fumonisin controls (see precautions for handling methanol solution)
4. 1 blue-labeled bottle of fumonisin-HRP conjugate solution
5. 1 green-labeled bottle of K-Blue Substrate solution
6. 1 red-labeled bottle of Red Stop solution
7. 28 dilution bottles prefilled with 7.9 mL of deionized water (dilution bottles not provided with product #8836)

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Extraction materials (items d through f available in kit form from Neogen, item 8052):
 - a. 70% methanol solution (Neogen item 8055, 8056) or ACS grade methanol
 - b. Distilled or deionized water
 - c. 250 mL graduated cylinder (Neogen item 9368)
 - d. Container with 125 mL capacity (Neogen item 9428)
 - e. Neogen filter syringes, Whatman #1 filter paper, or equivalent (Neogen item 9420, 9430)
 - f. Sample collection tubes (Neogen item 9421)
2. High-speed blender (Neogen item 9493, 9477)
3. Agri-Grind grinder or equivalent (Neogen item 9401, 9453)
4. Scale capable of weighing 2-25 g \pm 0.1 g (Neogen item 9427)
5. Microwell reader with a 650 nm filter (Neogen item 9303)
6. Pipettor, 12-channel (Neogen item 9273)
7. Pipettor, 100 μ L (Neogen item 9272, 9290)
8. Tips for 100 μ L and 12-channel pipettors (Neogen item 9410, 9407, 9417)
9. Paper towels or equivalent absorbent material
10. Plastic 1/2 gallon bucket for use as waste receptacle
11. Microwell holder (Neogen item 9402)
12. Timer (Neogen item 9426)
13. Waterproof marker
14. Wash bottle (Neogen item 9400)
15. 2 reagent boats for 12-channel pipettor (Neogen item 9450)
16. Distilled or deionized water

PRECAUTIONS

1. Methanol solution is highly flammable. Keep container tightly closed, and keep away from heat, sparks, open flame and those smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
2. Store test kit between 2-8°C (35-46°F) when not in use. Do not freeze.
3. Kits should be brought to room temperature 18-30°C (64-86°F) prior to use.
4. Avoid prolonged storage of kits at ambient temperatures.
5. Do not use kit components beyond expiration date.
6. Do not mix reagents from one kit serial with reagents from a different kit serial.
7. Do not run more than 24 wells per test.
8. Follow proper pipetting techniques, including priming of tips.
9. Use of incubation times other than those specified may give inaccurate results.
10. Treat all used liquids, including sample extract, and labware as if contaminated with fumonisin. Gloves and other protective apparel should be worn at all times.
11. To avoid cross-contamination, use clean pipette tips and glassware for each sample, and thoroughly wash all glassware between samples.
12. Commodities tested should have a pH of 6-8. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH contact your Neogen representative or Technical Services.

PROCEDURAL NOTES

1. **Substrate:** K-Blue Substrate is ready for use. The substrate should be clear — discard if it has turned blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until it is needed.
2. **Antibody wells:** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test is set to begin.

SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques. The sample should be ground and thoroughly mixed prior to proceeding with the extraction. Store samples at 2-8°C (35-46°F) until analyzed.

1. If not using Neogen's prepared solution, prepare a 70% methanol solution by mixing 7 parts ACS Grade methanol with 3 parts distilled or deionized water for each sample to be tested.
2. Obtain a representative sample. Grind the entire sample so that at least 75% of the ground material passes through a 20 mesh sieve, the particle size of a fine instant coffee.
3. Extract at a ratio of 1 part sample to 5 parts 70% methanol. Example: Blend 25 grams of ground sample with 125 mL of 70% methanol for 2 minutes in a high-speed blender.

Alternative method: Add 5 grams of ground sample to 25 mL of 70% methanol and shake vigorously for 3 minutes.

4. Filter the extract by pouring at least 5 mL through a Whatman #1 filter (or Neogen filter syringe) and collecting the filtrate as a sample.
5. Dilute sample by adding 100 µL of extract to a prefilled sample dilution bottle and mix well by swirling bottle. The sample is now ready for testing without further preparation. Repeat for each sample, making sure to label each bottle.

TEST PROCEDURE

Allow all reagents to warm to room temperature 18-30°C (64-86°F) before use.

1. Remove 1 red-marked mixing well for each sample to be tested plus 5 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant. Reseal the pack to protect the antibody. Mark one end of strip with a "1", and place strip in the well holder with the marked end on the left. Do not mark the inside or bottom of the wells.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Place 100 µL of conjugate from the blue-labeled bottle in each red-marked mixing well.
5. Using a new pipette tip for each, transfer 100 µL of controls and samples to the red-marked mixing wells as described below.

| | | | | | | | | | | | | |
|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|---------|
| 0 | 0.5 | 1 | 3 | 6 | S1 | S2 | S3 | S4 | S5 | S6 | S7 | Strip 1 |
| S8 | S9 | S10 | S11 | S12 | S13 | S14 | S15 | S16 | S17 | S18 | S19 | Strip 2 |

6. Using a 12-channel pipettor, mix the liquid in the wells by pipetting it up and down 3 times. Transfer 100 µL to the antibody-coated wells. Discard red-marked mixing wells.
7. Set timer for **5 minutes, mixing the wells for the first minute** of the room temperature incubation by sliding the microwell holder back and forth on a flat surface without splashing reagents from the wells.
8. The initial reaction is now complete. Shake out the contents of the antibody wells. Fill the wells with distilled or deionized water and dump them out. Repeat this step 5 times, then turn the wells upside-down and tap out on a paper towel until the remaining water has been removed.
9. Pour the needed volume of substrate from the green-labeled bottle into the green-labeled reagent boat.
10. With new tips on the 12-channel pipettor, prime and pipette 100 µL of substrate into the wells.
11. Set timer for **10 minutes, mixing the wells for the first minute** by sliding the microwell holder back and forth on a flat surface. Discard remaining substrate and rinse the reagent boat with water.
12. Pour Red Stop solution from the red-labeled bottle into the red-labeled reagent boat.
13. Eject the excess substrate from the 12-channel pipettor, prime the tips, and pipette 100 µL of Red Stop to each well. Mix by sliding back and forth on a flat surface. Discard the tips.
14. Wipe the bottom of the microwells with a dry cloth or towel and read in a microwell reader using a 650 nm filter. Air bubbles should be eliminated, as they could affect analytical results. Results should be read within 20 minutes after the addition of Red Stop.
15. Read and calculate results using Neogen's Stat Fax microwell reader, or equivalent. If using a strip/plate reader, calculate results using Neogen's Veratox for Windows software.

RETESTING

If positives occur in commodities not previously tested, confirm with an additional approved method prior to taking action.

PERFORMANCE CHARACTERISTICS

Limit of detection: 0.2 ppm (Determined by the mean average of 10 fumonisin free samples plus 2 standard deviations.)

Limit of quantitation: 0.5 ppm (Described as the lowest concentration point on the calibration curve that this test can reliably detect fumonisin.)

Range of quantitation: 0.5 - 6 ppm (For quantitating samples above 6 ppm contact Neogen Technical Services for dilution instructions.)

Validated matrixes: Barley, beet pulp*, corn, corn meal, corn germ meal*, corn gluten meal*, corn/soy blend, corn steep, DDGS*, DDGs wet cake*, millet, milo (grain sorghum), oats, oat hulls, oats (naked), pea fiber, pet food*, popcorn, potato (white), rice gluten, rice hulls, rough rice, rye, soy hydrolysate, soybeans, soybean meal, sunflower meal, wheat, and wheat bran.
*May require a pH adjustment.

NOTE: Neogen continues to validate new commodities. Please contact a representative for the latest commodity list.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Service can be reached between 8 a.m. and 6 p.m. Eastern time by calling 800/234-5333 or 517/372-9200 and asking for a Neogen sales representative or Technical Services. Assistance is available on a 24-hour basis by calling 800/867-0308. Training on this product, and all Neogen test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's Food Safety test kits, on Neogen's Web site at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TESTING KITS AVAILABLE FROM NEOGEN

Natural Toxins

- Aflatoxin, DON, Ochratoxin, Zearalenone, T-2/HT-2 Toxin, Fumonisin, Histamine

Foodborne Bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, Yeast and Mold, Total Plate Count, Generic *E. coli* and Total Coliforms, Protein Residues

Food Allergens

- Peanuts, Milk, Eggs, Almonds, Gluten, Soy, Hazelnut

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India

Neogen Food and Animal Security
+91 484 2306598, 2301582
info@neogenindia.com
www.neogenindia.com

©Neogen Corporation, 2015. Neogen, Veratox and K-Blue are registered trademarks of Neogen Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.

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