*Read instructions carefully before starting test*

# **ALERT<sup>®</sup>**

## **for Peanut Allergen**

### **Screening Test**

*REFRIGERATE AT 2–8°C (35–46°F) • DO NOT FREEZE*

#### **PEANUT ALLERGEN**

Food allergens are proteins in food that can create an immune response in sensitive individuals. Once ingested, food allergens can cause a number of reactions, ranging in severity from hives and itching to anaphylaxis. Anaphylaxis is a severe allergic reaction, involving vomiting, diarrhea, difficulty breathing, swelling of the mouth and tongue, and a rapid drop in blood pressure.

An estimated 3.5 to 4 percent of adults, and 6 to 8 percent of children, are sensitive in some degree to food allergens. More than 12 million people in the United States alone are known to have a food allergy, with an allergy to peanuts being one of the most prevalent.

Food manufacturers protect those with food allergies by clearly labeling their products with a list of ingredients. Testing for the presence of peanut protein residue components ensures food manufacturers that an unlabeled — and potentially dangerous — ingredient did not make its way into a food product.

#### **INTENDED USE**

Alert for Peanut Allergen is intended for the qualitative analysis of peanut protein residue in food products such as cookies, crackers, chocolate bars, ice cream, cereals, in clean-in-place rinses and on environmental surfaces.

#### **INTENDED USER**

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by peanuts or peanut products. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed the Neogen training.

## **ASSAY PRINCIPLES**

Alert for Peanut Allergen is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Peanut protein residue is extracted from samples by blending with a hot buffered salt solution (PBS). Extracted peanut residue is sampled and added to antibody-coated wells (capture antibody) where it binds to the antibody during an incubation. Any unbound peanut residue is washed away and a second antibody (detector antibody), which is enzyme labeled, is added. The detector antibody binds to the already bound peanut residue. After a second wash, the substrate is added. Color develops as a result of the presence of bound detector antibody. Red Stop reagent is added and the color of the resulting solution is observed.

## **STORAGE REQUIREMENTS**

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). Do not freeze.

## **MATERIALS PROVIDED**

1. 24 antibody-coated microwells
2. 1 yellow-labeled dropper bottle of 5 ppm peanut control
3. 1 blue-labeled dropper bottle of enzyme-labeled antibody conjugate
4. 1 green-labeled dropper bottle of K-Blue® Substrate solution
5. 1 red-labeled dropper bottle of Red Stop solution
6. 20 unlabeled empty sample dropper bottles with tips
7. 3 foil pouches of 10 mM PBS dry powder extraction solvent; each pouch contains enough powder to prepare 1 L of extraction solution
8. 1 wide-mouth bottle of 10 mM PBS-Tween wash buffer concentrate; a bottle contains enough concentrate to prepare 1 L of wash buffer
9. 2- 50 g cups of extraction additive
10. Plastic scoop to measure extraction additive

## **MATERIALS RECOMMENDED BUT NOT PROVIDED**

1. Allergen Extraction Kit (Neogen item #8429)
  - a. 20 disposable plastic extraction bottles
  - b. 20 sample collection tubes (12 x 75 mm) with caps
2. Allergen Environmental Swabbing Kit (Neogen item #8432S)
  - a. 100 sterile swabs
  - b. 100 dropper tips
3. High-speed blender with a 250 mL jar or 250 mL disposable plastic bottles (Neogen item #9493, #9477)
4. Scale capable of weighing  $5 \pm 0.1$  g (Neogen item #9427)
5. Shaker water bath, hot plate or equivalent heat source capable of maintaining  $60^{\circ}\text{C} \pm 1^{\circ}$
6. Thermometer
7. Whatman #4 filter paper or equivalent (Neogen item #9429)
8. Two 1 L bottles to prepare washing and extraction solutions (Neogen item #9472)
9. Paper towels or equivalent absorbent material
10. Microwell holder (Neogen item #9402)
11. Waterproof marker
12. Wash bottle (Neogen item #9400)
13. Distilled or deionized water
14. Timer (Neogen item #9426)
15. Graduated cylinder capable of measuring 125 mL (Neogen item #9368)

## PRECAUTIONS

1. Samples intended to be tested for peanut must be extracted separately from samples intended to be tested for other food allergens, such as egg and milk residues. The extraction additives for each type of test are designed specifically for the target food allergen.
2. Components of the Alert for Peanut Allergen test kit, such as controls and extraction reagents, may contain one or more of the following potentially allergic materials: milk, egg protein, peanut protein or soy protein. If allergic to any of these compounds, please use caution when using this product.
3. Concentrated food additives, colors and flavors may cause interferences on ELISA test methods. Contact Neogen's technical services for validation information.
4. Hydrolyzed and fermented proteins may not be detected using ELISA methods for allergen testing. Due to the nature of the proteins they may be undetectable in the assay, but there could still be active allergenic protein residue present.
5. The testing environment should be clean and dust-free.
6. Store test kit refrigerated, between 2–8°C (35–46°F), when not in use. Do not freeze test kits.
7. Bring kits to room temperature (18–30°C, 64–86°F) prior to use.
8. Avoid prolonged storage of kits at ambient temperatures.
9. Do not use kit components beyond expiration date.
10. Do not mix reagents from kit serial with those from a different serial.
11. Do not run more than 6 wells per test.
12. Use only incubation times specified. Others may give inaccurate results.

## PROCEDURAL NOTES

1. **Extraction solution.** Prepare extraction solution by adding a foil pouch of extraction solvent, 10 mM PBS, to 1 L distilled or deionized water. Swirl to mix thoroughly. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).
2. **Wash buffer.** Prepare the wash buffer solution by pouring all the wash buffer concentrate into an empty 1 L bottle. Rinse the wash buffer concentrate container with distilled or deionized water and pour into the 1 L bottle to ensure all the concentrate is used. Fill the 1 L bottle with additional distilled or deionized water, and swirl to assure thorough mixing. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).
3. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted and the test procedure is set to begin.

## SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques (see Neogen's Food Allergen Monitoring Handbook). The sample should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

1. Prepare the extraction solution as described in the procedural notes.
2. Preheat the extraction solution to 60°C (140°F).
3. Obtain a representative sample. If the sample is of a larger particle size, grind it to a very fine particle size.
4. Transfer 5 g of sample, or 5 mL of liquid sample, to a 250 mL disposable plastic bottle or a blender jar.
5. Add one level scoop of peanut extraction additive to the bottle/jar. Do not use the extraction additive from another allergen test kit.
6. Pour 125 mL of the 60°C (140°F) extraction solution into the bottle/jar and cap.
7. Extract by shaking (150 rpm) in a water bath at 60°C (140°F) for 15 minutes. Remove bottle from bath.

**Alternative:** Blend at high speed for 2 minutes.

8. Let the material settle for 5 minutes.  
**NOTE:** In rare cases, filtering may be necessary to achieve a supernatant free of suspended material. Filter the extract by pouring at least 5 mL through a Whatman #4 filter and collecting the filtrate as a sample, or centrifuging at 14,000 rpm for 5 minutes.
9. Using a disposable pipette, transfer approximately 1 mL of the supernatant (the top liquid portion of the extract) to a sample dropper bottle. Then, label and place a dropper tip on the bottle.
10. Begin the test procedure once the sample has cooled to room temperature (at least 15 minutes).

### ENVIRONMENTAL SWABBING PREPARATION AND EXTRACTION

1. Prepare the extraction solution as described in the procedural notes.
2. Prepare a working solution by combining 125 mL of the extraction solution to a separate container, and add one level scoop of extraction additive to a container. Prepare a fresh working solution daily. **NOTE:** Prepare 125 mL of working solution for every 25 environmental swabs to be tested.
3. Preheat the working solution to 60°C (140°F). Shake container to dissolve additive.
4. Gather the sample with a swab, using one of the following methods:  
**For dry surfaces:** Open a new swab and wet with extraction solution. Swab a 10 x 10 cm area by using a crosshatch technique. (Do not use working solution to moisten swabs.)  
**For wet surfaces:** Open a new swab and swab a 10 x 10 cm area by using a crosshatch technique. Do not moisten swab prior to use.
5. Return the swab to its original tube once sampling is complete. Remember to label each tube.
6. Remove the swab from its tube, and add 5 mL of working solution at 60°C to the tube. Mix by placing the swab back into the tube and shaking for 2 minutes by hand (inverting tube), or for 30 seconds with a vortex mixer.
7. Remove the swab from its tube.
8. Place a new sample dropper tip onto the tube. The solution in the tube now serves as the sample. **NOTE:** Use caution when inverting the tube, as some liquid may drip.

### TEST PROCEDURE

Allow the test kit and all reagents to warm to room temperature (18–30°C, 64–86°F) before using.

1. Remove 1 well for each sample to be tested plus 1 well for the control, and place into the well holder.
2. Mix each reagent by swirling its dropper bottle prior to use.
3. Add 3 drops from the yellow-labeled control dropper bottle to the first well. Add three drops from each sample extract to a respective well as indicated in the template below. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.

Control	S1	S2	S3	S4	S5
---------	----	----	----	----	----

4. Incubate microwells **10 minutes** at room temperature (18–30°C, 64–86°F).
5. Shake out the contents of the wells. Using a wash bottle filled with wash buffer, fill each well and shake out. Repeat 5 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
6. Add 3 drops from the blue-labeled conjugate dropper bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.
7. Incubate for **10 minutes** at room temperature (18–30°C, 64–86°F).
8. Shake out the contents of the wells. Using a wash bottle filled with wash buffer, fill each well and shake out. Repeat 5 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
9. Add 3 drops from the green-labeled substrate dropper bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.

10. Incubate for **10 minutes** at room temperature (18–30°C, 64–86°F).
11. Add 3 drops from the red-labeled Red Stop dropper bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface. The results are now ready to be interpreted.

### **INTERPRETATION OF RESULTS**

Visually compare the color of a sample well to the color of the control well. If the sample well has **more blue color** than the control well, the sample tests positive for peanut contamination of more than 5 ppm. If the sample well has **less blue color, or more red color**, than the control well, the sample contains **less than 5 ppm** of peanut contamination.

**Alternative:** Read wells (wipe bottom of wells first) in a microwell reader with a 650 nm filter. If the sample well has an optical density (OD) higher than the control well, the sample is positive for peanut contamination of more than 5 ppm. If the sample well has an OD lower than the control well, the sample contains less than 5 ppm of peanut contamination.

### **PERFORMANCE CHARACTERISTICS**

Allergen detection: This test detects peanut protein, and the results are expressed in terms of the detection of total peanut.

### **CUSTOMER SERVICE**

Neogen Customer Assistance and Technical Service can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

### **MSDS INFORMATION AVAILABLE**

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's Food Safety test kits, at [www.neogen.com](http://www.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

### **WARRANTY**

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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