Read instructions carefully before starting test

Reveal[®]

for

Ruminant in MBM

RUMINANT BY-PRODUCTS IN MEAL

The use of rendered animal products in animal feed is a common practice to supplement the feed's protein content. Concern about the consumption of ruminants (cattle, sheep) infected with bovine spongiform encephalopathy (BSE, a.k.a. "mad cow" disease) led the U.S. Food and Drug Administration to ban the feeding of rendered ruminant by-products, such as meat and bone meal (MBM), to other ruminants.

INTENDED USE

Reveal for Ruminant in MBM is an immunochromatographic lateral flow assay used for the qualitative analysis of beef and sheep MBM in nonruminant animal protein meal. The test can detect as little as 2% beef or sheep material in MBM. The test also will react with by-product material from goats. It will not detect material from other ruminants, such as deer or elk.

INTENDED USER

The test kit is designed for use by quality control personnel and others familiar with the possibility that meal containing ruminant by-products could enter the facility. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed Neogen training.

ASSAY PRINCIPLES

Reveal for Ruminant in MBM is a single-step lateral flow immunochromatographic assay. The extract is wicked through a reagent zone, which contains antibodies specific for heat stable beef and sheep muscle protein conjugated to colored particles. If beef or sheep by-product is present, it will be captured by the conjugated antibodies. The beef/sheep-antibody-particle complex then is wicked onto a membrane which contains a zone of antibody specific for beef and sheep muscle protein. This zone captures the complex allowing the particles to concentrate and form a visible line. If no beef or sheep MBM is present, no line will form. The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line always will form regardless of the presence of beef or sheep MBM, ensuring the strip is working properly.

STORAGE REQUIREMENTS

Store kit components at room temperature (18–30°C, 64–86°F) to ensure full shelf life.

MATERIALS PROVIDED

1. 25 test strips
2. 25 transfer pipettes
3. 25 sample tubes
4. 1 foil pouch of extraction solvent (each pouch is enough to prepare 1 gallon (3.8 L) of extraction solution in distilled or deionized water)
5. Peel-off label for labeling container of extraction solvent

MATERIALS REQUIRED BUT NOT PROVIDED

1. 12-ounce Mason-type jars (1 per sample)
2. Heat source to boil water
3. 1 gallon (3.8 L) container of distilled or deionized water
4. Scale capable of weighing 10 g (Neogen item 9427)
5. Sample tube rack (Neogen item 9466)
6. 100 mL graduated cylinder
7. Timer (Neogen item 9426)
8. Waterproof marker

PRECAUTIONS

1. Store test kit at room temperature (18–30°C, 64–86°F) when not in use.
2. Do not use kit components beyond expiration date.
3. Do not freeze test kits.
4. Treat all used liquids, including sample extract, and labware as if containing ruminant by-products.
5. To avoid cross-contamination, use clean labware for each sample, and thoroughly wash all labware between samples.

MEAL SAMPLE PREPARATION AND EXTRACTION

The samples to be tested should be collected according to accepted sampling techniques. Samples to be tested can be subsampled from representative samples collected from bulk animal protein meals for other tests. It is assumed that the samples collected are representative of the contents of the containers and are sufficiently homogenous with regard to sample contents. Extracted samples should be tested immediately.

1. Discard approximately 7 ounces (220 mL) from the 1 gallon (3.8 L) container of distilled or deionized water. Prepare extraction solution by adding the contents of the foil pouch to the container of water, and inverting the container approximately 20 times to completely dissolve.
2. Using the peel-off label supplied, label the extraction solution container indicating the solution's expiry date (6 months from date of preparation). Store any unused extraction solution at room temperature (18–30°C, 64–86°F). **NOTE:** Do not refrigerate.
3. Weigh 10 ± 0.1 g of sample and place in a Mason-type jar.
4. Add 100 ± 2 mL of extraction solution to the jar.
5. Heat in boiling water for **10 minutes**.
6. Remove the jar from the heat source. Swirl to mix. Do not allow the sample to cool prior to testing.

TEST PROCEDURE

1. Remove the appropriate number of sample tubes, label and place in tube rack.
2. Using a new transfer pipette for each sample, transfer approximately 0.5 mL of heated extract to a sample tube.
3. Place a strip with the sample end down into each sample tube.
4. Allow each strip to develop for **10 minutes** and read results.

INTERPRETATION OF RESULTS

Positive results: If a line forms in the test zone and another line forms in the control zone **within 10 minutes** (2 lines total), the sample is positive. A positive sample may be determined as soon as 2 lines are visible on the strip. If there is no line in the control zone, the test is invalid and the sample should be retested with another strip.

Negative results: If after the **full 10 minutes** there is no visible line in the test zone, but a visible line in the control zone, the sample is negative. If there is no line in the control zone, the test is invalid and the sample should be retested with another strip.

NOTE: Observations after 15 minutes may be inaccurate due to overdevelopment of the device. To avoid overdevelopment after 10 minutes, and to retain the test results, snip off the solid red bottom portion of the strip to stop the wicking of sample extract through the test and control zones.

PERFORMANCE CHARACTERISTICS

Reveal for Ruminant in MBM is designed to detect greater than or equal to 2% beef or sheep material in MBM. The intensity of the line in the test zone may vary by the amount of beef or sheep material detected.

CUSTOMER SERVICE


Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's Food Safety test kits, at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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
TESTING KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustacea, eggs, gliadin, hazelnut, lupine, milk, mustard, peanut, sesame, soy, walnut

Genetic modification

- CP4 (Roundup Ready®)

Ruminant by-products

- Meat and bone meal, feed

North America Neogen Headquarters

620 Leshler Place, Lansing, MI 48912 USA
800/234-5333 (USA/Canada) or 517/372-9200
Fax: 517/372-2006 • foodsafety@neogen.com
www.neogen.com

Europe, Middle East and Africa Neogen Europe

The Dairy School, Auchincruive, Ayr
KA6 5HW Scotland, UK
+ 44 (0) 1292 525 600 • Fax: + 44 (0) 1292 525 601
info_uk@neogeneurope.com
www.neogeneurope.com

Mexico Neogen Latinoamérica

Darwin No. 83, Col. Anzures, México, 11590 D.F.
(55) 5254-8235, (55) 5203-0111, (55) 5531-2837
Fax: (55) 5531-1647 • informacion@neogenlac.com
www.neogen.com

Brazil Neogen do Brasil

Rua: Alberto Guizo 760, Distrito Industrial João
Narezzi, Indaiatuba – SP Brasil, Cep: 13.347-402
Tel: +55 19 3935.3727
info@neogendobrasil.com.br • www.neogen.com