**Read instructions carefully before starting test**

# ALERT®

## Sulfites Detection Kit

### SULFITES

Sulfiting agents, sodium bisulfite and sodium metabisulfite, are used with shrimp, lobster and related crustaceans to prevent melanosis (a.k.a. "blackspot"). While sulfiting agents are very beneficial to the seafood industry, some consumers are highly allergic to sulfite residues in food. For this reason, FDA regulations state that seafood containing sulfite residues of 10 parts per million (ppm) or more must be labeled as containing sulfites. Monitoring and verification of sulfite residues are required as part of a Hazard Analysis and Critical Control Points (HACCP) program.

### INTENDED USE

This Alert for Sulfites Kit is a very simple and rapid semi-quantitative test for detecting sulfite residues in shrimp. This kit contains the materials necessary to perform 200 tests.

### TEST PRINCIPLE

Sulfite residues in a sample will produce color changes in a dye reagent provided with the kit. These color changes, and the rate of the changes, are indicative of the level of sulfite residues in the sample.

### MATERIALS PROVIDED

- 1 orange-capped bottle of Activator Solution
- 1 blue-capped bottle of Dye Reagent
- Color Chart

### STORAGE

This kit may be stored at room temperature (18–30°C, 64–86°F). Avoid high temperatures and direct sunlight. Do not allow kit to freeze.

### PRECAUTIONS

1. The Dye Reagent contains methanol. Keep bottle tightly closed. Methanol is highly flammable—keep away from heat, sparks and open flame. Methanol is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin and eyes. Dye Reagent will stain clothing. KEEP OUT OF THE REACH OF CHILDREN.
2. Discard sample after testing.

### PROCEDURAL NOTES

1. Activator Solution and Dye Reagent bottles are shipped in a foil pack.
2. When adding reagents, allow each drop to fall freely onto the sample, being careful not to touch the sample with the dropper tip. Hold the bottle perpendicular to the sample to ensure the correct size drop is dispensed.
3. Color changes are more noticeable if the reagents are added to a less-pigmented area of the shrimp. You may wish to test the whiter thorax area where the head was removed.
4. Pigmentation of different shrimp species may affect how reagent color compares to the color chart. Observe any color changes carefully.

### TEST PROCEDURE

1. Thaw frozen shrimp at room temperature. Remove the head of the shrimp and peel completely. For best results, test the less-pigmented muscle area.
2. Gently squeeze the orange-capped bottle to add 1 drop of Activator Solution to shrimp's exposed muscle.
3. Immediately add 1 drop of Dye Reagent from the blue-capped bottle to the same area as the Activator Solution. **Wait 1 minute.**
4. After 1 minute, note the color of the area where the dye was added. Use the following guidelines to determine the sample's sulfite level. The enclosed color chart can help determine color differences.
  - a. If the blue dye has not changed color, the shrimp was not treated with sulfiting agents.
  - b. If the blue dye has turned a violet color, the shrimp was treated with sulfiting agents. Sulfite levels are greater than 10 ppm, but less than 100 ppm.
  - c. If the blue dye turns clear, sulfite levels exceed 100 ppm.
5. Color changes after 1 minute should be disregarded.

### CUSTOMER SERVICE

Neogen Customer Assistance and Technical Service can be reached between 8 a.m. and 6 p.m. Eastern time by calling 800/234-5333 or 517/372-9200 and asking for a Neogen sales representative or Technical Services. Assistance is available on a 24-hour basis by calling 800/234-5333. Training on this product, and all Neogen test kits, is available.

### MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's Food Safety test kits, on Neogen's Web site at [www.neogen.com](http://www.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

### WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### TESTING KITS AVAILABLE FROM NEOGEN

#### Natural Toxins

- Aflatoxin, DON, Ochratoxin, Zearalenone, T-2/HT-2 Toxins, Fumonisin, Histamine

#### Foodborne Bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

#### Sanitation

- ATP, Yeast and Mold, Total Plate Count, Generic *E. coli* and Total Coliforms, Protein Residues

#### Food Allergens

- Almonds, Crustaceans, Eggs, Gliadin, Hazelnut, Lupine, Milk, Mustard, Peanuts, Sesame, Soy, Walnut

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