Read instructions carefully before using

Allergen Extraction Kit

*Designed for use with Neogen's Veratox[®], Alert[®] and Reveal[®]
allergen test kits*

INTENDED USE

This Allergen Extraction Kit contains the sample extraction materials to extract any unintended peanut, almond, crustacea, hazelnut, egg, lupine, milk, mustard, soy or gliadin residues in food.

STORAGE

This kit may be stored at room temperature (18–30°C, 64–86°F).

MATERIALS PROVIDED

1. 20 disposable extraction bottles
2. 20 transfer pipettes
3. Directions for use

MATERIALS REQUIRED BUT NOT PROVIDED

1. Shaker water bath adjusted to 60°C (140°F) with clamps to hold extraction bottles (room temperature for gliadin extraction and 30°C for crustacea extraction)
2. Graduated cylinder capable of measuring 125 mL
3. Scale capable of weighing 5 g (Neogen item 9427)
4. 1 L bottle (Neogen item 9472)
5. Timer (Neogen item 9426)
6. Distilled or deionized water
7. Laboratory grade ethanol, 190 proof (for gliadin extraction only)

PRECAUTIONS

1. The extraction additive provided with each allergen test kit is designed specifically for that target allergen. Separate extractions, using the appropriate additive, must be performed for each allergen tested for.
2. The extraction reagents of the allergen test kits may contain one or more of the following potentially allergic materials: milk, egg, peanut, soy, tree nuts or wheat. If allergic to any of these compounds, please exercise caution when using this product.
3. Good laboratory practices always should be followed.
4. Always use a new extraction bottle and sample collection tube for each sample to avoid cross-contamination.

SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques (see Neogen's Food Allergen Handbook). The sample should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

Veratox and Alert procedure for peanut, almond, crustacea, hazelnut, lupine, egg, milk, mustard or soy and Reveal procedure for peanut

1. Prepare the extraction solution (10 mM Phosphate Buffered Saline (PBS), pH 7.4) as described in the test kit insert's procedural notes section.
2. Preheat extraction solution to 60°C (140°F) by immersing the bottle containing the solution into the water bath and allowing it to reach 60°C (30°C for crustacea sample extraction).
3. Using an appropriate sampling and collection procedure, obtain a representative sample and grind it to a very fine particle size.
4. Transfer 5 g of sample, or 5 mL of liquid sample, into a 250 mL disposable plastic extraction bottle.
5. Add 1 level scoop of the appropriate extraction additive into the sample bottle.
6. Pour 125 mL of the 60°C (140°F) (or 30°C for crustacea) extraction solution to the sample bottle.
7. Cap the sample bottle to prevent contents from splashing during the extraction.
8. Extract by shaking (150 rpm) in a 60°C water bath for **15 minutes** (For crustacea, extract by shaking in a 30°C water bath for **30 minutes**). Remove the bottle from the bath.
9. Let the material settle for **5 minutes** before proceeding. Some samples may require a filtration to achieve a supernatant (the top liquid portion of the extract) free of suspended material. Filter the extract by pouring at least 5 mL through a Whatman no. 4 filter or filter syringe.
ALTERNATIVE: Centrifuge at 14,000 rpm for **5 minutes**. Use the clear supernatant as the sample. Do not filter samples for total milk analysis.
10. Remove approximately 1 mL of the supernatant and transfer to a sample collection tube. Begin the test procedure once the sample in the collection tube has cooled to room temperature (at least **15 minutes**).

Reveal procedure for total milk

1. Prepare the extraction solution (10 mM PBS, pH 7.4) as described in the test kit insert's procedural notes section.
2. Preheat extraction solution to 60°C (140°F) by immersing the bottle containing the solution into the water bath and allowing it to reach 60°C.
3. Using an appropriate sampling and collection procedure, obtain a representative sample and grind it to a very fine particle size.
4. Transfer 10 mL of liquid sample into a 250 mL disposable plastic extraction bottle.
5. Add 1 level scoop of the appropriate extraction additive into the sample bottle.
6. Pour 90 mL of the 60°C (140°F) extraction solution to the sample bottle.
7. Cap the sample bottle to prevent contents from splashing during the extraction.
8. Extract by shaking (150 rpm) in a 60°C water bath for **15 minutes**. Remove the bottle from the bath.

- Let the material settle for **5 minutes** before proceeding to the next step.
ALTERNATIVE: Centrifuge at 14,000 rpm for **5 minutes**. Use the clear supernatant as the sample. Do not filter samples for total milk analysis.
- Remove approximately 1 mL of the supernatant (the top liquid portion of the extract) and transfer to a sample collection tube. Begin the test procedure once the sample in the collection tube has cooled to room temperature (at least **15 minutes**).

Veratox and Alert procedure for gliadin or gliadin R5 (nonheat processed samples)

NOTE: Contact Neogen Technical Services for the gliadin cocktail procedure for heat processed samples.

For gliadin kits (Neogen items 8480 and 8481): These kits require the use of a 40% ethanol solution. To prepare, combine 4 parts of 190 proof ethanol with 6 parts distilled water.

For gliadin R5 kits (Neogen items 8510 and 8511): These kits require the use of a 60% ethanol solution. To prepare, combine 6 parts of 190 proof ethanol with 4 parts distilled water.

- Prepare appropriate ethanol solution (see above). Prepare sample extract dilution solution (10 mM PBS, pH 7.4) as detailed in the test kit insert's procedural section.
- Add 2 g ground sample, or 2 mL liquid sample, to a 250 mL clean extraction bottle.
- Add 1 scoop of extraction additive to the bottle.
- Add 20 mL (18 mL for liquid samples) of appropriate ethanol solution, cap the bottle tightly, then shake vigorously by hand for about **20 seconds** to ensure complete mixing.
- Extract by shaking (150 rpm) in a shaker for **15 minutes** at room temperature (a shaker water bath can work, but do not turn the heat on). Remove the bottle from shaker or bath.
- Let the bottle stand for about **10 minutes** to allow some of the sample to settle before withdrawing the clear extract.
- Remove the top from a new prefilled dilution dropper bottle and add 100 μ L of the upper layer of the extract to the dropper bottle. For Veratox for Gliadin, dilute 1:40 in PBS. For Veratox for Gliadin R5, dilute 1:50 in PBS.
- To mix, vortex the tube for **5 seconds**, or invert several times by hand.
- Test diluted samples within **2–3 hours** of extraction.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TESTING KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts

Genetic modification

- CP4 (Roundup Ready®)

Ruminant by-products

- Meat and bone meal, feed

Species identification

- Raw and cooked meat samples

North America

Neogen Headquarters

620 Leshar Place, Lansing, MI 48912 USA
800/234-5333 (USA/Canada) or 517/372-9200
Fax: 517/372-2006 • foodsafety@neogen.com
www.neogen.com

Europe, Middle East and Africa Neogen Europe

The Dairy School, Auchincruive, Ayr
KA6 5HW Scotland, UK
+ 44 (0) 1292 525 600
Fax: + 44 (0) 1292 525 601
info_uk@neogeneurope.com
www.neogeneurope.com

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