Read instructions carefully before starting test

Agri-Screen[®]

Ticket

Pesticide Detection Program

PESTICIDES

Without the use of modern pesticides, the high productivity demanded of modern agriculture to feed the increasing world population probably could not be achieved. The diverse nature of today's pesticides reflect the diverse nature of today's "pests." Pests include mammals and birds which eat crops; insects which either directly damage crops or cause crop-damaging diseases; fungi-causing diseases; human or animal parasites; and weeds which compete for the essential nutrients of desirable crops.

Pesticide manufacturers face the challenge of producing a chemical that is strong enough to protect a crop from pests, yet leaves the crop safe for human or animal consumption. Farmers face the challenge of applying the pesticide in precise quantities in ever-changing environmental conditions. The effects of pesticide residue on humans and animals vary as widely as the types of pesticides. Pesticide toxicity also varies considerably, depending on the route of absorption.

INTENDED USE

Agri-Screen Ticket is used to detect all major organophosphates, thiophosphates and carbamates in air, water, soil, produce, food, spills, surfaces and custom applications. Samples can be tested wet, dry, or with solvents.

ASSAY PRINCIPLES

Ticket is based on a simple biochemical principal. One side of Ticket has a disc saturated with cholinesterase, an enzyme in most living organisms, except plants, which controls muscle performance. Pesticides can inhibit an organism's ability to produce cholinesterase, and therefore, kill the organism. If enough pesticide is present in the sample, it will inactivate the cholinesterase and turn the disc white. A blue color indicates a negative result for the presence of pesticides.

STORAGE REQUIREMENTS

Ticket pesticide detectors and kit components can be stored at room temperature, 18–30°C (64–86°F).

MATERIALS PROVIDED

1. Pesticide detector Tickets
2. Activator ampoules*
3. 50 mL graduated beaker
4. Glass rod for breaking the activator ampoules and stirring the sample

*See precaution #5, Part #8920A and 8921A do not include activator ampoules.

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Water
2. Filter paper
3. Wide mouth jars
4. Leaf punch

PRECAUTIONS

1. Samples which are highly colored may obscure test results — use the combined sample method technique to produce a clear test sample.
2. Samples with pH above 8 or below 3 may produce a false negative test result. Above pH 8, pesticides will decompose. Below pH 3, the test will be inactive.
3. During development, the temperature of the Ticket should be maintained between 33–37°C (91.4–98.6°F). The test will work at lower temperatures, but sensitivity may be decreased.
4. Use of ethyl acetate as a sample extraction solvent requires removal of the reagent disc from the sensor as ethyl acetate will attack the plastic Ticket support.
5. You must use the bromine (activator ampoule) step to test for the presence of thiophosphate pesticides. Without the use of the bromine step, only organophosphates and carbamates will be detected. For no ampoule testing, follow the appropriate procedure skipping the “convert” step.
6. To avoid cross-contamination, thoroughly wash glassware between samples.

TEST PROCEDURES

WATER

1. **Collect:** Place 20 mL of sample in 50 mL beaker provided with kit.
2. **Convert:** Place activator ampoule in the 20 mL of sample and crush with glass rod. Allow to react for 3 minutes.
3. **Expose:** Remove Ticket from packet. Expose white disc only, and dip into beaker for 1 minute.
4. **Develop:** Remove Ticket from beaker, and expose the second disc. Fold Ticket at the perforation, and press discs together with finger and thumb for 3 minutes.
5. **Read results:** Open and read results immediately. White means pesticides were detected, and blue indicates they were not.

SOIL — WATER EXTRACTION (FIELD)

1. **Collect:** Collect representative soil sample.
2. **Extract:** Add equal parts of soil and water to a wide mouth jar, and shake vigorously for 1 minute.
3. **Filter:** Pour 20 mL of soil/water mixture through filter paper into the 50 mL beaker.
4. **Convert:** Place an activator ampoule in the filtered sample and crush ampoule with glass rod. Allow to react for 3 minutes.
5. **Expose:** Remove a Ticket from the packet, expose the white disc only, and dip the white disc into the beaker for 1 minute.
6. **Develop:** Expose the second disc. Fold Ticket at the perforation, and press discs together with finger and thumb for 3 minutes.
7. **Read results:** Open and read results immediately. White means pesticides were detected, and blue indicates they were not.

SURFACES

1. **Wipe/clean:** Wipe or clean the surface to be tested to remove dust and oil.
2. **Expose:** Remove Ticket and fold at the perforation. Wipe surface with white disc of the Ticket.
3. **Prepare:** Put 20 mL of clean water into the 50 mL graduated beaker and break the activator ampoule in the water with the glass rod.
4. **Convert:** Add 3 drops of the diluted activator solution to the white disc. Allow to react for 2 minutes.
5. **Develop:** Expose the second disc. Fold Ticket at the perforation, and press discs together with finger and thumb for 3 minutes.
6. **Read results:** Open and read results immediately. White means pesticides were detected, and blue indicates they were not.

SPILLS

1. **Collect:** Collect a spill sample of approximately 5 mL in the 50 mL graduated beaker.
2. **Dilute:** Add clean water to bring the volume of the sample and the water to 20 mL.
3. **Convert:** Place activator ampoule in the solution and crush ampoule with glass rod. Allow to react for 3 minutes.
4. **Expose:** Remove Ticket from packet, expose the white disc only, and dip white disc into the beaker solution for 1 minute.
5. **Develop:** Remove Ticket from the beaker and expose the second disc. Fold Ticket at the perforation, and press discs together with finger and thumb for 3 minutes.
6. **Read results:** Open and read results immediately. White means pesticides were detected, and blue indicates they were not.

AIR — PESTICIDE DRIFT AND WORKER EXPOSURE

1. **Expose:** Remove Ticket from the packet. Secure it to a tree, fence post or worker's clothing before spraying or handling pesticides.
2. **Collect:** After spraying, code and note the location of each Ticket, and remove to a clean environment.
3. **Prepare:** Break activator ampoule in 20 mL of clean water.
4. **Convert:** Add 3 drops of diluted activator solution to the white disc, and allow to react for 2 minutes.
5. **Develop:** Expose the second disc. Fold Ticket at the perforation, and press discs together with finger and thumb for 3 minutes.
6. **Read results:** Open and read results immediately. White means pesticides were detected, and blue indicates they were not.

AIR — WORKER RE-ENTRY

1. **Collect:** Collect 50 leaf punch samples.
2. **Extract:** Place the 50 leaf samples in a wide mouth jar, add 50 mL of clean water, and shake vigorously for 1 minute.
3. **Prepare:** Pour 20 mL of extract into the 50 mL beaker and place an activator ampoule into the beaker. Break the ampoule with a glass rod, and allow to react for 2 minutes.
4. **Expose:** Remove a Ticket from the packet, expose the white disc only, and dip the white disc into the beaker for 1 minute.
5. **Develop:** Expose the second disc. Fold Ticket at the perforation, and press discs together with finger and thumb for 3 minutes.
6. **Read results:** Open and read results immediately. White means pesticides were detected, and blue indicates they were not.

GRAINS AND PRODUCE — FIELD (for water soluble pesticides only)

Note: For the testing of non-water soluble pesticides, contact your Neogen representative.

1. **Collect:** Collect a representative sample.
2. **Extract:** Add equal weights of sample and water to a wide mouth jar. Shake for 1 minute.
3. **Prepare:** Pour 20 mL of water extract into the 50 mL beaker.
4. **Convert:** Break activator ampoule into the extract in the beaker, and allow to incubate for 2 minutes.
5. **Expose:** Remove a Ticket from the packet, expose the white disc only, and dip the white disc into the beaker for 1 minute.
6. **Develop:** Expose the second disc. Fold Ticket at the perforation, and press discs together with finger and thumb for 3 minutes.
7. **Read results:** Open and read results immediately. White means pesticides were detected, and blue indicates they were not.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

ORDERING INFORMATION

- 8920 AT-10 Kit – contains 10 tickets, 10 ampoules, beaker and rod
- 8920A AT-10 Kit – contains 10 tickets, beaker and rod
- 8921 AT-25 Kit – contains 25 tickets, 25 ampoules, beaker and rod
- 8921A AT-25 Kit – contains 25 tickets, beaker and rod
- 8950 Agri-Screen Tickets – contains 10 tickets
- 8951 Activator Ampoules – contains 10 ampoules
- 8956 Activator Ampoules – contains 5 ampoules
- 8953 Concentrator
- 8954 Beaker
- 8955 Rod, Glass

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts

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©Neogen Corporation, 2015. Neogen, Alert, K-Blue y Agri-Screen son marcas comerciales registradas de Neogen Corporation, Lansing, MI. Todas las otras marcas y nombres de productos son marcas comerciales o marcas comerciales registradas de sus respectivas compañías.