Read instructions carefully before starting test

ALERT[®] **for Gliadin**

Refrigerate at 2–8°C (35–46°F) • Do not freeze

GLIADIN

Gliadin is an alcohol-soluble protein found in wheat that belongs to a group of proteins called prolamins. Other prolamins include secalin, found in rye, and hordein, found in barley. Gluten consists of two groups of proteins (prolamins and glutelins) that are found in differing amounts in wheat, barley, rye and oats.

Gliadin and other prolamins have been identified as major causal agents in a number of disorders, including wheat allergy and gluten intolerance (celiac disease). Wheat allergy is a specific immune response to a number of wheat proteins, including gliadin, albumin, globulin, and glutenin. Celiac disease is a chronic reaction to gluten proteins that results in the poor absorption of nutrients in the small intestine.

Those with celiac disease must avoid gluten, and rely upon the correct labeling of food to make appropriate, safe food choices. Testing for the presence of gluten components ensures food manufacturers that an unlabeled—and potentially dangerous—ingredient did not make its way into a food product.

INTENDED USE

Alert for Gliadin is intended for the qualitative analysis of ingredients, clean-in-place solutions, finished food products, and environmental surfaces intended to be gluten-free for the presence of gliadins and prolamins found in wheat, rye and barley.

INTENDED USER

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by gluten. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed the Neogen training.

ALERT[®] for Gliadin

ASSAY PRINCIPLES

Alert for Gliadin is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Gliadin is extracted from samples with a 40% ethanol solution by shaking in a shaker or rotator. Extract is diluted in phosphate saline buffer and diluted samples are added to antibody-coated wells (capture antibody) where gluten will bind to the antibody during an incubation period. Any unbound gliadin is washed away and a second antibody, which is enzyme labeled (detector antibody) is added. The detector antibody binds to the gliadin during another incubation period. Unbound enzyme-labeled antibody is washed away and a one step substrate is added. Color develops as a result of the presence of bound-labeled antibody. A stopping reagent is added and the color of the solution is observed. Blue color indicates samples containing high levels of gliadin while purple or red samples contain little or no gliadin.

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F).

MATERIALS PROVIDED

1. 24 antibody-coated microwells
2. 1 yellow-labeled dropper bottle of 10 ppm gliadin control (25 ng/mL gliadin)
3. 1 blue-labeled dropper bottle of conjugate
4. 1 green-labeled dropper bottle of K-Blue[®] Substrate
5. 1 red-labeled dropper bottle of Red Stop solution
6. 20 prefilled dilution dropper bottles of 3.9 mL PBS
7. 1 foil pouch containing enough 10 mM PBS dry powder to prepare 1 L of dilution solvent
8. 1 wide-mouth bottle containing enough 10 mM PBS-Tween wash buffer concentrate to prepare 1 L in distilled or deionized water (pH 7.4)
9. 30 g of extraction additive in a specimen cup
10. Plastic 1 g scoop to measure extraction additive

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Orbital rotator, shaker or vortex to hold 50 cc centrifuge tubes for a 1 g sample or shaker or shaker water bath with clamps adjusted to hold 125 mL (4 oz) extraction bottles for a 2 g sample
2. Allergen Environmental Swabbing Kit (Neogen item 8432S)
3. Scale capable of weighing 1 + 0.1 g (Neogen item 9427), or scale capable of weighing 0.25 ± 0.01 g if using cocktail solution
4. 2 1 L bottles to prepare washing solution, extract solution (Neogen item 9472)
5. Timer (Neogen item 9426)
6. Microwell holder (Neogen item 9402)
7. Pipettor, 100 µL (Neogen item 9272, 9278, or 9276)
8. Pipette tips (Neogen item 9410)
9. Wash bottle (Neogen item 9400)
10. Paper towels or equivalent absorbent material
11. Waterproof marker
12. Distilled or deionized water
13. Laboratory grade ethanol (190 proof)
14. Special extraction additive for dark chocolate, cocoa and tannin (Neogen item 8482)
15. Cocktail solution for analysis of heat processed samples (Neogen item 8483)

PRECAUTIONS

1. Ethanol solution is highly flammable. Keep container tightly closed, and keep away from heat, sparks, open flame and those smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
2. Use a chemical hood when using the cocktail solution to analyze heat processed samples.
3. Components of Alert for Gliadin, such as controls and extraction additive, may contain one or more of the following potentially allergic materials: gluten, casein, almond protein and soy protein. If allergic to any of these compounds, use caution when using this product.
4. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze test kits, and avoid prolonged storage at ambient temperatures.
5. Bring kits to room temperature (18–30°C, 64–86°F) prior to use.
6. Do not use kit components beyond expiration date.
7. Do not mix reagents from one kit with reagents from a kit with a different serial number.
8. Do not run more than six wells per test.
9. Use only incubation times specified; others may give inaccurate results.

PROCEDURAL NOTES

1. **Dark chocolate, cocoa and tannin.** If testing commodities where dark chocolate, cocoa, and/or tannin are present, such as dark chocolate bars and cocoa powder, contact Neogen for a special extraction additive (Neogen item 8482). Add 1 scoop of this special extraction additive in addition to 1 scoop of the additive supplied with this test kit.
2. **Wash buffer.** Prepare the wash buffer solution by diluting the wash buffer concentrate in 960 mL of distilled or deionized water in an empty 1 L container. Ensure the transfer of all concentrate by rinsing the bottle several times with water, and swirl to ensure thorough mixing.
NOTE: Discard unused portions of extraction solution and wash buffer when the test kit has been used completely.
3. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.
4. **Sample extract dilution solution (PBS).** Prepare extract dilution solution by adding a foil pouch of dilution solvent, 10 mM PBS, to 1 L distilled or deionized water. Swirl to mix thoroughly.

SAMPLE PREPARATION AND EXTRACTION

For analyzing **heat processed samples**, follow extraction **procedure C**. For all commodities that were **not heat processed**, follow either extraction procedure **A or B**. Samples of an unknown origin should be extracted using extraction procedure **C**.

A. Extraction of non-heat processed samples with orbital shaker, rotator or vortex

1. Prepare 40% ethanol extraction solution by combining 4 parts ethanol with 6 parts distilled water.
2. Add 1 g ground sample, or 1 mL liquid sample, to a clean 50 cc tube.
3. Add 1 scoop of extraction additive to the tube.
4. Add 10 mL (9 mL for liquid samples) of 40% ethanol to the tube, cap tightly, then shake the tube vigorously by hand for about **20 seconds**, or vortex for **10 seconds**, to ensure complete mixing.
5. Extract by shaking (150 rpm) in an orbital shaker or rotator by laying down the tube on its side over the flat pad of the instrument, and holding it tightly using a rubber band or tape. Rotate or shake for **15 minutes** at room temperature. **ALTERNATIVE:** Vortex on high speed for **4 minutes**.
6. Remove the tube and let it stand in a rack for about **10 minutes** to enable the sample extract to settle before withdrawing the clear extract.
7. Dilute each sample 1:40 by withdrawing 100 µL of the upper layer of the extract and transferring it to a pre-filled dilution dropper bottle containing 3.9 mL of PBS. (Cap and tip must be removed first, then reinserted).
8. To mix, invert several times by hand.
9. Test diluted samples within **2–3 hours** of extraction.

B. Extraction of non-heat processed samples with shaker or shaker water bath

1. Prepare 40% ethanol extraction solution by combining 4 parts ethanol with 6 parts distilled water.
2. Add 2 g ground sample, or 2 mL liquid sample, to a 125 mL clean extraction bottle.
3. Add 1 scoop of extraction additive to the bottle.
4. Add 20 mL (18 mL for liquid samples) of 40% ethanol, cap the bottle tightly, then shake vigorously by hand for about **20 seconds** to ensure complete mixing.
5. Extract by shaking (150 rpm) in a shaker for **15 minutes** at room temperature (a shaker water bath can work, but do not turn the heat on). Remove the bottle from shaker or bath.
6. Let the bottle stand for about **10 minutes** to enable some of the sample to settle before withdrawing the clear extract.
7. Dilute each sample 1:40 by withdrawing 100 μ L of the upper layer of the extract and transferring it to a pre-filled dilution dropper bottle containing 3.9 mL of PBS. (Cap and tip must be removed first, then reinserted).
8. To mix, invert several times by hand.
9. Test diluted samples within **2–3 hours** of extraction.

C. Extraction of heat processed commodities

These commodities require a special pre-extraction cocktail solution that renature the heated gluten and allowed the proper determination of the gluten in the sample (Neogen item 8483). In analyzing heat processed samples, follow these instructions. Prefilled dilution bottles cannot be used for this extraction procedure.

1. Prepare 55% ethanol extraction solution by combining 55 parts ethanol with 45 parts distilled water.
2. Prepare sample extract dilution solution (PBS) as detailed in procedural note 4.
3. It is recommended to perform the extraction with the cocktail extraction additive under a chemical hood.
4. Weigh out 0.25 g sample to 50 cc screw cap centrifuge tube.
5. Add 2.5 mL of extraction cocktail solution (dilution factor 1:10).
6. Cap and vortex **10–20 seconds** to homogenize cocktail and sample.
7. Incubate **40 minutes** at 50°C (water bath or oven).
8. Remove samples and let cool for **5–10 minutes**.
9. Add 1 scoop of powdered additive.
10. Add 7.5 mL of 55% ethanol and Vortex again for **10–20 seconds** (the final concentration of ethanol will be 41% and the sample dilution to this point 1:40).
11. Shake (150–200 rpm) for 1 hour at room temperature on a rotator (tube on its side).
12. Centrifuge sample (if necessary) for **5 minutes** at 2500 rpm.
13. Dilute the sample 1:10 into 10 mM PBS, pH 7.4 (200 μ L sample into 1.8 mL PBS).
14. Samples are ready to run (1:400 final dilution).

TEST PROCEDURE

Allow the test kit and all reagents to warm to room temperature (18–30°C, 64–86°F) before using.

1. Remove 1 well for each sample to be tested plus 1 well for the control, and place into the well holder.
2. Mix each reagent by swirling its dropper bottle prior to use.
3. Add 3 drops from the yellow-labeled control dropper bottle to the first well. Add 3 drops from each diluted sample dropper bottle to a respective well as indicated in the template below. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.

Control S1 S2 S3 S4 S5

4. Incubate microwells **10 minutes** at room temperature (18–30°C, 64–86°F).
5. Shake out the contents of the wells. Using a wash bottle filled with wash buffer, fill each well and shake out. Repeat 5 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
6. Add 3 drops from the blue-labeled conjugate dropper bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.
7. Incubate for **10 minutes** at room temperature (18–30°C, 64–86°F).
8. Shake out the contents of the wells. Using a wash bottle filled with wash buffer, fill each well and shake out. Repeat 5 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
9. Add 3 drops from the green-labeled substrate dropper bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.
10. Incubate for **10 minutes** at room temperature (18–30°C, 64–86°F).
11. Add 3 drops from the red-labeled Red Stop dropper bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface. The results are now ready to be interpreted.

INTERPRETATION OF RESULTS

Visually compare the color of a sample well to the color of the control well. If the sample well has **more blue color** than the control well, the sample tests positive for gliadin contamination of **more than 10 ppm**. If the sample well has **less blue color or more red color** than the control well, the sample contains **less than 10 ppm** of gliadin contamination. **NOTE:** Standard controls were made from wheat gliadin and calculated as gliadin. Approximately 50% of the gluten is available as gliadin. Therefore, 10 ppm gliadin is equal to 20 ppm gluten.

ALTERNATIVE: Read wells (wipe bottom of wells first) in a microwell reader with a 650 nm filter. If the sample well has an optical density (OD) higher than the control well, the sample is positive for gliadin contamination of more than 10 ppm. If the sample well has an OD lower than the control well, the sample contains less than 10 ppm of gliadin contamination.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN**Natural toxins**

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts

Genetic modification

- CP4 (Roundup Ready®)

Ruminant by-products

- Meat and bone meal, feed

Species identification

- Raw and cooked meat samples

**North America****Neogen Headquarters**

620 Leshler Place, Lansing, MI 48912 USA
800/234-5333 (USA/Canada) or 517/372-9200
Fax: 517/372-2006 • foodsafety@neogen.com
www.neogen.com

Europe, Middle East and Africa**Neogen Europe**

The Dairy School, Auchincruive, Ayr
KA6 5HU Scotland, UK
+ 44 (0) 1292 525 600
Fax: + 44 (0) 1292 525 601
info_uk@neogeneurope.com
www.neogeneurope.com

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