*Read instructions carefully before starting test*

# Veratox<sup>®</sup> DON HS

## Quantitative High Sensitivity Test

*REFRIGERATE AT 2-8°C • DO NOT FREEZE*

### THE TOXIN

Deoxynivalenol (DON) is most commonly produced by the pink mold *Fusarium graminearum*. DON, a member of the trichothecene family, is produced by fungi living on cereal commodities such as wheat, corn, barley and ensilages. The toxicological effects attributed to DON include: nausea (vomiting), feed refusal, gastroenteritis, diarrhea, immunosuppression and blood disorders.

Accurate determination of the presence of the toxin is of major importance to those monitoring the quality of feed and food in which DON may occur. Testing these commodities for the toxin requires careful sampling, extraction, sanitation and quantitative analysis. The European Community (EC) is currently investigating maximum allowable levels for DON in foodstuffs. The Working Group on Agricultural Contaminants and the Standing Committee on the Food Chain and Animal Health have proposed the following maximum levels to be considered by the member states of the EC:

COMMODITIES	MAXIMUM LEVEL (ppb)
1. Unprocessed cereals other than durum wheat and corn	1,500
2. Unprocessed durum wheat and corn	2,000
3. Cereals for direct human consumption or use as a food ingredient, all products derived from cereals, except those specified in items 4, 5 and 6 below	500
4. Corn-based breakfast cereals and corn-based snacks	500
5. Wholemeal wheat flour, bran and pasta (dry)	750
6. Cereal-based food for infants and young children, and sold as such, and ingredients used for manufacture of these products	100

### INTENDED USE

Veratox for DON HS is intended for the quantitative analysis of DON in commodities like corn, wheat, barley, and processed cereal.

### INTENDED USER

The test kit is designed for use by quality control personnel and others familiar with food and feed possibly contaminated by DON. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed the Neogen training.

## ASSAY PRINCIPLES

Veratox for DON HS is a competitive direct enzyme-linked immunosorbent assay (CD-ELISA) which allows the user to obtain exact concentrations in parts per billion (ppb). Free DON in the samples and controls is allowed to compete with enzyme-labeled DON (conjugate) for the antibody binding sites. After a wash step, substrate is added, which reacts with the bound conjugate to produce blue color. More blue color means less DON. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form the standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of DON.

## STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2-8°C (35-46°F).

## MATERIALS PROVIDED

1. 48 antibody-coated microwells
2. 48 red-marked mixing wells
3. 5 yellow-labeled bottles of 0, 25, 50, 100 and 250 ppb DON controls
4. 1 blue-labeled bottle of DON-HRP conjugate solution
5. 1 green-labeled bottle of K-Blue® Substrate solution
6. 1 red-labeled bottle of Red Stop solution

## MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Extraction materials:
  - a. Distilled or deionized water
  - b. Graduated cylinder (Neogen item #9447/#9368)
  - c. Container with 125 mL capacity (Neogen item #9428)
  - d. Whatman #1 filter paper, Neogen filter syringe, or equivalent (Neogen item #9430/#9420)
  - e. Sample collection tubes (Neogen item #9421)
2. Agri-Grind grinder, or equivalent (Neogen item #9401)
3. Scale capable of weighing 5-50 grams (Neogen item #9427)
4. Microwell reader with a 650 nm filter (Neogen item #9302)
5. Pipettor, 12-channel (Neogen item #9273)
6. Pipettor, 100 µL (Neogen item #9272)
7. Tips for 100 µL and 12-channel pipettors (Neogen item #9410/#9407)
8. Paper towels or equivalent absorbent material
9. Plastic bucket for use as waste receptacle
10. Microwell holder (Neogen item #9402)
11. Timer (Neogen item #9426)
12. Waterproof marker
13. Wash bottle (Neogen item #9400)
14. 2 reagent boats for 12-channel pipettor (Neogen item #9435)
15. Distilled or deionized water

## PRECAUTIONS

1. Store test kit between 2-8°C (35-46°F) when not in use.
2. Do not use kit components beyond expiration date.
3. Do not mix reagents from one kit serial with reagents from a different kit serial.
4. Do not run more than 24 wells at one time.
5. Follow proper pipetting techniques, including the tip priming.
6. Use of incubation times other than those specified may give inaccurate results.
7. Kits should be brought to room temperature (18-30°C, 64-86°F) before use.
8. Avoid prolonged storage of kits at ambient temperatures.
9. Do not freeze test kits.
10. Commodity extracts should have a pH of 6-8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH contact your Neogen representative or Technical Services.
11. Treat all used liquids, including sample extract, and labware as if contaminated with DON. Gloves and other protective apparel should be worn at all times.
12. To avoid cross-contamination, use clean pipette tips and glassware for each sample, and thoroughly detoxify and wash all glassware between samples.

## PROCEDURAL NOTE

1. **Substrate:** K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until it is needed.
2. **Antibody wells:** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test is set to begin.

## SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques. The sample should be ground and thoroughly mixed prior to proceeding with the extraction. Store samples at 2-8°C (35-46°F) until analyzed.

1. Obtain a representative sample. Grind the entire sample so that at least 75% of the ground material passes through a 20 mesh sieve, the particle size of a fine instant coffee.
2. Using hand or mechanical means, vigorously shake 10 grams of ground sample in 50 mL of distilled or deionized water for 3 minutes.
3. Let material set for 2-3 minutes to enable some of the sample to settle before filtering extract.
4. Filter the extract by pouring at least 5 mL through a Whatman #1 filter (or Neogen filter syringe) and collecting the filtrate as a sample.
5. The sample is ready for testing.

## TEST PROCEDURE

Allow all reagents to warm to room temperature (18-30°C, 64-86°F) before use.

1. Remove 1 red-marked mixing well for each sample to be tested plus 5 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant and seal to protect the antibody. Mark one end of strip with a "1", and place strip in the well holder with the marked end on the left. Do not mark the inside or bottom of the wells.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Place 100 µL of conjugate from the blue-labeled bottle in each red-marked mixing well.
5. Using a new pipette tip for each, transfer 100 µL of controls and samples to the red-marked mixing wells as described below.

0	25	50	100	250	S1	S2	S3	S4	S5	S6	S7	Strip 1
S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19	Strip 2

6. Using a 12-channel pipettor, mix the liquid in the wells by pipetting it up and down 3 times. Transfer 100 µL to the antibody-coated wells. Mix by sliding the microwell holder back and forth on a flat surface for 10-20 seconds without splashing reagents from the wells. Incubate for **10 minutes** at room temperature (18-30°C, 64-86°F). Discard red-marked mixing wells.
7. Shake out the contents of the antibody wells. Fill the wells with distilled or deionized water and dump them out. Repeat this step 5 times, then turn the wells upside-down and tap out on a paper towel until the remaining water has been removed.
8. Pour the needed volume of substrate from the green-labeled bottle into the green-labeled reagent boat.
9. With new tips on the 12-channel pipettor, prime and pipette 100 µL of substrate into the wells and mix by sliding back and forth on a flat surface for 20 seconds.
10. Incubate **10 minutes**. Discard remaining substrate and rinse the reagent boat with water.
11. Pour Red Stop solution from the red-labeled bottle (same volume as the substrate) into the red-labeled reagent boat.
12. Eject the excess substrate from the 12-channel pipettor, prime the tips, and pipette 100 µL of Red Stop to each well. Mix by sliding back and forth on a flat surface. Discard the tips.
13. Wipe the bottom of the microwells with a dry cloth or towel and read in a microwell reader using a 650 nm filter. Air bubbles should be eliminated, as they could affect analytical results. Results should be read within 20 minutes after the addition of Red Stop.
14. Read and calculate results using Neogen's Stat-Fax microwell reader. If using an EL301 reader or other strip/plate reader, calculate results using Neogen's Veratox for Windows software.

## RETESTING

If positives occur in commodities not previously tested, confirm with an additional approved method prior to taking action. The test kit will not differentiate between DON and 3 acetyl DON.

## PERFORMANCE CHARACTERISTICS

Limit of quantitation: 25 ppb (0.025 ppm)

Range of quantitation: 25-250 ppb (For quantitating samples above 250 ppb, contact Neogen Technical Services for dilution instructions.)

Validation Matrixes: Barley, beet pulp\*, corn, corn bran\*\*, corn meal, corn germ meal\*, corn gluten feed\*, corn gluten meal\*, corn screenings, corn silage\*\*, corn/soy blend, corn starch, DDGS\*, flaxseed meal, hay\*\*, haylage\*\*, kamut, malted barley\*, milo, oats, petfood\*, popcorn, quinoa, rice, rice hulls, rye, soy flour, soy meal, sunflower meal, TMR\*\*, wheat, wheat bran, wheat flour and wheat midds.

\*May require pH adjustment.

\*\*Contact Neogen for special procedure.

## CUSTOMER SERVICE

Customer assistance and Technical Services for Veratox for DON HS can be reached by calling 800/234-5333 (USA/Canada) or 517/372-9200, 8 a.m. to 7 p.m. Eastern time. Assistance is available on a 24-hour basis by calling 800/867-0308. Training on this product, and all Neogen test kits, is available.

## MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's Food Safety test kits, at [www.neogen.com](http://www.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

## WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### TESTING KITS AVAILABLE FROM NEOGEN

#### Natural Toxins

- Aflatoxin, DON, Ochratoxin, Zearalenone, T-2/HT-2 Toxin, Fumonisin, Histamine

#### Foodborne Bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

#### Sanitation

- ATP, Yeast and Mold, Total Plate Count, Generic *E. coli* and Total Coliforms, Protein Residues

#### Food Allergens

- Peanuts, Milk, Eggs, Almonds, Gluten, Soy, Hazelnut

#### Genetic Modification

- CP4 (Roundup Ready®)

#### Ruminant By-products

- Meat and Bone Meal, Feed

620 Leshar Place • Lansing, MI 48912  
800/234-5333 (USA/Canada) or 517/372-9200 • fax: 517/372-2006  
e-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com) • [www.neogen.com](http://www.neogen.com)

©Neogen Corporation, 2009. Neogen, Veratox and K-Blue are registered trademarks of Neogen Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.

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