Read instructions carefully before starting test

Veratox[®] **for Mustard Allergen**

Quantitative Test

REFRIGERATE AT 2–8°C (35–46°F) • DO NOT FREEZE

MUSTARD ALLERGEN

Food allergens are proteins in food that can create an immune response in sensitive individuals. Once ingested, food allergens can cause a number of reactions, ranging in severity from hives and itching to anaphylaxis. Anaphylaxis is a severe allergic reaction, involving vomiting, diarrhea, difficulty breathing, swelling of the mouth and tongue, and a rapid drop in blood pressure.

An estimated 3.5% to 4% of adults, and 6% to 8% of children, are sensitive in some degree to food allergens. More than 12 million people in the United States alone are known to have a food allergy.

Food manufacturers protect those with food allergies by clearly labeling their products with a list of ingredients. Testing for the presence of mustard proteins assures food manufacturers that an unlabeled — and potentially dangerous — ingredient did not make its way into a food product.

INTENDED USE

Veratox for Mustard Allergen is intended for the quantitative analysis of mustard proteins in food products such as spices, sauces, dressings, meats and clean-in-place rinses.

INTENDED USER

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by mustard. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed the Neogen training.

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). Do not freeze.

ASSAY PRINCIPLES

The Veratox for Mustard Allergen test is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Mustard protein residue is extracted from samples with a Tris-EDTA solution by shaking in a heated water bath, followed by centrifugation or filtration. Extracted mustard protein is sampled and added to antibody-coated wells (capture antibody) where it binds to the antibody during an incubation. Any unbound mustard protein is washed away and a second antibody (detector antibody), which is enzyme labeled, is added. The detector antibody binds to the already bound mustard protein. After a second wash, substrate is added. Color develops as a result of the presence of bound detector antibody. Red Stop reagent is added and the color of the resulting solution is observed. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form a standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of mustard, expressed as ppm mustard.

MATERIALS PROVIDED

1. 48 antibody-coated microwells
2. 48 red-marked transfer wells
3. 5 yellow-labeled bottles of 0, 2.5, 5, 10, 25 ppm mustard controls
4. 2 blue-labeled bottles of enzyme-labeled antibody conjugate
5. 1 green-labeled bottle of K-Blue® Substrate
6. 1 red-labeled bottle of Red Stop solution
7. 5 pouches of Veratox mustard allergen extraction buffer. Each pouch is enough to prepare 1 L in distilled or deionized water.
8. 40 mL of 10 mM PBS-Tween washing reagent in a wide mouth bottle. Each bottle is enough to prepare 1 L in distilled or deionized water (pH 7.4)
9. 50 g of extraction additive in a specimen cup
10. Plastic scoop to measure extraction additive

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Allergen Extraction Kit (Neogen item 8429)
 - a. 20 disposable plastic extraction bottles
 - b. 20 sample collection tubes (12 x 75 mm) with caps
2. Allergen Environmental Swabbing Kit (Neogen item 8432S)
 - a. 100 sterile swabs
 - b. 100 dropper tips
3. Shaker water bath capable of maintaining $60 \pm 1^\circ\text{C}$ with holders for 250 mL extraction bottles (Neogen items 9298, 9299)
4. Whatman no. 4 filters or equivalent (Neogen item 9429)
5. Centrifuge (optional)
6. Filter funnel, 65 mm (Neogen item 9353)
7. Pipettor, 50–200 μL , adjustable (Neogen item 9276)
8. Pipettor, 12-channel (Neogen item 9273)
9. Pipette tips (Neogen item 9410, 9407, 9417)
10. Timer (Neogen items 9426, 9452)
11. Microwell reader with a 650 nm filter (Neogen item 9303)
12. Bottle, 1 L, to prepare washing solution (Neogen item 9472)
13. Bottle, 1 L heat safe, to prepare extract solution (Neogen item 9472)
14. Paper towels or equivalent absorbent material
15. Microwell holder (Neogen item 9402)
16. Waterproof marker
17. Wash bottle (Neogen item 9400)
18. Distilled or deionized water
19. 3 reagent boats for 12-channel pipettor (Neogen item 9435)

20. Graduated cylinder capable of measuring 125 mL (Neogen item 9368)
21. Scale capable of weighing 5 ± 0.1 g (Neogen item 9427)

PRECAUTIONS

1. Components of Veratox for Mustard Allergen, such as controls and extraction reagents, may contain one or more of the following potentially allergenic materials: mustard protein; casein, soy protein and/or tree nut protein. If allergic to any of these compounds, please use caution when handling this product.
2. Concentrated food additives, colors and flavors may cause interferences on ELISA test methods. Contact Neogen Technical Services for validation information.
3. Hydrolyzed and fermented proteins may not be detected using ELISA methods for allergen testing. Due to the breakdown of the proteins to small peptides or amino acids, proteins may become undetectable by this assay, but still could be allergenic and cause an allergic reaction.
4. Sponges should not be used for sample collection and allergen testing. Sample collection swabs other than Neogen swabs should be validated prior to use. General sponges and swabs may contain solutions or materials that may interfere with the test kit.
5. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze.
6. Bring kits to room temperature 18–30°C (64–86°F) prior to use.
7. Avoid prolonged storage of kits at ambient temperatures.
8. Do not use kit components beyond expiration date.
9. Do not mix reagents from one kit with reagents from a kit with a different serial number.
10. Do not run more than 24 wells per test.
11. Follow proper pipetting techniques (e.g., prime tips and use clean tips).
12. Use only incubation times specified. Others may give inaccurate results.
13. Use clean pipette tips and glassware for each sample to avoid cross-contamination. Thoroughly wash all glassware between samples.

PROCEDURAL NOTES

1. **Substrate:** K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until needed.
 2. **Conjugate:** The conjugate supplied with this kit is ready to use. One bottle is enough for 24 wells. Cover the reagent boat to keep the conjugate protected from direct light and contaminants.
 3. **Antibody wells:** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.
 4. **Mustard allergen extraction solution:** Prepare extraction solution by adding a pouch of Veratox Mustard Allergen extraction buffer, Tris/EDTA, to 1 L distilled or deionized water. Mix by swirling or shaking to thoroughly dissolve all buffer. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).
 5. **Wash buffer:** Prepare the wash buffer solution by pouring all the wash buffer concentrate into an empty 1 L container. Rinse the wash buffer concentrate bottle with distilled or deionized water and pour into the 1 L container to ensure all the concentrate is used. Fill the 1 L container with additional distilled or deionized water, and swirl to assure thorough mixing. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).
- NOTE:** Discard unused portions of extraction solution and wash buffer when the test kit has been used completely.

SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques (see Neogen's Food Allergen Handbook). The sample should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

1. Prepare the mustard allergen extraction solution as described in the procedural notes.
2. Preheat extraction solution to 60°C (140°F) by immersing the bottle containing the solution into the water bath and allowing it to reach 60°C.
3. Using an appropriate sampling and collection procedure, obtain a representative sample and grind it to a very fine particle size.
4. Transfer 5 g of sample or 5 mL of liquid sample into a 250 mL disposable extraction bottle.
5. Add 1 level scoop of the extraction additive to the sample bottle.
6. Pour 125 mL of the 60°C (140°F) extraction solution to the sample bottle.
7. Cap the sample bottle to prevent contents from splashing during the extraction.
8. Extract by shaking (150 rpm) in a water bath at 60°C (140°F) for **15 minutes**. Remove the bottle from the bath.
9. Let material settle for **5 minutes** to enable some of the sample to settle before proceeding to the next step.
10. Filter the extract by pouring at least 5 mL through a Whatman no. 4 filter and collecting the filtrate as a sample.
ALTERNATIVE: Centrifuge at 14,000 rpm for **5 minutes (20 minutes** for lower speeds). Use the clear supernatant as a sample.
11. Allow extracts to cool to room temperature before beginning analysis.
12. Discard extracts after completion of analysis.

ENVIRONMENTAL SWABBING PREPARATION AND EXTRACTION

1. Prepare the extraction solution as described in the procedural notes.
2. Prepare a working solution by combining 125 mL of the extraction solution to a separate container, and add 1 level scoop of extraction additive to a container. Prepare a fresh working solution daily.
NOTE: Prepare 125 mL of working solution for every 25 environmental swabs to be tested.
3. Preheat the working solution to 60°C (140°F). Shake container to dissolve additive.
4. Gather the sample with a swab, using one of the following methods:
For dry surfaces: Open a new swab and wet with extraction solution. Swab a 10 x 10 cm area by using a crosshatch technique. (Do not use working solution to moisten swabs.)
For wet surfaces: Open a new swab and swab a 10 x 10 cm area by using a crosshatch technique. Do not moisten swab prior to use.
5. Return the swab to its original tube once sampling is complete. Remember to label each tube.
6. Remove the swab from its tube, and add 5 mL of working solution at 60°C to the tube. Mix by placing the swab back into the tube and shaking for **2 minutes** by hand (inverting tube), or for **30 seconds** with a vortex mixer.
7. Remove the swab from its tube.
8. Place a new sample dropper tip onto the tube. The solution in the tube now serves as the sample.
NOTE: Use caution when inverting the tube, as some liquid may drip.

TEST PROCEDURE FOR QUANTITATION

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

1. Remove 1 red-marked mixing well for each sample to be tested plus 5 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of the strip with a "1", and place strip in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Using a new pipette tip for each, transfer 150 µL of controls and sample extracts to the red-marked transfer wells as shown in the template below. Only run up to two 12-well strips at a time.

0	2.5	5	10	25	S1	S2	S3	S4	S5	S6	S7	Strip 1
S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19	Strip 2

5. Place tips on the 12-channel pipettor and transfer 100 µL of the controls and sample extracts to the antibody-coated wells. Mix for **20 seconds** by sliding the well holder back and forth on a flat surface.
6. Incubate microwells **10 minutes** at room temperature (18–30°C, 64–86°F). Discard the red-marked transfer wells.
7. Empty the contents of the wells into a sink. With a wash bottle fill each antibody well with the wash buffer solution and dump out. Repeat the washing 5 times, then turn the wells upside down and tap out on a paper towel until the remaining washing solution is removed.
8. Pour the needed volume of conjugate from the blue-labeled bottle into a clean reagent boat.
9. Using the 12-channel pipettor and new tips, transfer 100 µL of the conjugate into all the wells and mix for **20 seconds** by sliding the well holder back and forth on a flat surface.
10. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
11. Wash all wells with the wash buffer solution as described in step 7.
12. Pour the needed volume of substrate solution from the green-labeled bottle into a clean reagent boat.
13. Place new tips on the 12-channel pipettor and transfer 100 µL of substrate into each well and mix for **20 seconds**. Do not eject tips.
14. Incubate for **10 minutes** at room temperature (18–30°C, 64–86°F).
15. Pour the needed volume of Red Stop solution from the red-labeled bottle into a clean reagent boat.
16. With the same tips used to dispense the substrate, transfer 100 µL of Red Stop into each well and mix for **20 seconds**.
17. Wipe the bottom of the microwells and read in a microwell reader with a 650 nm filter. Results should be read within **20 minutes** after the addition of Red Stop solution.
18. Interpret the test's results using Neogen's Stat Fax microwell reader, or an equivalent strip reader. If using a strip reader, calculate the results using Neogen's Veratox software.

TEST PROCEDURE FOR SCREENING

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

1. Remove 1 well for each sample to be tested plus 1 well for the control, and place into the well holder.
2. Choose one yellow-labeled control bottle to serve as the screening level for the test. 5 ppm is recommended.
3. Mix each reagent by swirling its bottle prior to use.
4. Add 100 µL from the yellow-labeled control bottle to the first well. Add 100 µL from each sample extract to a respective well as indicated in the template below. For environmental swabs, add 3 drops from swab tube with dropper tip. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.

Control	S1	S2	S3	S4	S5
---------	----	----	----	----	----

5. Incubate microwells **10 minutes** at room temperature (18–30°C, 64–86°F).
6. Shake out the contents of the wells. Using a wash bottle filled with wash buffer solution, fill each well and shake out. Repeat 5 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
7. Add 100 µL from the blue-labeled conjugate bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.
8. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
9. Shake out the contents of the wells. Using a wash bottle filled with wash buffer solution, fill each well and shake out. Repeat 5 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
10. Add 100 µL from the green-labeled substrate bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.
11. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
12. Add 100 µL from the red-labeled Red Stop bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface. The results are now ready to be interpreted.
13. Visually compare the color of a sample well to the color of the control well. If the sample well has more blue color than the control well, the sample tests positive for mustard contamination of more than the control used. If the sample well has less blue color, or more red color, than the control well, the sample contains less than the control used of mustard contamination.

ALTERNATIVE: Read wells (wipe bottom of wells first) in a microwell reader with a 650 nm filter. If the sample well has an optical density (OD) higher than the control well, the sample is positive for mustard contamination of more than the control used. If the sample well has an OD lower than the control well, the sample contains less than the control used of mustard contamination.

PERFORMANCE CHARACTERISTICS

Limit of quantitation: 2.5 ppm mustard (Described as the lowest concentration point on the calibration curve that this test can reliably detect mustard seed protein.)

Range of quantitation: 2.5–25 ppm mustard (For quantitating samples above 25 ppm, contact a Neogen representative for dilution instructions.)

Allergen detection: This test detects mustard proteins and the results are expressed as ppm of mustard.

Cross-reactivity: This test cross-reacts with seeds from the *Brassicaceae* family (e.g., rapeseed and broccoli seed) due to their similarity in structure to mustard protein.

Protein conversion: Test kit results are expressed as mustard. To express results as protein, multiply the mustard result by 0.26 (e.g., 2.5 ppm mustard x 0.26 = 0.65 ppm mustard protein. *On average, mustard contains 26% protein.

*Source: USDA National Nutrient Database for Standard Reference Release 28, #02024 - Spices, mustard seed, ground.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's Food Safety test kits, at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/en/terms-and-conditions.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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