*Read instructions carefully before starting test*

# **ALERT<sup>®</sup>**

## **for Soy Allergen**

### **Screening Test**

*Refrigerate at 2–8°C (35–46°F) • Do not freeze.*

#### **SOY ALLERGEN**

Food allergens are proteins in food that can create an immune response in sensitive individuals. Once ingested, food allergens can cause a number of reactions, ranging in severity from hives and itching to anaphylaxis. Anaphylaxis is a severe allergic reaction, involving vomiting, diarrhea, difficulty breathing, swelling of the mouth and tongue, and a rapid drop in blood pressure.

An estimated 3.5 to 4 percent of adults, and 6 to 8 percent of children, are sensitive in some degree to food allergens. More than 12 million people in the United States alone are known to have a food allergy.

Food manufacturers protect those with food allergies by clearly labeling their products with a list of ingredients. Testing for the presence of soy ensures food manufacturers that an unlabeled — and potentially dangerous — ingredient did not make its way into a food product.

#### **INTENDED USE**

Alert for Soy is intended for the qualitative analysis of soy protein residue in products such as cookies, crackers, meat, cereals, granola/power bars, ice cream, salad dressings, cereals, juices, milk, in clean-in-place rinses and on environmental surfaces.

#### **INTENDED USER**

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by soy or soy products. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed the Neogen training.

## ASSAY PRINCIPLES

Alert for Soy is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Soy protein residue is extracted from samples by blending with a hot buffered salt solution (PBS). Extracted soy protein is sampled and added to antibody-coated wells (capture antibody) where it binds to the antibody during an incubation. Any unbound soy protein is washed away and a second antibody (detector antibody), which is enzyme labeled, is added. The detector antibody binds to the already bound soy protein. After a second wash, the substrate is added. Color develops as a result of the presence of bound detector antibody. Red Stop reagent is added and the color of the resulting solution is observed.

## STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). Do not freeze.

## MATERIALS PROVIDED

1. 24 antibody-coated microwells
2. 1 yellow-labeled dropper bottle of 2.5 ppm soy
3. 1 blue-labeled dropper bottle of enzyme-labeled antibody conjugate
4. 1 green-labeled dropper bottle of K-Blue® Substrate solution
5. 1 red-labeled dropper bottle of Red Stop solution
6. 20 unlabeled empty sample droppers with tips
7. 3 foil pouches of 10 mM PBS dry powder extraction solvent; each pouch contains enough powder to prepare 1 L of extraction solution
8. 1 wide-mouth bottle of 10 mM PBS-Tween wash buffer concentrate; a bottle contains enough concentrate to prepare 1 L of wash buffer
9. Extraction additive in a specimen cup
10. Plastic scoop to measure extraction additive

## MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Allergen Extraction Kit (Neogen item 8429)
  - a. 20 disposable extraction bottles
  - b. 20 disposable transfer pipettes
2. Allergen Environmental Swabbing Kit (Neogen item 8432s)
  - a. 100 sterile swabs
  - b. 100 disposable dropper tips
3. High-speed blender with a 250 mL jar (Neogen items 9493/9477)
4. Scale capable of weighing  $5 \pm 0.1$  g (Neogen item 9427)
5. Shaker water bath, hot plate or equivalent heat source capable of maintaining  $60 \pm 1^\circ\text{C}$
6. Thermometer
7. Whatman no. 4 filter paper or equivalent (Neogen item 9429)
8. Two 1 L bottles to prepare washing and extraction solutions (Neogen item 9472)
9. Paper towels or equivalent absorbent material
10. Microwell holder (Neogen item 9402)
11. Waterproof marker
12. Wash bottle (Neogen items 9400, 9399)
13. Distilled or deionized water
14. Timer (Neogen item 9426)
15. Graduated cylinder capable of measuring 125 mL (Neogen item 9368)

## PRECAUTIONS

1. Samples intended to be tested for soy must be extracted separately from samples intended to be tested for other food allergens. The extraction additives for each type of test are designed specifically for the target food allergen.
2. Components of the Alert for Soy Allergen test kit, such as controls and extraction reagents, may contain one or more of the following potentially allergenic materials: milk, egg protein, peanut protein or soy protein. If allergic to any of these compounds, please use caution when using this product.
3. Concentrated food additives, colors and flavors may cause interferences on ELISA test methods. Contact Neogen Technical Services for validation information.
4. Hydrolyzed and fermented proteins may not be detected using ELISA methods for allergen testing. Due to the nature of the proteins they may be undetectable in the assay, but there could still be active allergenic protein residue present.
5. Sponges should not be used for sample collection and allergen testing. Sample collection swabs other than Neogen swabs should be validated prior to use. General sponges and swabs may contain solutions or materials that may interfere with the test kit.
6. The testing environment should be clean and dust-free.
7. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze.
8. Bring kits to room temperature 18–30°C (64–86°F) prior to use.
9. Avoid prolonged storage of kits at ambient temperatures.
10. Do not use kit components beyond expiration date.
11. Do not mix reagents from a kit with one lot number with reagents from a kit with a different lot number.
12. Do not run more than 6 wells per test.
13. Use only incubation times specified. Others may give inaccurate results.

## PROCEDURAL NOTES

1. **Extraction solution.** Prepare extraction solution by adding one foil pouch containing 10 mM Phosphate Buffered Saline (PBS) to 1 L distilled or deionized water. Swirl to mix thoroughly. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).
2. **Wash buffer.** Prepare the wash buffer solution by pouring all the wash buffer concentrate into an empty 1 L bottle. Rinse the wash buffer concentrate container with distilled or deionized water and pour into the 1 L bottle to ensure all the concentrate is used. Fill the 1 L bottle with additional distilled or deionized water, and swirl to assure thorough mixing. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).  
**NOTE:** Discard unused portions of extraction solution and wash buffer when the test kit has been used completely.
3. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.

## SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques (see Neogen's Food Allergen Handbook). The sample should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

1. Prepare the extraction solution as described in the procedural notes.
2. Preheat the extraction solution to 60°C (140°F).
3. Obtain a representative sample. If the sample is of a larger particle size, grind it to a very fine particle size.

4. Transfer 5 g of sample, or 5 mL of liquid sample, to a 250 mL plastic bottle or blender jar.
5. Add one level scoop of soy extraction additive to the plastic bottle or blender jar. Do not use the extraction additive from another allergen test kit.
6. Pour 125 mL of the 60°C (140°F) extraction solution into the plastic bottle or blender jar and cap to prevent contents from splashing during extraction.
7. Extract by shaking (150 rpm) in a water bath at 60°C (140°F) for **15 minutes**. Remove bottle from bath.  
**ALTERNATIVE:** Blend at high speed for **2 minutes**.
8. Let the material settle for **5 minutes**.  
**NOTE:** In rare cases, filtering may be necessary to achieve a supernatant free of suspended material. Filter the extract by pouring at least 5 mL through a Whatman no. 4 filter and collecting the filtrate as a sample, or centrifuging at 14,000 rpm for **5 minutes**.
9. Using a disposable pipette, transfer approximately 1 mL of the supernatant (the top liquid portion of the extract) to a sample dropper bottle. Label and place a dropper tip on the bottle.
10. Begin the test procedure once the sample has cooled to room temperature (at least **15 minutes**).

### ENVIRONMENTAL SWABBING PREPARATION AND EXTRACTION

1. Prepare the extraction solution as described in the procedural notes.
2. Prepare a working solution by combining 125 mL of the extraction solution to a separate container, and add one level scoop of extraction additive to a container. Prepare a fresh working solution daily.  
**NOTE:** Prepare 125 mL of working solution for every 25 environmental swabs to be tested.
3. Preheat the working solution to 60°C (140°F). Shake container to dissolve additive.
4. Gather the sample with a swab, using one of the following methods:  
**For dry surfaces:** Open a new swab and wet with extraction solution. Swab a 10 x 10 cm area by using a crosshatch technique. (Do not use working solution to moisten swabs.)  
**For wet surfaces:** Open a new swab and swab a 10 x 10 cm area by using a crosshatch technique. Do not moisten swab prior to use.
5. Return the swab to its original tube once sampling is complete. Remember to label each tube.
6. Remove the swab from its tube, and add 5 mL of working solution at 60°C to the tube. Mix by placing the swab back into the tube and shaking for **2 minutes** by hand (inverting tube), or for **30 seconds** with a vortex mixer.
7. Remove the swab from its tube.
8. Place a new sample dropper tip onto the tube. The solution in the tube now serves as the sample.  
**NOTE:** Use caution when inverting the tube, as some liquid may drip.

### TEST PROCEDURE

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

1. Remove 1 well for each sample to be tested plus 1 well for the control, and place into the well holder.
2. Mix each reagent by swirling its dropper bottle prior to use.
3. Add 3 drops from the yellow-labeled control dropper bottle to the first well. Add three drops from each sample extract to a respective well as indicated in the template below. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.

Control	S1	S2	S3	S4	S5
---------	----	----	----	----	----

4. Incubate microwells **10 minutes** at room temperature 18–30°C (64–86°F).
5. Shake out the contents of the wells. Using a wash bottle filled with wash buffer, fill each well and shake out. Repeat 10 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.

6. Add 3 drops from the blue-labeled conjugate dropper bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.
7. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
8. Shake out the contents of the wells. Using a wash bottle filled with wash buffer, fill each well and shake out. Repeat 10 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
9. Add 3 drops from the green-labeled substrate dropper bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.
10. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
11. Add 3 drops from the red-labeled Red Stop dropper bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface. The results are now ready to be interpreted.

### INTERPRETATION OF RESULTS

Visually compare the color of a sample well to the color of the control well. If the sample well has **more blue color** than the control well, the sample tests positive for soy contamination of **more than 2.5 ppm**. If the sample well has **less blue color**, than the control well, **or more red color**, the sample contains **less than 2.5 ppm** of soy contamination.

**ALTERNATIVE:** Read wells (wipe bottom of wells first) in a microwell reader with a 650 nm filter. If the sample well has an optical density (OD) higher than the control well, the sample is positive for soy contamination of more than 2.5 ppm. If the sample well has an OD lower than the control well, the sample contains less than 2.5 ppm of soy contamination.

### PERFORMANCE CHARACTERISTICS

Allergen detection: This test detects soy proteins and the results are expressed as ppm of soy flour. 2.5 ppm of soy flour is equivalent to 10 ppm soy protein isolate (SPI).

**NOTE:** Beginning November 2012, the kit control of 10 ppm SPI was revised to be expressed as 2.5 ppm of soy flour. Reagents and test procedures are unchanged. Only the expression of the control value changed. Prior results of greater than or less than 10 ppm SPI will now be greater than or less than 2.5 ppm soy flour.

### CUSTOMER SERVICE

Neogen Customer Assistance and Technical Service can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

### SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's Food Safety test kits, at [www.neogen.com](http://www.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

### TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit [www.neogen.com/Corporate/termsconditions.html](http://www.neogen.com/Corporate/termsconditions.html).

### WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## TEST KITS AVAILABLE FROM NEOGEN

### Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

### Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts

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