*Read instructions carefully before starting test*

# **Reveal<sup>®</sup>** **for** **Histamine**

*STORE AT ROOM TEMPERATURE (18–30°C)*

## **HISTAMINE**

High levels of histamine may develop in a variety of fish species as they decompose. These species include tuna, mahi-mahi, marlin, bluefish, sardines, anchovy, bonito, herring and mackerel. Ingestion of histamine may cause scombroid poisoning in humans, which may lead to a variety of symptoms, including rash, nausea, vomiting, diarrhea, hypotension, heart palpitations and muscle weakness. Paralysis and death have also been reported in cases of scombroid poisoning.

Because of its potential for human illness, the U.S. Food and Drug Administration (FDA) has ruled that extensive refrigeration records and/or histamine testing must be included in a hazard analysis and critical control point (HACCP) program for relevant fish species. The FDA has set an action level of 50 parts per million (ppm) for histamine in domestic and imported fish.

## **INTENDED USE**

Reveal<sup>®</sup> for Histamine is intended for the qualitative analysis of histamine in scombroid species of fish, such as tuna and mahi-mahi.

## **INTENDED USER**

The test kit is designed for use by quality control personnel and others familiar with histamine analysis in fish. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed Neogen training.

## **ASSAY PRINCIPLES**

Reveal for Histamine is a single-step lateral flow assay based on a competitive immunoassay format. The extract is wicked through a reagent zone, which contains histamine-specific antibodies conjugated to colored particles. If histamine is present, it will be captured by the particle-antibody complex. The histamine-antibody-particle complex then is wicked onto a membrane, which contains a zone of histamine conjugated to a protein carrier. This zone captures only uncomplexed antibody-particle conjugate, allowing the colored particles to concentrate and form a visible line. The membrane also contains a control zone where the antibody-particle complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of histamine, ensuring the strip is functioning properly. In addition, the control line also is used as a reference for comparison to the test zone for visual interpretation of the test results. If histamine is present above a pre-determined limit, a line lighter than the control line or no line will form. Therefore, a positive sample will result in a light line or absence of a line in the test zone. Alternatively, a negative sample will form a line of the same or darker intensity as the control line in the test zone.

## **STORAGE REQUIREMENTS**

Store kit components at room temperature (18–30°C, 64–86°F) to ensure full shelf life. Test strips should remain capped in their original sample tubes until used to ensure optimal performance.

## **MATERIALS PROVIDED**

1. 25 histamine test strips
2. 25 sample cups
3. 25 bottles of 7 mL sample diluent

## **MATERIALS RECOMMENDED BUT NOT PROVIDED**

1. Graduated cylinder (Neogen item 9368)
2. Sample extraction bottles, 250 mL (Neogen item 9399)
3. Sample collection tubes (Neogen item 9421)
4. Blender (Neogen item 9493, 9477)
5. Scale capable of weighing 10–50 grams (Neogen item 9427)
6. Pipettor, 100 µL (Neogen items 9272, 9278, 9276)
7. Pipette tips (Neogen items 9410, 9407, 9417)
8. Paper towels or equivalent absorbent material
9. Rack for Reveal sample cups (Neogen item 9475)
10. Timer (Neogen item 9426)
11. Waterproof marker
12. Distilled or deionized water
13. AccuScan® reader system (Neogen item 9565) (Optional)

## **PRECAUTIONS**

1. The test strips must remain inside the stay-dry tube before use.
2. Store test kit at room temperature (18–30°C, 64–86°F) when not in use.
3. Do not use kit components beyond expiration date.
4. Treat all used liquids, including sample extract, and labware as if contaminated with histamine. Gloves and other protective apparel should be worn at all times.
5. To avoid cross-contamination, use clean equipment for each sample, and thoroughly wash all equipment between samples.
6. Glassware should not be used for extraction purposes. Histamine may adhere to glass, therefore using glassware may affect test results.

## **SAMPLE PREPARATION**

### **A. Canned tuna: AOAC 937.07b**

Place entire contents of can/pouch (meat and liquid) into a blender. Blend until homogenous. Store samples at 2–8°C (35–46°F) until analyzed.

### **B. Fresh or thawed frozen raw fish: AOAC 937.07a**

Clean and eviscerate three fish. Cut three cross-sectional pieces 2.5 cm (1 inch) thick, from back of the pectoral fin, halfway to the vent and one posterior to the vent. Debone slices and blend or grind combined samples until homogenous. Store samples at 2–8°C (35–46°F) until analyzed.

## **SAMPLE EXTRACTION**

1. Add 10 g of the homogenous mixture to a clean disposable extraction bottle containing 190 mL of distilled or deionized water.
2. Tightly cap and vigorously shake the bottle for **15–20 seconds** to suspend the fish tissue in the water.
3. Wait approximately **5 minutes**, and then shake the bottle for **15–20 seconds** to resuspend the fish tissue.
4. Wait an additional **5 minutes**, and again shake the bottle for **15–20 seconds** to resuspend the fish tissue. Allow the tissue to settle to the bottom of the bottle for about **30 seconds**. The sample is now ready for extract dilution.

## SAMPLE EXTRACT DILUTION

Because of the sensitivity of the assay, the sample extract must be diluted before performing the test.


1. Using a clean pipette tip, add 100  $\mu\text{L}$  of the extract to a new bottle of sample diluent provided with the kit. Gently swirl to mix. The sample is now ready to test.
2. Repeat for all samples, using a new sample diluent bottle for each sample.

## TEST PROCEDURE

1. Remove the appropriate number of sample cups and place in sample cup rack.
2. Using a pipettor and new tip, add 200  $\mu\text{L}$  of diluted sample to a sample cup. Repeat for each sample using a new tip and cup.
3. Place a new histamine test strip with the sample end down into the sample cup.
4. Allow the strip to develop in the sample cup for **5 minutes**.
5. Remove strip and read results visually against a white background or by using one of Neogen's AccuScan readers.

## TEST INTERPRETATION

1. If no control line appears, the test is not considered valid.
2. If the test line has an intensity the same as or darker than the control line, the sample is interpreted as negative (less than 50 ppm histamine).
3. If the test line has an intensity lighter than the control line, the sample is interpreted as positive (greater than 50 ppm).
4. If the test line is absent, the sample is interpreted as highly positive.

**Notes:**

1. Observations should be made immediately following the 5 minute development. Observations after 5 minutes may be inaccurate due to over-development of the device.
2. The strips can be read using Neogen's AccuScan readers\*. This eliminates subjectivity by interpreting and storing sample results.
3. Positives can be confirmed by Neogen's quantitative Veratox® for Histamine\* kit or HPLC.
4. Methods are available for alternative detection limits. Please contact Neogen Technical Services for further instructions.

*\*Contact Neogen or visit [www.neogen.com](http://www.neogen.com) for more information on these products.*

**VALIDATED MATRICES**

Tuna and mahi-mahi.

**CUSTOMER SERVICE**

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

**MSDS INFORMATION AVAILABLE**

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at [www.neogen.com](http://www.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

**TERMS AND CONDITIONS**

For Neogen's full terms and conditions, please visit [www.neogen.com/Corporate/termsconditions.html](http://www.neogen.com/Corporate/termsconditions.html).

**WARRANTY**

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## TEST KITS AVAILABLE FROM NEOGEN

### Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

### Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts

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