Read instructions carefully before starting test

Veratox[®] for Chloramphenicol

Quantitative Test

REFRIGERATE AT 2–8°C (35–46°F) • DO NOT FREEZE

CHLORAMPHENICOL

Chloramphenicol was among the first broad-spectrum antibiotics. It was isolated from *Streptomyces venezuelae* in 1947, and its chemical structure was elucidated shortly thereafter. The relatively simple structure allows the compound to be produced cheaply by completely synthetic methods. It is bacteriostatic against a range of gram-positive and gram-negative bacteria. However, serious side-effects related to suppression of bone marrow activity are observed in patients undergoing treatment with chloramphenicol. Some of the side effects are dose related and reversible upon ending treatment. The most serious side-effect, potentially fatal aplastic anemia, may develop after ending treatment and is unrelated to dose.

Because of the side-effects, chloramphenicol is generally not prescribed in developed countries where alternative antibiotics are readily available. Chloramphenicol is banned for use in production of food animals in the United States, Canada, and the European Union. However, it continues to be used in some countries due to its broad-specificity, availability, and low cost.

INTENDED USE

This test is intended for quantitative determination of chloramphenicol in shrimp or honey. Positive samples should be confirmed using an instrumental method such as mass spectrometry.

INTENDED USER

This kit is intended for use by laboratory personnel familiar with the analysis of chloramphenicol in shrimp or honey. The method requires extraction of the shrimp tissue or honey with an organic solvent. The user should be experienced in the use of immunosorbent assays (ELISA). Good technique is important for obtaining accurate results.

ASSAY PRINCIPLES

Veratox for Chloramphenicol is an enzyme-linked immunosorbent assay (ELISA). The basis of this test is a competition between a chloramphenicol-enzyme conjugate and any free chloramphenicol from the sample for a limited number of binding sites on the antibody-coated microplate wells.

Samples are extracted, and the extracts (and solutions of chloramphenicol standards) are added to the microplate wells. A solution of the chloramphenicol-enzyme conjugate is also added to each of the wells. At this point, time is allowed for any free chloramphenicol and the chloramphenicol-enzyme conjugate to compete for the binding sites on the plate. Next, the plate is washed to remove unbound material and a substrate solution is added to each well. During an ensuing reaction period, any bound enzyme conjugate reacts with the substrate to produce a blue-colored solution. The color-forming reaction is terminated by adding Stop Solution to the wells of the plate, turning the blue color to yellow. The absorbance of each well is then measured using a microplate reader set at 450 nm.

The extent of color formation is inversely proportional to the amount of chloramphenicol present in that well. Samples that contain little or no chloramphenicol will have a relatively large absorbance while a sample that contains more chloramphenicol will have a smaller absorbance. Quantitative results may be obtained by comparing the absorbance of the samples to a standard curve generated from the wells containing the chloramphenicol standards.

STORAGE REQUIREMENTS

This kit can be used until the expiration date on the label when stored at 2–8°C (35–46°F). Do not freeze.

MATERIALS PROVIDED

1. 1 plate (96 wells) of antibody-coated microwells
2. 1 vial of chloramphenicol 0.0 ng/mL standard, 2 mL
3. 1 vial of chloramphenicol 0.02 ng/mL standard, 2 mL
4. 1 vial of chloramphenicol 0.05 ng/mL standard, 2 mL
5. 1 vial of chloramphenicol 0.1 ng/mL standard, 2 mL
6. 1 vial of chloramphenicol 0.5 ng/mL standard, 2 mL
7. 1 vial of chloramphenicol 2 ng/mL standard, 2 mL
8. 1 bottle of ready-to-use chloramphenicol enzyme conjugate, 20 mL
9. 1 bottle of sample dilution buffer, 100 mL
10. 1 bottle of wash buffer concentrate, 100 mL (dilute 1:10 with deionized water before use)
11. 1 bottle of substrate solution, 28 mL
12. 1 bottle of stop solution, 28 mL

MATERIALS REQUIRED BUT NOT PROVIDED

1. Pipettor, 100 µL (Neogen item #9272)
2. Pipettor, multichannel, 100 µL
3. Pipette tips for 100 µL and multichannel pipettors (Neogen item #9410)
4. 3 reagent boats for multichannel pipettor
5. Paper towels or equivalent absorbent material
6. Pipettor capable of measuring 1.0 mL
7. Pipettor capable of measuring 0.5 mL (for isooctane/chloroform)
8. Graduated or volumetric pipettor capable of measuring 6.0 mL of solvent
9. Disposable pipettes capable of measuring 4.0 mL of solvent
10. Glass test tubes capable of being centrifuged at 2,000 g (16 x 100 mm preferred)
11. Caps for glass test tubes
12. Ethyl Acetate (ACS reagent grade or better)
13. Chloroform
14. Isooctane (2,2,4-trimethylpentane)
15. Deionized or distilled water
16. Blender or food processor
17. Centrifuge capable of 2,000 g, swinging bucket rotor preferred
18. Heating block
19. Compressed air or nitrogen and necessary tubing
20. Water bath
21. Vortex

22. Automated plate washer
23. Microplate reader with 450 nm filter (650 nm reference filter is optional)
24. Balance capable of measuring 3.0 g
25. Timer (Neogen item #9426)
26. Graduated cylinder to dilute wash buffer
27. Waterproof marker
28. Pipette aid for dispensing solvents
29. Stainless steel spatulas for transferring shrimp into test tubes
30. Thermometer

PRECAUTIONS

1. Do not use components beyond expiration date.
2. Store test kit at 2–8°C (35–46°F) when not in use. Do not freeze test kits.
3. Avoid prolonged storage of kits at ambient temperatures.
4. Prior to performing the test, allow ELISA kit components, standards, and samples to reach room temperature (20–23°C, 68–74°F). This usually takes about 60 minutes.
5. Do not mix reagents from one kit lot with reagents from a different kit lot.
6. Use aseptic technique when opening and removing reagents from vials and bottles.
7. Always pour substrate out of the bottle into a clean reagent boat. Do not pipette out of the bottle. Unclean tips can contaminate the substrate. Do not put excess substrate back into the bottle.
8. Do not pipette reagents by mouth. Always use clean pipette tips for the buffer, enzyme conjugate, standards, and samples.
9. All specimens should be considered potentially hazardous. Exercise proper handling precautions.
10. The plate wells are for single use only. Each well can accommodate only one sample or standard. Do not reuse wells.
11. Use of incubation times other than those specified may give inaccurate results.
12. Most of the reagents contain anti-microbial preservatives. Avoid skin contact. If contact occurs, flush immediately with water.
13. Avoid contacting skin or eyes with stop solution. In case of contact, rinse with large amounts of water.
14. Do not smoke, eat, or drink in areas where specimens or reagents are being handled.
15. Dispose of reagents in accordance with all applicable governmental regulations.

PROCEDURAL NOTES

1. Allow the strips to warm to room temperature (20–23°C, 68–74°F) before removing them from the foil pouch.
2. Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after the samples are extracted and the test procedure is set to begin.
3. Desiccant bag must remain sealed in the foil pouch with unused strips. Remove excess air before sealing.
4. Dispense only the amounts of enzyme conjugate and substrate needed for immediate testing into clean vessels. From these vessels, pipette to the wells of the plate. **Do not put excess reagents back into original bottles.** Contamination of these reagents may cause enzyme conjugate degradation or premature development of blue color in the substrate.
5. Before opening the standards vials, tap vials in an upright position to remove any liquid in the cap.
6. The standards provided with the kit are ready to use. Do not dilute the standards in the same manner as sample extracts.
7. Before pipetting a reagent, prime the pipette tip once with that reagent (i.e. fill the tip with the desired amount of reagent and dispense back into the same vial and repeat).
8. When pipetting into the wells, do not allow the pipette tip to touch the inside of the well or any of the reagents already in the well.
9. Timing and speed of reagent addition can influence performance of the test. This is especially true for the addition of substrate and stop solution. Substrate and stop solution should be applied to the plate wells as quickly as possible, preferably with a multichannel pipettor.
10. A good wash step is critical to achieving good well-to-well reproducibility. All liquid must be removed from the wells before applying substrate. After washing, the plate should be tamped on absorbent paper to remove any remaining residue.

EXTRACTION OF SHRIMP TISSUE

1. Remove the heads and shells from the shrimp to be tested. Grind at least 5 shrimp (≥ 50 g) in a blender or food processor for approximately 1 minute or until the material appears to be a uniform paste.
2. Weigh 3.0 g of the blended tissue into a glass tube.
Optional: Centrifuge the tube briefly (approximately 200 g for 2 minutes to transfer the tissue into the bottom of the tube).
3. Add 6 mL of ethyl acetate to the tube. Using a stainless steel spatula, scrape the tissue off the walls of the tube and break up any large aggregates. All of the tissue should be suspended in the solvent.
4. Cap the tube, and vortex the tube vigorously for 30 seconds.
5. Centrifuge the tube for 10 minutes at approximately 700 g.
6. Remove 4 mL of the ethyl acetate layer, and transfer it to a clean glass tube.
7. Evaporate the ethyl acetate to just dryness at 80°C using compressed air or nitrogen and a heating block.
8. Add 0.5 mL of 2:3 isooctane/chloroform (v/v), and vortex briefly to ensure all of the residue is lifted off the walls of the tube.
9. Add 1.0 mL sample dilution buffer and vortex for 10 seconds. A longer vortex time promotes the formation of emulsions.
10. Centrifuge the tube for 10 minutes at 2,000 g.
NOTE: After step 10, if there is no distinct separation of layers or the upper layer remains an emulsion, heat the tube briefly (≤ 5 min) in a water bath at 80°C. (Loosely cap the tube to avoid evaporating the contents.) Then centrifuge again for 10 min at 2000 g and proceed with step 11.
11. The upper layer is the aqueous solution. Carefully pipette 100 μ L portions of this layer into the wells of the plate for the test. It is important that none of the solvent or emulsion is transferred during this process.

EXTRACTION OF HONEY

1. Weigh 3.0 g of honey and place into a tube.
2. Add 3.0 mL of water. Cap the tube and vortex for 30 seconds, or until the honey is dissolved to form a uniform solution.
3. Add 6.0 mL ethyl acetate to the tube.
4. Place the tube on a rocker or similar mixer and mix gently for 10 minutes.
5. Centrifuge the tube for 5 minutes at approximately 700 g.
6. Remove 2.0 mL of the ethyl acetate phase (upper layer) and transfer it to a clean glass tube.
7. Evaporate the ethyl acetate until just before dry at 80°C using compressed air or nitrogen and a heating block.
8. Add 1.0 mL sample dilution buffer and vortex for 30 seconds. The residue will not completely dissolve, but it must be suspended and dispersed in the buffer. If necessary, allow the residue to soak for a few minutes and vortex until the residue is broken up and suspended. If the honey is crude or may contain pieces of the comb, see note below.
9. Carefully pipette 100 μ L of the solution into the wells of the plate for the test.
NOTE: Raw honey may contain a significant amount of fatty compounds from the comb. In this case, after evaporating the ethyl acetate (Step 7) add 0.5 mL of 2:3 isooctane/chloroform (v/v) to the residue and vortex for 10 seconds. Then add 1.0 mL sample dilution buffer and vortex for 10 seconds. Centrifuge the tube for 10 minutes at 2000 g. Carefully pipette 100 μ L of the aqueous (upper) layer into the wells of the plate for the test. This treatment is only necessary for fatty samples, but can be done for any sample.

TEST PROCEDURE

1. Dilute wash buffer concentrate ten-fold with deionized water (i.e. 20 mL of wash buffer concentrate plus 180 mL of deionized water). Mix thoroughly.
2. Determine the number of wells to be used. Remove unneeded strips from the plate holder, and return them to the foil pouch.
3. Add 100 μ L of sample or standards to the appropriate wells in duplicate.
NOTE: Do not dilute the provided standards. **Wait 10 minutes before proceeding to Step 4.**
4. Add 100 μ L of diluted enzyme conjugate to each well. Use an multi-channel pipettor for rapid addition.
5. Incubate at room temperature for 10 minutes.
NOTE: Keep plate away from drafts and temperature fluctuations.

6. After the incubation is complete, dump the liquid from the wells. Tamp the plate on a clean lint-free towel to remove any liquid remaining in the wells.
7. Wash each well with 300 μ L of diluted wash buffer. For manual wash procedures repeat for a total of 3 washings, invert and tap dry between each washing. If an automated plate washer is used, wash the plate 3 times with 300 μ L of diluted wash buffer. After completing the last wash step, tamp the plate on a clean lint-free towel to remove any liquid remaining in the wells, and wipe the bottom of the wells with a lint-free towel to remove any liquid on the outside of the wells.
8. Add 100 μ L of the substrate solution to each well. Use a multi-channel pipettor for best results.
9. Incubate at room temperature for 10 minutes. Gently shake plate periodically.
10. Add 100 μ L of stop solution to each well to stop the reaction. Use a multi-channel pipettor for best results. Mix by shaking plate gently.
11. Measure the absorbance of each well at 450 nm using a microplate reader.

INTERPRETATION OF RESULTS

Calculate the average of the duplicate absorbance values for each standard and sample. Express the mean absorbance (B) of each sample and standard as a fraction relative to the mean absorbance of the 0.0 ng/mL standard (B_0). Do this by dividing each value by the absorbance of the 0.0 ng/mL standard (B/B_0). The 0.0 ng/mL standard will become 1.0.

Prepare a calibration curve by plotting $\log(B/B_0)$ versus $\log(\text{chloramphenicol concentration})$ for each of the standards. Using the value of $\log(B/B_0)$ for each sample, the concentration of chloramphenicol in each sample as ng/mL can be read or calculated from the curve. To account for the extraction procedure, the result must be multiplied by 500 to express the result as pg/g tissue (ppt).

SPECIFICITY AND SENSITIVITY

The specificity was determined by testing the cross-reactivity of related compounds.

Compound	Cross reactivity
Chloramphenicol	100%
Chloramphenicol succinate	130%
Chloramphenicol glucuronide	26%
Chloramphenicol base	< 0.01%
Florfenicol	< 0.01%
Thiamphenicol	< 0.01%

The lower limit of detection is approximately 0.02 ng/mL. This is the lowest concentration that can be determined to be significantly different from the zero standard under ideal conditions. Allowing for the extraction, this is equivalent to 10 ppt (pg/g) in shrimp tissue. Care should be taken in selecting an appropriate cutoff value.

Typical recovery rates using the recommended extraction procedure are $\geq 80\%$.

CUSTOMER SERVICE/TECHNICAL SUPPORT

Neogen Customer Assistance and Technical Service can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's Food Safety test kits, at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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