Read instructions carefully before starting test

ALERT[®]

for Gliadin R5

Refrigerate at 2–8°C (35–46°F) • Do not freeze

GLIADIN

Gliadin is an alcohol-soluble protein found in wheat that belongs to a group of proteins called prolamins. Other prolamins include secalin, found in rye, and hordein, found in barley. Gluten consists of two groups of proteins (prolamins and glutelins) that are found in differing amounts in wheat, barley, rye and oats.

Gliadin and other prolamins have been identified as major causal agents in a number of disorders, including wheat allergy and gluten intolerance (celiac disease). Wheat allergy is a specific immune response to a number of wheat proteins, including gliadin, albumin, globulin, and glutenin. Celiac disease is a chronic reaction to gluten proteins that results in the poor absorption of nutrients in the small intestine.

Those who must avoid gluten rely upon the correct labeling of food to make appropriate, safe food choices. Testing for the presence of gluten components ensures food manufacturers that an unlabeled—and potentially dangerous—ingredient did not make its way into a food product.

INTENDED USE

Alert for Gliadin R5 is intended for the qualitative analysis of ingredients, clean-in-place solutions, finished food products, and environmental surfaces intended to be gluten-free for the presence of gliadin and prolamins found in wheat, rye and barley.

INTENDED USER

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by gluten. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed the Neogen training.

ASSAY PRINCIPLES

Alert for Gliadin R5 is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Gliadin is extracted from samples with a 60% ethanol solution by shaking in a shaker or rotator. Extract is diluted in Phosphate Buffered Saline (PBS) and diluted samples are added to antibody-coated wells (capture antibody) where gliadin will bind to the antibody during an incubation period. Any unbound gliadin is washed away and a second antibody, which is R5 enzyme labeled (detector antibody) is added. The detector antibody binds to the gliadin during another incubation period. Unbound enzyme-labeled antibody is washed away and a one step substrate is added. Color develops as a result of the presence of bound-labeled antibody. A stopping reagent is added and the color of the solution is observed. Blue color indicates samples containing high levels of gliadin while purple or red samples contain little or no gliadin.

ALERT[®] for Gliadin R5

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F).

MATERIALS PROVIDED

1. 24 antibody-coated microwells
2. 1 yellow-labeled dropper bottle of 10 ppm gliadin control (20 ng/mL gliadin)
3. 1 blue-labeled dropper bottle of conjugate
4. 1 green-labeled dropper bottle of K-Blue[®] Substrate
5. 1 red-labeled dropper bottle of Red Stop solution
6. 20 prefilled dilution dropper bottles of 4.9 mL PBS
7. 1 foil pouch containing enough 10 mM PBS dry powder to prepare 1 L of dilution solvent
8. 1 widemouth bottle containing enough 10 mM PBS-Tween wash buffer concentrate to prepare 1 L in distilled or deionized water (pH 7.4)
9. 30 g of extraction additive in a specimen cup
10. Plastic 1 g scoop to measure extraction additive

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Orbital rotator or shaker to hold 50 cc centrifuge tubes for a 1 g sample or shaker or shaker water bath with clamps adjusted to hold 125 mL (4 oz) extraction bottles for a 2 g sample
2. Allergen Environmental Swabbing Kit (Neogen item 8432S)
3. Scale capable of weighing 1 ± 0.1 g (Neogen item 9427), or scale capable of weighing 0.25 ± 0.01 g if using renaturing cocktail solution
4. Gliadin renaturing cocktail solution for analysis of heat processed samples (Neogen item 8515, 8515S, or 8515B)
5. 2 1 L bottles to prepare washing solution, extract solution (Neogen item 9472)
6. Timer (Neogen item 9426)
7. Microwell holder (Neogen item 9402)
8. Pipettor, 100 μ L (Neogen item 9272, 9278 or 9276)
9. Pipette tips (Neogen item 9410)
10. Wash bottle (Neogen item 9400)
11. Paper towels or equivalent absorbent material
12. Waterproof marker
13. Distilled or deionized water
14. Laboratory grade ethanol (190 proof)
15. Centrifuge (optional)

PRECAUTIONS

1. Ethanol solution is highly flammable. Keep container tightly closed, and keep away from heat, sparks, open flame and those smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
2. Components of Alert for Gliadin R5, such as controls and extraction additive, may contain one or more of the following potentially allergic materials: gluten, casein, almond protein and soy protein. If allergic to any of these compounds, use caution when using this product.
3. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze test kits, and avoid prolonged storage at ambient temperatures.
4. Bring kits to room temperature (18–30°C, 64–86°F) prior to use.
5. Do not use kit components beyond expiration date.
6. Do not mix reagents from one kit with reagents from a kit with a different serial number.
7. Do not run more than 6 wells per test.
8. Use only incubation times specified; others may give inaccurate results.

PROCEDURAL NOTES

1. **Sample extract dilution solution (PBS).** Prepare extract dilution solution by adding a foil pouch of dilution solvent, 10 mM PBS, to 1 L distilled or deionized water. Swirl to mix thoroughly.
2. **Wash buffer.** Prepare the wash buffer solution by pouring the wash buffer concentrate into a 1 L container. Add 960 mL of distilled or deionized water. Swirl to assure thorough mixing.
NOTE: Discard unused portions of extraction solution and wash buffer when the test kit has been used completely.
3. **Extraction additive.** Use extraction additive with all the samples utilizing procedure A or B (non-heat processed samples). Heat processed samples (procedure C) should **only** use the extraction additive if samples contain buckwheat, chestnut flour or tannins/phenolic compounds such as chocolate, coffee, cocoa, wine, herbs or fruits.
4. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.

SAMPLE PREPARATION AND EXTRACTION

For analyzing **heat processed samples**, follow extraction **procedure C**. For all commodities that were **not heat processed**, follow either extraction procedure **A or B**. Samples of an unknown origin should be extracted using extraction procedure **C**.

A. Extraction of non-heat processed samples with orbital shaker or rotator

1. Prepare 60% ethanol extraction solution by combining 6 parts ethanol with 4 parts distilled water.
2. Add 1 g ground sample, or 1 mL liquid sample, to a clean 50 cc tube.
3. Add 1 level scoop of extraction additive to the tube.
4. Add 10 mL (9 mL for liquid samples) of 60% ethanol to the tube, cap tightly, then shake the tube vigorously by hand for about **1 minute**, or vortex for **30 seconds**, to ensure complete mixing.
5. Extract by shaking (150 rpm) in an orbital shaker or rotator by laying down the tube on its side over the flat pad of the instrument, and holding it tightly using a rubber band or tape. Rotate or shake for **10 minutes** at room temperature.
6. Remove the tube and let it stand in a rack for about **10 minutes** to enable the sample extract to settle before withdrawing the clear extract. If necessary, centrifuge sample for **10 minutes** at $\geq 2,500$ g at room temperature.
7. Dilute each sample 1:50 by withdrawing 100 μ L of the upper layer of the extract and transferring it to a prefilled dilution dropper bottle containing 4.9 mL of PBS. (Cap and tip must be removed first, then reinserted).
8. To mix, invert several times by hand or vortex for **5 seconds**.
9. Test diluted samples within **2–3 hours** of extraction.

B. Extraction of non-heat processed samples with shaker or shaker water bath

1. Prepare 60% ethanol extraction solution by combining 6 parts ethanol with 4 parts distilled water.
2. Add 2 g ground sample, or 2 mL liquid sample, to a 125 mL clean extraction bottle.
3. Add 1 level scoop of extraction additive to the bottle.
4. Add 20 mL (18 mL for liquid samples) of 60% ethanol, cap the bottle tightly, then shake vigorously by hand for about **20 seconds** to ensure complete mixing.
5. Extract by shaking (150 rpm) in a shaker for **10 minutes** at room temperature (a shaker water bath can work, but do not turn the heat on). Remove the bottle from shaker or bath.
6. Let the bottle stand for about **10 minutes** to enable some of the sample to settle before withdrawing the clear extract. If necessary, centrifuge sample for **10 minutes** at $\geq 2,500$ g at room temperature.
7. Dilute each sample 1:50 by withdrawing 100 μ L of the upper layer of the extract and transferring it to a prefilled dilution dropper bottle containing 4.9 mL of PBS. (Cap and tip must be removed first, then reinserted.)
8. To mix, invert several times by hand or vortex for **5 seconds**.
9. Test diluted samples within **2–3 hours** of extraction.

C. Extraction of heat processed commodities

Heat processed commodities require the gliadin renaturing cocktail solution (Neogen item 8515) that renatures the heated sample and allows the proper determination of the gliadin in the sample. Prefilled dilution bottles **cannot** be used for this extraction procedure.

1. Prepare 80% ethanol extraction solution by combining 8 parts ethanol with 2 parts distilled water.
2. Prepare sample extract dilution solution (PBS) as detailed in procedural note 1.
3. Weigh out 0.25 g sample to 50 cc screw cap centrifuge tube.
4. Add 2.5 mL of renaturing cocktail solution.
5. If samples contain buckwheat, chestnut flour or tannins/phenolic compounds such as chocolate, coffee, cocoa, wine, herbs or fruits, add 1 level scoop of extraction additive. For all other commodity types, do not add extraction additive.
6. Cap and vortex **30 seconds** to homogenize cocktail and sample.
7. Incubate **40 minutes** at 50°C (water bath or oven).
8. Remove samples and let cool for **5–10 minutes**.
9. Add 7.5 mL of 80% ethanol and vortex again for **10–20 seconds**.
10. Shake (150–200 rpm) for **1 hour** at room temperature on a rotator (tube on its side).
11. Centrifuge sample (if necessary) for **10 minutes** at $\geq 2,500$ g at room temperature.
12. Dilute the sample 1:12.5 into 10 mM PBS, pH 7.4 (200 μ L sample into 2.3 mL PBS)
13. Samples are ready to run.

TEST PROCEDURE

Allow the test kit and all reagents to warm to room temperature (18–30°C, 64–86°F) before using.

1. Remove 1 well for each sample to be tested plus 1 well for the control, and place into the well holder.
2. Mix each reagent by swirling its dropper bottle prior to use.
3. Add 3 drops from the yellow-labeled control dropper bottle to the first well. Add 3 drops from each diluted sample dropper bottle to a respective well as indicated in the template below. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.

Control	S1	S2	S3	S4	S5
---------	----	----	----	----	----

4. Incubate microwells **10 minutes** at room temperature (18–30°C, 64–86°F).
5. Shake out the contents of the wells. Using a wash bottle filled with wash buffer, fill each well and shake out. Repeat 5 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
6. Add 3 drops from the blue-labeled conjugate dropper bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.
7. Incubate for **10 minutes** at room temperature (18–30°C, 64–86°F).
8. Shake out the contents of the wells. Using a wash bottle filled with wash buffer, fill each well and shake out. Repeat 5 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
9. Add 3 drops from the green-labeled substrate dropper bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.
10. Incubate for **10 minutes** at room temperature (18–30°C, 64–86°F).
11. Add 3 drops from the red-labeled Red Stop dropper bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface. The results are now ready to be interpreted.

INTERPRETATION OF RESULTS

Visually compare the color of a sample well to the color of the control well. If the sample well has **more blue color** than the control well, the sample tests positive for gliadin contamination of **more than 10 ppm**. If the sample well has **less blue color or more red color** than the control well, the sample contains **less than 10 ppm** of gliadin contamination. **NOTE:** Standard controls were made from wheat gliadin and calculated as gliadin. Approximately 50% of the gluten is available as gliadin. Therefore, 10 ppm gliadin is equal to 20 ppm gluten.

ALTERNATIVE: Read wells (wipe bottom of wells first) in a microwell reader with a 650 nm filter. If the sample well has an optical density (OD) higher than the control well, the sample is positive for gliadin contamination of more than 10 ppm. If the sample well has an OD lower than the control well, the sample contains less than 10 ppm of gliadin contamination.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN**Natural toxins**

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts

Genetic modification

- CP4 (Roundup Ready®)

Ruminant by-products

- Meat and bone meal, feed

Species identification

- Raw and cooked meat samples

**North America****Neogen Headquarters**

620 Leshler Place, Lansing, MI 48912 USA
800/234-5333 (USA/Canada) or 517/372-9200
Fax: 517/372-2006 • foodsafety@neogen.com
www.neogen.com

Europe, Middle East and Africa**Neogen Europe**

The Dairy School, Auchincruive, Ayr
KA6 5HW Scotland, UK
+ 44 (0) 1292 525 600
Fax: + 44 (0) 1292 525 601
info_uk@neogeneurope.com
www.neogeneurope.com

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