Read instructions carefully before starting test

Reveal[®] Q+

for Fumonisin

FGIS 2018-116

THE TOXIN

Discovered in 1989, fumonisins are a family of mycotoxins produced by different species of the mold *Fusarium*. These molds commonly infect corn (in fact, they are considered ubiquitous in corn) and rice, hence the potential for fumonisins to be found in feed and foodstuffs is high. Fumonisin affect various animals differently and have been linked to esophageal cancer in humans. The Environmental Protection Agency classifies fumonisins as category II-B carcinogens.

Horses are extremely sensitive to low amounts of fumonisin, which can cause leukoencephalomalacia (liquefaction of the brain). In swine, research has shown fumonisin attacks the cardiopulmonary system causing pulmonary edema, as well as liver and pancreatic lesions.

The FDA has issued advisory levels for fumonisins as follows:

For human consumption (total fumonisins)

Commodity	Level
Degermed dry milled corn products	2 ppm
Whole/partially dry milled corn products, dry milled corn bran, cleaned corn for mass production	4 ppm
Cleaned corn for popcorn	3 ppm

Animal feeds (corn/corn by-products, total fumonisins)

For	Level
Equids and rabbits	5 ppm, ≤ 20% of diet
Swine and catfish	20 ppm, ≤ 50% of diet
Breeding ruminants, breeding poultry and breeding mink	30 ppm, ≤ 50% of diet
Ruminants ≥ 3 months old being raised for slaughter and mink being raised for pelt production	60 ppm, ≤ 50% of diet
Poultry being raised for slaughter	100 ppm, ≤ 50% of diet
All other species or classes of livestock and pet animals	10 ppm, ≤ 50% of diet

European Union regulations for fumonisin (Total B₁+B₂)

Foodstuffs	Maximum levels
Unprocessed corn	2 ppm
Corn flour, meal, grits, germ, and oil	1 ppm
Corn-based foods intended for direct human consumption unless otherwise listed	1 ppm
Corn-based baby foods intended for infants	0.2 ppm

INTENDED USE/USER

Reveal Q+ for Fumonisin is intended for the quantitative analysis of corn for fumonisin. The test kit is designed for use by quality control personnel and others familiar with commodities possibly contaminated by fumonisin.

ASSAY PRINCIPLES

Reveal Q+ for Fumonisin is a single-step lateral flow immunochromatographic assay based on a competitive immunoassay format. The extract is wicked through a reagent zone, which contains antibodies specific for fumonisin conjugated to colloidal gold particles. If fumonisin is present, it will be captured by the particle-antibody complex. The fumonisin-antibody-particle complex then is wicked onto a membrane, which contains a zone of fumonisin conjugated to a protein carrier. This zone captures any uncomplexed fumonisin antibody, allowing the particles to concentrate and form a visible line. As the level of fumonisin in a sample increases, free fumonisin will form a complex with the antibody-gold particles. This allows less antibody-gold to be captured in the test zone. Therefore, as the concentration of fumonisin in the sample increases, the test line density decreases. Algorithms programmed into the readers convert these line densities into a quantitative result displayed in parts per million (ppm). The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line always will form regardless of the presence of fumonisin, ensuring the strip is functioning properly.

STORAGE REQUIREMENTS

Store kit components at room temperature (18–30°C, 64–86°F) to ensure full shelf life. Test strips should remain capped in their original tubes until used to ensure optimal performance.

MATERIALS PROVIDED

1. 25 Reveal Q+ for Fumonisin test strips
2. 25 red sample dilution cups
3. 25 clear sample cups
4. 1 bottles of sample diluent
5. Instructions for use

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Extraction materials
 - a. 65% ethanol solution (NEOGEN item 8073, 8074)
 - b. Sample collection cups with lids – 125 mL (NEOGEN item 9428, 9428B)
 - c. Sample collection tubes with caps (NEOGEN item 9421, 9421B)
 - d. Filter syringes (NEOGEN item 9272, 9860)
- OR**
- e. Centrifuge, mini (NEOGEN item 9330)
- f. Microcentrifuge tubes (NEOGEN item 9172)
- g. Dispensing pump or graduated cylinder (NEOGEN item 9448, 9447)
2. Agri-Grind grinder or equivalent (NEOGEN item 9401, 9453)
3. Scale capable of weighing 5–50 g ± 0.1 g (NEOGEN item 9427)
4. Timer (NEOGEN item 9426, 9452)
5. Reveal sample cup rack (NEOGEN item 9475)
6. AccuScan Pro, AccuScan Gold or Raptor Integrated Analysis Platform (NEOGEN item 9680 or 9696)
7. Pipettor, 100 µL (NEOGEN item 9272, 9860)
8. Pipette tips, 100 µL (NEOGEN item 9407, 9410, 9417)
9. Pipettor, 400 µL fixed (NEOGEN item 9693)
10. Pipette tips, 200–1,000 µL (NEOGEN item 9464, 9487, 9293, 9292)
11. Raptor cartridges (NEOGEN item 9681)

PRECAUTIONS

1. The test strips must remain inside the stay-dry tube before use.
2. Ethanol is highly flammable. Keep container tightly closed and away from heat, sparks, open flame and those who are smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
3. Store test kit at room temperature (18–30°C, 64–86°F) when not in use. Do not freeze.
4. Do not use kit components beyond expiration date.
5. Treat all used liquids, including sample extract, and labware as if contaminated with fumonisin. Gloves and other protective apparel should be worn at all times.
6. To avoid cross-contamination, use clean glassware for each sample, and thoroughly wash all glassware between samples.
7. Ensure the device lot number and the curve details match the lot ID number selected on the reader. Failure to update the lot-specific QR code within the AccuScan Pro and AccuScan Gold reader will cause inaccurate results.

ACCUSCAN READER SET UP

AccuScan Pro reader

1. Enter the lot-specific QR code by selecting the QR code icon on the reader. Place the QR code into the cartridge and insert the cartridge into the reader. **NOTE:** For instructions on manually entering sample IDs, see the AccuScan Pro user manual.
2. Return to the home screen and select the test strip icon. Touch the **Mycotoxin** category, then select the **Fumonisin Q+** test type.

AccuScan Gold reader

1. Enter the lot specific QR code by selecting **Scan QR** from the main screen. Place the lot specific QR code into the white cartridge adapter labeled Cal/QR and place the cartridge into the reader.
2. The valid code will be scanned by the reader and provide information on the lot number and expiry date. Verify this information is correct and then add the lot ID to the reader by pressing **Add Lot ID**.

NOTE: The lot ID for the current lot will now be stored with the test ID (e.g., fumonisin) and can be selected when running a test.

SAMPLE PREPARATION

The sample to be tested should be collected according to accepted sampling techniques (see FGIS sampling protocol or contact your NEOGEN representative). Obtain a representative sample (minimum 100 g). Grind the sample so at least 95% of the ground material passes through a 20 mesh sieve (about the particle size of fine espresso).

SAMPLE EXTRACTION

1. If not using NEOGEN's prepared solution, prepare a 65% ethanol solution by mixing 6.5 parts ethanol with 3.5 parts distilled or deionized water for each sample.
2. Weigh out 10 g \pm 0.1 of sample into an extraction cup.
3. Add 50 mL of 65% ethanol to the extraction cup.
4. Vigorously shake, using hand or mechanical means, for 3 minutes, or blend for 1 minute.
5. Allow the sample to settle, then filter with a filter syringe or Whatman #4 filter paper to collect a minimum of 3 mL filtrate into a sample collection tube. Or, you may also pipette 1 mL of sample into a 2.0 mL microcentrifuge tube, and centrifuge for 30 seconds using a microcentrifuge (approx. 2,000 x g).

ALTERNATE – FGIS method

1. Combine 50 g of ground sample with 250 mL of 65% ethanol and shake vigorously for **3 minutes** or blend for **1 minute**. Allow sample to settle for 1–2 minutes.
2. Filter the extract by pouring 5 mL through a filter syringe and collecting the filtrate.

VALIDATED/VERIFIED MATRICES

1:5 extraction: (1 part sample to 5 parts 65% ethanol), broken rice, canola meal, corn grits, cornmeal, corn, green lentils, oatmeal, soybean meal, wheat flour, wheat.

1:4 extraction: (1 part sample to 4 parts 65% ethanol), barley, corn gluten meal, DDGS, DDG wet cake, DDG syrup, hominy, milo, oats, rough rice, wheat bran, wheat midds. Please contact NEOGEN for additional extraction information.

1:3 extraction: (1 part sample to 3 parts 65% ethanol), brewers rice, corn protein concentrate, lentil fiber, navy beans, rice bran, rice hulls, sorghum, sunflower meal. Please contact NEOGEN for additional extraction information.

NOTE: NEOGEN continues to validate new commodities. Please contact a representative for the latest validated commodity list.

TEST PROCEDURE

1. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary.
2. Add 100 μL of sample extract to a red dilution cup.
3. Add 200 μL of sample diluent to the red dilution cup with the sample. Mix by pipetting up and down 5 times.
4. Transfer 100 μL of diluted sample extract into a clear sample cup.
5. Place a new Reveal Q+ for Fumonisin test strip with the sample end down into the sample cup and set timer for **6 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
6. Remove the strip from the sample cup after it has developed for **6 minutes** and read immediately (within 30 seconds).
7. For the AccuScan Gold Reader: Select Category: Mycotoxin Q Plus
8. **Test Name: Q+ Fumonisin (results in ppm)**

For the AccuScan Pro Reader: **Q+ Fumonisin NE (results in ppb)**
Select Category: Mycotoxin
Test Name: Fumonisin Q+

DILUTION PROCEDURE

Samples greater than 6 ppm must be diluted and re-tested.

1. Add 200 μL of sample filtrate to a sample collection tube. Add 200 μL of 65% ethanol. Mix well by pipetting up and down 5 times.
2. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary.
3. Add 100 μL of diluted sample filtrate (from step one) to a new red sample dilution cup.
4. Add 200 μL of sample diluent to the red dilution cup with the sample. Mix by pipetting up and down 5 times.
5. Transfer 100 μL of diluted sample extract (from step 4) into a clear sample cup.
6. Place a new Reveal Q+ for Fumonisin test strip with the sample end down into the sample cup and set timer for 6 minutes. Ensure the test **strip** comes into contact with liquid and begins to wick.
7. Remove the strip from the sample cup after it has developed for 6 minutes and read immediately (within 30 seconds).

- 8. For the AccuScan Gold Reader:**
Select Category: Mycotoxin Q Plus **Test Name: Q+ Fumonisin**
- For the AccuScan Pro Reader:**
Select Category: Mycotoxin **Test Name: Fumonisin Q+**

NOTE: the AccuScan reader will not calculate your dilution. Final result displayed will need to be multiplied by 2.

READING TEST RESULTS

Test strips should be read within **30 seconds** of completion of the **6 minute** incubation. Refer to **AccuScan Reader Set Up** for test selection and set up information.

1. Select the assay type (e.g., fumonisin) from the menu and ensure the device lot number matches the lot ID number selected on the reader. **NOTE:** Failure to update the lot-specific QR code will cause inaccurate results.
2. Fully insert the Reveal Q+ test strip into the R-labeled cartridge adapter with the sample end first and results facing out.

3. Insert the cartridge with test strip upside-down into the AccuScan Gold reader (the test lines will face downward into the reader) or test-strip side up for the AccuScan Pro. The reader will automatically begin analyzing the cartridge. **CAUTION:** Removing cartridge prior to completion can result in invalid readings.
4. The AccuScan reader will analyze the test strip and results will be displayed and stored in the reader.


NOTES

1. Ensure device is fully inserted into cartridge.
2. Readings should be made between **6–6.5 minutes**. Readings after 6.5 minutes may be inaccurate due to overdevelopment of the device.

TEST PROCEDURE – RAPTOR INTEGRATED ANALYSIS PLATFORM

1. Place the appropriate number of red sample dilution cups into a sample cup rack. Label cups if necessary.
2. Add 200 μL of sample filtrate to each red sample dilution cup.
3. Add 400 μL of sample diluent to each red sample dilution cup. Mix by pipetting up and down five times.
4. Fully insert a Reveal Q+ for Fumonisin test strip into a Raptor cartridge.
5. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor Integrated Analysis Platform reader.
 - a. The bar code on the test strip will be read – the Raptor reader identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader on the front of the Raptor will turn on automatically.
 - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the system.
6. Enter a sample ID if desired. Hit the Accept button to go to the Add Sample screen.
7. Add 400 μL of sample from the red sample dilution cup (from step 3) to the Raptor cartridge.
 - a. The Raptor reader will start automatically.
 - b. Additional samples can be started in the other ports while the first sample is processing
8. Results will be displayed on the Raptor screen after the 6-minute testing is complete.

DILUTION PROCEDURE – RAPTOR INTEGRATED ANALYSIS PLATFORM

Samples greater than 6 ppm must be diluted and re-tested.

1. Place the appropriate number of red sample dilution cups into a sample cup rack. Label cups if necessary.
2. Add 200 µL of sample filtrate to a sample collection tube.
3. Add 200 µL of 65% ethanol. Mix well by pipetting up and down 5 times.
4. Transfer 200 µL of diluted sample (from step 3) to a new red sample dilution cup.
5. Add 400 µL of sample diluent (pink-labeled bottle) to the red sample dilution cup. Mix by pipetting up and down five times.
6. Fully insert a Reveal Q+ for Fumonisin test strip into a Raptor cartridge.
7. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor Integrated Analysis Platform reader.
 - a. The bar code on the test strip will be read – the Raptor reader identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader on the front of the Raptor will turn on automatically.
 - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the system.
8. Enter a sample ID if desired. Hit the Accept button to go to the add sample screen.
9. Add 400 µL of sample extract from the red sample dilution cup (from step 5) to the Raptor cartridge.
10. The Raptor reader will start automatically
11. Additional samples can be started in the other ports while the first sample is processing
12. Results will be displayed on the Raptor screen after the 6-minute testing is complete. The final result will need to be multiplied by 2.

PERFORMANCE CHARACTERISTICS

Reveal Q+ for Fumonisin is designed for the quantitative analysis of fumonisin.

1. Limit of detection: 0.3 ppm
2. Range of detection: 0.3–6 ppm

NOTE: Samples greater than 6 ppm must be diluted and re-tested. Results below the range of detection should be reported as less than 0.3 ppm.

CUSTOMER SERVICE

NEOGEN Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all NEOGEN test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of NEOGEN's test kits, on NEOGEN's Web site at NEOGEN.com, or by calling NEOGEN at 800.234.5333 or 517.372.9200.

TERMS AND CONDITIONS

For NEOGEN's full terms and conditions, please visit neogen.com/terms-and-conditions/

WARRANTY

NEOGEN Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ergot alkaloids, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, coconut crustaceans, eggs, gliadin, hazelnut, milk, multi-treenut, mustard, peanuts, sesame, soy, walnuts

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