*Read instructions carefully before starting test*

# for Coconut Allergen

## Quantitative Test

REFRIGERATE AT 2–8°C (35–46°F) • DO NOT FREEZE

### **COCONUT ALLERGEN**

Food allergens are proteins in food that can create an immune response in sensitive individuals. Once ingested, food allergens can cause a number of reactions, ranging in severity from hives and itching to anaphylaxis. Anaphylaxis is a severe allergic reaction, involving vomiting, diarrhea, difficulty breathing, swelling of the mouth and tongue, and a rapid drop in blood pressure.

An estimated 32 million people in the United States are sensitive in some degree to food allergens (about 10% of the population). Food manufacturers protect those with food allergies by clearly labeling their products with a list of ingredients. Testing for the presence of coconut components ensures food manufacturers that an unlabeled — and potentially dangerous — ingredient did not make its way into a food product.

### **INTENDED USE**

Veratox® for Coconut Allergen is intended for the full quantitative analysis or simple screening, of coconut protein residues in food products such as cookies, chocolate bars, and beverages.

### **INTENDED USER**

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by coconut or coconut products. Since technique is very important, operators should be trained by a NEOGEN representative or someone who has completed the NEOGEN training.

### **ASSAY PRINCIPLES**

The Veratox for Coconut Allergen test is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Coconut protein is extracted from samples with a phosphate buffered salt solution (PBS) by shaking in a heated water bath and clarified by centrifugation or filtration. The extracted coconut protein is added to antibody-coated wells (capture antibody) where it binds to the antibody during an incubation. Next, a washing step is performed to remove unbound proteins, and then an enzyme-labeled antibody (detector antibody), is added to the microwells. The detector antibody binds coconut protein captured in the microwells. After a second washing step, the enzyme substrate is added. A blue color will develop when the coconut protein is present as a result of the bound detector antibody. Lastly, Red Stop

reagent is added to stop the color development reaction. The optical density of color development is measured in a microwell reader. The coconut protein concentration for unknown samples is determined by plotting the samples' optical densities against those generated from a standard curve.

### **STORAGE REQUIREMENTS**

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). Do not freeze.

### **MATERIALS PROVIDED**

1. 48 antibody-coated microwells
2. 48 red-marked transfer wells
3. 6 yellow-labeled bottles of 0, 1, 2.5, 5, 10, and 25 ppm coconut controls
4. 2 blue-labeled bottles of enzyme-labeled antibody conjugate
5. 1 green-labeled bottle of K-Blue® Substrate
6. 1 red-labeled bottle of Red Stop solution
7. 5 foil pouches of phosphate-buffered saline (PBS) powder. Each pouch contains PBS powder sufficient to prepare 1.0 L of 10 mM PBS buffer at pH 7.4.
8. 40 mL of 10 mM PBS-Tween washing reagent concentrate in a wide mouth bottle. Each bottle contains wash reagent concentrate sufficient to prepare 1 L in distilled or deionized water (pH 7.4)
9. 50 g of extraction additive in a specimen cup
10. Plastic scoop to measure extraction additive

### **MATERIALS RECOMMENDED BUT NOT PROVIDED**

1. Allergen Extraction Kit (NEOGEN item 8429)
  - a. 20 disposable plastic extraction bottles
  - b. 20 sample collection tubes (12 x 75 mm) with caps
2. Shaker water bath capable of maintaining  $60 \pm 1^\circ\text{C}$  with holders for 250 mL disposable plastic bottles (NEOGEN items 9298, 9299)
3. Whatman No. 4 filters or equivalent (NEOGEN item 9429)
4. Centrifuge (optional)
5. Pipettor, 50–200  $\mu\text{L}$ , adjustable (NEOGEN item 9276)
6. Pipettor, 12-channel (NEOGEN item 9273)
7. Pipette tips (NEOGEN item 9410, 9407, 9417)
8. Timer (NEOGEN items 9426, 9452)
9. Microwell reader with a 650 nm filter (NEOGEN item 9303)
10. 1 L bottle to prepare washing solution (NEOGEN item 9472)
11. 1 L heat safe bottle to prepare extract solution (NEOGEN item 9472)
12. Paper towels or equivalent absorbent material
13. Microwell holder (NEOGEN item 9402)
14. Waterproof marker
15. Wash bottle (NEOGEN item 9400)
16. Distilled or deionized water
17. 3 reagent boats for 12-channel pipettor (NEOGEN item 9435)
18. Graduated cylinder capable of measuring 125 mL (NEOGEN item 9368)
19. Scale capable of weighing  $5 \pm 0.1$  g (NEOGEN item 9427)

## PRECAUTIONS

1. Components of Veratox for Coconut Allergen, including controls, contain coconut. If allergic, please use caution when handling this product.
2. Concentrated food additives, colors and flavors may cause interferences on ELISA test methods. Contact NEOGEN's technical services for validation information.
3. Hydrolyzed and fermented proteins may not be detected using ELISA methods for allergen testing. Due to the breakdown of proteins to small peptides or amino acids, they may become undetectable by this assay, but still could be allergenic and cause an allergic reaction.
4. Sponges should not be used for sample collection and allergen testing. Sample collection swabs other than NEOGEN swabs should be validated prior to use. General sponges and swabs may contain solutions or materials that may interfere with the test kit.
5. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze test kits.
6. Do not use kit components beyond expiration date.
7. Do not mix reagents from one kit serial with reagents from a different kit serial.
8. Do not run more than 24 wells per test.
9. Follow proper pipetting techniques (e.g., prime tips and use clean tips).
10. Use only incubation times specified. Others may give inaccurate results.
11. Bring kits to room temperature 18–30°C (64–86°F) prior to use.
12. Avoid prolonged storage of kits at ambient temperatures.
13. Use clean pipette tips and glassware for each sample to avoid cross-contamination. Thoroughly wash all glassware between samples.

## PROCEDURAL NOTES

1. Substrate. K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. Do not return unused substrate to the bottle. Cover the reagent boat to keep the substrate protected from light until needed.
2. Conjugate. The conjugate supplied with this kit is ready to use. One bottle is enough for 24 wells. Cover the reagent boat to keep the conjugate protected from direct light and contaminants.
3. Antibody wells. Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.
4. Extraction solution. Prepare extraction solution by adding a foil pouch of PBS powder to 1 L distilled or deionized water. Swirl to mix thoroughly. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).
5. Wash buffer. Prepare the wash buffer solution by pouring all the wash buffer concentrate into an empty 1 L container. Rinse the wash buffer concentrate bottle with distilled or deionized water and pour into the 1 L container to ensure all the concentrate is used. Fill the 1 L container with additional distilled or deionized water, and swirl to assure thorough mixing. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).

**NOTE:** Discard unused portions of extraction solution and wash buffer when the test kit has been used completely.

## SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques (see NEOGEN's Food Allergen Handbook). The sample should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

1. Prepare the extraction solution as described in the procedural notes.
2. Preheat extraction solution to 60°C (140°F) by immersing the bottle containing the solution into the water bath and allowing it to reach 60°C.
3. Using your sampling and collection procedure, obtain a representative sample and grind it to a very fine particle size.
4. Transfer 5 g of sample or 5 mL of liquid sample into a 250 mL disposable plastic bottle.
5. Add one level scoop of the extraction additive to the sample bottle.
6. Pour 125 mL of the 60°C (140°F) extraction solution to the sample plastic bottle.
7. Cap the bottle to prevent contents from splashing during the extraction.
8. Extract by shaking (150 rpm) in a water bath at 60°C (140°F) for 15 minutes. Remove the bottle from the bath.
9. Let material settle for 5 minutes to enable some of the sample to settle before proceeding to the next step.
10. Filter the extract by pouring at least 5 mL through a Whatman No. 4 filter and collecting the filtrate as a sample.  
Alternative: Centrifuge at 14,000 rpm for 5 minutes (20 minutes for lower speeds). Use the clear supernatant as a sample.
11. Allow extracts to cool to room temperature before beginning analysis.
12. Discard extracts after completion of analysis.

## TEST PROCEDURE FOR QUANTITATION

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

1. Remove 1 red-marked mixing well for each sample to be tested plus 6 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of the strip with a "1," and place strip in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Using a new pipette tip for each, transfer 150 µL of controls and sample extracts to the red-marked transfer wells as shown in the template below. Only run up to two 12-well strips at a time.

0	1	2.5	5	10	25	S1	S2	S3	S4	S5	S6	Strip 1
S7	S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	Strip 2

5. Place tips on the 12-channel pipettor and transfer 100 µL of the controls and sample extracts to the antibody-coated wells. Mix for 20 seconds by sliding the well holder back and forth on a flat surface.
6. Incubate microwells 10 minutes at room temperature 18–30°C (64–86°F). Discard the red-marked transfer wells.

7. Empty the contents of the wells into a sink. With a wash bottle fill each antibody well with the wash buffer solution and dump out. Repeat the washing five times, then turn the wells upside down and tap out on a paper towel until the remaining washing solution is removed.
8. Pour the needed volume of conjugate from the blue-labeled bottle into a clean reagent boat.
9. Using the 12-channel pipettor and new tips, transfer 100  $\mu$ L of the conjugate into all the wells and mix for 20 seconds by sliding the well holder back and forth on a flat surface.
10. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
11. Wash all wells with the wash buffer solution as described in step 7.
12. Pour the needed volume of substrate solution from the green-labeled bottle into a clean reagent boat.
13. Place new tips on the 12-channel pipettor and transfer 100  $\mu$ L of substrate into each well and mix for 20 seconds. Do not eject tips.
14. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
15. Pour the needed volume of Red Stop solution from the red-labeled bottle into a clean reagent boat.
16. With the same tips used to dispense the substrate, transfer 100  $\mu$ L of Red Stop into each well and mix for 20 seconds.
17. Wipe the bottom of the microwells and read in a microwell reader with a 650 nm filter. Results should be read within 20 minutes after the addition of Red Stop solution.
18. Interpret the test's results using NEOGEN's Stat-Fax microwell reader, or an equivalent strip reader. If using a strip reader, calculate the results using NEOGEN's Veratox software.

### TEST PROCEDURE FOR SCREENING

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

1. Remove 1 well for each sample to be tested plus 1 well for the control, and place into the well holder.
2. Choose one yellow-labeled control bottle to serve as your screening level for the test. 2.5 ppm is recommended.
3. Mix each reagent by swirling its bottle prior to use.
4. Add 100  $\mu$ L from the yellow-labeled control bottle to the first well. Add 100  $\mu$ L from each sample extract to a respective well as indicated in the template below. For environmental swabs, add three drops from swab tube with dropper tip. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.

Control	S1	S2	S3	S4	S5
---------	----	----	----	----	----

5. Incubate microwells 10 minutes at room temperature 18–30°C (64–86°F).
6. Shake out the contents of the wells. Using a wash bottle filled with wash buffer solution, fill each well and shake out. Repeat five times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
7. Add 100  $\mu$ L from the blue-labeled conjugate bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.

8. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
9. Shake out the contents of the wells. Using a wash bottle filled with wash buffer solution, fill each well and shake out. Repeat five times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
10. Add 100 µL from the green-labeled substrate bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.
11. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
12. Add 100 µL from the red-labeled Red Stop bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface. The results are now ready to be interpreted.
13. Visually compare the color of a sample well to the color of the control well. If the sample well has more blue color than the control well, the sample tests positive for coconut protein contamination of more than the control used. If the sample well has less blue color, or more red color, than the control well, the sample contains less than the control used of coconut protein contamination.  
Alternative: Read wells (wipe bottom of wells first) in a microwell reader with a 650 nm filter. If the sample well has an optical density (OD) higher than the control well, the sample is positive for coconut protein contamination of more than the control used. If the sample well has an OD lower than the control well, the sample contains less than the control used of coconut protein contamination.

#### **PERFORMANCE CHARACTERISTICS**

Limit of quantitation: 1.0 ppm coconut protein (Described as the lowest concentration point on the calibration curve that this test can reliably detect coconut allergen.)

Range of quantitation: 1.0–25 ppm coconut protein (For quantitating samples above 25 ppm, contact a NEOGEN representative for dilution instructions.)

Allergen detection: This test detects coconut protein, and the results are expressed as ppm of coconut protein.

Cross-reactivity: Cross-reacts with: Brazil nut, oat, cherry pit, peach pit, cashew, poppy seed, and hazelnut. Please see Veratox for Coconut validation report for additional information.

**CUSTOMER SERVICE**

NEOGEN Customer Assistance and Technical Service can be reached by using the contact information on the back of this booklet. Training on this product, and all NEOGEN test kits, is available.

**SDS INFORMATION AVAILABLE**

Safety data sheets (SDS) are available for this test kit, and all of NEOGEN's Food Safety test kits, at [foodsafety.NEOGEN.com](http://foodsafety.NEOGEN.com), or by calling NEOGEN at 800.234.5333 or 517.372.9200.

**TERMS AND CONDITIONS**

For NEOGEN's full terms and conditions, please visit [NEOGEN.com/terms-and-conditions/](http://NEOGEN.com/terms-and-conditions/)

**WARRANTY**

NEOGEN Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## TESTING KITS AVAILABLE FROM NEOGEN

### Natural toxins

- Aflatoxin, DON, ergot alkaloids, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

### Food allergens

- Almonds, coconut, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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