Read instructions carefully before starting test

Reveal[®] Q+

for Zearalenone

Quantitative Test

FGIS 2018-115

THE TOXIN

Zearalenone primarily is produced by the mold *Fusarium graminearum*, which also commonly produces deoxynivalenol (DON). Hence, there is evidence if zearalenone is detected, there is a high probability other fusarial mycotoxins may be present. Zearalenone is classified as an estrogenic mycotoxin because it frequently causes estrogenic responses in animals.

When zearalenone-contaminated feed or grain is eaten by livestock, it can cause a wide variety of reproductive problems. In swine, it causes vulvovaginitis, low birth weights, fetal reabsorption, aborted pregnancies, reduced litter sizes, abnormal estrus, and feminization of immature males. Zearalenone can delay the breeding process and cost the producer significant economic and physical losses.

Livestock producers increasingly are becoming aware of zearalenone problems, and have looked for ways to reduce risks related to contaminated feed.

The best protection against mycotoxins is monitoring for their presence in feeds and foods, which means testing all along the pathway from initial harvest of grains to the finished product.

The FDA has not issued any advisory levels for zearalenone.

European Union regulations for zearalenone are as follows:

Foodstuffs	Maximum levels ppb
Unprocessed corn	200
Other unprocessed cereals and corn for direct human consumption	100
Other cereals for direct human consumption	75
Bread, pastries, breakfast cereals, snacks, corn-based snacks and corn-based breakfast cereals	50
Baby foods and cereals intended for infants	20

INTENDED USE/USER

Reveal Q+ for Zearalenone is intended for the quantitative analysis of corn and wheat products for zearalenone. The test kit is designed for use by quality control personnel and others familiar with commodities possibly contaminated by zearalenone.

ASSAY PRINCIPLES

Reveal Q+ for Zearalenone is a single-step lateral flow immunochromatographic assay based on a competitive immunoassay format. The extract is wicked through a reagent zone, which contains antibodies specific for zearalenone conjugated to colloidal gold particles. If zearalenone is present, it will be captured by the particle-antibody complex. The zearalenone-antibody-particle complex then is wicked onto a membrane, which contains a zone of zearalenone conjugated to a protein carrier. This zone captures any uncomplexed zearalenone antibody, allowing the particles to concentrate and form a visible line. As the level of zearalenone in a sample increases, free zearalenone will complex with the antibody-gold particles. This allows less antibody-gold to be captured in the test zone. Therefore, as the concentration of zearalenone in the sample increases, the test line density decreases. Algorithms programmed into the readers convert these line densities into a quantitative result displayed in parts per billion (ppb). The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line always will form regardless of the presence of zearalenone, ensuring the strip is functioning properly.

STORAGE REQUIREMENTS

Store kit components at room temperature (18–30°C, 64–86°F) to ensure full shelf life. Test strips should remain capped in their original tubes until used to ensure optimal performance.

MATERIALS PROVIDED

1. 25 Reveal Q+ for Zearalenone test strips
2. 25 red sample dilution cups
3. 25 clear sample cups
4. 1 bottle of sample diluent
5. Instructions for use

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. 65% ethanol solution (NEOGEN item 8073, 8074)
2. Sample collection cups with lids (NEOGEN item 9428)
3. Agri-Grind grinder or equivalent (NEOGEN item 9401, 9453)
4. Scale capable of weighing 5–50 g \pm 0.1 g (NEOGEN item 9427)
5. Timer (NEOGEN item 9426)
6. Reveal sample cup rack (NEOGEN item 9475)
7. AccuScan reader: AccuScan Pro or AccuScan Gold
8. Dispensing pump or graduated cylinder (NEOGEN item 9448, 9447)
9. Filter syringe (NEOGEN item 9420)
10. Sample collection tubes with caps (NEOGEN item 9421, 9421B)
11. Pipettor, 100 μ L (NEOGEN item 9272, 9278)
12. Pipettor, 200 μ L fixed (NEOGEN item 9488)
13. Pipette tips, 20–200 μ L (NEOGEN item 9407, 9410, 9417)

PRECAUTIONS

1. The test strips must remain inside the stay-dry tube before use.
2. Ethanol is highly flammable. Keep container tightly closed, and away from heat, sparks, open flame and those who are smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
3. Store test kit at room temperature (18–30°C, 64–86°F) when not in use. Do not freeze.
4. Do not use kit components beyond expiration date.
5. Treat all used liquids, including sample extract, and labware as if contaminated with zearalenone. Gloves and other protective apparel should be worn at all times.
6. To avoid cross-contamination, use clean glassware for each sample, and thoroughly wash all glassware between samples.
7. Ensure the device lot number and the curve details match the lot ID number selected on the reader. Failure to update the lot-specific QR code within the AccuScan Pro and AccuScan Gold will cause inaccurate results.

ACCUSCAN READER SET UP

AccuScan Pro reader

1. Enter the lot-specific and commodity-specific QR code by selecting the QR code icon on the reader. Place the QR code into the cartridge and insert the cartridge into the reader. **NOTE:** For instructions on manually entering sample IDs, see the AccuScan Pro user manual.
2. Return to the home screen and select the test strip icon. Touch the mycotoxin category, then select the **Zear Q+ Corn** or **Zear Q+ Wheat** test type.

AccuScan Gold reader

1. Enter the lot-specific and the commodity-specific QR code by selecting **Scan QR** from the main screen. Place the specific QR code into the white cartridge adapter labeled Cal/QR and place the cartridge into the reader.
2. The valid code will be scanned by the reader and provide information on the lot number and expiry date. Verify this information is correct and then add the lot ID to the reader by pressing **Add Lot ID**. **NOTE:** The lot ID for the current lot will now be stored with the test ID (e.g., **Q+ Zearalenone Corn** or **Q+ Zearalenone Wheat**) and can be selected when running a test.

3. Return to the main screen and select the **Mycotoxin Q+** category then select the **Q+ Zear** test type.
4. The reader will prompt the user to select the lot ID of the kit being tested. A sample ID can be added at this time.

SAMPLE PREPARATION

The sample to be tested should be collected according to accepted sampling techniques (see FGIS sampling protocol or contact your NEOGEN representative). Obtain a representative sample (minimum 100 g). Grind the sample so at least 95% of the ground material passes through a 20- mesh sieve (about the particle size of fine espresso). If not using NEOGEN's prepared solution, prepare a 65% ethanol solution by mixing 6.5 parts ethanol with 3.5 parts distilled or deionized water for each sample.

PROTOCOL FOR CORN CURVE SET

Sample Extraction

1. Extract at a ratio of 1 part sample to 3 parts 65% ethanol. For example, combine 10 g of ground sample with 30 mL of 65% ethanol.
2. Vigorously shake, using hand or mechanical means (250 rpm) for **3 minutes**, or blend for **1 minute**.
3. Allow the sample to settle, then filter at least 4 mL with a filter syringe, or Whatman No. 1 filter paper. Alternatively, pipette sample into a 2.0 mL microcentrifuge tube and centrifuge for 30 seconds.
4. The sample is now ready for testing.

FGIS Method Sample Extraction

1. Combine 50 g of ground sample with 150 mL of 65% ethanol.
2. Vigorously shake, using hand or mechanical means (250 rpm) for **3 minutes**.
3. Allow the sample to settle for 1 minute, then filter at least 3-5 mL with a filter syringe
4. The sample is now ready for testing.

Test Procedure

1. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary
2. Add 100 µL of sample extract to the red sample cup.
3. Add 200 µL of sample diluent to the red dilution cup with the sample extract. Mix by pipetting up and down 5 times.
4. Transfer 100 µL of diluted sample extract into a new clear sample cup.
5. Place a new Reveal Q+ for Zearalenone test strip with the sample end down into the sample cup and set timer for **6 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
6. Remove the strip from the sample cup after it has developed for 6 minutes and read immediately (within 30 seconds).
7. **For the AccuScan Gold Reader:** **Select Category: Mycotoxin Plus**
 Test Name: Q+ Zear Corn
For the AccuScan Pro Reader: **Select Category: Mycotoxin**
 Test Name: Zear Q+ Corn

Dilution Procedure

Samples greater than 1200 ppb must be diluted and re-tested.

1. Add 100 µL sample filtrate to a sample collection tube.
2. Add 200 µL 65 % ethanol solution to the sample collection tube. Mix well.
3. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary
4. Add 100 µL of diluted sample extract (from step 2) to the red sample cup.
5. Add 200 µL of sample diluent to the red dilution cup with the sample extract. Mix by pipetting up and down 5 times.
6. Transfer 100 µL of diluted sample extract into a new clear sample cup.
7. Place a new Reveal Q+ for Zearalenone test strip with the sample end down into the sample cup and set timer for **6 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
8. Remove the strip from the sample cup after it has developed for 6 minutes and read immediately (within 30 seconds).
9. **For the AccuScan Gold Reader: Select Category: Mycotoxin Plus
Test Name: Q+ Zear Corn**
**For the AccuScan Pro Reader: Select Category: Mycotoxin
Test Name: Zear Q+ Corn**
Final result displayed will need to be multiplied by 3.

PROTOCOL FOR WHEAT CURVE SET

Sample Extraction

1. Extract at a ratio of 1 part sample to 5 parts 65% ethanol. For example, combine 10 g of ground sample with 50 mL of 65% ethanol.
2. Vigorously shake, using hand or mechanical means (250 rpm) for **3 minutes**, or blend for **1 minute**.
3. Allow the sample to settle, then filter at least 4 mL with a filter syringe, or Whatman No. 1 filter paper. Alternatively, pipette sample into a 2.0 mL microcentrifuge tube and centrifuge for 30 seconds.
4. The sample is now ready for testing.

FGIS Method Sample Extraction

1. Combine 50 g of ground sample with 250 mL of 65% ethanol.
2. Vigorously shake, using hand or mechanical means (250 rpm) for **3 minutes**.
3. Allow the sample to settle for 1 minute, then filter at least 3–5 mL with a filter syringe
4. The sample is now ready for testing. **NOTE: when the test is complete, multiply the results by 1.66 for wheat.**

Test Procedure

1. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary
2. Add 100 µL of sample extract to the red sample cup.
3. Add 200 µL of sample diluent to the red dilution cup with the sample extract. Mix by pipetting up and down 5 times.
4. Transfer 100 µL of diluted sample extract into a new clear sample cup.
5. Place a new Reveal Q+ for Zearalenone test strip with the sample end down into the sample cup and set timer for **6 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
6. Remove the strip from the sample cup after it has developed for 6 minutes and read immediately (within 30 seconds).
7. **For the AccuScan Gold Reader: Select Category: Mycotoxin Q+ Plus
Test Name: Q+ Zear Wheat**
**For the AccuScan Pro Reader: Select Category: Mycotoxin
Test Name: Zear Q+ Wheat**

Dilution Procedure

Samples greater than 1200 ppb must be diluted and re-tested.

1. Add 100 µL sample filtrate to a sample collection tube.
 2. Add 200 µL 65 % ethanol solution to the sample collection tube. Mix well.
 3. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary
 4. Add 100 µL of diluted sample extract (from step 2) to the red sample cup.
 5. Add 200 µL of sample diluent to the red dilution cup with the sample extract. Mix by pipetting up and down 5 times.
 6. Transfer 100 µL of diluted sample extract into a new clear sample cup.
 7. Place a new reveal Q+ for Zearalenone test strip with the sample end down into the sample cup and set timer for **6 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
 8. Remove the strip from the sample cup after it has developed for 6 minutes and read immediately (within 30 seconds).
 9. **For the AccuScan Gold Reader: Select Category: Mycotoxin Q+ Plus
Test Name: Q+ Zear Wheat**
**For the AccuScan Pro Reader: Select Category: Mycotoxin
Test Name: Zear Q+ Wheat**
- Final result displayed will need to be multiplied by 3.

READING TEST RESULTS

Test strips should be read within **30 seconds** of the 6 minute incubation. Refer to

AccuScan Reader Set Up for test selection and set up information.

1. Select the assay type (e.g., zearalenone) from the menu and ensure the device lot number matches the lot ID number selected on the reader. **NOTE:** Failure to update the lot-specific QR code will cause inaccurate results.
2. Fully insert the Reveal Q+ test strip into the black R-labeled cartridge adapter with the sample end first and results facing out.

3. Insert the cartridge with test strip upside-down into the AccuScan Gold reader (the test lines will face downward into the reader) or test-strip side up for the AccuScan Pro. The reader will automatically begin analyzing the cartridge.

CAUTION: Removing cartridge prior to completion can result in invalid readings.

4. The AccuScan reader will analyze the test strip and results will be displayed and stored in the reader.

NOTES

1. Ensure device is fully inserted into cartridge.
2. Readings should be made within 30 seconds of the 6-minute incubation time. Readings after 6.5 minutes may be inaccurate due to over-development of the device.
3. The strips must be read using NEOGEN's AccuScan Pro or AccuScan Gold readers.

PERFORMANCE CHARACTERISTICS

Corn Curve Commodities

Limit of Detection: 50 ppb

Range of Detection: 50 – 1200 ppb

Wheat Curve Commodities

Limit of Detection: 25 ppb

Range of Detection: 25 – 1200 ppb

Samples greater than 1200 ppb must be diluted and re-tested following the dilution protocols.

VALIDATED MATRICES

Curve set 1: Corn

Curve set 2: Wheat

NOTE: NEOGEN continues to validate new commodities. Please contact a representative for the latest validated commodity list.

CUSTOMER SERVICE

NEOGEN Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all NEOGEN test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of NEOGEN's test kits, on NEOGEN's website at NEOGEN.com, or by calling NEOGEN at 800.234.5333 or 517.372.9200.

TERMS AND CONDITIONS

For NEOGEN's full terms and conditions, please visit www.NEOGEN.com/Corporate/termsconditions.html.

WARRANTY

NEOGEN Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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©NEOGEN Corporation, 2020. NEOGEN, AccuScan y Reveal son marcas comerciales registradas de NEOGEN Corporation. Todas las demás marcas y nombres de producto son marcas comerciales o marcas comerciales registradas de sus respectivas compañías.