*Read instructions carefully before starting test*

# Veratox<sup>®</sup> MAX

## for Total Aflatoxin

### Quantitative Test with Aqueous Extraction

REFRIGERATE AT 2–8°C • DO NOT FREEZE

**FGIS – 2016-088**

#### THE TOXIN

Aflatoxin is a toxic and carcinogenic substance produced by certain strains of the molds *Aspergillus flavus* and *A. parasiticus*. The four principle types of aflatoxin are: B<sub>1</sub>, B<sub>2</sub>, G<sub>1</sub> and G<sub>2</sub>. Aflatoxin B<sub>1</sub> is the most frequently encountered of the group and the most toxic. The commodities most affected by aflatoxin are corn, peanuts, cottonseed, milo, and the majority of tree nuts.

The effects in animals of ingesting excessive amounts of the toxin range from chronic health and performance problems to death. Aflatoxin has been shown to cause liver damage or cancer, decreased milk and egg production, immune suppression and interference with reproductive efficiency.

The FDA has set maximum allowable levels of aflatoxin in food and feed. Therefore, accurate determination of the presence of the toxin is of major importance to those monitoring the quality of food and feed in which aflatoxin may occur. Testing these commodities for the toxin requires careful sampling, chemical extraction, sanitation and quantitative analysis.

The FDA has issued regulatory levels for aflatoxin as follows:

| For                                             | Level   | Commodities                  |
|-------------------------------------------------|---------|------------------------------|
| Humans                                          | 20 ppb  | All food except milk         |
| All animal species                              | 20 ppb  | All feed (exceptions below)  |
| <b>Exceptions:</b>                              |         |                              |
| Breeding cattle, breeding swine, mature poultry | 100 ppb | Corn                         |
| Finishing swine (>100 lbs.)                     | 200 ppb | Corn                         |
| Finishing beef cattle                           | 300 ppb | Corn                         |
| All animal species                              | 300 ppb | Cottonseed meal used in feed |

#### INTENDED USE

Veratox MAX for Total Aflatoxin is intended for the quantitative analysis of aflatoxin in corn and corn gluten meal. This kit is designed to detect the four principle types of aflatoxin: B<sub>1</sub>, B<sub>2</sub>, G<sub>1</sub>, and G<sub>2</sub> with superior cross-reactivity when compared to other available kits. Veratox for Total Aflatoxin is the first ELISA microwell test to utilize an aqueous extraction procedure.

## INTENDED USER

This test kit is designed for use by quality control personnel and others familiar with food and feed possibly contaminated by aflatoxin. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed the Neogen training.

## ASSAY PRINCIPLES

Veratox MAX for Total Aflatoxin is a direct competitive ELISA in a microwell format which allows the user to obtain exact concentrations in parts per billion (ppb). Free aflatoxin in the samples and controls is allowed to compete with enzyme-labeled aflatoxin (conjugate) for the antibody binding sites. After a wash step, substrate is added, which reacts with the bound conjugate to produce blue color. More blue color means less aflatoxin. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form the standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of aflatoxin.

## STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). The MAX 2 packets can be used until the expiration date listed on the label when stored at room temperature 18–30°C (64–86°F).

## MATERIALS PROVIDED

1. 48 antibody-coated microwells
2. 48 red-marked mixing wells
3. 4 yellow-labeled bottles of 0, 5, 15 and 50 ppb aflatoxin controls
4. 1 blue-labeled bottle of aflatoxin-HRP conjugate solution
5. 1 green-labeled bottle of K-Blue® Substrate solution
6. 1 red-labeled bottle of Red Stop solution
7. 25 MAX 2 Aqueous Extraction Packets
8. 1 pink-labeled bottle of dilution diluent

## MATERIALS REQUIRED BUT NOT PROVIDED

1. Extraction materials
    - a. MAX 2 Aqueous Extraction Packets (Neogen item 8036)
    - b. MAX 2 G-50 Aqueous Extraction Packets (Neogen item 8036G)
    - c. 250 mL graduated cylinder or dispensing pump (Neogen item 9368, 9448)
    - d. Container with 125 mL capacity (Neogen item 9428)
    - e. Filter syringes, Whatman #4 filter paper, or equivalent (Neogen item 9420, 9519, 9429)
    - f. Sample collection tubes (Neogen item 9421)
  - or**
  - g. Centrifuge, mini (Neogen item 9330)
  - h. Micro-centrifuge tubes (Neogen item 9172)
  2. Agri-Grind grinder or equivalent (Neogen item 9401, 9453)
  3. Scale capable of weighing 5-50 grams (Neogen item 9427)
  4. Microwell reader with a 650 nm filter (Neogen item 9303)
  5. Pipettor, 12-channel (Neogen item 9273)
  6. Pipettor, 100 µL (Neogen item 9272, 9290)
  7. Pipette tips for 100 µL and 12-channel pipettors (Neogen item 9410, 9407, 9417)
  8. Paper towels or equivalent absorbent material
  9. Plastic bucket for use as waste receptacle
  10. Microwell holder (Neogen item 9402)
  11. Timer (Neogen item 9426)
  12. Waterproof marker
  13. Wash bottle (Neogen item 9400)
  14. 2 reagent boats for 12-channel pipettor (Neogen item 9450)
  15. Distilled or deionized water
  16. Shaker, Rockit (Neogen item 9390)
  17. Blender (Neogen item 9493)
  18. Blender jar, mini (Neogen item 9477)
-

## PRECAUTIONS

1. Store test kit refrigerated, between 2–8°C (35–46°F) when not in use. Do not freeze.
2. Kits should be brought to room temperature 18–30°C (64–86°F) for at least 2 hours prior to use.
3. Avoid prolonged storage of kits at ambient temperatures.
4. Do not use kit components beyond expiration date.
5. Do not mix reagents from one kit serial with reagents from a different kit serial.
6. Do not run more than 24 wells at a time.
7. Follow proper pipetting techniques, including priming of tips.
8. Use of incubation times other than those specified may give inaccurate results.
9. Treat all used liquids, including sample extract, and labware as if contaminated with aflatoxin.
10. To avoid cross-contamination, use clean pipette tips and glassware for each sample, and thoroughly detoxify and wash all glassware between samples.
11. Commodities tested should have a pH of 6–8. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH contact your Neogen representative or Technical Services.

## PROCEDURAL NOTE

1. K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until it is needed.
2. **Antibody wells:** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test is set to begin.

## SAMPLE PREPARATION

The sample to be tested should be collected according to accepted sampling techniques (see FGIS sampling protocol or contact your Neogen representative). Obtain a representative sample (minimum 100 g). Grind the entire sample so that at least 95% of the ground material passes through a 20 mesh sieve (about the particle size of a fine instant espresso).

## SAMPLE EXTRACTION

1. Weigh out 10 g  $\pm$  0.1 g of sample into container cup. Pour contents of 1 MAX 2 Aqueous Extraction Packet into the container cup.
2. Add 50 mL of distilled or deionized water over the sample and extraction powder.
3. Vigorously shake, using hand or mechanical means, for 3 minutes. Let sample settle for 10 minutes.
4. Filter the extract with a filter syringe or through a Whatman #4 filter paper to collect a minimum of 3 mL sample filtrate into a sample collection tube. You may also pipette 1 mL of sample into a 2.0 mL microcentrifuge tube, and centrifuge for 30 seconds using a micro-centrifuge.
5. The sample is now ready for testing.

## FGIS METHOD SAMPLE EXTRACTION

1. Weigh out 50 g  $\pm$  0.1 g of sample into container cup. Pour contents of 1 MAX2 G-50 Aqueous Extraction Packets into the container cup.
2. Add 250 mL of distilled or deionized water over the sample and extraction powder.
3. Vigorously shake using mechanical means for 3 minutes or blend for 1 minute. Let sample settle for 10 minutes.
4. Filter the extract by pouring at least 5 mL through a Whatman #4 filter, or use a filter syringe.
5. The sample is now ready for testing.

## TEST PROCEDURE

Allow all reagents to warm to room temperature 18–30°C (64–86°F) for at least two hours prior to use.

1. Remove 1 red-marked mixing well for each sample to be tested plus 4 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of strip with a “1”, and place strip in the well holder with the marked end on the right. Do not mark the inside or bottom of the wells.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Place 100 µL of conjugate from the blue-labeled bottle in each red-marked mixing well.
5. Using a new pipette tip for each, transfer 100 µL of controls and samples to the red-marked mixing wells as described below.

|    |     |     |     |     |     |     |     |     |     |     |     |         |
|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|---------|
| 0  | 5   | 15  | 50  | S1  | S2  | S3  | S4  | S5  | S6  | S7  | S8  | Strip 1 |
| S9 | S10 | S11 | S12 | S13 | S14 | S15 | S16 | S17 | S18 | S19 | S20 | Strip 2 |

6. Using a 12-channel pipettor, mix the liquid in the wells by pipetting it up and down 5 times. Transfer 100 µL to the antibody-coated wells. Discard the red-marked mixing wells.
7. Set timer for 5 minutes, mixing the wells for the first 30 seconds of the room temperature incubations by sliding the microwell holder back and forth on a flat surface without splashing reagents from wells.
8. Shake out the contents of the antibody wells. Fill the wells with distilled or deionized water and dump them out. Repeat this step 5 times, then turn the wells upside-down and tap out on a paper towel until the remaining water has been removed.
9. Pour the needed volume of substrate from the green-labeled bottle into the green-labeled reagent boat.
10. With new tips on the 12-channel pipettor, prime and pipette 100 µL of substrate into the wells.
11. Set timer for 5 minutes, mixing the wells for the first 30 seconds by sliding back and forth on a flat surface. Discard remaining substrate and rinse the reagent boat with water.
12. Pour Red Stop solution from the red-labeled bottle into the red-labeled reagent boat.
13. Pipette 100 µL of Red Stop to each well. Mix by sliding back and forth on a flat surface. Discard the tips.
14. Wipe the bottom of the microwells with a dry cloth or towel and read in a microwell reader using a 650 nm filter. Air bubbles should be eliminated, as they could affect analytical results. Results must be read within 20 minutes after the addition of Red Stop.
15. Read and calculate results using Neogen's Stat Fax microwell reader, or equivalent. If using a strip/plate reader, calculate results using Neogen's Veratox software.

## RETESTING

If positives occur in commodities not previously tested, confirm with an additional approved method prior to taking action.

## PERFORMANCE CHARACTERISTICS

Limit of detection: 2.5 ppb

Limit of quantitation: 5 ppb (Described as the lowest concentration point on the calibration curve that this test can reliably detect aflatoxin.)

Range of quantitation: 5–50 ppb.

10–100 ppb quantitation:

Dilute the filtered sample 1:1 with supplied diluent (pink-labeled bottle) using a sample collection cup (i.e. 200 µL filtered sample with 200 µL dilution diluent). Mix well. Run diluted sample following the test procedure. Multiply your result by 2.

30–300 ppb quantitation:

Dilute the filtered sample 1:5 with supplied diluent (pink-labeled bottle) using a sample collection cup (i.e. 100 µL filtered sample with 500 µL dilution diluent). Mix well. Run diluted sample following the test procedure. Multiply your result by 6.

Neogen validated matrixes: Corn and corn gluten meal.

*\*May require a pH adjustment or special extraction procedure.*

**NOTE:** Neogen continues to validate new commodities. Please contact a representative for the latest validated commodity list.

## CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

## SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's Food Safety test kits, at [foodsafety.neogen.com](http://foodsafety.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

## TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit [www.neogen.com/en/terms-and-conditions](http://www.neogen.com/en/terms-and-conditions).

## WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## TEST KITS AVAILABLE FROM NEOGEN

### Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

### Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, milk, multi-treenut, mustard, peanuts, sesame, soy, walnuts

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©Neogen Corporation, 2018. Neogen, Veratox y K-Blue son marcas registradas de Neogen Corporation. Todas las otras marcas y nombres de productos son marcas o marcas registradas de sus respectivas compañías.