Read instructions carefully before starting test

Reveal[®] Q+ MAX

for DON

Quantitative Test

With Aqueous Extraction

THE TOXIN

Deoxynivalenol (DON) most commonly is produced by the pink mold *Fusarium graminearum*. DON, a member of the trichothecene family, is produced by fungi living on cereal commodities such as wheat, corn, barley and ensilages. The toxicological effects attributed to DON include nausea (vomiting), feed refusal, gastroenteritis, diarrhea, immunosuppression and blood disorders.

Pigs have been shown to be highly sensitive to DON. They will refuse to eat feeds when DON levels of ≥ 1 parts per million (ppm) are present. DON and its analogs cause toxic effects in other species as well, with varying degrees of sensitivity. DON has been implicated in causing problems in processed food, including off-flavor in ready-to-eat cereals and adverse effects on dough quality. Accurate determination of the presence of the toxin is of major importance to those monitoring the quality of feed and food in which DON may occur. Testing these commodities for the toxin requires careful sampling, extraction, sanitation and quantitative analysis.

The FDA has issued advisory levels for DON as follows:

For	Level	Commodities
Humans	1 ppm	Finished wheat products (flour, bran and germ)
Ruminating beef, feedlot cattle, and chickens	10 ppm in <50% of diet (5 ppm total diet)	All grains, grain by-products
Swine	5 ppm in <20% of diet (1 ppm total diet)	All grains, grain by-products
All other animals	5 ppm in <40% of diet (2 ppm total diet)	All grains, grain by-products

European Union regulations for DON as follows:

Foodstuffs	Maximum levels (ppb)
Unprocessed cereals other than durum wheat, oats and maize	1250
Unprocessed durum wheat and oats	1750
Cereal flour, including maize flour, maize grits and maize meal	750
Bread, pastries, biscuits, cereal snacks and breakfast cereals	500
Processed cereal-based foods for infants and young children and baby food	200

INTENDED USE/USER

Reveal Q+ MAX for DON is intended for the quantitative analysis of DON levels in grain and grain products, such as wheat, barley, corn and oats. The test kit is designed for use by quality control personnel and others familiar with commodities possibly contaminated by DON.

ASSAY PRINCIPLES

Reveal Q+ MAX for DON is a single-step lateral flow immunochromatographic assay based on a competitive immunoassay format. The extract is wicked through a reagent zone, which contains antibodies specific for DON conjugated to colloidal gold particles. If DON is present, it will be captured by the particle-antibody complex. The DON-antibody-particle complex then is wicked onto a membrane, which contains a zone of DON conjugated to a protein carrier. This zone captures any uncomplexed DON antibody, allowing the particles to concentrate and form a visible line. As the level of DON in a sample increases, free DON will complex with the antibody-gold particles. This, in turn, allows less antibody-gold to be captured in the test zone. Therefore, as the concentration of DON in the sample increases, the test line density decreases. Algorithms programmed into the AccuScan® reader convert these line densities into a quantitative result displayed in parts per million (ppm). The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of DON, ensuring the strip is functioning properly.

STORAGE REQUIREMENTS

Store kit components at room temperature between 18–30°C (64–86°F) to ensure full shelf life. Test strips should remain capped in their original tubes until used to ensure optimal performance.

MATERIALS PROVIDED

1. 25 Reveal Q+ MAX for DON test strips
2. 25 red sample dilution cups
3. 25 clear sample cups
4. 2 bottles of blue-labeled sample diluent
5. 25 MAX 1 Aqueous Extraction Packets
6. Instructions for use

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Extraction Materials
 - a. Sample collection cups with lids 125 mL (Neogen item 9428, 9428B)
 - b. Sample collection tubes with caps (Neogen item 9421, 9421B)
 - c. Filter syringes, Whatman No. 4 filter paper or equivalent (Neogen item 9420, 9429, 9519)**OR**
 - d. Centrifuge, mini (Neogen item 9330)
 - e. Microcentrifuge tubes (Neogen item 9372)
 - f. Dispensing pump or graduated cylinder (Neogen item 9448, 9447)
 - g. Distilled or deionized water
2. Agri-Grind grinder or equivalent (Neogen item 9401,9453)
3. Scale capable of weighing 5–50 ± 0.1 g (Neogen item 9427)
4. Timer (Neogen item 9426, 9452)

5. Reveal sample cup rack (Neogen item 9475)
6. Pipettor, 100 μ L (Neogen item 9860)
7. Pipette tips, 1–200 μ L (Neogen items 9407, 9410, 9417)
8. Pipettor, 100 μ L–1,000 μ L (Neogen item 9463)
9. Pipette tips, 100–1,000 μ L (Neogen items 9464, 9487)
10. AccuScan Gold reader (Neogen item 9595) or AccuScan Pro reader
11. Raptor Integrated Analysis Platform (Neogen item 9680)
 - a. Raptor Cartridges (Neogen item 9681)
 - b. Raptor Exact Volume Transfer Pipettes (Neogen item 9682)
12. MAX 1 Aqueous Extraction packets (Neogen item 8089)
13. MAX 1–G50 Aqueous Extraction packets (Neogen item 8089G)

PRECAUTIONS

1. The test strips must remain inside the stay-dry tube before use.
2. Store test kit at room temperature 18–30°C (64–86°F) when not in use. Do not freeze.
3. Do not use kit components beyond expiration date.
4. Treat all used liquids, including sample extract, and labware as if contaminated with DON.
5. Ensure the device lot number and the curve details match the lot ID number selected on the reader. Failure to update the lot-specific QR code within the AccuScan Gold reader will cause inaccurate results.
6. Commodity extracts should have a pH of 6–8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH contact a Neogen representative or Technical Services.

SAMPLE PREPARATION

The sample to be tested should be collected according to accepted sampling techniques (see GIPSA sampling protocol or contact your Neogen representative). Obtain a representative sample (minimum 100 g). Grind the sample so at least 95% of the ground material passes through a 20 mesh sieve (about the particle size of fine espresso).

SAMPLE EXTRACTION

1. Weigh out 10 ± 0.1 g of sample into extraction cup.
2. Add contents of one (1) MAX 1 aqueous extraction packet to the extraction cup.
3. Add 50 mL distilled or deionized water to the extraction cup.
4. Vigorously shake, using hand or mechanical means, for 3 minutes, or blend for 1 minute.
5. Allow the sample to settle, then filter with a filter syringe or Whatman No. 4 filter paper.
6. You may also pipette 1 mL of sample into a 1.5 mL micro-centrifuge tube, and centrifuge for 30 seconds using a micro-centrifuge.
7. The sample is now ready for testing.

SAMPLE EXTRACTION (GIPSA METHOD)

1. Weigh out $50 \text{ g} \pm 0.1$ g of sample into extraction container.
2. Add contents of one (1) MAX 1-G50 aqueous extraction packet to the extraction cup.
3. Add 250 mL distilled or deionized water to the extraction container.
4. Vigorously shake, using hand or mechanical means, for 3 minutes, or blend for 1 minute.
5. Allow the sample to settle, then filter with a filter syringe or Whatman #4 filter paper or you may also pipette 1 ml of sample into a 1.5 mL micro-centrifuge tube, and centrifuge for 30 seconds using a micro-centrifuge.

ACCUSCAN GOLD READER SET UP

1. Enter the lot specific QR code by selecting **Scan QR** from the main screen. Place the lot specific QR code into the white cartridge adapter labeled Cal/QR and place the cartridge into the reader.
2. The valid code will be scanned by the reader and provide information on the lot number and expiry date. Verify this information is correct and then add the lot ID to the reader by pressing **Add Lot ID**.
NOTE: The lot ID for the current lot will now be stored with the test ID (e.g., aflatoxin, DON) and can be selected when running a test.

ACCUSCAN PRO READER SET UP

1. Enter the lot specific QR code by selecting the QR code icon on the reader. Place the QR code into the cartridge and insert the cartridge into the reader.
NOTE: For instructions on manually entering sample IDs, see the AccuScan Pro user manual.

TEST PROCEDURE – ACCUSCAN READER

1. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary.
2. Add 100 µL of sample extract to each red sample dilution cup.
3. Add 1,000 µL of sample diluent (blue labeled bottle) to the red dilution cup with sample extract. Mix by pipetting up and down 5 times.
4. Using a new pipette tip, transfer 100 µL of diluted sample extract into a new clear sample cup.
5. Place a new Reveal Q+ MAX for DON test strip with the sample end down into the sample cup and set timer for **3 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
6. Remove the strip from the sample cup after it has developed for **3 minutes**.
7. For the AccuScan Gold reader: Select CATEGORY: Mycotoxin Q+MAX
 TEST NAME: Q+M DON ppm

For the AccuScan Pro reader: Select CATEGORY: Mycotoxin Q+MAX
 TEST NAME: DON Q+MAX

Dilution Procedure

Samples **greater than 6 ppm** will need to be diluted and re-tested.

1. Add 100 µL sample filtrate to a sample collection tube.
2. Add 500 µL distilled or deionized water to the sample collection tube. Mix well by pipetting up and down 5 times, or cap the tube and shake.
3. Add 100 µL of diluted sample extract (from step 2) to each red sample dilution cup.
4. Add 1,000 µL of sample diluent (blue labeled bottle) to the red dilution cup with sample extract. Mix by pipetting up and down 5 times.
5. Using a new pipette tip, transfer 100 µL of diluted sample extract into a new clear sample cup.
6. Place a new Reveal Q+ MAX for DON test strip with the sample end down into the clear sample cup and set timer for **3 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
7. Remove the strip from the sample cup after it has developed for **3 minutes**.
8. For the AccuScan Gold reader: Select CATEGORY: Mycotoxin Q+MAX
 TEST NAME: Q+M DON ppm

For the AccuScan Pro reader: Select CATEGORY: Mycotoxin Q+MAX
 TEST NAME: DON Q+MAX

NOTE: The AccuScan reader will not calculate your dilution. Final result displayed will need to be multiplied by 6.


READING TEST RESULTS

NOTE: Test strips should be read within **30 seconds** of completion of the 3 minute incubation. Refer to **AccuScan Reader Set Up** for test selection and set up information.

1. Select the assay type (e.g., DON) from the menu and ensure the device lot number matches the lot ID number selected on the reader.

NOTE: Failure to update the lot-specific QR code will cause inaccurate results.

2. Fully insert the Reveal Q+ test strip into the R-labeled cartridge adapter with the sample end first and results facing out.

3. Insert the cartridge with test strip upside-down into the AccuScan Gold reader (the test lines will face downward into the reader). If using the AccuScan Pro reader, insert the cartridge right side up. The reader will automatically begin analyzing the cartridge.

CAUTION: Removing cartridge prior to completion can result in invalid readings.

4. The AccuScan reader will analyze the test strip and results will be displayed and stored in the reader.


NOTES

1. Ensure device is fully inserted into cartridge.
2. Readings should be made between **3–3.5 minutes**. Readings after 3.5 minutes may be inaccurate due to overdevelopment of the device.

TEST PROCEDURE – RAPTOR INTEGRATED ANALYSIS PLATFORM

1. Place the appropriate number of red sample dilution cups into a sample cup rack. Label cups if necessary.
2. Add 100 μL of sample extract to each red sample dilution cup.
3. Add 1,000 μL of sample diluent to each red sample dilution cup. Mix by pipetting up and down five times.
4. Fully insert a Reveal for DON Q+ MAX test strip into a Raptor cartridge.
5. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor Integrated Analysis Platform reader.
 - a. The bar code on the test strip will be read — the Raptor reader identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader in the front of the Raptor will turn on automatically.
 - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the system.
6. Enter Sample ID if desired.
7. Add 400 μL of sample extract from the red sample dilution cup to the Raptor cartridge.
 - a. The Raptor reader will start automatically.
 - b. Additional samples can be started in the other ports while the first sample is processing.
8. Results will be displayed on the Raptor screen after the 3 minute testing is complete.

DILUTION PROCEDURE – RAPTOR INTEGRATED ANALYSIS PLATFORM

Samples **greater than 6 ppm** will need to be diluted and re-tested.

1. Add 100 μL sample filtrate to a sample collection tube.
2. Add 500 μL distilled or deionized water to the sample collection tube. Mix well by pipetting up and down five times, or cap the tube and shake.
3. Add 100 μL of diluted sample extract (from step 2) to each red sample dilution cup.
4. Add 1,000 μL of sample diluent (blue labeled bottle) to the red dilution cup with sample extract. Mix by pipetting up and down five times.
5. Fully insert a Reveal DON Q+ MAX test strip into a Raptor cartridge.
6. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor reader.
 - a. The bar code on the test strip will be read — the Raptor reader identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader in the front of the Raptor will turn on automatically.
 - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the system.
7. Enter Sample ID if desired.
8. Add 400 μL of sample extract from the red sample dilution cup to the Raptor cartridge.
 - a. The Raptor reader will start automatically.
 - b. Additional samples can be started in the other ports while the first sample is processing.
9. Results will be displayed on the Raptor reader's screen after the three-minute testing is complete.

NOTE: The Raptor reader will not calculate your dilution. Final result displayed will need to be multiplied by 6.

PERFORMANCE CHARACTERISTICS

1. Limit of detection: 0.3 ppm
2. Range of quantitation: 0.3–30 ppm (with dilution)

NOTE: All samples over 6 ppm must be diluted and retested. See protocol on page 4. Samples below the limit of quantitation should be reported as less than 0.3 ppm.

VALIDATED MATRICES

Alfalfa	Malted barley ¹	Sweet potato, powder
Barley ¹	Millet	Sweet potato, pellets
Beet pulp	Milo (grain sorghum)	Triticale
Brewer's rice	Navy bean	Wheat midds ¹
Broken rice	Oatmeal	Wheat bran
Brown rice	Oats ¹	Wheat ¹
Buckwheat	Oat fiber	Wheat flour ¹
Buckwheat groats*	Oat groats	Wheat germ
Canary seed	Pearled barley ¹	Wheat gluten
Canned pet food	Peas	Wheat, waxy
Chickpea	Pea flower	White sesame
Corn gluten feed*	Pea protein	Whole rye
Corn gluten meal ¹	Peas, yellow	
Corn oil	Rapeseed	
Corn silage*	Quinoa	
Corn ¹	Potato protein	
DDGS ¹ *	Red spring wheat	
DDG wet cake*	Potato with skins powder*	
DDG syrup*	Red dog flour	
Defatted rice bran	Rolled oats	
Dried blueberries	Rough rice ¹	
Dried peas	Rice bran	
Flaxseed*	Rice flour	
Grass silage*	Safflower	
Green lentils	Soybean meal	
Haylage*	Soybean	
Kamut (Khorasan wheat)	Soy flour	
Kiwicha	Soy hulls	
Lentil fiber	Spelt	
Lentils	Sunflower nugget	
Linseed, brown	Sunflower meal	
Linseed, gold	Sunflower seeds	

¹ = Validated by USDA-GIPSA

*pH adjustment step may be necessary

06-2017

NOTE: Neogen continues to validate new commodities. Please contact a representative for the latest validated commodity list.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/en/terms-and-conditions.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TESTING KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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Neogen Food and Animal Security

+91 484 2306598, 2301582
info@neogenindia.com
www.neogenindia.com

©Neogen Corporation, 2018. Neogen, AccuScan, Raptor and Reveal are registered trademarks of Neogen Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.

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