Read instructions carefully before starting test

Reveal[®] Q+ MAX

for Aflatoxin Quantitative Test With Aqueous Extraction

FGIS 2016-087

THE TOXIN

Aflatoxin is a toxic and carcinogenic substance produced by certain strains of the molds *Aspergillus flavus* and *A. parasiticus*. There are four principle types of aflatoxin: B₁, B₂, G₁ and G₂. Aflatoxin B₁ is the most frequently encountered of the group and the most toxic. The commodities most affected by aflatoxin are corn, peanuts, cottonseed, milo, and the majority of tree nuts.

The effects on animals ingesting excessive amounts of aflatoxin range from chronic health and performance problems to death. Aflatoxin has been shown to cause liver damage or cancer, decreased milk and egg production, immune suppression, and interference with reproductive efficiency.

Many countries have set maximum allowable levels of aflatoxin in food and feed. Therefore, accurate determination of the toxin's presence is of major importance to those monitoring the quality of food and feed in which aflatoxin may occur. Testing these commodities for the toxin requires careful sampling, chemical extraction, quality assurance, and quantitative analysis.

The FDA has issued regulatory levels for aflatoxin as follows:

For	Level	Commodities
Humans	20 ppb	All food except milk
All animal species	20 ppb	All feed (exceptions below)
Exceptions:		
Breeding cattle, breeding swine, mature poultry	100 ppb	Corn
Finishing swine (> 100 lbs.)	200 ppb	Corn
Finishing beef cattle	300 ppb	Corn
Finishing beef cattle, swine, poultry	300 ppb	Cottonseed meal

The European Union has issued regulations for aflatoxin (total B₁, B₂, G₁, and G₂) as follows (as referenced in the Commission Regulation (EC) No. 1881/2006):

Foodstuffs	Levels
Groundnuts (peanuts) and other oilseeds, subject to sorting and treatment before human consumption or use as an ingredient in foodstuffs; Almonds, pistachio, hazelnuts, Brazil nuts and apricot kernels subject to sorting and treatment before human consumption or use as an ingredient in foodstuffs	15 ppb
Tree nuts, other than those listed above subject to sorting and treatment before human consumption or use as an ingredient in foodstuffs; Almonds, pistachio, hazelnuts, Brazil nuts and apricot kernels for direct human consumption or use as an ingredient in foodstuffs; Maize, rice or dried fruits subject to sorting and treatment before human consumption or as an ingredient in foodstuffs	10 ppb
Groundnuts (peanuts) and other oilseeds for direct human consumption or use as an ingredient in foodstuffs; Dried fruit and processed products thereof, for direct human consumption or use as an ingredient in foodstuffs	4 ppb

INTENDED USE/USER

Reveal Q+ MAX for Aflatoxin is intended for the quantitative analysis of corn and corn products for aflatoxin. Testing can be complete in 6 minutes after extraction and can detect a wide range of aflatoxin. This test utilizes an aqueous extraction thus eliminating the need for the disposal of hazardous waste. The test kit is designed for use by quality control personnel and others familiar with commodities possibly contaminated with aflatoxin.

Corn, peanuts, pistachios, almonds, and walnuts were validated by AOAC Research Institute Performance Tested Method. All other matrices were validated internally by Neogen.

ASSAY PRINCIPLES

Reveal Q+ MAX for Aflatoxin is a single-step lateral flow immunochromatographic assay based on a competitive immunoassay format. The extract is wicked through a reagent zone, which contains antibodies specific for aflatoxin conjugated to colloidal gold particles. If aflatoxin is present, it will be captured by the particle-antibody complex. The aflatoxin-antibody-particle complex then is wicked onto a membrane, which contains a zone of aflatoxin conjugated to a protein carrier. This zone captures any uncomplexed aflatoxin antibody, allowing the particles to concentrate and form a visible line. As the level of aflatoxin in a sample increases, free aflatoxin will complex with the antibody-gold particles. This allows less antibody-gold to be captured in the test zone. Therefore, as the concentration of aflatoxin in the sample increases, the test line density decreases. Algorithms programmed into the AccuScan® readers convert these line densities into a quantitative result displayed in parts per billion (ppb). The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of aflatoxin, ensuring the strip is functioning properly.

STORAGE REQUIREMENTS

Store kit components at room temperature 18–30°C (64–86°F) to ensure full shelf life. Test strips should remain capped in their original tubes until used to ensure optimal performance.

MATERIALS PROVIDED

1. 25 Reveal Q+ MAX for Aflatoxin test strips
2. 25 red sample dilution cups
3. 25 clear sample cups
4. 1 bottle of sample diluent
5. 25 MAX 1 Aqueous Extraction packets
6. Kit instructions for use

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Extraction Materials
 - a. Sample collection cups with lids 125 mL (Neogen item 9428, 9428B)
 - b. Sample collection tubes with caps (Neogen item 9421, 9421B)
 - c. Filter syringes, Whatman #4 filter paper or equivalent (Neogen item 9420, 9429, 9519)
OR
 - d. Centrifuge, mini (Neogen item 9330)
 - e. Microcentrifuge tubes (Neogen item 9172)
 - f. Dispensing pump or graduated cylinder (Neogen item 9448, 9447)
2. Agri-Grind grinder or equivalent (Neogen item 9401,9453)
3. Scale capable of weighing 5–50 grams \pm 0.1 g (Neogen item 9427)
4. Timer (Neogen item 9426,9452)
5. Reveal sample cup rack (Neogen item 9475)
6. Pipettor, 100 μ L (Neogen item 9860, 9290)
7. Pipette tips, 1-200 μ L (Neogen item 9407, 9410, 9417)
8. Distilled or deionized water
9. AccuScan Pro (Neogen item 9565), or AccuScan Gold (Neogen item 9595) or Raptor Integrated Analysis Platform (Neogen item 9680)
10. Raptor cartridges (Neogen item 9681)
11. Pipettor, 250 μ L (Neogen item 9692)
12. Pipette tips, 100 μ L–1 mL (Neogen item 9464,9487)
13. Pipettor, 400 μ L fixed (Neogen item 9693)
14. Pipette tips, 200–1000 μ L (Neogen item 9293,9292)
15. Raptor 400 μ L exact Volume pipettes (Neogen item 9682)
16. MAX 1 Aqueous Extraction packets (Neogen item 8089)
17. MAX 1–G50 Aqueous Extraction packets (Neogen item 8089G)

PRECAUTIONS

1. The test strips must remain inside the stay-dry tube before use.
2. Store test kit at room temperature 18–30°C (64–86°F) when not in use. Do not freeze.
3. Do not use kit components beyond expiration date.
4. Treat all used liquids, including sample extract, and labware as if contaminated with aflatoxin. Consult the appropriate safety data sheets (SDS) for proper handling of aflatoxin.
5. Ensure the device lot number and the curve details match the lot ID number selected on the reader. Failure to update the lot-specific QR code within the AccuScan reader will cause inaccurate results.
6. Commodity extracts should have a pH of 6–8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH contact a Neogen representative or Technical Services.

ACCUSCAN GOLD READER SETUP

1. Enter the lot specific QR code by selecting **Scan QR** from the main screen. Place the lot specific QR code into the white cartridge adapter labeled Cal/QR and place the cartridge into the reader.
2. The valid code will be scanned by the reader and provide information on the lot number and expiry date. Verify this information is correct and then add the lot ID to the reader by pressing **Add Lot ID**.

NOTE: The lot ID for the current lot will now be stored with the test ID (e.g., aflatoxin) and can be selected when running a test.

ACCUSCAN PRO READER SETUP

1. Enter the lot specific QR code by selecting the QR code icon on the reader. Place the QR code into the cartridge and insert the cartridge into the reader.

NOTE: For instructions on manually entering sample IDs, see the AccuScan Pro user manual.

SAMPLE PREPARATION

The sample to be tested should be collected according to accepted sampling techniques (see FGIS sampling protocol or contact your Neogen representative). Obtain a representative sample (minimum 100 gm). Grind the sample so at least 95% of the ground material passes through a 20 mesh sieve (about the particle size of fine espresso).

SAMPLE EXTRACTION

1. Weigh out 10 g \pm 0.1 g of sample into extraction cup.
2. Add contents of one (1) MAX 1 aqueous extraction packet to the extraction cup.
3. Add 50 mL distilled or deionized water to the extraction cup.
4. Vigorously shake, using hand or mechanical means, for 3 minutes, or blend for 1 minute.
5. Allow the sample to settle, then filter with a filter syringe or Whatman #4 filter paper to collect a minimum of 3 mL filtrate into a sample collection tube. Or, you may also pipette 1 mL of sample into a 2.0 mL microcentrifuge tube, and centrifuge for 1 minute using a microcentrifuge (approx. 2,000 x g).

SAMPLE EXTRACTION (FGIS METHOD)

1. Weigh out 50 g \pm 0.1 g of sample into extraction container.
2. Add contents of one (1) MAX 1-G50 aqueous extraction packet to the extraction cup.
3. Add 250 mL distilled or deionized water to the extraction container.
4. Vigorously shake, using hand or mechanical means, for 3 minutes.
5. Allow the sample to settle, then filter with a filter syringe or Whatman #4 filter paper to collect a minimum of 3 mL filtrate into a sample collection tube. Or, you may also pipette 1 mL of sample into a 2.0 mL micro-centrifuge tube, and centrifuge for 1 minute using a micro-centrifuge (approx. 2,000 x g).

TEST PROCEDURE

1. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary.
2. Add 100 μ L of sample diluent to each **red** sample dilution cup.
3. Add 100 μ L of sample extract to a red dilution cup with sample diluent. Mix by pipetting up and down 5 times.
4. Transfer 100 μ L of diluted sample extract into a new **clear** sample cup.
5. Place a new Reveal Q+ MAX for Aflatoxin test strip with the sample end down into the clear sample cup and set timer for **6 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
6. Remove the strip from the sample cup after it has developed for **6 minutes**.
7. **For the AccuScan Gold reader: Select CATEGORY: Mycotoxin Q+ MAX
TEST NAME: Q+M Afla**

**For the AccuScan Pro reader: Select CATEGORY: Mycotoxin Q+ MAX
TEST NAME: Afla Q+ MAX**

DILUTION PROCEDURE

Samples greater than 50 ppb will need to be diluted and retested.

1. Add 100 μ L sample filtrate to a sample collection tube. Add 500 μ L distilled or deionized water. Mix well.
2. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary.
3. Add 100 μ L of sample diluent to each **red** sample dilution cup.
4. Add 100 μ L of sample (from step 1) to a red dilution cup with sample diluent. Mix by pipetting up and down 5 times.
5. Transfer 100 μ L of sample extract from the red dilution cup into a **clear** sample cup.
6. Place a new Reveal Q+ MAX for Aflatoxin test strip with the sample end down into the clear sample cup and set timer for **6 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
7. Remove the strip from the sample cup after it has developed for **6 minutes**.
8. **For the AccuScan Gold reader: Select CATEGORY: Mycotoxin Q+ MAX
TEST NAME: Q+M Afla**

Final result displayed will need to be multiplied by 6.

**For the AccuScan Pro reader: Select CATEGORY: Mycotoxin Q+ MAX
TEST NAME: Afla Q+ MAX**

Final result displayed will need to be multiplied by 6.

READING TEST RESULTS

NOTE: Test strips should be read within **30 seconds** of completion of the 6 minute incubation.

1. Select the assay type (e.g., aflatoxin) from the menu and ensure the device lot number matches the lot ID number selected on the reader. **NOTE:** Failure to update the lot-specific QR code will cause inaccurate results.
2. Fully insert the Reveal Q+ test strip into the R-labeled cartridge adapter with the sample end first and results facing out.

3. Insert the cartridge with test strip upside-down into the AccuScan Gold reader (the test lines will face downward into the reader). If using the AccuScan Pro reader, insert the cartridge right side up. The reader will automatically begin analyzing the cartridge.

CAUTION: Removing cartridge prior to completion can result in invalid readings.

4. The AccuScan reader will analyze the test strip and results will be displayed and stored in the reader.


NOTES

1. Ensure device is fully inserted into cartridge.
2. Readings should be made between **6 and 6.5 minutes**. Readings after 6.5 minutes may be inaccurate due to overdevelopment of the device.
3. The strips must be read using Neogen's AccuScan Gold or AccuScan Pro readers.

TEST PROCEDURE — RAPTOR INTEGRATED ANALYSIS PLATFORM

1. Place the appropriate number of red sample dilution cups into a sample cup rack. Label cups if necessary.
2. Add 250 μL of sample extract to each red sample dilution cup.
3. Add 250 μL of sample diluent to red sample dilution cup. Mix by pipetting up and down five times.
4. Fully insert a Reveal Q+ MAX for Aflatoxin test strip into a Raptor cartridge.
5. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor.
 - a. The bar code on the test strip will be read. The system identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader in the front of the Raptor will turn on automatically.
 - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the reader.
6. Enter Sample ID if desired.
7. Add 400 μL of sample extract from the red sample dilution cup to the Raptor cartridge.
 - a. The Raptor system will start automatically.
 - b. Additional samples can be started in the other ports while the first sample is processing.
8. Results will be displayed on the Raptor screen after the 6-minute testing period is complete.

DILUTION PROCEDURE - RAPTOR INTEGRATED ANALYSIS PLATFORM

Samples **greater than 50 ppb** will need to be diluted and re-tested.

1. Add 100 μL sample filtrate to a sample collection tube.
2. Add 500 μL distilled or deionized water to the sample collection tube. Mix well by pipetting up and down 5 times, or cap the tube and shake.
3. Add 250 μL of diluted sample extract (from step 2) to each red sample dilution cup.
4. Add 250 μL of sample diluent (blue labeled bottle) to the red dilution cup with sample extract. Mix by pipetting up and down five times.
5. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor.
 - a. The bar code on the test strip will be read. The system identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader in the front of the Raptor will turn on automatically.
 - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the reader.
6. Enter Sample ID if desired.
7. Add 400 μL of sample extract from the red sample dilution cup to the Raptor cartridge.
 - a. The Raptor system will start automatically.
 - b. Additional samples can be started in the other ports while the first sample is processing.
8. Results will be displayed on the Raptor screen after the 6-minute testing period is complete.

NOTE: The Raptor reader will not calculate your dilution. **Final result displayed will need to be multiplied by 6.**

PERFORMANCE CHARACTERISTICS

1. Limit of detection: 2 ppb
2. Range of detection: 3–300 ppb (with dilution).
Samples greater than 50 ppb need to be diluted and re-tested following the dilution procedure on page 5.

VALIDATED MATRICES

Barley, brown rice, cashews, corn (including dent or field corn, corn flour, cracked corn, corn grits or polenta and corn screenings), cottonseed, pistachio, popcorn, sorghum, wheat, walnuts.

1:6 extraction – almond, peanuts, peanut butter (Example: 10 g matrix and 60 g distilled or deionized water. Add one MAX 1 packet.)

1:7 extraction – white rice (Example: 10 g matrix and 70 g distilled or deionized water. Add one MAX 1 packet.)

NOTE: Neogen continues to validate new commodities. Please contact a representative for the latest validated commodity list.

AOAC Performance Tested Method (PTM) Validated Matrices: corn, almonds, pistachios, walnuts, peanuts.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/en/terms-and-conditions.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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