Read instructions carefully before starting test

Reveal[®] Q+ MAX **for Ergot Alkaloids**

Quantitative Test with Aqueous Extraction

Refrigerate at 2–8°C • DO NOT FREEZE

THE TOXIN

Ergot alkaloids are toxic and carcinogenic substances. The six ergot alkaloids studied (ergometrine, ergotamine, ergosine, ergocristine, ergocryptine and ergocornine) have been defined by the European Food Safety Authority (EFSA) as being of major importance to human and animal health.

Throughout history, there have been several deadly episodes of food poisoning by ergot alkaloids, known as St. Anthony's fire or ergotism. Ergot alkaloids are secondary metabolites produced by fungi of the *Claviceps* genus (chiefly *Claviceps purpurea*), which are common pathogens of cereals and pasture grasses.

The best protection against these mycotoxins is monitoring for their presence in feeds and foods. That means testing all along the pathway from initial harvest of grains to the finished product.

INTENDED USE/USER

Reveal[®] Q+ MAX for Ergot Alkaloids is intended for the quantitative analysis of coarse ground grains of wheat and rye for ergot alkaloids. Testing can be completed in 8 minutes after extraction and can detect within two ranges, 50–500 parts per billion (ppb) and 500–5,000 ppb, of ergot alkaloids depending on regional regulations. This test uses an aqueous extraction, thus eliminating the need for the disposal of hazardous waste. The test kit is designed for use by quality control personnel and others familiar with commodities possibly contaminated with ergot alkaloids.

ASSAY PRINCIPLES

Reveal Q+ MAX for Ergot Alkaloids is a single-step lateral flow immunochromatographic assay based on a competitive immunoassay format. The extract is wicked through a reagent zone, which contains antibodies specific for ergot alkaloids conjugated to colloidal gold particles. If ergot alkaloids are present, they will be captured by the particle-antibody complex. The ergot alkaloid -antibody-particle complex is then wicked onto a membrane, which contains a zone of ergot alkaloids conjugated to a protein carrier. This zone captures any uncomplexed ergot alkaloid antibody, allowing the particles to concentrate and form a visible line. As the level of ergot alkaloids in a sample increases, free ergot alkaloids will complex with the antibody-gold particles. This allows less antibody-gold to be captured in the test zone. Therefore, as the concentration of ergot alkaloids in the sample increases, the test line density decreases.

Algorithms programmed into NEOGEN's Raptor® platform convert these line densities into a quantitative result displayed in parts per billion (ppb). The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form, regardless of the presence of ergot alkaloids, ensuring the strip is functioning properly.

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). Test strips should remain in the original tube until used to ensure optimal performance.

MATERIALS PROVIDED

1. 25 Reveal Q+ MAX for Ergot Alkaloids test strips
2. 50 red sample dilution cups
3. 2 bottles of sample diluent
4. 25 MAX 1 Aqueous Extraction packets
5. Kit instructions for use

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Extraction materials
 - a. Sample collection cups with lids, 125 mL (NEOGEN item 9428, 9428B)
 - b. Sample collection tubes with caps (NEOGEN item 9421, 9421B)
 - c. Centrifuge, mini (NEOGEN item 9330)
 - d. Microcentrifuge tubes (NEOGEN item 9172)
 - e. Dispensing pump or graduated cylinder (NEOGEN item 9448, 9447)
2. Agri-Grind grinder or equivalent (NEOGEN item 9401, 9453)
3. Scale capable of weighing 5–50 grams ± 0.1 g (NEOGEN item 9427)
4. Timer (NEOGEN item 9426, 9452)
5. Reveal sample cup rack (NEOGEN item 9475)
6. Pipettor, adjustable 100–1000 µL
7. Pipette tips, 100–1000 µL (NEOGEN item 9464, 9487)
8. MAX 1 Aqueous Extraction packets (NEOGEN item 8089)
9. Raptor Integrated Analysis Platform (NEOGEN item 9680)
10. Raptor cartridges (NEOGEN item 9681)

PRECAUTIONS

1. The test strips must remain inside the stay-dry tube before use.
2. Store test kit at refrigerated temperature 2–8°C (35–46°F)
3. Kit should be brought to room temperature 18–30°C (64–86°F) for at least **30 minutes** prior to use.
4. Do not use kit components beyond expiration date.
5. Commodity extracts should have a pH of 6–8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH, contact a NEOGEN representative or Technical Services.
6. Treat all used liquids, including sample extract and labware, as if contaminated with ergot alkaloids.
7. MAX 1 packets may be stored at room temperature 18–30°C (64–86°F) if desired.

SAMPLE PREPARATION


The sample to be tested should be collected according to accepted sampling techniques (see FGIS sampling protocol or contact your NEOGEN representative). Obtain a representative sample (minimum 100 g). Grind the sample so at least 95% of the ground material passes through a 20-mesh sieve.

SAMPLE EXTRACTION

1. Weigh out 10 g ± 0.1 g of sample into an extraction cup.
2. Add contents of one (1) MAX 1 aqueous extraction packet to the extraction cup.
3. Add 50 mL distilled or deionized water to the extraction cup.
4. Vigorously shake, using hand or mechanical means, for **3 minutes**.
5. Allow the sample to settle for a minimum of 2 minutes.
6. Pipette 1 mL of sample into a 2 mL microcentrifuge tube, and centrifuge for **1 minute**.

TEST PROCEDURE RAPTOR INTEGRATED ANALYSIS PLATFORM – CURVE 1 (50–500 PPB)

1. Place the appropriate number of red sample dilution cups into a sample cup rack. Label cups if necessary.
2. Add 100 µL of sample extract to a red sample dilution cup.
3. Add 700 µL of sample diluent to the red dilution cup with sample extract. Mix by pipetting up and down 5 times.
4. Fully insert a Reveal Q+ MAX for Ergot Alkaloid test strip into a Raptor cartridge.
5. Insert the Raptor cartridge containing the test strip into the Raptor port.
 - a. The bar code on the test strip will be read. The system identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader in the front of the Raptor will turn on automatically.
 - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the reader.
6. Enter a Sample ID if desired.
7. Select CURVE 1 for 50–500 ppb range.

8. Add 400 μL of sample extract from the red sample dilution cup to the Raptor cartridge.
 - a. The Raptor system will start automatically.
 - b. Additional samples can be started in the other ports while the first sample is processing.
9. Results will be displayed on the Raptor screen after the 8-minute testing period is complete.

NOTE: Samples greater than 500 ppb (run on Curve 1) should be re-tested following test protocol for Curve 2.

TEST PROCEDURE RAPTOR INTEGRATED ANALYSIS PLATFORM – CURVE 2 (500–5,000 PPB)

1. Place the appropriate number of red sample dilution cups into a sample cup rack. Label cups if necessary.
2. Add 100 μL of sample extract to a red sample dilution cup.
3. Add 700 μL of sample diluent to the red dilution cup with sample extract. Mix by pipetting up and down 5 times.
4. In a new red dilution cup, add 200 μL of diluted sample (from step 3). then add 800 μL of sample diluent to the cup and mix by pipetting up and down 5 times.
5. Fully insert a Reveal Q+ MAX for Ergot Alkaloid test strip into a Raptor cartridge.
6. Insert the Raptor cartridge containing the test strip into the Raptor port.
 - a. The bar code on the test strip will be read. The system identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader in the front of the Raptor will turn on automatically.
 - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the reader.
7. Enter a Sample ID if desired.
8. Select CURVE 2 for 500–5,000 ppb range.
9. Add 400 μL of sample extract (from Step 4) from the red sample dilution cup to the Raptor cartridge.
 - a. The Raptor system will start automatically.
 - b. Additional samples can be started in the other ports while the first sample is processing.
10. Results will be displayed on the Raptor screen after the 8-minute testing period is complete.

NOTE: Samples greater than 5000 ppb (run on Curve 2) should be diluted and re-tested following test protocol for Curve 2.

PERFORMANCE CHARACTERISTICS

Limit of detection: 40 ppb

Range of quantitation: Curve 1: 50–500 ppb
Curve 2: 500–5,000 ppb

NOTE: Samples greater than 500 ppb (run on Curve 1) should be re-tested following test protocol for Curve 2. Samples greater than 5000 ppb (run on Curve 2) should be diluted and re-tested following test protocol for Curve 2.

VALIDATED MATRICES

Ground grains of wheat and rye

CUSTOMER SERVICE

NEOGEN Customer Assistance and Technical Service can be reached by using the contact information on the back of this booklet. Training on this product, and all NEOGEN test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of NEOGEN's test kits, on NEOGEN's Web site at foodsafety.neogen.com, or by calling NEOGEN at 800.234.5333 or 517.372.9200.

TERMS AND CONDITIONS

Please visit www.NEOGEN.com/en/terms-and-conditions for NEOGEN's full terms and conditions.

WARRANTY

NEOGEN Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TESTING KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, coconut, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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