BetaStar[®]

Advanced

for Tetracyclines

INTRODUCTION

BetaStar[®] Advanced for Tetracyclines is an assay for the rapid detection of tetracycline, chlortetracycline and oxytetracycline in raw, commingled cow milk at or below established tolerance and/or safe levels.

REACTION MECHANISM

BetaStar Advanced for Tetracyclines is a single step lateral flow immunochromatographic assay based on a competitive immunoassay format. The milk is wicked through a reagent zone, which contains antibodies and receptors conjugated to colloidal gold particles. If tetracycline, chlortetracycline or oxytetracycline are present, they will be captured by their particle-antibody complex. The drug-antibody-particle complex then is wicked onto a membrane. The membrane contains one test line, consisting of a drug conjugated to a protein carrier designed to capture the tetracyclines. This test line captures any antibody gold particle complexes (in absence of antibiotics) to concentrate and form a visible line.

As the level of a particular drug in the sample increases, free drug molecules will complex with the antibody gold particles. This allows less antibody-gold to be captured into the test lines. Therefore, as the concentration of a drug in the sample increases, the test line density decreases. The membrane also contains a control line where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of the drug, ensuring the strip is functioning properly. This assay is meant to be performed using the Raptor[®] Integrated Analysis System.

MATERIALS REQUIRED

The BetaStar Advanced for Tetracyclines kit contains:

1. 25 test strips for detecting tetracyclines stored in tubes with a desiccant in the lid
2. 25 disposable 400 μ L poly-pipettes

MATERIALS REQUIRED BUT NOT PROVIDED

1. Raptor Integrated Analysis System (Neogen item number: 9680)
2. Package of 50 Raptor cartridges for test strips (Neogen item number: 9681)
3. Raptor Verification Standards (cartridges containing test strip standards) (Neogen item number: 9685)
4. OPTIONAL: Seiko Instruments Inc. (SII), Smart Label Printer 620 (Neogen item number: 9566). Print-out required for NCIMS testing.
5. Tetracycline Positive Control (Oxytetracycline at 300 ± 30 ppb) (Neogen item number: TET010)

TEST PREPARATION

1. BetaStar Advanced for Tetracyclines is designed for storage under normal ambient environmental conditions: 18–30°C (64–86°F).
2. One test strip will need to be loaded into a single cartridge prior to performing the assay.
3. The test strip may be loaded into any of the three cartridge positions.
4. Remove only the number of test strips and cartridges you intend to use.
5. Test strips and cartridges that have been removed from the kit must be kept clean and dry for optimum performance.

POSITIVE AND NEGATIVE CONTROLS FOR NCIMS TESTING IN THE U.S.

The performance of the test kit must be validated with positive and negative controls each day that the test kit is used and with each new lot of reagents. The controls must be run before milk samples are evaluated.

Negative control: For raw commingled cow milk or previously frozen milk that tested at a 1.15 or higher ratio reading with BetaStar Advanced for Tetracyclines. Must be used within 72 hours of testing, and temperature maintained at 0–4.4°C.

Positive control: Use the oxytetracycline positive control available from Neogen (item #TET010). Follow the insert instructions for preparation of the material. Test for suitability each time prepared must produce 0.85 or lower ratio reading. Store prepared controls at 0–4.4°C up to 48 hours.

PERFORMANCE INFORMATION – SENSITIVITY

Dose response information

BetaStar Advanced for Tetracyclines Performance Information

| ppb | Oxytetracycline | Chlortetracycline | Tetracycline |
|-----|-----------------|-------------------|--------------|
| 0 | 0% | 0% | 0% |
| 50 | 0% | 0% | 0% |
| 100 | 0% | 0% | 0% |
| 150 | 40% | 0% | 17% |
| 200 | 100% | 13% | 93% |
| 250 | 100% | 87% | 100% |
| 300 | 100% | 100% | 100% |

| Tetracyclines Detected | Tetracycline | Chlortetracycline | Oxytetracycline |
|------------------------|--------------|-------------------|-----------------|
| Tolerance (ppb) * | 300 | 300 | 300 |
| Calculated 90/95 ** | 245 ppb | 254 ppb | 190 ppb |

* The tolerance level for each of the tetracycline antibiotics has been set by the FDA as the sum of all tetracyclines.

** The 90/95 sensitivity level was calculated from data collected during validation testing by a third party laboratory. Sensitivity was based on testing multiple levels of each tetracycline antibiotic at 30 samples per concentration.

PERFORMANCE INFORMATION — SELECTIVITY

The BetaStar Advanced for Tetracyclines test does not cross-react with the following drugs at levels of 100 ppb: amoxicillin, ampicillin, ceftiofur, cephapirin, cloxacillin, penicillin, sulfadiazine, sulfanilamide, sulfathiazole, sulfamethazine, sulfapyridine, sulfadimethoxine, gentamicin, neomycin, streptomycin, ivermectin, erythromycin, pirlimycin, tilmicosin, novobiocin, furosemide, trichlormethiazide, thiabendazole, chlorothiazide, oxytocin, phenylbutazone, dexamethasone, dipyrone, flunixin, enrofloxacin, para-aminobenzoic acid (PABA) and florphenicol.

TRAINING

The BetaStar Advanced for Tetracyclines test is recommended for use by personnel who have received training by a Neogen representative. In the United States, under NCIMS recommendations, trained individuals need to maintain proficiency by regular use and/or state sponsored training/certification programs. Individuals needing training should contact Neogen Technical Services at 800/234-5333.

OPERATING INSTRUCTIONS

For all NCIMS testing, refer to the current 2400 form for BetaStar Advanced for Tetracyclines (Appendix N Bulk Milk Tanker Screening Form).

Sample Information

- Raw, commingled cow milk must be used between 0°C and 4.4°C.
- Samples must be refrigerated and testing must be completed within 72 hours of sampling.
- Thoroughly mix the sample before testing.
- For NCIMS purposes, milk samples cannot be frozen at any time. The testing must be completed within 72 hours of sampling. Completion includes the time from the initial test, any necessary retesting of the initial positive, confirmation and/or producer traceback testing.

Precautions and Warnings

Do not use kits past their expiration date.

DAILY PERFORMANCE AND OPERATION CHECK

1. Reader calibration occurs automatically when the Raptor Integrated Analysis Platform is powered on. If the calibration is unsuccessful, the reader will not operate. Calibration must be performed on each day of use.
2. Reader calibration checks should be done daily with the calibration check cartridges provided by Neogen. The positive cartridge must have an average ratio value of less than 0.85 and negative cartridge must have an average ratio value of greater than 1.15; if they do not meet specifications, contact Neogen.
3. Run a positive control and a negative control before use on each new lot of kits. The results must be appropriate for the control. The positive control must have a ratio value of less than 0.85 and negative control must have a ratio value greater than 1.15; if they do not meet specifications, contact Neogen.
4. Run a positive control and a negative control daily to assure appropriate results are achieved. If results are incorrect, the control samples should be re-run. If the problem persists, discontinue testing and contact Neogen for assistance.

TEST PROCEDURE

1. Insert the cartridge containing the test strip into any of the three ports within the Raptor Integrated Analysis System.
 - a. The bar code on the test device will be read. If the QR code for the lot of devices has not been entered into the Raptor, the bar code reader in the front of the Raptor will turn on automatically. Scan the QR code found on the container storing the test devices.
 - b. The assay temperature for the BetaStar Advanced for Tetracyclines is 65°C. This temperature is programmed in the Raptor Integrated Analysis System.
 - c. When a cartridge is inserted into one of the 3 ports, the port will automatically begin to adjust to the proper temperature.
 - d. The user will not be able to progress to the next step until the incubator reaches 65°C.
2. Enter the sample ID by either scanning or using the on-screen keyboard.
3. Mix the milk sample prior to adding the sample into the cartridge.
4. Pipette a 400 µL milk sample into the sample port at the back of the cartridge. Press “Next.” The analysis will occur automatically.
5. Insert the used pipette into the sample port located in the back of the cartridge. This will prevent double loading the same sample or loading a second sample into the same cartridge.
6. After the incubation, the results will appear on the reader screen. In addition, results may be printed using SII Smart Label Printer 620 connected to the Raptor.

TEST INTERPRETATION

If the ratio of test to control line intensities as determined by the reader is ≥ 1.0 , the test is **negative**. If the ratio is < 1.0 , the test is **positive**. If the test is invalid, the reader will indicate “invalid.” Re-test the control/sample. *NOTE: Visual interpretation of the result is not allowed under the NCIMS Milk Regulatory Programs, Appendix N, PMO.*

The reader printout is formatted to provide test information as follows:

Reader ID: Reader serial number

User Name: Operator of the test

Time & Date: Time and date of test performed

Sample ID: Unique identifier for sample

Temperature: 65.0°C

Port: One of three possible ports the test can be run within

Test A: Test description of the test in the first position within the cartridge (ex: BetaStar Adv. for Tet)

Lot: Six digit lot number for the test

General Result: Positive or negative result for the sample. If at least one of up to three devices is positive, the general result will be positive.

Ratio: Numerical value which determines positive or negative results

Result: Positive or negative result for specific test or test line

FOR NCIMS TESTING IN THE U.S.

- The original test result, if **negative**, is reported as **Not Found**.
- The original test result, if **positive**, is considered an **Initial Positive**.

RE-TEST OF INITIAL POSITIVE

1. For milk samples yielding an Initial Positive result, promptly re-test the same sample in duplicate, along with positive and negative controls.
2. Providing controls are valid, if both duplicate assays are negative, this indicates a negative result (reported as **Not Found**).
3. Providing both controls are valid, if one or both duplicate samples test positive, the sample is a “**Tetracycline Presumptive Positive Test.**” In accordance with Appendix N of the PMO, report all presumptive positives to the appropriate regulatory agency.
4. If either control is invalid, repeat the testing. If either control is invalid again, contact Neogen for technical assistance at 800/234-5333.

NCIMS TESTING

All NCIMS testing must be performed in accordance with Appendix N of the current PMO. All requirements of the PMO must be followed for all additional testing. Milk testing must be completed within 72 hours of sampling, including the time from the initial test, any necessary retesting of an initial positive, confirmation and/or producer trace-back testing.

REMARKS

1. If performing only one test, the cartridge may be inserted into any one of the three ports. Similarly, the test strip may be placed into any of the three slots found on the cartridge.
2. A single milk sample can be tested for both tetracycline and beta-lactam antibiotics at the same time using the Raptor Integrated Analysis System. In order to test for tetracyclines alongside beta-lactam antibiotics, the tetracycline assay test strips provided in this kit could be inserted into either of the two remaining cartridge positions available in the test strip cartridge provided in the BetaStar Advanced for Beta-lactams test kit (purchased separately).
3. If the test sample does not migrate on the strip, the test is invalid. This situation will occur when the test is performed on abnormal milk, such as clotted milk, or if the procedure has not been performed properly.
4. Any test strips removed from the storage tube that remain unused at the end of the testing day must be discarded.

PRECAUTIONS

1. When handling the BetaStar Advanced for Tetracyclines test strips, ensure hands are clean and dry. This will protect against contamination of the test strip.
2. Confirm that the test strip is correctly inserted into the cartridge prior to inserting cartridge into the ports in the Raptor.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

Please visit <http://www.neogen.com/en/terms-and-conditions> for Neogen's full terms and conditions.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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Neogen Food and Animal Security

+91 484 2306598, 2301582
info@neogenindia.com
www.neogenindia.com

©Neogen Corporation, 2018. Neogen, BetaStar and Raptor are registered trademarks of Neogen Corp., Lansing, MI. Patent: <http://www.neogen.com/Corporate/patents.html>

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©Neogen Corporation, 2018. Neogen, BetaStar y Raptor son marcas registradas de Neogen Corp., Lansing, MI.