*Read instructions carefully before starting test*

# BetaStar<sup>®</sup> S

## Combo

### assay for beta-lactam and tetracycline antibiotics

#### **INTRODUCTION**

BetaStar<sup>®</sup> S Combo is a rapid detection assay for beta-lactams, tetracyclines, and desfuoylceftiofur (the metabolite of ceftiofur). The desfuoylceftiofur line detects the sum of the parent drug ceftiofur and its metabolite desfuoylceftiofur.

#### **REACTION MECHANISM**

BetaStar S Combo is a single-step lateral flow immunochromatographic assay based on a competitive immunoassay format. The sample is wicked through a reagent zone, which contains antibodies and receptors conjugated to colloidal gold particles. If beta-lactams, tetracyclines, or desfuoylceftiofur are present, they will be captured by their particle-antibody complex. The drug-antibody-particle complex then is wicked onto a membrane, which contains a separate zone for each drug conjugated to a protein carrier. These independent zones capture any uncomplexed drug-antibody, allowing the particles to concentrate and form a visible line. As the level of a particular drug in the sample increases, free drug molecules will complex with the antibody-gold particles. This allows less antibody-gold to be captured in the test zones. Therefore, as the concentration of a drug in the sample increases, the test line density decreases. The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of the drug, ensuring the strip is functioning properly.

## **MATERIALS REQUIRED**

### **BetaStar S Combo 50 kit (BSC050) contains:**

1. 50 plastic vials
2. 2 containers with 25 test strips
3. 50 disposable pipettes

### **BetaStar S Combo Raptor 100 kit (BSCR100) contains:**

1. 100 Raptor cartridges
2. 4 containers with 25 test strips
3. 100 disposable pipettes

## **MATERIALS REQUIRED FOR VISUAL OR SOLO END-POINT**

1. Heater block/incubator capable of maintaining a temperature of  $47.5 \pm 1^\circ\text{C}$

## **MATERIALS REQUIRED FOR RAPTOR OR RAPTOR SOLO**

1. Raptor Integrated Analysis Platform (Neogen item 9680)
2. Raptor Solo Platform (Neogen item 9696)
3. Raptor cartridge included in Raptor kit (Neogen item 9681)

## **TEST PREPARATION WITH INCUBATOR**

1. A daily temperature check of the heater block is recommended. Ensure the heater block has been turned on and preheated and the temperature is maintained at  $47.5 \pm 1^\circ\text{C}$ .
2. BetaStar S Combo is designed for use under normal ambient environmental conditions ( $15\text{--}30^\circ\text{C}$ ). Remove the kit from the refrigerator and leave the test strip container at room temperature ( $15\text{--}30^\circ\text{C}$ ) for 10–15 minutes prior to opening to prevent condensation.
3. Test strips that have been removed from the test strip container must be kept clean and dry.

## **TEST PROCEDURE WITH INCUBATOR**

1. Mix milk sample.
2. Pipette 0.4 mL milk sample into the bottom of the vial. This is achieved by pressing the pipette tip to the bottom of the vial to release the sample.
3. Place the vial into the heater block.
4. Place the test strip into the vial that is in the heater block. The arrows on the test strip must be oriented downward in the vial. Incubate the test strip in the vial for 5 minutes at  $47.5 \pm 1^\circ\text{C}$ .
5. After the 5-minute incubation, remove the test strip from the vial.

## **VISUAL TEST INTERPRETATION**

At the completion of the 5-minute incubation, remove the test strip from the vial. Immediately compare the intensities of the antibiotic test lines (lines 2, 3 and 4) to the control line (line 1). If the intensity of the test line is greater than or equal to the control line, the milk is negative for the presence of the antibiotic. If the intensity of the test line is less than the intensity of the control line, the milk sample is positive for the antibiotic.

## **TEST INTERPRETATION USING THE RAPTOR SOLO IN END-POINT MODE (OPTIONAL)**

At the completion of the five-minute incubation, remove the test strip from the vial. Place the test strip into the Raptor Solo end-point test strip holder. Insert the holder into the Raptor Solo reader. Select the Dairy Antibiotics test category and select the BetaStar S Combo test. Press the Run Test icon to analyze the test strip. Read the data in the analysis windows. If the ratio of the test line intensity to the control line intensity is  $\geq 1.0$ , the test is negative. If the ratio is  $<1.0$ , the test is positive.

## **RAPTOR® INTEGRATED ANALYSIS PLATFORM AND RAPTOR SOLO**

BetaStar S Combo can be used with the Raptor® Integrated Analysis Platform (Neogen item 9680) and Raptor Solo (Neogen item 9696). A Raptor cartridge is required (Neogen item 9681).

## **TEST PROCEDURE FOR RAPTOR AND RAPTOR SOLO**

1. Insert the test strip into the Raptor cartridge in any of the available slots.
2. Insert the cartridge containing the test strip into any of the 3 ports within the Raptor Integrated Analysis Platform or single port of the Raptor Solo.
  - a. The bar code on the test device will be read. If the QR code for the lot of devices has not been entered into the Raptor, the bar code reader in the front of the Raptor will turn on automatically. Scan the QR code found on the container storing the test devices.
  - b. The assay temperature for the BetaStar S Combo test is 65°C. This temperature is programmed in the Raptor.
  - c. When the cartridge is inserted into a port, the port will automatically begin to adjust to the proper temperature.
  - d. The user will not be able to proceed until the incubator temperature reaches 65°C.
3. Enter the sample ID by either scanning or using the on-screen keyboard.
4. Mix the milk sample prior to adding the sample into the cartridge.
5. Pipette slowly 0.4 mL milk sample into the sample port located in the back of the cartridge. Do not seat the pipette directly into the back of the cartridge. Press “Next.” The analysis will take place automatically.
6. Insert the used pipette into the sample port located in the back of the cartridge. This will prevent double loading the same sample or loading a second sample into the same cartridge.
7. After a 5-minute incubation, the results will appear on the Raptor screen.

## **REMARKS**

1. Milk sample should be kept refrigerated between 2–8° just prior to testing.
2. All reagents must be kept refrigerated between 2–8°C. Before opening the test strip container, it should be equilibrated to room temperature for at least 10 minutes.
3. If the test sample does not migrate on the strip, the test is invalid. This situation will occur when the test is performed on abnormal milk such as clotted milk, or if the procedure has not been performed properly.
4. When using the Raptor Solo in End-point mode the End-point track must be installed.

## **PRECAUTIONS**

When handling the BetaStar S Combo test strip, ensure hands are clean and dry. This will protect against contamination of the test strip.

## **CUSTOMER SERVICE**

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

## **SDS INFORMATION AVAILABLE**

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at [foodsafety.neogen.com](https://foodsafety.neogen.com), or by calling Neogen at 800.234.5333 or 517.372.9200.

## **TERMS AND CONDITIONS**

Please visit [www.neogen.com/en/terms-and-conditions](https://www.neogen.com/en/terms-and-conditions) for Neogen's full terms and conditions.

## **WARRANTY**

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## TEST KITS AVAILABLE FROM NEOGEN

### Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

### Food allergens

- Almonds, crustaceans, coconut, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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