Read instructions carefully before starting test

BetaStar® 4D

Assay for beta-lactams, tetracyclines, chloramphenicol and streptomycin residues in raw, commingled cow milk

INTRODUCTION

BetaStar® S 4D is a rapid detection assay for beta-lactams, tetracyclines, chloramphenicol and streptomycin residues.

REACTION MECHANISM

BetaStar S 4D is a single-step lateral flow immunochromatographic assay based on a competitive immunoassay format. The milk is wicked through a reagent zone, which contains antibodies and receptors conjugated to colloidal gold particles. If beta-lactams, tetracyclines, chloramphenicol or streptomycin are present, they will be captured by their particle-antibody complex. The drug-antibody-particle complex then is wicked onto a membrane, which contains a separate zone for each drug conjugated to a protein carrier. These independent zones capture any uncomplexed drug antibody, allowing the particles to concentrate and form a visible line. As the level of a particular drug in the sample increases, free drug molecules will complex with the antibody-gold particles. This allows less antibody-gold to be captured in the test zones. Therefore, as the concentration of a drug in the sample increases, the test line density decreases. Algorithms programmed into the AccuScan® Pro reader convert these line densities into positive or negative results for each drug. The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of the drug, ensuring the strip is functioning properly.

MATERIALS REQUIRED

BetaStar S 4D 50 kit (BS4D050) contains:

1. 50 plastic vials
2. 2 containers with 25 test strips
3. 50 disposable pipettes
4. 1 positive/negative control pack

BetaStar S 4D Raptor 100 kit (BS4DR100) contains:

1. 100 disposable Raptor cartridges
2. 4 containers with 25 test strips
3. 100 disposable pipettes
4. 1 positive/negative control pack

MATERIALS REQUIRED FOR ACCUSCAN PRO OR GOLD

1. Heater block/incubator capable of maintaining a temperature of $47.5 \pm 1^\circ\text{C}$
2. AccuScan Pro or AccuScan Gold reader

MATERIALS REQUIRED FOR RAPTOR OR RAPTOR SOLO

1. Raptor Integrated Analysis Platform (NEOGEN item 9680)
2. Raptor Solo Platform (NEOGEN item 9696)
3. Raptor cartridge included in Raptor kit (NEOGEN item 9681)

TEST PREPARATION ACCUSCAN READER

1. A daily temperature check of the heater block is recommended. Ensure the heater block has been turned on and preheated and the temperature is maintained at $47.5 \pm 1^\circ\text{C}$.
2. BetaStar S 4D is designed for use under normal ambient environmental conditions ($15\text{--}30^\circ\text{C}$). Remove the kit from the refrigerator and leave the test strip container at room temperature ($15\text{--}30^\circ\text{C}$) for **10–15 minutes** prior to opening to prevent condensation.
3. Test strips that have been removed from the test strip container must be kept clean and dry.

ACCUSCAN PRO READER SET UP

1. Enter the lot-specific QR code by selecting the QR code icon on the reader.
2. Place the QR code into the code device cartridge and insert the cartridge into the reader.
3. Return to the home screen and select the test strip icon. Touch the **Dairy Residues** category, then select BetaStar S 4D test type.

ACCUSCAN GOLD READER SET UP

1. Enter the lot specific QR code by selecting the Scan QR button on the main screen.
2. Place the QR code into the Cal/QR device cartridge and insert the cartridge into the reader.
3. Return to the main screen and select the **Dairy Antibiotics** category and **BetaStar S 4D** assay type.
4. Select the lot number for the kit in use from the drop-down menu available.

PRECAUTIONS

1. Milk sample should be kept refrigerated between $2\text{--}8^\circ\text{C}$
2. Milk sample must be from at least four animals and must be commingled for at least 24 hours prior to testing.

TEST PROCEDURE

1. Mix the refrigerated milk sample.
2. Pipette 0.4 mL milk sample into the vial.
3. Place the vial into the heater block.
4. Immediately place the test strip into the vial that is in the heater block. The word “sample” on the test strip must be oriented downward in the vial. Incubate the test strip in the vial for 10 minutes at $47.5 \pm 1^\circ\text{C}$.
5. After the 10-minute incubation, remove the test strip from the vial.
6. Read the test strip using the AccuScan Pro or Accuscan Gold Reader.
NOTE: BetaStar S 4D is not a visual test. Results must be determined using the AccuScan Pro or AccuScan Gold Reader.

TEST INTERPRETATION USING THE ACCUSCAN READER

At the completion of the 10-minute incubation, remove the test strip from the vial. Place the test strip onto the AccuScan test strip holder. Insert the holder into the AccuScan reader. The results will be displayed on the touch screen. For each residue, if the ratio of the test line intensity to the control line intensity is greater than the predetermined lot specific value, the test is negative. If the ratio is less than or equal to this value, the test is positive.

NOTE: The image to the left only is meant to indicate line position. BetaStar S 4D is not a visual test. Results must be determined using the AccuScan reader.

RAPTOR TEST PROCEDURE

1. Insert the cartridge containing the test strip into any of the three ports within the Raptor Integrated Analysis Platform.
 - a. The bar code on the test device will be read. If the QR code for the lot of devices has not been entered into the Raptor, the bar code reader in the front of the Raptor will turn on automatically. Scan the QR code found on the container storing the test devices.
 - b. The assay temperature for the BetaStar S 4D test is 65°C . This temperature is programmed in the Raptor.
 - c. When a cartridge is inserted into one of the three ports, the port will automatically begin to adjust to the proper temperature.
 - d. The user will not be able to proceed until the incubator temperature reaches 65°C .
2. Enter the sample ID by either scanning or using the on-screen keyboard.
3. Mix the refrigerated milk sample prior to adding the sample into the cartridge.
4. Make sure the port has reached 65°C . Follow the prompt on-screen.
5. Pipette slowly 0.4 mL milk sample into the sample port located in the back of the cartridge. Do not seat the pipette directly into the back of the cartridge during sample addition. Press “Next.” The analysis will take place automatically.
6. Insert the used pipette into the sample port located in the back of the cartridge. This will prevent double-loading the same sample or loading a second sample into the same cartridge.
7. After a 10-minute incubation, the results will appear on the Raptor screen.

NEGATIVE CONTROL PREPARATION

1. Mix buffer bottle for 10 seconds.

POSITIVE CONTROL PREPARATION

1. Add 1.2 mL (3 x 0.4 mL) of negative buffer by dispensing 0.4 mL three times with a new pipette into positive control vial.
2. Mix for 15 seconds or until pellet is completely dissolved.

REMARKS

1. All reagents must be kept refrigerated between 2–8°C. Before opening the test strip container, it should be equilibrated to room temperature for at least **10 minutes**.
2. If the test sample does not migrate on the strip, the test is invalid. This situation will occur when the test is performed on abnormal milk such as clotted milk, or if the procedure has not been performed properly.
3. The positive control is designed for the 4D test, which detects penicillin G at 4 ppb, tetracycline at 10 ppb, streptomycin at 200 ppb and chloramphenicol at 0.3 ppb.
4. When handling the BetaStar S 4D test strip, ensure hands are clean and dry. This will protect against contamination of the test strip.

NOTE: The concentration of the antibiotics is established to accommodate the freeze-dry process and degradation over the shelf life of the product. The level of antibiotics in the control ensure a positive response over the shelf life of the product.

CUSTOMER SERVICE

NEOGEN Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all NEOGEN test kits, is available.

Control line
Tetracycline line
Chloramphenicol line
Streptomycin line
Beta-lactam line

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of NEOGEN's test kits, on NEOGEN's website at foodsafety.neogen.com, or by calling NEOGEN at 800.234.5333 or 517.372.9200.

TERMS AND CONDITIONS

Please visit www.neogen.com/en/terms-and-conditions for NEOGEN's full terms and conditions.

WARRANTY

NEOGEN Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ergot alkaloids, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, coconut, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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