Read instructions carefully before starting test

Veratox[®]

for Clenbuterol

CLENBUTEROL

Clenbuterol is a beta agonist with some structural and pharmacological similarities to ractopamine and zilpaterol, but its effects are more potent and longer lasting. Clenbuterol increases the rate at which body fat is metabolized while increasing the body's basal metabolic rate (BMR). It is commonly used for smooth muscle-relaxant properties as a bronchodilator and tocolytic.

Clenbuterol is not approved for use in food producing animals in many countries like the U.S., Canada and Mexico. FDA lists clenbuterol as a prescription only, restricted-use drug for the treatment of horses. In addition, clenbuterol is a therapeutic drug for asthma, approved for human use in some countries in Europe (Bulgaria and Russia) and Asia (China).

INTENDED USE

Veratox[®] for Clenbuterol is a competitive enzyme immunoassay for the quantitative analysis of clenbuterol in urine, muscle, liver, kidney and serum. Veratox for Clenbuterol enables international and government regulatory agencies, food manufacturers and processors, as well as quality assurance organizations, to detect clenbuterol in urine, muscle, liver, kidney and serum in response to customer concerns about food safety. Veratox for Clenbuterol is an accurate, rapid, cost-effective and easy-to-use diagnostic tool compared to instrument based diagnostics like LC-MS-MS. The unique features of the kit are:

- Rapid sample prep without a column
- High sensitivity and specificity
- A quick ELISA assay (just over 1.5 hours)

ASSAY PRINCIPLES

The method is based on a competitive colorimetric ELISA assay. The clenbuterol antibody has been coated in the plate wells. During the analysis, sample is added along with the clenbuterol-horseradish peroxidase (Clenbuterol-HRP) conjugate. If the clenbuterol residue is present in the sample, it will compete for the clenbuterol antibody, thereby preventing the Clenbuterol-HRP from binding to the antibody attached to the well. The resulting color intensity, after addition of the HRP substrate (TMB), has an inverse relationship with the clenbuterol residue concentration in the sample.

Veratox for Clenbuterol has the capacity for 96 determinations or testing of 36 samples in duplicate (assuming 24 wells for standards). Return any unused microwells to the foil bag and reseal them with the desiccant provided in the original package.

STORAGE REQUIREMENTS

Store kit at 2–8°C (35–46°F). Some components should be stored at frozen temperatures if not used within 1 month. The shelf life of the kit is 12 months when properly stored.

MATERIALS PROVIDED

Kit Contents	Amount	Storage
Clenbuterol Ab-coated Microtiter Plate	1 x 96-well plate (8 wells x 12 strips)	2–8°C
Clenbuterol Standards: Negative control (white cap tube) 0.04 ng/mL (yellow cap tube) 0.08 ng/mL (orange cap tube) 0.16 ng/mL (pink cap tube) 0.32 ng/mL (purple cap tube) 0.64 ng/mL (blue cap tube) 10 ng/mL (spiking, red cap tube)	1.8 mL 1.8 mL 1.8 mL 1.8 mL 1.8 mL 1.8 mL 1.8 mL	2–8°C
100X Clenbuterol-HRP Conjugate	175 µL	2–8°C*
Clenbuterol-HRP Diluent	20 mL	2–8°C
20X Wash Solution**	28 mL	2–8°C
Stop Buffer**	14 mL	2–8°C
TMB Substrate**	12 mL	2–8°C
Balance Buffer	12 mL	2–8°C
20X Sample Extraction Buffer	10 mL	2–8°C
Sample Balance Buffer I	10 mL	2–8°C
Sample Balance Buffer II (powdered)	62 g	2–8°C
10X PBS	25 mL	2–8°C
Liver Extraction Buffer	10 mL	2–8°C
Serum Extraction Buffer A	12 mL	2–8°C
Serum Balance Buffer B	25 mL	2–8°C

*If the kit will be unused for over 1 month, store 100X Clenbuterol-HRP Conjugate at -20°C or in a freezer.

**These components (within their expiry) can be used interchangeably with other Veratox for Clenbuterol kits.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Microtiter plate reader with a 450 nm filter (Neogen item 9303)
2. 20–200 µL pipettes (Neogen item 9276)
3. Multichannel pipette: 50–300 µL (optional) (Neogen item 9385)
4. Centrifuge 4,000 x g
5. Timer (Neogen item 9426)
6. Wash bottle (Neogen item 9400)
7. Paper towels or equivalent absorbent material
8. Veratox software (Neogen item 9305)
9. Lab station (Neogen item 9481)
10. Distilled or deionized water
11. Blender or food processor to homogenize sample
12. Rotary evaporator or other means to evaporate down a sample
13. Acetonitrile
14. HPLC grade hexane

PRECAUTIONS

Neogen strongly recommends that you read the following warnings and precautions to ensure your full awareness of ELISA techniques and other details you should pay close attention to when running the assays. Periodically, optimizations and revisions are made to the kit and insert. Therefore, it is important to follow the protocol included with the kit. If you need further assistance, you may contact your local distributor or Neogen at foodsafety@neogen.com.

1. The standards contain clenbuterol. Handle with particular care.
2. Do not use the kit past the expiration date.
3. Do not intermix reagents from different kits or lots except for components with the same lot numbers within their expiration dates. CLENBUTEROL-HRP CONJUGATES AND PLATES ARE KIT AND LOT SPECIFIC.
4. Make sure HRP-Conjugate and diluent are mixed in correct volumes.
5. Try to maintain a laboratory temperature of 20–25°C (68–77°F). Avoid running assays under or near air vents, as this may cause excessive cooling, heating and/or evaporation. Also, do not run assays in direct sunlight, as this may cause excessive heat and evaporation. Cold bench tops should be avoided by placing several layers of paper towel or some other insulation material under the assay plates during incubation.
6. Make sure you are using only distilled or deionized water since water quality is very important.
7. Follow proper pipetting techniques, including priming tips by filling and dispensing solution once before use.
8. Incubations of assay plates should be timed as precisely as possible. Be consistent when adding standards to the assay plate. Add your standards first and then your samples.
9. Add standards to plate only in the order from low concentration to high concentration as this will minimize the risk of compromising the standard curve.
10. Always refrigerate plates in sealed bags with a desiccant to maintain stability. Prevent condensation from forming on plates by allowing them to equilibrate to room temperature 20–25°C (68–77°F) while in the sealed bag.

SAMPLE PREPARATION AND EXTRACTION

Be sure samples are properly stored. In general, samples should be refrigerated at 2–4°C for no more than 1–2 days. Freeze samples to a minimum of -20°C if they need to be stored for a longer period. Frozen samples can be thawed at room temperature 20–25°C (68–77°F) or in a refrigerator before use. Preparation protocols for samples other than below can be made available upon request. Please contact your local distributor or Neogen at foodsafety@neogen.com.

1. Preparation of 1X Sample Extraction Buffer

Mix 1 volume of 20X Sample Extraction Buffer with 19 volumes of distilled water.

2. Preparation of 1X PBS

Mix 1 volume of 10X PBS Buffer with 9 volumes of distilled water.

3. Preparation of 1X Sample Balance Buffer

Add 10 mL of the Sample Balance Buffer I and 62 g of Sample Balance Buffer II (powdered) to a 250 mL plastic bottle. Add 190 mL of distilled water to the bottle. Mix well to obtain a homogeneous solution.

Urine

1. Centrifuge urine at 4000 x *g* for 5 minutes.
2. Follow step 1 of Procedure A in the ELISA Testing Protocol.

NOTE: The total in each well is 200 µL for the assay.

Dilution factor: 2.

Muscle/Liver/Kidney

1. Remove fat from muscle/liver/kidney. Homogenized the sample with a suitable mixer.
2. Add 2 mL of 1X Sample Extraction Buffer to 1 g of the homogenized sample, vortex for 5 minutes at maximum speed.
3. Add 3 mL of acetonitrile and 2.5 mL of 1X Sample Balance Buffer, 100 µL of Liver Extraction Buffer to the sample. Vortex for 1 minute at maximum speed.
4. Centrifuge for 5 minutes at 4,000 x *g* at room temperature 20–25°C (68–77°F), transfer 900 µL of the supernatant to a new tube.
5. Use a rotary evaporator to dry the sample in a 60°C water bath under reduced pressure. Alternatively, the sample can be dried by blowing nitrogen gas in a 60°C water bath.
6. Add 1 mL of hexane to dissolve the sample and then add 0.6 mL of 1X PBS, vortex for 1 minute at maximum speed.
7. Centrifuge the sample for 5 minutes at 4,000 x *g*, aspirate off the upper hexane layer.
8. Use 100 µL of the lower aqueous layer for the assay.
9. Follow step 1 in Procedure A of the ELISA Testing Protocol.
NOTE: The total in each well is 200 µL for the assay.
Final Dilution factor: 4 (when accounting for ELISA dilution).

Serum

1. To 500 µL of serum in a 1.5-mL micro-centrifuge tube, add 125.0 µL of Serum Extraction Buffer A.
2. Vortex samples for 3 minutes.
3. Centrifuge sample for 5 minutes at maximum speed.
4. Transfer 150 µL of supernatant to a new 1.5-mL micro-centrifuge tube containing 450 µL of Serum Balance Buffer B.
5. Vortex sample for 2 minutes.
6. Centrifuge sample for 5 minutes at maximum speed.
7. Use 200 µL of supernatant per well in the assay.
8. Follow step 1 in Procedure B of the ELISA testing protocol.
NOTE: Dilution factor: 5.25

TEST PROCEDURE

Reagent Preparation

IMPORTANT: All reagents should be brought up to room temperature before use (1–2 hours at 20–25°C / 68–77°F); Make sure you read “Precautions” section on page 3. Solutions should be prepared just prior to ELISA test. All reagents should be mixed by gently inverting or swirling prior to use. Prepare volumes that are needed for the number of wells being run. Do not return the reagents to the original stock tubes/bottles. Using disposable reservoirs when handling reagents can minimize the risk of contamination and is recommended.

1. **Preparation of 1X Wash Solution**
Mix 1 volume of the 20X Wash Solution with 19 volumes of distilled water.
2. **Preparation of 1X Clenbuterol-HRP**
Mix 1 volume of the 100X Clenbuterol-HRP with 99 volumes of the Clenbuterol-HRP Diluent. Prepare just prior to use on the ELISA.
NOTE: The 1X Clenbuterol-HRP should be used within 5 minutes after preparation.

ELISA TESTING PROTOCOL

Label the individual strips that will be used and aliquot reagents as the following example:

Component	Volume per Reaction	24 Reactions
1X Clenbuterol-HRP Conjugate	100 µL	2.4 mL
1X Wash Solution	2.0 mL	48 mL
Stop Buffer	100 µL	2.4 mL
TMB Substrate	100 µL	2.4 mL

PROCEDURE A FOR URINE, MUSCLE, LIVER AND KIDNEY SAMPLES

1. Add 200 µL of each Clenbuterol Standards in duplicate into different sample wells.
NOTE: Add standards to plate only in the order from low concentration to high concentration.
2. Add 100 µL of Balance Buffer to each well corresponding with a sample. Add 100 µL of the sample to the well. Mix by sliding back and forth on a flat surface for **1 minute**.
NOTE: The total volume in each well is 200 µL.
3. Incubate the plate for **60 minutes** at room temperature 20–25°C (68–77°F).
4. Aspirate all fluid from each well after the incubation, invert the plate and repeatedly tap the plate on paper towels to dry. Dry the plate with paper towels only.
NOTE: Perform the next step immediately after drying the plate. Do not allow the plate to air dry between working steps.
5. Add 100 µL of fresh prepared 1X Clenbuterol-HRP Conjugate to each well and mix by sliding back and forth on a flat surface for **1 minute**.
6. Incubate the plate for **15 minutes** at room temperature 20–25°C (68–77°F).
7. Wash the wells 3 times with 250 µL of 1X Wash Solution per well. After the last wash, invert the plate and gently tap the plate dry on paper towels.
NOTE: Perform the next step immediately after drying the plate. Do not allow the plate to air dry between working steps.
8. Add 100 µL of the TMB Substrate to each well. Start the 15 minute incubation immediately after adding the substrate. Mix by sliding back and forth on a flat surface for **1 minute** while incubating. Total incubation time for TMB Substrate is 15 minutes.
NOTE: Do not put any substrate back into the original container to avoid any potential contamination. Any substrate solution exhibiting coloration is indicative of deterioration and should be discarded. Covering the microtiter plate while incubating is recommended.
9. Add 100 µL of Stop Buffer to each well to stop the enzyme reaction.
10. Read the plate as soon as possible following the addition of Stop Buffer on a plate reader with a 450 nm wavelength.
NOTE: Before reading, use a lint-free wipe on the bottom of the plate to ensure no moisture or fingerprints interfere with the readings.

PROCEDURE B FOR SERUM SAMPLES

1. Add 200 μL of each Clenbuterol Standard in duplicate into different sample wells.
NOTE: Add standards to plate only in the order from low concentration to high concentration.
2. Add 200 μL of the sample to the corresponding sample well. Mix by sliding back and forth on a flat surface for 1 minute.
NOTE: The total volume in each well is 200 μL .
3. Incubate the plate for **60 minutes** at room temperature 20–25°C (68–77°F).
4. Wash the wells 2 times with 250 μL of 1X Wash Solution per well after incubation. Invert the plate and repeatedly tap the plate on paper towels to dry. Dry the plate with paper towels only.
NOTE: Perform the next step immediately after drying the plate. Do not allow the plate to air dry between working steps.
5. Add 100 μL of fresh prepared 1X Clenbuterol-HRP Conjugate to each well and mix by sliding back and forth on a flat surface for **1 minute**.
6. Incubate the plate for **25 minutes** at room temperature 20–25°C (68–77°F).
7. Wash the wells 3 times with 250 μL of 1X Wash Solution per well. After the last wash, invert the plate and gently tap the plate dry on paper towels.
NOTE: Perform the next step immediately after drying the plate. Do not allow the plate to air dry between working steps.
8. Add 100 μL of the TMB Substrate to each well. Start the 15 minute incubation immediately after adding the substrate. Mix by sliding back and forth on a flat surface for **1 minute** while incubating. Total incubation time for TMB Substrate is 15 minutes.
NOTE: Do not put any substrate back into the original container to avoid any potential contamination. Any substrate solution exhibiting coloration is indicative of deterioration and should be discarded. Covering the microtiter plate while incubating is recommended.
9. Add 100 μL of Stop Buffer to each well to stop the enzyme reaction.
10. Read the plate as soon as possible following the addition of Stop Buffer on a plate reader with a 450 nm wavelength.
NOTE: Before reading, use a lint-free wipe on the bottom of the plate to ensure no moisture or fingerprints interfere with the readings.


CLENBUTEROL CONCENTRATION CALCULATIONS

A standard curve can be constructed by plotting the mean relative absorbance (%) obtained from each reference standard against its concentration in ng/mL on a logarithmic curve.

$$\text{Relative absorbance (\%)} = \frac{\text{absorbance standard (or sample)} \times 100}{\text{absorbance zero standard}}$$

Use the mean relative absorbance values for each sample to determine the corresponding concentration of the tested drug in ng/mL from the standard curve. Neogen's Veratox Software for Windows is available to calculate test results. Please contact your local distributor or foodsafety@neogen.com for further information.

The following figure is a typical clenbuterol standard curve.

PERFORMANCE CHARACTERISTICS

Sensitivity (Detection Limit)

Sample Type	Detection Limit (ng/g or ppb)
Urine	0.2
Liver/Kidney	0.3
Muscle	0.15
Serum	0.5

Specificity (Cross-Reactivity)

Analytes	Cross-Reactivity (%)
Clenbuterol	100
Mabuterol	36
Brombuterol	33
Clenpenterol	28
Cimbuterol	28
Mapenterol	23
Salbutamol	6
Carbuterol	5
Terbutaline	4
Cimaterol	3
Fenoterol	< 1%
Ractopamine	< 1%
Zilpaterol	< 1%

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/en/terms-and-conditions.

WARRANTY

Neogen makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

OTHER RESIDUE DIAGNOSTICS KITS

RT-96-TT-AMS	Alert for Ractopamine — qualitative microwell assay, 96 wells
9551	Veratox for Chloramphenicol — range 10–1,000 ppt, 96 wells
8416	Veratox for Malachite Green — range 1–4 ppb, 48 wells
DR073	Veratox for Avermectins — range 6.4–300 ppb, 96 wells
DR107	Veratox for Florfenicol — range 0.15–100 ppm, 96 wells
DR081	Veratox for Oxytetracycline — range 3.6–324 ppb, 96 wells

North America
Neogen Headquarters
800/234-5333 (USA/Canada)
foodsafety@neogen.com
foodsafety.neogen.com

Europe, Middle East and Africa
Neogen Europe
+ 44 (0) 1292 525 600
info_uk@neogeneurope.com
www.neogeneurope.com

Mexico
Neogen Latinoamerica
+52 (55) 5254-8235
informacion@neogenlac.com
www.neogenlac.com

Brazil
Neogen do Brasil
+55 19 3935.3727
info@neogendobrasil.com.br
www.neogendobrasil.com.br

China
Neogen Bio-Scientific Technology
+86 21 6271 7013
info@neogenchina.com.cn
www.neogenchina.com.cn

India
Neogen Food and Animal Security
+91 484 2306598, 2301582
info@neogenindia.com
www.neogenindia.com

© Neogen Corporation, 2018. Neogen, Veratox and K-Blue are registered trademarks of Neogen Corporation, Lansing, MI. All other brand and product names are trademarks or registered trademarks of their respective companies. Patent: <http://www.neogen.com/en/patents>

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