Read instructions carefully before starting test

Veratox[®] for Tylosin

Refrigerate at 2–8°C (35–46°F) • Do not freeze

TYLOSIN

Tylosin is a macrolide antibiotic used in treating illnesses found in animals produced for food. Tylosin can be administered to individual animals through injection, but is more commonly used in feed to treat groups of animals.

INTENDED USE

Veratox[®] for Tylosin is a competitive enzyme immunoassay for the quantitative analysis of tylosin in muscle.

Veratox for Tylosin enables international and government regulatory agencies, food manufacturers and processors, as well as quality assurance organizations, to detect tylosin in muscle and to satisfy customer concerns about food safety.

ASSAY PRINCIPLES

The method is based on a competitive colorimetric ELISA assay. The drug of interest has been coated in the plate wells. During the analysis, sample is added along with the primary antibody specific for the target drug. If the target is present in the sample, it will compete for the antibody, thereby preventing the antibody from binding to the drug attached to the well. The secondary antibody, tagged with a peroxidase enzyme, targets the primary antibody that is complexed to the drug coated on the plate wells. The resulting color intensity, after addition of substrate, has an inverse relationship with the target concentration in the sample.

Veratox for Tylosin has the capacity for 96 determinations or testing of 36 samples in duplicate (assuming 24 wells for standards). Return any unused microwells to the foil bag and reseal them with the desiccant provided in the original package.

STORAGE REQUIREMENTS

Store kit at 2–8°C (35–46°F). Some components should be stored at frozen temperatures if not used within 1 month. The shelf life of the kit is 12 months when properly stored.

MATERIALS PROVIDED

Kit Contents	Amount	Storage
Tylosin-coated Plate	1 x 96-well plate (8 wells x 12 strips)	2–8°C
Tylosin Standards: 0 ng/mL (white cap tube) 1 ng/mL (yellow cap tube) 2.5 ng/mL (orange cap tube) 5 ng/mL (pink cap tube) 10 ng/mL (purple cap tube) 25 ng/mL (blue cap tube) Tylosin Spiking Stock 1000 ng/mL (red cap tube)	0.8 mL each standard	2–8°C
Tylosin Antibody #1*	12 mL	2–8°C*
100X HRP-Conjugated Antibody #2	250 µL	2–8°C*
Antibody #2 Diluent	20 mL	2–8°C
Tylosin Resuspension Buffer	2 x 25 mL	2–8°C
20X Wash Solution	28 mL	2–8°C
Stop Buffer	14 mL	2–8°C
TMB Substrate**	12 mL	2–8°C

* If the kit will be unused for over 1 month, store Tylosin Antibody #1 and 100X HRP-Conjugated Antibody #2 at -20°C or in a freezer.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Microtiter plate reader with a 450 nm filter (Neogen item 9303)
2. 20–200 µL pipettes (Neogen item 9276)
3. Multichannel pipette: 50–300 µL (optional) (Neogen item 9385)
4. 1000 µL pipettor (Neogen 9337)
5. Vortex Mixer (Neogen 9494)
6. Centrifuge 4,000 x *g*
7. Timer (Neogen item 9426)
8. Wash bottle (Neogen item 9400)
9. Paper towels or equivalent absorbent material
10. Veratox software (Neogen item 9305)
11. Lab station (Neogen item 9481)
12. Distilled or deionized water
13. Blender or food processor to homogenize sample
14. Nitrogen gas evaporator
15. Methanol
16. N-Hexane

PRECAUTIONS

Neogen strongly recommends that you read the following precautions to ensure your full awareness of ELISA techniques and other details you should pay close attention to when running the assays. More information can also be found in the Technical Solutions & Troubleshooting section. Periodically, optimizations and revisions are made to the kit and manual. Therefore, it is important to follow the protocol included with the kit. If you need further assistance, you may contact Neogen at 800/234-5333 or 517/372-9200.

1. The standards contain tylosin. Handle with particular care.
2. Do not use the kit past the expiration date.
3. Do not intermix reagents from different kits or lots except for components with the same part numbers within their expiration dates. ANTIBODIES AND PLATES ARE KIT- AND LOT-SPECIFIC. Make sure that antibody #2 and diluent are mixed in correct volumes.
4. Try to maintain a laboratory temperature of 20–25°C (68–77°F). Avoid running assays under or near air vents, as this may cause excessive cooling, heating and/or evaporation. Also, do not run assays in direct sunlight, as this may cause excessive heat and evaporation. Cold bench tops should be avoided by placing several layers of paper towel under the assay plates during incubation.
5. Make sure you are using only distilled or deionized water since water quality is very important.
6. Incubations of assay plates should be timed as precisely as possible. Be consistent when adding standards to the assay plate. Add your standards first and then your samples.
7. Follow proper pipetting techniques, including priming tips by filling and dispensing solution once before use.
8. Add standards to plate only in the order from low concentration to high concentration as this will minimize the risk of compromising the standard curve.
9. Always refrigerate plates in sealed bags with a desiccant to maintain stability. Prevent condensation from forming on plates by allowing them to equilibrate to room temperature (20–25°C / 68–77°F) while in the sealed bag.

SAMPLE PREPARATION

Be sure samples are properly stored. In general, samples should be refrigerated at 2–4°C for no more than 48 hours. Freeze samples to a minimum of -20°C if they need to be stored for a longer period. Frozen samples can be thawed at room temperature (20–25°C / 68–77°F) or in a refrigerator before use.

SAMPLE EXTRACTION

Muscle

1. Add 4 mL of 100% methanol and 1 mL (\pm 0.1 mL) of n-hexane to 1.0 g (\pm 0.05 g) of a room-temperature homogenized sample in a 15 mL plastic conical tube.
2. Vortex sample manually for **1 minute**.
NOTE: Tilt the tube at 45° or on its side to aid in the breaking up of the pellet.
3. Open and close the cap to release any built-up gas. Incubate the sample for **15 minutes** at 45°C, then vortex sample for an additional **1 minute**. Open and close the cap once again to release any built-up gas.
4. Centrifuge sample for **10 minutes** at 4,000 x g.
5. Transfer 0.5 mL of the methanol layer to a new 15 mL tube.
NOTE: To ensure that neither hexane nor the interphase is transferred to the new tube, perform this step carefully by slowly puncturing the interphase with the pipette tip and transferring only from the middle of the methanol layer.
6. Use a nitrogen-blowing apparatus in a 60°C water bath to evaporate the methanol and concentrate the sample to a dried residue (about 40 minutes at 8 PSI).
7. Add 1 mL of Tylosin Resuspension Buffer to the dried residue. Scrape each tube's bottom back-and-forth vigorously against a micro-centrifuge rack at least 10 times to disrupt the pellet and aid in its dissolution.
8. Vortex sample well for **1 minute**, followed by a **10 minute** incubation at room temperature (20–25°C / 68–77°F). Then, vortex the sample for an additional **1 minute**.
9. Use 50 μ L per well for the assay.

NOTE: Dilution factor: 10.

REAGENT PREPARATION

IMPORTANT: All reagents should be brought up to room temperature before use (1–2 hours at 20–25°C / 68–77°F). Make sure you read the Precautions section on page 3. Solutions should be prepared prior to running the ELISA test. All reagents should be mixed by gently inverting or swirling prior to use. Prepare volumes that are needed for the number of wells being run. Do not return the reagents to the original stock tubes/bottles. Using disposable reservoirs when handling reagents can minimize the risk of contamination and is recommended.

A. Preparation of 1X HRP-Conjugated Antibody #2

Mix 1 volume of the 100X HRP-Conjugated Antibody #2 with 99 volumes of Antibody #2 Diluent.

B. Preparation of 1X Wash Solution

Mix 1 volume of the 20X Wash Solution with 19 volumes of distilled water.

TEST PROCEDURE

Label the individual strips that will be used and prepare reagent aliquots based on the number of reactions:

Component	Volume Per Reaction	24 Reactions
Tylosin Antibody #1	100 μ L	2.4 mL
1X HRP-Conjugated Antibody #2	150 μ L	3.6 mL
1X Wash Solution	2.0 mL	48 mL
Stop Buffer	100 μ L	2.4 mL
TMB Substrate	100 μ L	2.4 mL


1. Add 50 μ L of each Tylosin Standard in duplicate into different wells.
NOTE: Add standards to plate in order from low concentration to high concentration (white, yellow, orange, pink, purple, blue) and close each standard vial tightly.
2. Add 50 μ L of each sample in duplicate into different sample wells.
3. Add 100 μ L of Antibody #1 and mix by sliding back and forth on a flat surface for **1 minute**.
4. Incubate the plate for **30 minutes** at room temperature (20–25°C / 68–77°F).
5. Shake out the contents of the wells into a waste container. Wash each well with 250 μ L of 1X Wash Solution and dump out, repeat the wash step 3 times. After the last wash, invert the plate and gently tap the plate dry on paper towels.
NOTE: Perform the next step immediately (within 1 minute) after plate washings. Do not allow the plate to air dry between working steps.
6. Add 150 μ L of the 1X Antibody #2 solution. Incubate the plate for **30 minutes** at room temperature (20–25°C / 68–77°F).
NOTE: Avoid direct sunlight and cold bench tops during the incubation. Covering the microtiter plate while incubating is recommended.
7. Shake out the contents of the wells into a waste container. Wash each well with 250 μ L of 1X Wash Solution and dump out, repeat the wash step 3 times. After the last wash, invert the plate and gently tap the plate dry on paper towels.
NOTE: Perform the next step immediately (within 1 minute) after plate washings. Do not allow the plate to air dry between working steps.
8. Add 100 μ L of TMB Substrate to each well. Incubate the plate for **15 minutes** at room temperature (20–25°C / 68–77°F) in the dark. Time the reaction immediately after adding the substrate. Mix by sliding back and forth on a flat surface for **1 minute** while incubating.
NOTE: Do not put any substrate back into the original container to avoid any potential contamination. Covering the microtiter plate while incubating is recommended.
9. After incubation, add 100 μ L of Stop Buffer to stop the enzyme reaction.
10. Read the plate as soon as possible following the addition of Stop Buffer on a plate reader with a 450 nm wavelength.
NOTE: Before reading, use a lint-free wipe on the bottom of the plate to ensure no moisture or fingerprints interfere with the readings.

INTERPRETATION OF RESULTS

A standard curve can be constructed by plotting the mean relative absorbance (%) obtained from each reference standard against its concentration in ng/mL on a logarithmic curve.

$$\text{Relative absorbance (\%)} = \frac{\text{absorbance standard (or sample)} \times 100}{\text{absorbance zero standard}}$$

Use the mean relative absorbance values for each sample to determine the corresponding concentration of the tested drug in ng/mL from the standard curve. Veratox software is available upon request to evaluate the results. Please contact Neogen for further information. The following figure is a typical Tylosin standard curve:

Sensitivity (Assay Range)

Sample Type	Assay Range (ng/g or ppb)
Muscle	10 – 250

Specificity (Cross-Reactivity)

Analytes	Cross-Reactivity (%)
Tylosin	100

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

OTHER RESIDUE DIAGNOSTICS KITS

RT-96-TT-AMS	ALERT for Ractopamine — qualitative microwell assay, 96 wells
9551	Veratox for Chloramphenicol — range 10–1,000 ppt, 96 wells
8416	Veratox for Malachite Green — range 1–4 ppb, 48 wells
DR021	Veratox for Clenbuterol — range 0.08–1.28 ppb, 96 wells
DR107	Veratox for Florfenicol — range 0.15–100 ppm, 96 wells
DR073	Veratox for Avermectins — range 6.4–300 ppb, 96 wells

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