Read instructions carefully before starting test

Veratox[®] for Sulfonamide

Refrigerate at 2–8°C (35–46°F) • Do not freeze

SULFONAMIDE

Sulfonamides are a class of antibiotics used in both feed and water for the treatment of intestinal infections and other diseases in food producing animals. Sulfonamide residues can occur in food of animal origin. A maximum residue limit (MRL) for all substances of the sulfonamide class is 100 ppb in muscle, fat, liver and kidney.

INTENDED USE

Veratox[®] for Sulfonamide is a competitive enzyme immunoassay for the quantitative analysis of sulfonamide in kidney or muscle.

Veratox for Sulfonamide enables international and government regulatory agencies, food manufacturers and processors, as well as quality assurance organizations, to detect sulfonamide in kidney or muscle in response to customer concerns about food safety. Veratox for Sulfonamide is an accurate, rapid, cost-effective and easy-to-use diagnostic tool compared to instrument based diagnostics like LC-MS-MS.

ASSAY PRINCIPLES

The method is based on a competitive colorimetric ELISA assay. The sulfonamide antibody has been coated in the plate wells. During the analysis, sample is added along with the primary antibody specific for the target drug. If sulfonamide residue is present in the sample, it will compete for the sulfonamide antibody, thereby preventing the antibody from binding to the antibody attached to the well. The secondary antibody, tagged with a peroxidase enzyme, targets the primary antibody that is complexed to the drug coated on the plate wells. The resulting color intensity, after addition of the HRP substrate (TMB), has an inverse relationship with the sulfonamide residue concentration in the sample.

Veratox for Sulfonamide has the capacity for 96 determinations or testing of 36 samples in duplicate (assuming 24 wells for standards). Return any unused microwells to the foil bag and reseal them with the desiccant provided in the original package.

STORAGE REQUIREMENTS

Store kit at 2–8°C (35–46°F). Some components should be stored at frozen temperatures if not used within 1 month. The shelf life of the kit is 12 months when properly stored.

MATERIALS PROVIDED

Kit Contents	Amount	Storage
Sulfonamide Ab-coated Microtiter Plate	1 x 96-well plate (8 wells x 12 strips)	2–8°C
SMX Standards: 0 ng/mL (white cap tube) 0.50 ng/mL (yellow cap tube) 1.5 ng/mL (orange cap tube) 5.0 ng/mL (pink cap tube) 15.0 ng/mL (purple cap tube) 50.0 ng/mL (blue cap tube) 1000 ng/mL (spiking, red cap tube)	0.8 mL 0.8 mL 0.8 mL 0.8 mL 0.8 mL 0.8 mL 0.8 mL	2–8°C
Sulfonamide Antibody #1	12 mL	2–8°C*
100X HRP-Conjugated Antibody #2	250 µL	2–8°C*
Sulfonamide Antibody #2 Diluent	20 mL	2–8°C
10X Sulfonamide Extraction Buffer A	25 mL	2–8°C
20X Wash Solution**	28 mL	2–8°C
Stop Buffer**	14 mL	2–8°C
TMB Substrate**	12 mL	2–8°C
3X Kidney Extraction Buffer	20 mL	2–8°C
10X Sample Extraction Buffer E	25 mL	2–8°C
Concentrate of Tissue Extraction Buffer	62 g	2–8°C
Tissue Clean Up Mix	12 g	2–8°C

*If the kit will be unused for over 1 month, store Sulfonamide Antibody #1 and 100X HRP-Conjugated Antibody #2 at -20°C or in a freezer.

**These components (within their expiry) can be used interchangeably with other Veratox for Sulfonamide kits.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Microtiter plate reader with a 450 nm filter (Neogen item 9303)
2. Incubator
3. Tissue Mixer
4. Acetonitrile
5. Ethyl Acetate
6. 20–200 µL pipettes (Neogen item 9276)
7. Multichannel pipette: 50–300 µL (optional) (Neogen item 9385)
8. Centrifuge 4,000 x *g*
9. Timer (Neogen item 9426)
10. Wash bottle (Neogen item 9400)
11. Paper towels or equivalent absorbent material
12. Veratox software (Neogen item 9305)
13. Lab station (Neogen item 9481)
14. Distilled or deionized water
15. Blender or food processor to homogenize sample
16. Rotary evaporator or other means to evaporate down a sample with nitrogen

PRECAUTIONS

Neogen strongly recommends that you read the following warnings and precautions to ensure your full awareness of ELISA techniques and other details you should pay close attention to when running the assays. Periodically, optimizations and revisions are made to the kit and insert. Therefore, it is important to follow the protocol coming with the kit. If you need further assistance, you may contact your local distributor or Neogen at foodsafety@neogen.com.

1. The standards contain sulfonamide. Handle with particular care.
2. Do not use the kit past the expiration date.
3. Do not intermix reagents from different kits or lots except for components with the same lot numbers within their expiration dates. SULFONAMIDE-HRP CONJUGATES AND PLATES ARE KIT AND LOT SPECIFIC.
4. Make sure HRP-Conjugate and diluent are mixed in correct volumes.
5. Try to maintain a laboratory temperature of 20–25°C (68–77°F). Avoid running assays under or near air vents, as this may cause excessive cooling, heating and/or evaporation. Also, do not run assays in direct sunlight, as this may cause excessive heat and evaporation. Cold bench tops should be avoided by placing several layers of paper towel or some other insulation material under the assay plates during incubation.
6. Make sure you are using only distilled or deionized water since water quality is very important.
7. Follow proper pipetting techniques, including priming tips by filling and dispensing solution once before use.
8. Incubations of assay plates should be timed as precisely as possible. Be consistent when adding standards to the assay plate. Add your standards first and then your samples.
9. Add standards to plate only in the order from low concentration to high concentration as this will minimize the risk of compromising the standard curve.
10. Always refrigerate plates in sealed bags with a desiccant to maintain stability. Prevent condensation from forming on plates by allowing them to equilibrate to room temperature (20–25°C / 68–77°F) while in the sealed bag.

SAMPLE PREPARATION AND EXTRACTION

Be sure samples are properly stored. In general, samples should be refrigerated at 2–4°C for no more than 48 hours. Freeze samples to a minimum of -20°C if they need to be stored for a longer period. Frozen samples can be thawed at room temperature (20–25°C / 68–77°F) or in a refrigerator before use. Preparation protocols for samples other than below can be made available upon request. Please contact your local distributor or Neogen at foodsafety@neogen.com.

1. **Preparation of 1X Kidney Extraction Buffer**
Mix 1 part of 3X Kidney Extraction Buffer with 2 parts of deionized water. Vortex very well.
2. **Preparation of 1X Sulfonamide Extraction Buffer A**
Mix 1 volume of 10X Sulfonamide Extraction Buffer A with 9 volumes of deionized water. Vortex very well.
3. **Preparation of 1X Sample Extraction Buffer E**
Mix 1 volume of 10X Sample Extraction Buffer E with 9 volumes of deionized water. Vortex very well.
4. **Preparation of 1X Tissue Extraction Buffer**
Add the entire bag of powder from the Concentrate of Tissue Extraction Buffer bag to a 250 mL plastic bottle. Dissolve the powder by adding 180 mL of deionized water and shaking well for 5 minutes. Let the solution sit for 10 minutes, and then shake again to dissolve the remaining powder. It is normal for a small amount of particles to remain at the bottom of the bottle. Add increments of 1 mL of deionized water and shake until there are no particles remaining.

Kidney/Muscle

1. Homogenize 10 g of tissue with a suitable mixer
2. Weigh 1 g (\pm 0.05 g) of sample into a 15 mL plastic conical tube. Force the entire sample to the bottom of the tube.
3. Incubate each sample for **15 minutes** at 37°C.
NOTE: If performing an analyte-recovery experiment, add the target spiking concentration to each tube, and then incubate the sample.
4. Add exactly 2.5 mL of acetonitrile and 0.5 mL of deionized water to each sample.
NOTE: Tightly close each cap so as to not allow any solvents to escape during vortexing.
5. Vortex each sample thoroughly for **30 seconds** to ensure all the solvent is able to fully penetrate the tissue. During vortexing, tilt the tube at an angle to forcibly agitate the mixture. Allow the mixture to settle to the bottom of the tube.
6. Slowly open each cap to allow any built-up gas to escape. Tightly close each cap. Centrifuge each sample for **10 minutes** at 4,000 x *g* at room temperature (20–25°C / 68–77°F).
7. Slowly transfer 1 mL of the supernatant to a new 15 mL tube containing 0.5 mL of 1X Kidney Extraction Buffer.
8. Add 2.5 mL of ethyl acetate with a glass pipette. Vortex each sample manually or using a multi-tube vortexer for **5 minutes**.
9. Slowly open each cap to allow any built-up gas to escape. Tightly close each cap. Centrifuge each sample for **10 minutes** at 4,000 x *g* at room temperature (20–25°C / 68–77°F).
10. Slowly transfer 0.5 mL of the upper layer to a new 15 mL plastic conical tube.
11. Use a nitrogen-blowing apparatus in a 60°C water bath to evaporate the solvent and concentrate the sample to a dried residue.
12. To this residue, slowly add exactly 1.0 mL of 1X Sulfonamide Extraction Buffer A.
NOTE: Use a new pipette tip for each sample.
13. Tightly close each tube's cap. Scrape each tube's bottom back-and-forth vigorously against a micro-centrifuge rack at least 10 times to disrupt the pellet and aid in its dissolution.
NOTE: The solution should become white and foamy depending on the force used during this step.
14. Vortex each sample for **1 minute**.
15. Sonicate each sample for **1 minute** in a 37°C water bath.
NOTE: Tap the tube against the edge or bottom of the sonicator while sonicating. Make sure that any solid particles are off of the side of the wall of the tube.
16. Vortex for an additional **1 minute**, then incubate at room temperature (20–25°C / 68–77°F) for 10 minutes before continuing to the next step.
NOTE: If any solid particles are still apparent in the solution after incubation, repeat from Step 13 above. Due to the nature of certain samples along with how well each sample was extracted, it may take more or less time to dissolve the dried residue in the final resuspension. If the residue is not fully dissolved, it will give inaccurate and inconsistent results during the ELISA.
17. Vortex each sample for **15 seconds** before using 50 μ L per well for the ELISA.
NOTE: Dilution factor: 20.

TEST PROCEDURE

Reagent Preparation

IMPORTANT: All reagents should be brought up to room temperature before use (1–2 hours at 20–25°C / 68–77°F); Make sure you read “Precautions” section on page 3. Solutions should be prepared just prior to ELISA test. All reagents should be mixed by gently inverting or swirling prior to use. Prepare volumes that are needed for the number of wells being run. Do not return the reagents to the original stock tubes/bottles. Using disposable reservoirs when handling reagents can minimize the risk of contamination and is recommended.

1. **Preparation of 1X Wash Solution**
Mix 1 volume of the 20X Wash Solution with 19 volumes of distilled water.
2. **Preparation of 1X HRP-Conjugated Antibody #2**
Mix 1 volume of the 100X HRP-Conjugated Antibody #2 with 99 volumes of the Sulfonamide Antibody #2 Diluent.

ELISA Testing Protocol

Label the individual strips that will be used and aliquot reagents as the following example:

Component	Volume per Reaction	24 Reactions
Sulfonamide Antibody #1	100 µL	2.4 mL
1X HRP-Conjugated Antibody #2	150 µL	3.6 mL
1X Wash Solution	2.0 mL	48 mL
Stop Buffer	100 µL	2.4 mL
TMB Substrate	100 µL	2.4 mL

1. Add 50 µL of each SMX Standard in duplicate into different wells.
NOTE: Add standards to plate only in the order from low concentration to high concentration (white, yellow, orange, pink, purple, blue).
2. Add 50 µL of each sample in duplicate into different sample wells.
3. Add 100 µL of Sulfonamide Antibody # 1 to each well. Mix by sliding back and forth on a flat surface tapping against a hard force point for **1 minute**.
NOTE: A hard force point can be created by placing finger tips firmly against the flat surface.
4. Incubate the plate for 30 minutes at room temperature (20–25°C / 68–77°F).
NOTE: Keep plate covered and in the dark.
5. Shake out the contents of the wells into a waste container. Wash each well with 250 µL of 1X Wash Solution and dump out, repeat the wash step 3 times. After the last wash, invert the plate and tap the plate dry on paper towels.
NOTE: After the addition of wash buffer to each well, let the plate incubate at room temperature (20–25°C / 68–77°F) for 20 seconds to allow equilibration of the buffer.
6. Add 150 µL of freshly prepared 1X HRP-Conjugated Antibody # 2 to each well. Mix by sliding back and forth on a flat surface for **1 minute**. Incubate the plate for **30 minutes** at room temperature (20–25°C / 68–77°F) in the dark.
7. Shake out the contents of the wells into a waste container. Wash each well with 250 µL of 1X Wash Solution and dump out, repeat the wash step 3 times. After the last wash, invert the plate and tap the plate dry on paper towels.
8. Add 100 µL of TMB substrate to each well. Incubate the plate for **15 minutes** at room temperature (20–25°C / 68–77°F) in the dark. Mix by sliding back and forth on a flat surface for **1 minute** while incubating. Start the 15 minute incubation immediately after adding the substrate. Total incubation time for TMB Substrate is 15 minutes.
NOTE: Do not put any substrate back into the original container to avoid any potential contamination. Any substrate solution exhibiting coloration is indicative of deterioration and should be discarded. Covering the microtiter plate while incubating is recommended.
9. After incubation, add 100 µL of Stop Buffer to each well to stop the enzyme reaction.
10. Read the plate as soon as possible following the addition of Stop Buffer on a plate reader with a 450 nm wavelength.
NOTE: Before reading, use a lint-free wipe on the bottom of the plate to ensure no moisture or fingerprints interfere with the readings.


SULFONAMIDE CONCENTRATION CALCULATIONS

A standard curve can be constructed by plotting the mean relative absorbance (%) obtained from each reference standard against its concentration in ng/mL on a logarithmic curve.

$$\text{Relative absorbance (\%)} = \frac{\text{absorbance standard (or sample)} \times 100}{\text{absorbance zero standard}}$$

Use the mean relative absorbance values for each sample to determine the corresponding concentration of the tested drug in ng/mL from the standard curve. Neogen's Veratox Software for Windows is available to calculate test results. Please contact your local distributor or foodsafety@neogen.com for further information.

The following figure is a typical sulfonamide standard curve.

PERFORMANCE CHARACTERISTICS**Sensitivity (Assay Range)**

Sample Type	Assay Range (ng/g or ppb)
Kidney/muscle	10 – 1,000

Specificity (Cross-Reactivity)

Analytes	Cross-Reactivity (%)
Sulfadimethoxine	100
Sulfapyridine	>100
Sulfamerazine	>100
Sulfamethazine	>100
Sulfamethoxyipyridazine	>100
Sulfamethizole	>100
Sulfadiazine	100
Sulfamethoxazole	>100
Sulfachlorpyridazine	>100
Sulfaquinoxaline	71
Sulfathiazole	71

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

OTHER RESIDUE DIAGNOSTICS KITS

RT-96-TT-AMS	ALERT for Ractopamine — qualitative microwell assay, 96 wells
DR021	Veratox for Clenbuterol — range 0.08 – 1.28 ppb, 96 wells
DR107	Veratox for Florfenicol — range 0.15 – 100 ppm, 96 wells
DR024	Veratox for Fluoroquinolone — range 0.5 – 8 ppb, 96 wells
DR026	Veratox for Tylosin — range 10 – 250 ppb, 96 wells
DR081	Veratox for Oxytetracycline — range 3.6–324 ppb, 96 wells
DR073	Veratox for Avermectins — range 6.4 – 300 ppb, 96 wells

North America
Neogen Headquarters
800/234-5333 (USA/Canada)
foodsafety@neogen.com
foodsafety.neogen.com

Europe, Middle East and Africa
Neogen Europe
+ 44 (0) 1292 525 600
info_uk@neogeneurope.com
www.neogeneurope.com

Mexico
Neogen Latinoamerica
+52 (55) 5254-8235
informacion@neogenlac.com
www.neogenlac.com

Brazil
Neogen do Brasil
+55 19 3935.3727
info@neogendobrasil.com.br
www.neogendobrasil.com.br

China
Neogen Bio-Scientific Technology
+86 21 6271 7013
info@neogenchina.com.cn
www.neogenchina.com.cn

India
Neogen Food and Animal Security
+91 484 2306598, 2301582
info@neogenindia.com
www.neogenindia.com

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