

*Read instructions carefully before starting test*

# **Veratox**<sup>®</sup> **for Oxytetracycline**

*Refrigerate at 2–8°C (35–46°F) • Do not freeze*

## **OXYTETRACYCLINE**

Oxytetracycline is a broad spectrum antibiotic which is part of the Tetracycline class of antibiotics. Oxytetracycline is FDA approved for use in food producing animals and is frequently employed in animal production for its excellent antibacterial and pharmacokinetic properties.

## **INTENDED USE**

Veratox<sup>®</sup> for Oxytetracycline is a competitive enzyme immunoassay for the quantitative analysis of Oxytetracycline in muscle.

Veratox for Oxytetracycline enables international and government regulatory agencies, food manufacturers and processors, as well as quality assurance organizations, to detect Oxytetracycline in muscle and to satisfy customer concerns about food safety. Veratox for Oxytetracycline is an accurate, rapid, cost-effective and easy-to-use diagnostic tool compared to instrument based methods like HPLC and GC.

## **ASSAY PRINCIPLES**

The method is based on a competitive colorimetric ELISA assay. The drug of interest has been coated in the plate wells. During the analysis, sample is added along with the primary antibody specific for the target drug. If the target is present in the sample, it will compete for the antibody, thereby preventing the antibody from binding to the drug attached to the well. The secondary antibody, tagged with a peroxidase enzyme, targets the primary antibody that is complexed to the drug coated on the plate wells. The resulting color intensity, after addition of substrate, has an inverse relationship with the target concentration in the sample.

Veratox for Oxytetracycline has the capacity for 96 determinations or testing of 36 samples in duplicate (assuming 24 wells for standards). Return any unused microwells to the foil bag and reseal them with the desiccant provided in the original package.

## **STORAGE REQUIREMENTS**

Store kit at 2–8°C (35–46°F). Some components should be stored at frozen temperatures if not used within 1 month. The shelf life of the kit is 12 months when properly stored.

## MATERIALS PROVIDED

| Kit Contents  | Amount   | Storage                   |
|---|--|---------------------------|
| Oxytetracycline-coated Plate  | 1 x 96-well plate (8 wells x 12 strips)  | 2–8°C                     |
| <b>Empty tubes for Standards:</b><br>Negative control (white cap tube)<br>0.15 ng/mL (yellow cap tube)<br>0.75 ng/mL (orange cap tube)<br>1.5 ng/mL (pink cap tube)<br>4.5 ng/mL (purple cap tube)<br>13.5 ng/mL (blue cap tube)<br><b>Oxytetracycline Spiking Stock</b><br>50,000 ng/mL (red cap tube) | 8 empty dark plastic tubes total (2 extra included for dilution). Reuse tubes for standard preparation. See page 4 (reagent preparation, part C) for further directions.<br><br>0.8 mL | 2–8°C<br><br><b>-20°C</b> |
| Oxytetracycline Antibody #1   | 12 mL  | 2–8°C                     |
| 100X HRP-Conjugated Antibody #2   | 300 µL   | 2–8°C*                    |
| Antibody #2 Diluent**   | 20 mL  | 2–8°C                     |
| 20X Wash Solution**   | 28 mL  | 2–8°C                     |
| Stop Buffer**   | 14 mL  | 2–8°C                     |
| TMB Substrate**   | 12 mL  | 2–8°C                     |
| 5X OXYTET Extraction Buffer   | 2 x 25 mL  | 2–8°C                     |
| 10X TET Sample Diluent  | 28 mL  | 2–8°C                     |
| TET Balance Buffer Concentrate  | 5 mL   | 2–8°C                     |

\* If the kit will be unused for over 1 month, store 100X HRP-Conjugated Antibody #2 at -20°C or in a freezer.

\*\*These components (within their expiry) can be used interchangeably with other Veratox for Oxytetracycline kits.

## MATERIALS REQUIRED BUT NOT PROVIDED

1. Microtiter plate reader with a 450 nm filter (Neogen item 9303)
2. 20–200 µL pipettes (Neogen item 9276)
3. Multichannel pipette: 50–300 µL (optional) (Neogen item 9385)
4. 1000 µL pipettor (Neogen 9337)
5. Vortex Mixer (Neogen 9494)
6. Centrifuge 4,000 x g
7. Timer (Neogen item 9426)
8. Wash bottle (Neogen item 9400)
9. Paper towels or equivalent absorbent material
10. Veratox software (Neogen item 9305)
11. Lab station (Neogen item 9481)
12. Distilled or deionized water
13. Blender or food processor to homogenize sample

## PRECAUTIONS

Neogen strongly recommends that you read the following precautions to ensure your full awareness of ELISA techniques and other details you should pay close attention to when running the assays. More information can also be found in the Technical Solutions & Troubleshooting section. Periodically, optimizations and revisions are made to the kit and manual. Therefore, it is important to follow the protocol included with the kit. If you need further assistance, you may contact Neogen at 800/234-5333 or 517/372-9200.

1. The standards contain Oxytetracycline. Handle with particular care.
2. Do not use the kit past the expiration date.
3. Do not intermix reagents from different kits or lots except for components with the same part numbers within their expiration dates. ANTIBODIES AND PLATES ARE KIT- AND LOT-SPECIFIC.
4. Try to maintain a laboratory temperature of 20–25°C (68–77°F). Avoid running assays under or near air vents, as this may cause excessive cooling, heating and/or evaporation. Also, do not run assays in direct sunlight, as this may cause excessive heat and evaporation. Cold bench tops should be avoided by placing several layers of paper towel under the assay plates during incubation.
5. Make sure you are using only distilled or deionized water since water quality is very important.
6. Incubations of assay plates should be timed as precisely as possible. Be consistent when adding standards to the assay plate. Add your standards first and then your samples.
7. Follow proper pipetting techniques, including priming tips by filling and dispensing solution once before use.
8. Add standards to plate only in the order from low concentration to high concentration as this will minimize the risk of compromising the standard curve.
9. Always refrigerate plates in sealed bags with a desiccant to maintain stability. Prevent condensation from forming on plates by allowing them to equilibrate to room temperature (20–25°C / 68–77°F) while in the sealed bag.

### Tips for Optimal ELISA Performance

1. Allow the entire kit to equilibrate at room temperature for at least two hours before starting any ELISA assay.
2. Avoid light as much as possible during sample preparation and ELISA assay.
3. For plate washing steps: After addition of 250 µL wash buffer to the wells, incubate the plate for 20–30 seconds; shake the plate gently before pouring out the wash buffer. Repeat this procedure for each of the three washes.
4. Incubate the plate for 55 minutes at room temperature (20–25°C / 68–77°F).

## SAMPLE PREPARATION

Be sure samples are properly stored. In general, samples should be refrigerated at 2–4°C for no more than 48 hours. Freeze samples to a minimum of -20°C if they need to be stored for a longer period. Frozen samples can be thawed at room temperature (20–25°C / 68–77°F) or in a refrigerator before use.

### A. Preparation of 1X OXYTET Extraction Buffer

Mix 1 volume of 5X OXYTET Sample Extraction Buffer with 4 volumes of distilled water.

### B. Preparation of 1X TET Sample Diluent

Mix 1 volume of 10X TET Sample Diluent with 9 volumes of distilled water.

### C. Preparation of 1X OXYTET Sample Balance Buffer

Combine 1.5 mL of TET Balance Buffer Concentrate and 1 mL of distilled water.

## SAMPLE EXTRACTION

### Muscle

1. To 1 g of a room temperature homogenized sample in a 15 mL plastic conical tube, add 5.0 mL of 1X OXYTET Extraction Buffer.
2. Vortex sample manually for **3 minutes**. Tilt the tube at 45° or on its side to aid in the breaking up of the pellet. Incubate sample for 5 minutes at 60°C, then vortex sample for an additional **1 minute**.
3. Centrifuge sample for **10 minutes** at 4,000 x g.
4. Transfer 100 µL of the supernatant to 75 µL of 1X OXYTET Sample Balance Buffer in a 2 mL tube. Briefly vortex the tube, then add 325 µL of 1X TET Sample Diluent.
5. Vortex sample well for **1 minute**. Incubate sample for 10 minutes at room temperature (20–25°C / 68–77°F), then vortex sample for an additional **1 minute**.
6. Centrifuge sample for **5 minutes** at maximum speed (~ 12,000 x g) in a micro-centrifuge.
7. Use 75 µL per well in the assay.

**NOTE:** Dilution factor: 24.

## REAGENT PREPARATION

**IMPORTANT:** All reagents should be brought up to room temperature before use (1–2 hours at 20–25°C / 68–77°F). Make sure you read the Precautions section on page 3. Solutions should be prepared prior to running the ELISA test. All reagents should be mixed by gently inverting or swirling prior to use. Prepare volumes that are needed for the number of wells being run. Do not return the reagents to the original stock tubes/bottles. Using disposable reservoirs when handling reagents can minimize the risk of contamination and is recommended.

### A. Preparation of 1X Wash Solution

Mix 1 volume of the 20X Wash Solution with 19 volumes of distilled water.

### B. Preparation of 1X HRP-Conjugated Antibody #2

Mix 1 volume of the 100X HRP-Conjugated Antibody #2 with 99 volumes of Antibody #2 Diluent.

### C. Preparation of Oxytetracycline Working Standards

Bring the 50,000 ng/mL spiking solution to room temperature then vortex well for 1 minute. The Oxytetracycline standards must be prepared fresh weekly, or each time an ELISA is performed. To make the working standards, serially dilute the 50,000 ng/mL spiking stock with 1X TET Sample Diluent. Vortex each standard for at least 1 minute before continuing to the next dilution. The working standards remain stable for up to 1 week when stored at 4°C. If standards have been stored for more than one week, dump out contents and keep the dark plastic tubes supplied with kit for preparing new standards in the future.

| Oxytetracycline Source       | Volume of Oxytetracycline Source | Volume of 1X TET Standard Diluent | Working Standards             |
|------------------------------|----------------------------------|-----------------------------------|-------------------------------|
| 50,000 ppb<br>(red cap tube) | 20 µL                            | 980 µL                            | 1,000 ppb                     |
| 1,000 ppb                    | 100 µL                           | 900 µL                            | 100 ppb                       |
| 100 ppb                      | 135 µL                           | 865 µL                            | 13.5 ppb<br>(blue cap tube)   |
| 13.5 ppb                     | 350 µL                           | 700 µL                            | 4.5 ppb<br>(purple cap tube)  |
| 4.5 ppb                      | 350 µL                           | 700 µL                            | 1.5 ppb<br>(pink cap tube)    |
| 1.5 ppb                      | 500 µL                           | 500 µL                            | 0.75 ppb<br>(orange cap tube) |
| 0.75 ppb                     | 200 µL                           | 800 µL                            | 0.15 ppb<br>(yellow cap tube) |
| N/A                          | N/A                              | 500 µL                            | Negative Control vial         |

## TEST PROCEDURE

Label the individual strips that will be used and prepare reagent aliquots based on the number of reactions:

| Component                     | Volume Per Reaction | 24 Reactions |
|-------------------------------|---------------------|--------------|
| Oxytetracycline Antibody #1   | 100 $\mu$ L         | 2.4 mL       |
| 1X HRP-Conjugated Antibody #2 | 150 $\mu$ L         | 3.6 mL       |
| 1X Wash Solution              | 2.5 mL              | 60 mL        |
| Stop Buffer                   | 100 $\mu$ L         | 2.4 mL       |
| TMB Substrate                 | 100 $\mu$ L         | 2.4 mL       |

1. Add 75  $\mu$ L of each Oxytetracycline Standard.  
**NOTE:** Add standards to plate in order from low concentration to high concentration, starting with the negative control (white, yellow, orange, pink, purple, blue, red.)
2. Add 75  $\mu$ L of each sample in duplicate into different sample wells. Because of the nature of some samples, a precipitate may form after sample preparation. If a precipitate has formed, vortex the sample thoroughly before adding to wells.
3. Add 100  $\mu$ L of Antibody #1 and mix by sliding back and forth on a flat surface for **1 minute**.
4. Incubate the plate for **55 minutes** at room temperature (20–25°C / 68–77°F).
5. Wash the plate 3 times with 250  $\mu$ L of 1X Wash Solution. After the last wash, invert the plate and gently tap the plate dry on paper towels.  
**NOTE:** Perform the next step immediately after plate washings. Do not allow the plate to air dry between working steps.
6. Add 150  $\mu$ L of the 1X Antibody #2 solution and mix by sliding back and forth on a flat surface for **1 minute**. During sliding, tap well holder edge against a firm item—such as a hand/finger on the flat surface—be careful not to spill liquid out of the wells. Incubate the plate for **20 minutes** at room temperature (20–25°C / 68–77°F).  
**NOTE:** Avoid direct sunlight and cold bench tops during the incubation. Covering the microtiter plate while incubating is recommended.
7. Wash the plate 3 times with 250  $\mu$ L of 1X Wash Solution. After the last wash, invert the plate and gently tap the plate dry on paper towels.  
**NOTE:** Perform the next step immediately (within 1 minute) after plate washings. Do not allow the plate to air dry between working steps.
8. Add 100  $\mu$ L of TMB Substrate to each well. Incubate the plate for **15 minutes** at room temperature (20–25°C / 68–77°F). Time the reaction immediately after adding the substrate. Mix by sliding back and forth on a flat surface for **1 minute** while incubating.  
**NOTE:** Do not put any substrate back into the original container to avoid any potential contamination. Covering the microtiter plate while incubating is recommended.
9. After incubation, add 100  $\mu$ L of Stop Buffer to stop the enzyme reaction.
10. Read the plate as soon as possible following the addition of Stop Buffer on a plate reader with 450 nm primary filter and 630 nm differential filter wavelengths. (Before reading, use a lint-free wipe on the bottom of the plate to ensure no moisture or fingerprints interfere with the readings.)


## INTERPRETATION OF RESULTS

A standard curve can be constructed by plotting the mean relative absorbance (%) obtained from each reference standard against its concentration in ng/mL on a logarithmic curve.

$$\text{Relative absorbance (\%)} = \frac{\text{absorbance standard (or sample)} \times 100}{\text{absorbance zero standard}}$$

Use the mean relative absorbance values for each sample to determine the corresponding concentration of the tested drug in ng/mL from the standard curve. Veratox software is available upon request to evaluate the results. Please contact Neogen for further information. The following figure is a typical Oxytetracycline standard curve:

**Oxytetracycline Standard Curve**



### Sensitivity (Detection Limit)

| Sample Type | Detection Limit (ng/g or ppb) |
|-------------|-------------------------------|
| Muscle      | 99                            |

### Specificity (Cross-Reactivity)

| Analytes          | Cross-Reactivity (%) |
|-------------------|----------------------|
| Oxytetracycline   | 100                  |
| Doxycycline       | 100                  |
| Tetracycline      | >100                 |
| Chlortetracycline | >100                 |



## CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

## SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at [foodsafety.neogen.com](http://foodsafety.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

## TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit [www.neogen.com/Corporate/termsconditions.html](http://www.neogen.com/Corporate/termsconditions.html).

## WARRANTY

Neogen makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### OTHER RESIDUE DIAGNOSTICS KITS

|                     |  |
|---------------------|--|
| <b>RT-96-TT-AMS</b> | <b>ALERT for Ractopamine</b> — qualitative microwell assay, 96 wells |
| <b>9551</b>         | <b>Veratox for Chloramphenicol</b> — range 10–1,000 ppt, 96 wells    |
| <b>8416</b>         | <b>Veratox for Malachite Green</b> — range 1–4 ppb, 48 wells         |
| <b>DR021</b>        | <b>Veratox for Clenbuterol</b> — range 0.08–1.28 ppb, 96 wells       |
| <b>DR107</b>        | <b>Veratox for Florfenicol</b> — range 0.15–100 ppm, 96 wells        |
| <b>DR073</b>        | <b>Veratox for Avermectins</b> — range 6.4–300 ppb, 96 wells         |



#### North America

##### Neogen Headquarters

620 Leshar Place, Lansing, MI 48912 USA  
800/234-5333 (USA/Canada) or  
517/372-9200

Fax: 517/372-2006  
[foodsafety@neogen.com](mailto:foodsafety@neogen.com)  
[foodsafety.neogen.com](http://foodsafety.neogen.com)

#### Europe, Middle East and Africa

##### Neogen Europe

The Dairy School, Auchincruive, Ayr  
KA6 5HU Scotland, UK  
+ 44 (0) 1292 525 600

Fax: + 44 (0) 1292 525 601  
[info\\_uk@neogeneurope.com](mailto:info_uk@neogeneurope.com)  
[foodsafety.neogen.com](http://foodsafety.neogen.com)

#### Mexico

##### Neogen Latinoamérica

Prolongación 5 de Mayo #27  
Col. Parque Industrial Naucalpan  
Naucalpan, Estado de Mexico C.P. 53489  
+52 (55) 5254-8235, +52 (55) 5203-0111,  
+52 (55) 5531-2837

Fax: +52 (55) 5531-1647  
[informacion@neogenlac.com](mailto:informacion@neogenlac.com)  
[foodsafety.neogen.com](http://foodsafety.neogen.com)

#### Brazil

##### Neogen do Brasil

Rua: Alberto Guizo 760, Distrito Industrial João  
Narezzi, Indaiatuba – SP Brasil, Cep: 13.347-402  
Tel: +55 19 3935.3727

[info@neogendobrasil.com.br](mailto:info@neogendobrasil.com.br)  
[foodsafety.neogen.com](http://foodsafety.neogen.com)

© Neogen Corporation, 2016. Neogen, Veratox and K-Blue are registered trademarks of Neogen Corporation, Lansing, MI. All other brand and product names are trademarks or registered trademarks of their respective companies. Patent: <http://www.neogen.com/Corporate/patents.html>

DR081B

V-Oxytetracycline\_0816