Read instructions carefully before starting test

Veratox[®]

for Lincomycin

*Refrigerate at 2–8°C (35–46°F)
For extended storage, review kit instructions.*

LINCOMYCIN

Lincomycin is a lincosamide antibiotic that comes from the actinomycete *Streptomyces lincolnensis*. Although similar in structure, antibacterial spectrum, and mechanism of action to macrolides, lincomycin is also effective against other organisms including actinomycetes, mycoplasma, and some species of *Plasmodium*.

INTENDED USE

Veratox[®] for Lincomycin is a competitive enzyme immunoassay for the quantitative analysis of lincomycin in muscle. Veratox for Lincomycin enables international and government regulatory agencies, food manufacturers and processors, as well as quality assurance organizations, to detect lincomycin in muscle in response to customer concerns about food safety.

ASSAY PRINCIPLES

The method is based on a competitive colorimetric ELISA assay. Lincomycin, the drug of interest, has been coated in the plate wells. During the analysis, sample is added along with the lincomycin-horseradish peroxidase (Lincomycin-HRP) conjugate. If the target is present in the sample, it will compete for the antibody, thereby preventing the antibody from binding to the drug attached to the well. The resulting color intensity, after addition of the HRP substrate (TMB), has an inverse relationship with the target concentration in the sample.

Veratox for Lincomycin has the capacity for 96 determinations or testing of 36 samples in duplicate (assuming 24 wells for standards). Return any unused microwells to the foil bag and reseal them with the desiccant provided in the original package.

STORAGE REQUIREMENTS

Store kit at 2–8°C (35–46°F). Some components should be stored at frozen temperatures if not used within 1 month. The shelf life of the kit is 12 months when properly stored.

MATERIALS PROVIDED

Kit Contents	Amount	Storage
Lincomycin Ab-coated Microtiter Plate	1 x 96-well plate (8 wells x 12 strips)	2–8°C
Lincomycin Standards: Negative control (white cap tube) 0.1 ng/mL (yellow cap tube) 0.25 ng/mL (orange cap tube) 0.6 ng/mL (pink cap tube) 1.5 ng/mL (purple cap tube) 3.75 ng/mL (blue cap tube) 1,000 ng/mL (spiking, red cap tube)	1.8 mL 1.8 mL 1.8 mL 1.8 mL 1.8 mL 1.8 mL 1.8 mL	2–8°C
Lincomycin-HRP Conjugate	6 mL	2–8°C*
200X Sample Extraction Buffer	1 mL	2–8°C
20X Wash Solution	28 mL	2–8°C
Stop Buffer	14 mL	2–8°C
TMB Substrate	12 mL	2–8°C
Sample Balance Buffer	0.5 mL	2–8°C
10X Lincomycin Dilution Buffer	6 mL	2–8°C

*If the kit will be unused for over 1 month, store Lincomycin-HRP Conjugate at -20°C or in a freezer.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Microtiter plate reader with a 450 nm filter (Neogen item 9303)
2. 20–200 µL pipettes (Neogen item 9276)
3. Multichannel pipette: 50–300 µL (optional) (Neogen item 9385)
4. 1000 µL pipettor (Neogen 9337)
5. Vortex Mixer (Neogen 9494)
6. Centrifuge 4,000 x *g*
7. Timer (Neogen item 9426)
8. Wash bottle (Neogen item 9400)
9. Paper towels or equivalent absorbent material
10. Veratox software (Neogen item 9305)
11. Lab station (Neogen item 9481)
12. Distilled or deionized water
13. Blender or food processor to homogenize sample
14. Heater block or incubator capable of at least 70°C

PRECAUTIONS

Neogen strongly recommends that you read the following warnings and precautions to ensure your full awareness of ELISA techniques and other details you should pay close attention to when running the assays. Periodically, optimizations and revisions are made to the kit and insert. Therefore, it is important to follow the protocol included with the kit. If you need further assistance, you may contact your local distributor or Neogen at foodsafety@neogen.com.

1. The standards contain lincomycin. Handle with particular care.
2. Do not use the kit past the expiration date.
3. Do not intermix reagents from different kits or lots except for components with the same lot numbers within their expiration dates. ANTIBODIES AND PLATES ARE KIT-AND-LOT SPECIFIC.
4. Try to maintain a laboratory temperature of 20–25°C (68–77°F). Avoid running assays under or near air vents, as this may cause excessive cooling, heating and/or evaporation. Also, do not run assays in direct sunlight, as this may cause excessive heat and evaporation. Cold bench tops should be avoided by placing several layers of paper towel or some other insulation material under the assay plates during incubation.
5. Make sure you are using only distilled or deionized water since water quality is very important.
6. Follow proper pipetting techniques, including priming tips by filling and dispensing solution once before use.
7. Incubations of assay plates should be timed as precisely as possible. Be consistent when adding standards to the assay plate. Add your standards first and then your samples.
8. Add standards to plate only in the order from low concentration to high concentration as this will minimize the risk of compromising the standard curve.
9. Always refrigerate plates in sealed bags with a desiccant to maintain stability. Prevent condensation from forming on plates by allowing them to equilibrate to room temperature (20–25°C / 68–77°F) while in the packaging.

SAMPLE PREPARATION AND EXTRACTION

Be sure samples are properly stored. In general, samples should be refrigerated at 2–4°C for no more than 48 hours. Freeze samples to a minimum of -20°C if they need to be stored for a longer period. Frozen samples can be thawed at room temperature (20–25°C / 68–77°F) or in a refrigerator before use. Preparation protocols for samples other than below can be made available upon request. Please contact your local distributor or Neogen at foodsafety@neogen.com.

1. Preparation of 1X Sample Extraction Buffer

Mix 1 volume of 10X Sample Extraction Buffer with 9 volumes of distilled water.

2. Preparation of 1X Lincomycin Dilution Buffer

Mix 1 volume of the 10X Lincomycin Dilution Buffer with 9 volumes of distilled or deionized water.

Muscle

1. Remove any fat from sample. Homogenize sample with a suitable mixer.
2. To 0.5 g of the homogenized sample in a 15 mL plastic conical tube, add 4.5 mL of 1X Sample Extraction Buffer.
3. Vortex sample manually for 3 minutes at maximum speed.
4. Centrifuge sample for 5 minutes at 4,000 x *g* at room temperature (20–25°C / 68–77°F).
5. Transfer 500 µL of the supernatant to a new 1.5 mL tube. Avoid any fat on the surface of the sample.
6. Incubate sample for 20 minutes at 60°C either in a heater block or incubator.
7. Vortex sample for 2 minutes.
8. Incubate sample for 20 minutes at -20 °C.
9. Centrifuge sample for 5 minutes at 4,000 x *g* at room temperature (20–25°C / 68–77°F).
10. Transfer 200 µL of the supernatant to a new 1.5 mL tube containing 5 µL of Sample Balance Buffer, then add 1.2 mL of 1X Lincomycin Dilution Buffer (2 x 600 µL portions).
11. Vortex sample for 2 minutes. Incubate sample for 10 minutes. Vortex for an additional 30 seconds.
12. Use 100 µL per well for the assay.

NOTE: Dilution factor: 70.

TEST PROCEDURE

Reagent Preparation

IMPORTANT: All reagents should be brought up to room temperature before use (1–2 hours at 20–25°C / 68–77°F); Make sure you read “Precautions” section on page 3. Solutions should be prepared just prior to ELISA test. All reagents should be mixed by gently inverting or swirling prior to use. Prepare volumes that are needed for the number of wells being run. Do not return the reagents to the original stock tubes/bottles. Using disposable reservoirs when handling reagents can minimize the risk of contamination and is recommended.

1. Preparation of 1X Wash Solution

Mix 1 volume of the 20X Wash Solution with 19 volumes of distilled water.

ELISA Testing Protocol

Label the individual strips that will be used and aliquot reagents as the following example:

Component	Volume per Reaction	24 Reactions
Lincomycin-HRP	50 µL	1.2 mL
1X Wash Solution	2.0 mL	48 mL
Stop Buffer	100 µL	2.4 mL
TMB Substrate	100 µL	2.4 mL

1. Add 100 µL of each Lincomycin Standard in duplicate into different wells. Add standards to plate only in the order from low concentration to high concentration.
2. Add 100 µL of each sample in duplicate into different sample wells.
3. Add 50 µL of Lincomycin-HRP Conjugate to each well. Mix by sliding back and forth on a flat surface for **1 minute**.
4. Incubate the plate for **30 minutes** at room temperature (20–25°C / 68–77°F).
5. Wash the plate 3 times with 250 µL 1X Wash Solution. After the last wash, invert the plate and gently tap the plate dry on paper towels.
NOTE: Perform the next step immediately after drying the plate. Do not allow the plate to air dry between working steps.
6. Add 100 µL of TMB substrate to each well. Mix well by sliding back and forth on a flat surface manually for 1 minute. Incubate the plate for 20 minutes at room temperature (20–25°C / 68–77°F) in the dark. Time the reaction immediately after adding the substrate.
NOTE: Do not put any substrate back to the original container to avoid any potential contamination. Any substrate solution exhibiting coloration is indicative of deterioration and should be discarded.
7. After incubation, add 100 µL of Stop Buffer to stop the enzyme reaction.
8. Read the plate as soon as possible following the addition of Stop Buffer on a plate reader with a 450 nm wavelength filter.
NOTE: Before reading, use a lint-free wipe on the bottom of the plate to ensure no moisture or fingerprints interfere with the readings.


LINCOMYCIN CONCENTRATION CALCULATIONS

A standard curve can be constructed by plotting the mean relative absorbance (%) obtained from each reference standard against its concentration in ng/mL on a logarithmic curve.

$$\text{Relative absorbance (\%)} = \frac{\text{absorbance standard (or sample)} \times 100}{\text{absorbance zero standard}}$$

Use the mean relative absorbance values for each sample to determine the corresponding concentration of the tested drug in ng/mL from the standard curve. Neogen's Veratox Software for Windows is available to calculate test results. Please contact your local distributor or foodsafety@neogen.com for further information.

The following figure is a typical lincomycin standard curve.

PERFORMANCE CHARACTERISTICS

Sensitivity (Detection Limit)

Sample Type	Assay Range (ng/g or ppb)
Muscle	7 – 262.5

Specificity (Cross-Reactivity)

Analytes	Cross-Reactivity (%)
Lincomycin	100

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

OTHER RESIDUE DIAGNOSTICS KITS

RT-96-TT-AMS	ALERT for Ractopamine — qualitative microwell assay, 96 wells
9551	Veratox for Chloramphenicol — range 10–1,000 ppt, 96 wells
8416	Veratox for Malachite Green — range 1–4 ppb, 48 wells
DR073	Veratox for Avermectins — range 6.4–300 ppb, 96 wells
DR021	Veratox for Clenbuterol — range 0.08–1.28 ppb, 96 wells
DR081	Veratox for Oxytetracycline — range 3.6–324 ppb, 96 wells

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Neogen Food and Animal Security

+91 484 2306598, 2301582
info@neogenindia.com
www.neogenindia.com

© Neogen Corporation, 2016. Neogen, Veratox and K-Blue are registered trademarks of Neogen Corporation, Lansing, MI. All other brand and product names are trademarks or registered trademarks of their respective companies.

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