Read instructions carefully before starting test

Reveal[®] **for** **Histamine**

STORE AT ROOM TEMPERATURE (18–30°C)

HISTAMINE

High levels of histamine may develop in a variety of fish species as they decompose. These species include tuna, mahi-mahi, marlin, bluefish, sardines, anchovy, bonito, herring and mackerel. Ingestion of histamine may cause scombroid poisoning in humans, which may lead to a variety of symptoms, including rash, nausea, vomiting, diarrhea, hypotension, heart palpitations and muscle weakness. Paralysis and death have also been reported in cases of scombroid poisoning.

Because of its potential for human illness, the U.S. Food and Drug Administration (FDA) has ruled that extensive refrigeration records and/or histamine testing must be included in a hazard analysis and critical control point (HACCP) program for relevant fish species. The FDA has set an action level of 50 parts per million (ppm) for histamine in domestic and imported fish.

INTENDED USE

Reveal[®] for Histamine is intended for the qualitative analysis of histamine in scombroid species of fish, such as tuna and mahi-mahi.

INTENDED USER

The test kit is designed for use by quality control personnel and others familiar with histamine analysis in fish. Since technique is very important, operators should be trained by a NEOGEN representative or someone who has completed NEOGEN training.

ASSAY PRINCIPLES

Reveal for Histamine is a single-step lateral flow assay based on a competitive immunoassay format. The extract is wicked through a reagent zone, which contains histamine-specific antibodies conjugated to colored particles. If histamine is present, it will be captured by the particle-antibody complex. The histamine-antibody-particle complex then is wicked onto a membrane, which contains a zone of histamine conjugated to a protein carrier. This zone captures only uncomplexed antibody-particle conjugate, allowing the colored particles to concentrate and form a visible line. The membrane also contains a control zone where the antibody-particle complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of histamine, ensuring the strip is functioning properly. In addition, the control line also is used as a reference for comparison to the test zone for visual interpretation of the test results. If histamine is present above a pre-determined limit, a line lighter than the control line or no line will form. Therefore, a positive sample will result in a light line or absence of a line in the test zone. Alternatively, a negative sample will form a line of the same or darker intensity as the control line in the test zone.

STORAGE REQUIREMENTS

Store kit components at room temperature (18–30°C, 64–86°F) to ensure full shelf life. Test strips should remain capped in their original sample tubes until used to ensure optimal performance.

MATERIALS PROVIDED

1. 25 histamine test strips
2. 25 sample cups
3. 1 bottle sample diluent
4. 25 red dilution cups

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Graduated cylinder (NEOGEN item 9368)
2. Sample extraction bottles, 250 mL (NEOGEN item 9399)
3. Sample collection tubes (NEOGEN item 9421)
4. Blender (NEOGEN item 9493, 9477)
5. Scale capable of weighing 10–50 grams (NEOGEN item 9427)
6. Pipettor, 100 µL (NEOGEN items 9272, 9278, 9276)
7. Pipette tips (NEOGEN items 9410, 9407, 9417)
8. Pipettor, 1000 µL, (NEOGEN item 9337)
9. Pipette tips (NEOGEN items 9464, 9487)
10. Paper towels or equivalent absorbent material
11. Rack for Reveal sample cups (NEOGEN item 9475)
12. Timer (NEOGEN item 9426)
13. Waterproof marker
14. Distilled or deionized water

PRECAUTIONS

1. The test strips must remain inside the stay-dry tube before use.
2. Store test kit at room temperature (18–30°C, 64–86°F) when not in use.
3. Do not use kit components beyond expiration date.
4. Treat all used liquids, including sample extract, and labware as if contaminated with histamine. Gloves and other protective apparel should be worn at all times.
5. To avoid cross-contamination, use clean equipment for each sample, and thoroughly wash all equipment between samples.
6. Glassware should not be used for extraction purposes. Histamine may adhere to glass, therefore using glassware may affect test results.

SAMPLE PREPARATION

A. **Canned tuna: AOAC 937.07b**

Place entire contents of can/pouch (meat and liquid) into a blender. Blend until homogenous. Store samples at 2–8°C (35–46°F) until analyzed.

B. **Fresh or thawed frozen raw fish: AOAC 937.07a**

Clean and eviscerate three fish. Cut three cross-sectional pieces 2.5 cm (1 inch) thick, from back of the pectoral fin, halfway to the vent and one posterior to the vent. Debone slices and blend or grind combined samples until homogenous. Store samples at 2–8°C (35–46°F) until analyzed.

SAMPLE EXTRACTION

1. Add 10 g of the homogenous mixture to a clean disposable extraction bottle containing 190 mL of distilled or deionized water.
2. Tightly cap and vigorously shake the bottle for **15–20 seconds** to suspend the fish tissue in the water.
3. Wait approximately **5 minutes**, and then shake the bottle for **15–20 seconds** to resuspend the fish tissue.
4. Wait an additional **5 minutes**, and again shake the bottle for **15–20 seconds** to resuspend the fish tissue. Allow the tissue to settle to the bottom of the bottle for about **30 seconds**. The sample is now ready for extract dilution.

TEST PROCEDURE

1. Remove the appropriate number of clear sample cups and red dilution cups and place in a sample cup rack.
2. Add 100 µL of sample extract to each red dilution cup.
3. Add 1000 µL of sample diluent to the red dilution cup. Mix well by pipetting up and down 5 times.
4. Transfer 100 µL of diluted sample extract to a clear sample cup.
5. Place a new Reveal for Histamine test strip with the sample end down into the sample cup and set timer for 5 minutes. Ensure the test strip comes into contact with liquid and begins to wick.
6. Allow the test strip to develop in the sample cup for 5 minutes.
7. Remove the strip and read results visually by holding the strip up to line of sight. Strip must be read within 30 seconds after the 5 minute incubation.

TEST INTERPRETATION

Negative result:

If the test line intensity is greater than or equal to the control line, the sample is below the detection limit (less than 50 ppm histamine).

NOTE: The presence of the control line, regardless of intensity, indicates a valid test.


Positive result:

If the test line intensity is less than the control line, the sample is greater than the detection limit (greater than 50 ppm histamine).

If the test line is absent, and the control line is present, the sample is interpreted as highly positive.

Invalid result:

If the control line is absent, the test is invalid.

NOTES:

1. Observations should be made immediately following the 5 minute development. Observations after 5 minutes may be inaccurate due to over-development of the device.
2. Positives can be confirmed by NEOGEN's Reveal® Q+ for Histamine, Veratox for Histamine or HPLC.

VALIDATED MATRICES

Tuna and mahi-mahi.

CUSTOMER SERVICE

NEOGEN Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all NEOGEN test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of NEOGEN's test kits, on NEOGEN's website at neogen.com, or by calling NEOGEN at 800.234.5333 or 517.372.9200.

TERMS AND CONDITIONS

For NEOGEN's full terms and conditions, please visit neogen.com/terms-and-conditions/.

WARRANTY

NEOGEN Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ergot alkaloids, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, coconut, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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