Read instructions carefully before starting test

 **Veratox[®]**
for Zearalenone
Quantitative Test

REFRIGERATE AT 2-8°C • DO NOT FREEZE

THE TOXIN

Zearalenone is primarily produced by the mold *Fusarium graminearum*, which also commonly produces deoxynivalenol. Hence, there is evidence that if zearalenone is detected, there is a high probability that other fusarial mycotoxins may be present. Zearalenone is classified as an estrogenic mycotoxin because it frequently causes estrogenic responses in animals.

When zearalenone-contaminated feed or grain is eaten by livestock, it can cause a wide variety of reproductive problems. In swine, it causes vulvovaginitis, low birth weights, fetal reabsorption, aborted pregnancies, reduced litter sizes, abnormal estrus and feminization of immature males. Zearalenone can delay the breeding process and cost the producer significant economic and physical losses.

Livestock producers are becoming increasingly aware of zearalenone problems, and have looked for ways to reduce risks related to contaminated feed.

The best protection against mycotoxins is monitoring for their presence in feeds and foods. That means testing all along the pathway from initial harvest of grains to the finished product.

INTENDED USE

Veratox for Zearalenone is intended for the quantitative analysis of zearalenone in such commodities as corn, wheat, oats, barley, pet food and DDGS.

INTENDED USER

The test kit is designed for use by quality control personnel and others familiar with food and feed possibly contaminated by zearalenone. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed the Neogen training.

ASSAY PRINCIPLES

Veratox for Zearalenone is a competitive direct ELISA in a microwell format which allows the user to obtain exact concentrations in parts per billion (ppb). Free zearalenone in the samples and controls is allowed to compete with enzyme-labeled zearalenone (conjugate) for the antibody binding sites. After a wash step substrate is added, which reacts with the bound conjugate to produce blue color. More blue color means less zearalenone. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form the standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of zearalenone.

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2-8°C (35-46°F).

MATERIALS PROVIDED

1. 48 antibody-coated microwells
2. 48 red-marked mixing wells
3. 5 yellow-labeled bottles of 0, 25, 75, 150 and 500 ppb zearalenone controls (see precautions for handling of methanol solution)
4. 1 blue-labeled bottle of zearalenone-HRP conjugate solution
5. 1 green-labeled bottle of K-Blue® Substrate solution
6. 1 red-labeled bottle of Red Stop solution
7. Directions for use

MATERIALS REQUIRED BUT NOT PROVIDED

1. Extraction materials (items c through e are available in kit form from Neogen, item 8052):
 - a. 70% methanol (ACS Grade, Neogen item 8055, 8056)
 - b. 250 mL graduated cylinder (Neogen item 9368)
 - c. Container with 125 mL capacity (Neogen item 9428)
 - d. Neogen filter syringes, Whatman #1 filter paper, or equivalent (Neogen item 9420, 9430)
 - e. Sample collection tubes (Neogen item 9421)
 - f. Sample dilution tubes 5-10mL (Neogen item 9468)
2. Agri-Grind grinder or equivalent (Neogen item 9401, 9453)
3. Scale capable of weighing 2-25 grams (Neogen item 9427)
4. Microwell reader with a 650 nm filter (Neogen item 9303)
5. Pipettor, 12-channel (Neogen item 9273)
6. Pipettor, 100 μ L (Neogen item 9290)
7. Pipette tips for 100 μ L and 12-channel pipettors (Neogen item 9410, 9407, 9417)
8. Paper towels or equivalent absorbent material
9. Plastic 1/2 gallon bucket for use as waste receptacle
10. Microwell holder (Neogen item 9402)
11. Timer (Neogen item 9426)
12. Waterproof marker
13. Wash bottle (Neogen item 9400)
14. 2 reagent boats for 12-channel pipettor (Neogen item 9450)
15. Distilled or deionized water

PRECAUTIONS

1. Methanol solution is highly flammable. Keep container tightly closed, and keep away from heat, sparks, open flame and those smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
2. Store test kit between 2-8°C (35-46°F) when not in use. Do not freeze.
3. Do not use kit components beyond expiration date.
4. Do not mix reagents from one kit serial with reagents from a different kit serial.
5. Do not run more than 24 wells at one time.
6. Follow proper pipetting techniques, including the proper priming of tips.
7. Incubation times other than those specified may give inaccurate results.
8. Kits should be at room temperature 18-30°C (64-86°F) prior to use.
9. Avoid prolonged storage of kits at ambient temperatures.
10. Treat all used liquids, including sample extract, and labware as if contaminated with zearalenone. Gloves and other protective apparel should be worn at all times.
11. To avoid cross-contamination, use clean pipette tips and glassware for each sample, and thoroughly detoxify and wash all glassware between samples.
12. Commodity extracts should have a pH of 6-8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH contact your Neogen representative or Technical Services.

PROCEDURAL NOTES

1. **Substrate:** K-Blue Substrate is ready for use. The substrate should be clear to a very light blue—discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until it is needed.
2. **Antibody wells:** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted and the test procedure is set to begin.

SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques. The sample should be ground and thoroughly mixed prior to proceeding with the extraction. Store samples at 2-8°C (35-46°F) until analyzed. **Note:** If you are using Neogen's Mycotoxin Extraction Kit, follow the instructions in that kit for the extraction procedure. If you are preparing your own extraction solution, continue with the instructions that follow.

1. If not using Neogen's prepared solution, prepare a 70% methanol solution by mixing 7 parts ACS Grade methanol with 3 parts distilled or deionized water for each sample to be tested.
2. Obtain a representative sample. Grind the entire sample so that at least 75% of the ground material passes through a 20 mesh sieve, the particle size of a fine instant coffee.
3. Extract at a ratio of 1 part sample with 5 parts 70% methanol. Example: Add 5 grams of ground sample to 25 mL of 70% methanol/water and shake vigorously for 3 minutes, or blend for 2 minutes.
4. Filter the extract by pouring at least 5 mL through a Whatman #1 filter (or Neogen filter syringe) and collecting the filtrate as a sample.
5. Dilute sample extract 1:5 in water by mixing 1 mL of sample extract with 4 mL of distilled water in a clean sample dilution tube.
6. The sample is now ready for testing.

TEST PROCEDURE – FOR QUANTITATION

Allow reagents to warm to room temperature 18-30°C (64-86°F) prior to use.

1. Remove 1 red-marked mixing well for each sample to be tested plus 5 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant and reseal the pack to protect the antibody. Mark one end of strip with a "1", and place strip in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Place 100 µL of conjugate from the blue-labeled bottle in each red-marked mixing well.
5. Using a new pipette tip for each, transfer 100 µL of controls and diluted samples to the red-marked mixing wells as shown below.

0	25	75	150	500	S1	S2	S3	S4	S5	S6	S7	Strip 1
S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19	Strip 2

6. Using a 12-channel pipettor, mix the liquid in the wells by pipetting it up and down 3 times. Transfer 100 µL to the antibody-coated wells. Mix by sliding the microwell holder back and forth on a flat surface for 10 seconds without splashing the reagents from the wells. Incubate for **5 minutes** at room temperature 18-30°C (64-86°F). Discard the red-marked mixing wells.
7. The initial reaction is now complete. Shake out the contents of the antibody wells.
8. Fill each antibody well with distilled or deionized water and dump them out. Repeat this step 5 times, then turn the wells upside down and tap out on a paper towel until the remaining water has been removed.
9. Pour the needed volume of substrate from the green-labeled bottle into the green-labeled reagent boat and, with new tips, pipette 100 µL of substrate into the wells and mix by sliding back and forth on a flat surface for 10 seconds. Incubate **5 minutes**. Discard remaining substrate and rinse the reagent boat with water.
10. Pour Red Stop solution from the red-labeled bottle (same volume as prepared for substrate) into the red-labeled reagent boat. Using the same pipette tips as were used to dispense substrate, add 100 µL Red Stop to each well and mix by sliding back and forth on a flat surface. Discard the tips.

11. Wipe the bottom of the microwells with a dry cloth or towel and read in a microwell reader using a 650 nm filter. Air bubbles should be eliminated, as they could affect analytical results. Results should be read within 20 minutes after the addition of Red Stop.
12. Read and calculate results using Neogen's Stat Fax microwell reader, or equivalent. If using an EL301 reader or other strip/plate reader, calculate results using Neogen's Veratox for Windows software.

TEST PROCEDURE – FOR SCREENING

Allow reagents to warm to room temperature, 18-30°C (64-86°F), prior to use.

1. Remove 1 red-marked mixing well for each sample to be tested plus 1 red-marked well for the control, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant and reseal the pack to protect the antibody. Mark one end of strip with a "1", and place strip in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Place 100 µL of conjugate from the blue-labeled bottle in each red-marked mixing well.
5. **Choose one control** level to screen samples against for each test performed. Using a new pipette tip for each, transfer 100 µL of the diluted samples and chosen control to the red-marked mixing wells as shown below. Thoroughly mix the sample and control with the conjugate as you add them to the well by depressing the plunger 5 times. Do not use more than 6 wells at one time.

Control	S1	S2	S3	S4	S5
---------	----	----	----	----	----

6. Using a new tip for each, transfer 100 µL from each mixing well to the corresponding antibody-coated well. Mix by sliding the microwell holder back and forth on a flat surface for 10 seconds without splashing the reagents from the wells. Incubate for **5 minutes** at room temperature, 18-30°C (64-86°F).
7. The initial reaction is now complete. Shake out the contents of the antibody wells.
8. Fill each antibody well with distilled or deionized water and dump them out. Repeat this step 5 times, then turn the wells upside down and tap out on a paper towel until the remaining water has been removed.
9. Pipette 100 µL of substrate into each well and mix by sliding back and forth on a flat surface for 10 seconds. Incubate **5 minutes**.
10. Pipette 100 µL Red Stop to each well and mix by sliding back and forth on a flat surface. Discard the tip.
11. Wipe the bottom of the microwells with a dry cloth or towel.
Microwells may be read visually or using a 650 nm filter. If the sample well is blue or more blue than the control well, the sample contains less toxin than the control. If the sample well shows less blue color (more red color) than the control, the sample contains more toxin than the control. For optimum observation of color differences, place the wells on a white surface and read looking down through the solution.

RETESTING

If positives occur in commodities not previously tested, confirm with an additional approved method prior to taking action.

PERFORMANCE CHARACTERISTICS

Limit of detection: 10 ppb (Determined by the mean average of 10 zearalenone free samples plus 2 standard deviations.)

Limit of quantitation: 25 ppb (Described as the lowest concentration point on the calibration curve that this test can reliably detect zearalenone.)

Range of quantitation: 25-500 ppb (For quantitating samples above 500 ppb contact Neogen Technical Services for dilution instructions.)

Validated matrixes: Barley, corn, corn gluten meal*, corn silage, DDGs*, oats, oat flour, oat hulls*, oats (naked), pea fiber, pet food*, popcorn, potato, rice, rice (brown), rice flour (white), rice hulls, rye, soybean meal, tapioca, wheat, wheat bran*.

*Commonly acidic and may require a pH adjustment.

NOTE: Neogen continues to validate new commodities. Please contact a representative for the latest validated commodity list.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's Food Safety test kits, on Neogen's Web site at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TESTING KITS AVAILABLE FROM NEOGEN

Natural Toxins

- Aflatoxin, DON, Ochratoxin, Zearalenone, T-2/HT-2 Toxin, Fumonisin, Histamine

Foodborne Bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, Yeast and Mold, Total Plate Count, Generic *E. coli* and Total Coliforms, Protein Residues

Food Allergens

- Peanuts, Milk, Eggs, Almonds, Gluten, Soy, Hazelnut

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©Neogen Corporation, 2015. Neogen, Agri-Screen, Veratox and K-Blue are registered trademarks of Neogen Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.

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