Read instructions carefully before starting test

Agri-Screen®

Fumonisin Screening Test

REFRIGERATE AT 2-8°C • DO NOT FREEZE

THE TOXIN

Discovered in 1989, fumonisins are a family of mycotoxins produced by different species of the mold *Fusarium*. These molds commonly infect corn (in fact, they are considered ubiquitous in corn) and rice, hence the potential for fumonisins to be found in feed and foodstuffs is high. Fumonisin affect various animals differently and have been linked to esophageal cancer in humans. The Environmental Protection Agency classifies fumonisins as Category II-B carcinogens.

Horses are extremely sensitive to low amounts of fumonisin, which can cause leukoencephalomalacia (liquefaction of the brain). In swine, research has shown fumonisin attacks the cardiopulmonary system causing pulmonary edema, as well as liver and pancreatic lesions.

The FDA has issued a final guidance for total fumonisins ($FB_1+FB_2+FB_3$) in food and animal feeds:

HUMAN FOODS (total fumonisins)

Degermed dry milled corn products.....	2 ppm
Whole/partially dry milled corn products, dry milled corn bran, cleaned corn for masa production	4 ppm
Cleaned corn for popcorn	3 ppm

ANIMAL FEEDS (corn/corn by-products, total fumonisins)

Equids and rabbits.....	5 ppm, < 20% of diet
Swine and catfish.....	20 ppm, < 50% of diet
Breeding ruminants, breeding poultry and breeding mink.....	30 ppm, < 50% of diet
Ruminants > 3 months old being raised for slaughter and mink being raised for pelt production	60 ppm, < 50% of diet
Poultry being raised for slaughter.....	100 ppm, < 50% of diet
All other species or classes of livestock and pet animals.....	10 ppm, < 50% of diet

INTENDED USE

Agri-Screen for Fumonisin is intended for the screening for fumonisin in such commodities as corn, cornmeal and rice.

INTENDED USER

The test kit is designed for use by quality control personnel and others familiar with food and feed possibly contaminated by fumonisin. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed the Neogen training.

ASSAY PRINCIPLES

Agri-Screen for Fumonisin is a competitive direct enzyme-linked immunosorbent assay (CD-ELISA) which allows the user to determine how the concentration of fumonisin in a sample compares to the concentration of fumonisin in a supplied control of 5 parts per million (ppm). A methanol/water solution is used to extract any existing fumonisin from a ground sample. After a dilution, extracted fumonisin in the filtrate is mixed with enzyme-labeled fumonisin (conjugate). The mixed solution is transferred to antibody-coated wells, where free fumonisin and conjugate compete for antibody binding sites. After a wash step, substrate is added. Color develops as a result of the presence of bound conjugate. Red stopping reagent is added and the color of the resulting solution is observed. Blue color indicates negative samples. Red indicates strong positives.

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2-8°C (35-46°F).

MATERIALS PROVIDED

1. 24 antibody-coated microwells
2. 24 red-marked mixing wells
3. 1 yellow-labeled bottle of 5 ppm fumonisin control (see precautions for handling of methanol solution)
4. 1 blue-labeled bottle of fumonisin-HRP conjugate solution
5. 1 green-labeled bottle of K-Blue® Substrate solution
6. 1 red-labeled bottle of Red Stop solution
7. 14 pink-labeled bottles for sample dilution, prefilled with 7.9 mL of deionized water (*Dilution bottles not provided with Neogen item #8811*)
8. Directions for use

MATERIALS REQUIRED BUT NOT PROVIDED

1. Extraction materials (items d through f available in kit form from Neogen, item #8052):
 - a. 70% methanol solution (Neogen item #8055)
 - b. Distilled or deionized water
 - c. 250 mL graduated cylinder (Neogen item #9368)
 - d. Container with 125 mL capacity (Neogen item #9428)
 - e. Neogen filter syringes, Whatman #1 filter paper, or equivalent (Neogen item #9420/#9430)
 - f. Sample collection tubes (Neogen item #9421)
2. 50 mL graduated cylinder (Neogen item #9447)
3. High-speed blender (Neogen item #9493/#9477)
4. Agri-Grind grinder or equivalent (Neogen item #9401/#9453)
5. Scale capable of weighing 5-25 grams (Neogen item #9427)
6. Pipettor 100 µL (Neogen item #9272/#9278)
7. Tips for 100 µL pipettor (Neogen item #9410/#9407)
8. Paper towels or equivalent absorbent material

9. Plastic bucket for use as waste receptacle
10. Microwell holder (Neogen item #9402)
11. Timer (Neogen item #9426)
12. Waterproof marker
13. Wash bottle (Neogen item #9400)
14. Distilled or deionized water

PRECAUTIONS

1. Methanol solution is highly flammable. Keep container tightly closed, and keep away from heat, sparks, open flame and those smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
2. Store test kit between 2-8°C (35-46°F) when not in use. Do not freeze test kits.
3. Kits should be at room temperature (18-30°C, 64-86°F) prior to use.
4. Avoid prolonged storage of kits at ambient temperatures.
5. Do not use kit components beyond expiration date.
6. Do not mix reagents from one kit serial with reagents from a different kit serial.
7. Do not run more than 6 wells per test, unless using a multichannel pipettor.
8. Follow proper pipetting techniques, including priming pipette tips by filling and dispensing solution once before use.
9. Use of incubation times other than those specified may give inaccurate results.
10. Treat all used liquids, including sample extract, and labware as if contaminated with fumonisin. Gloves and other protective apparel should be worn at all times.
11. To avoid cross-contamination, use clean pipette tips and glassware for each sample, and thoroughly wash all glassware between samples.
12. Commodities tested should have a pH of 6-8. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH contact your Neogen representative or Technical Services.
13. Do not use substrate that has turned blue prior to use.

CHANGING THIS TEST KIT'S SCREENING LEVEL

This test kit can be used to screen samples for fumonisin at levels other than the 5 ppm supplied control. Contact Neogen for more information on screening samples at a level other than 5 ppm.

SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques. Store samples at 2-8°C (35-46°F) until analyzed.

1. If not using Neogen's prepared solution, prepare a 70% methanol solution by mixing 7 parts ACS Grade methanol with 3 parts distilled or deionized water for each sample to be tested.
2. Obtain a representative sample. Grind the sample so at least 75% of the ground material passes through a 20 mesh sieve, about the particle size of fine instant coffee.
3. Extract at a ratio of 1 part sample to 5 parts 70% methanol. Blend 25 grams of ground sample with 125 mL of 70% methanol for 2 minutes in a high-speed blender. **Alternative method:** Add 10 grams of ground sample to 50 mL of 70% methanol and shake vigorously for 3 minutes.
4. Filter at least 5 mL of the extract by using a Neogen filter syringe or pouring through a Whatman #1 filter. Collect the filtrate as the sample.
5. Dilute the sample by adding 100 µL of extract to a prefilled sample dilution bottle and mix well by swirling bottle. The sample is now ready for testing. Repeat for each sample, making sure to label each bottle.

TEST PROCEDURE

Allow all reagents to warm to room temperature (18-30°C, 64-86°F) prior to use.

1. Remove 1 red-marked mixing well from the foil pack for each sample to be tested, and one for the control, and place in the well holder. **Note:** Do not run more than 6 wells at a time unless you are using a multichannel pipettor. Contact Neogen for more information.
2. Remove an equal number of antibody-coated wells. Return wells which will not be used immediately to the foil pack and reseal to protect the antibody. Mark one end of the strip with a "1", and place in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Using a new pipette tip, add 100 µL of conjugate from the blue-labeled bottle to each red-marked mixing well. Discard the tip.
5. Using a new tip, add 100 µL of the control from the yellow-labeled bottle to the first well of the red-marked strip. Thoroughly mix by inserting the tip in the liquid and pipetting up and down 5 times. Discard the tip.
6. Using a new tip, add 100 µL of the first sample to the second red-marked well. Thoroughly mix by inserting the tip in the liquid and pipetting up and down 5 times. Discard the tip. Repeat the process for each additional sample in a following red-marked well.
7. Using a new tip for each, transfer 100 µL from each red-marked well to the corresponding antibody-coated well. Discard the red-marked wells.
8. Set timer for **5 minutes, mixing the wells for the first minute** of the room temperature incubation by sliding the microwell holder back and forth on a flat surface without splashing reagents from the wells.
9. The initial reaction is now completed. Shake out the contents of the wells.
10. Fill each well with distilled or deionized water and shake out. Repeat 5 times. Remove all water droplets by turning wells upside down and vigorously tapping on an absorbent paper.
11. Using a new tip, add 100 µL of substrate from the green-labeled bottle to each well. Discard the tip.
12. Set timer for **10 minutes, mixing the wells for the first minute** by sliding the microwell holder back and forth on a flat surface.
13. Using a new tip, add 100 µL of Red Stop from the red-labeled bottle to each well. Discard the tip. Mix as before. Visually check to confirm thorough mixing.

INTERPRETATION OF RESULTS

If a sample well is **as blue** or **darker blue** than the control well, the sample contains **less than 5 ppm** of fumonisin. If a sample well shows **less blue** color, or **more red** color, than the control, the sample contains **more than 5 ppm** of fumonisin. For optimum observation of color differences, place the wells on a white surface and read looking down through the solution.

Alternative: Microwells may be read in a microwell reader. Wipe the bottom of wells, blank reader on air using a 650 nm filter, and compare sample readings to the 5 ppm control reading.

RETESTING

If positives occur in commodities not previously tested, confirm with an additional approved method prior to taking action.

VALIDATED MATRICES

Barley, beet pulp*, corn, corn meal, corn germ*, corn gluten meal*, corn/soy blend, DDGS*, milo, pea fiber, petfood*, popcorn, rice, rice hulls, rough rice, soybeans, soybean meal, soyflower meal, wheat and wheat bran.

**May require a pH adjustment.*

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Service can be reached between 8 a.m. and 7 p.m. Eastern time by calling 800/234-5333 or 517/372-9200 and asking for a Neogen sales representative or Technical Services. Assistance is available on a 24-hour basis by calling 800/234-5333. Training on this product, and all Neogen test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's test kits, on Neogen's Web site at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TESTING KITS AVAILABLE FROM NEOGEN

Natural Toxins

- Aflatoxin, DON, Ochratoxin, Zearalenone, T-2/HT-2 Toxins, Fumonisin, Histamine

Foodborne Bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, Yeast and Mold, Total Plate Count, Generic *E. coli* and Total Coliforms, Protein Residues

Food Allergens

- Almonds, Crustacea, Eggs, Gliadin, Hazelnut, Lupine, Milk, Mustard, Peanuts, Sesame, Soy, Walnut

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©Neogen Corporation, 2011. Neogen, Veratox y K-Blue son marcas comerciales registradas de Neogen Corporation. Todas las otras marcas y nombres de productos son marcas comerciales o marcas comerciales registradas de sus respectivas compañías.