Read instructions carefully before starting test

 **Veratox**[®]
for T-2/HT-2
Quantitative Test

Refrigerate at 2–8°C • Do not freeze

THE TOXINS

T-2/HT-2 toxins are trichothecene mycotoxins produced by several species of *Fusarium* molds. As T-2 toxin is readily metabolized to HT-2 toxin, and the toxins have been shown to produce numerous adverse effects on many animals, these two mycotoxins are frequently evaluated together.

Animals affected by the toxins include swine, dairy cattle, poultry, dogs, cats and horses. Effects of the toxins include digestive disorders, hemorrhage, edema, oral lesions, dermatitis, and blood disorders. Damage caused by the toxins to the digestive track is irreversible. In the most severe cases, these toxins will cause death. T-2 toxin is the principal causal toxin in the human disease alimentary toxic aleukia.

Poultry studies have shown T-2 intoxication has led to a reduction in weight gain and other problems such as beak lesions, poor feathering, motor function impairment and increased susceptibility to *Salmonella* spp.

The best protection against these mycotoxins is monitoring for their presence in feeds and foods. That means testing all along the pathway from initial harvest of grains to the finished product.

INTENDED USE

The Veratox[®] for T-2/HT-2 test kit is a competitive direct enzyme-linked immunosorbent assay (CD-ELISA) for the quantitative analysis of T-2/HT-2 toxins in such commodities as corn, barley, oats, rice, rye, soy and wheat.

INTENDED USER

The test kit is designed for use by quality control personnel and others familiar with food and feed possibly contaminated by T-2 and/or HT-2 toxins. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed the Neogen training.

ASSAY PRINCIPLES

Veratox for T-2/HT-2 is a competitive direct ELISA in a microwell format that allows the user to obtain exact concentrations in parts per billion (ppb) of either T-2, HT-2 or a combination of both. Free T-2 or HT-2 toxins in the samples and controls are allowed to compete with enzyme-labeled HT-2 toxin (conjugate) for the antibody binding sites. After a wash step substrate is added, which reacts with the bound conjugate to produce blue color. More blue color means less T-2/HT-2 toxins. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form the standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of T-2/HT-2 toxins.

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F).

MATERIALS PROVIDED

1. 48 antibody-coated microwells
2. 48 red-marked mixing wells
3. 5 yellow-labeled bottles of 0, 25, 50, 100 and 250 ppb T-2 toxin controls (see precautions for handling of methanol solution)
4. 1 blue-labeled bottle of HT-2 toxin-HRP conjugate solution
5. 1 green-labeled bottle of K-Blue[®] Substrate solution
6. 1 red-labeled bottle of Red Stop solution
7. Directions for use

MATERIALS REQUIRED BUT NOT PROVIDED

1. Extraction materials (items c through e are available in kit form, Neogen item 8052):
 - a. 70% methanol (ACS Grade) (Neogen item 8055, 8056)
 - b. 250 mL graduated cylinder (Neogen item 9368)
 - c. Container with 125 mL capacity (Neogen item 9428)
 - d. Neogen filter syringes, Whatman no. 1 filter paper, or equivalent (Neogen items 9420, 9430)
 - e. Sample collection tubes (Neogen item 9421)
2. Agri-Grind grinder or equivalent (Neogen item 9401, 9453)
3. Scale capable of weighing 2–25 grams (Neogen item 9427)
4. Microwell reader with a 650 nm filter (Neogen item 9303)
5. Pipettor, 12-channel (Neogen item 9273)
6. Pipettor, 100 μ L (Neogen item 9272, 9290)
7. Pipette tips for 100 μ L and 12-channel pipettors (Neogen items 9410, 9407, 9417)
8. Paper towels or equivalent absorbent material
9. Plastic bucket for use as waste receptacle
10. Microwell holder (Neogen item 9402)
11. Timer (Neogen item 9426)
12. Waterproof marker
13. Wash bottle (Neogen item 9400)
14. 2 reagent boats for 12-channel pipettor (Neogen item 9450)
15. Distilled or deionized water

PRECAUTIONS

1. Methanol solution is highly flammable. Keep container tightly closed, and keep away from heat, sparks, open flame and those smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
2. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze.
3. Do not use kit components beyond expiration date.
4. Do not mix reagents from a kit with one serial number with reagents from a kit with a different serial number.
5. Do not run more than 24 wells at one time.
6. Follow proper pipetting techniques, including the proper priming of tips.
7. Incubation times other than those specified may give inaccurate results.
8. Kits should be at room temperature 18–30°C (64–86°F) prior to use.
9. Avoid prolonged storage of kits at ambient temperatures.
10. Treat all used liquids, including sample extract, and labware as if contaminated with T-2 or HT-2 toxins. Gloves and other protective apparel should be worn at all times.
11. To avoid cross-contamination, use clean pipette tips and glassware for each sample, and thoroughly detoxify and wash all glassware between samples.
12. Commodity extracts should have a pH of 6–8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH contact a Neogen representative or Technical Services.

PROCEDURAL NOTES

1. **Substrate:** K-Blue Substrate is ready for use. The substrate should be clear to a very light blue—discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until it is needed.
2. **Antibody wells:** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted and the test procedure is set to begin.

SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques. The sample should be ground and thoroughly mixed prior to proceeding with the extraction. Store samples at 2–8°C (35–46°F) until analyzed. **NOTE:** If using Neogen's Mycotoxin Extraction Kit, follow the instructions in that kit for the extraction procedure. If preparing your own extraction solution, continue with the instructions that follow.

1. If not using Neogen's prepared solution, prepare a 70% methanol solution by mixing 7 parts ACS grade methanol with 3 parts deionized or distilled water for each sample to be tested.
2. Obtain a representative sample. Grind the entire sample so at least 75% of the ground material passes through a 20 mesh sieve, the particle size of a fine instant coffee.
3. Extract at a ratio of 1 part sample with 5 parts 70% methanol. **EXAMPLE:** Add 5 g of ground sample to 25 mL of 70% methanol/water and shake vigorously for **3 minutes**, or blend for **2 minutes**.
4. Filter the extract by pouring at least 5 mL through a Whatman no. 1 filter (or Neogen filter syringe) and collecting the filtrate as a sample.
5. Dilute the filtrate 1:1 (e.g., 1 mL in 1 mL) with deionized or distilled water and mix.
6. The sample is now ready for testing without further preparation.

TEST PROCEDURE (FOR QUANTITATION)

Allow reagents to warm to room temperature 18–30°C (64–86°F) prior to use.

1. Remove 1 red-marked mixing well for each sample to be tested plus 5 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells that will not be used immediately to the foil pack with desiccant and reseal the pack to protect the antibody. Mark one end of strip with a "1", and place strip in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Place 100 µL of conjugate from the blue-labeled bottle in each red-marked mixing well.
5. Using a new pipette tip for each, transfer 100 µL of controls and samples to the red-marked mixing wells as shown below:

0	25	50	100	250	S1	S2	S3	S4	S5	S6	S7	Strip 1
S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19	Strip 2

6. Using a 12-channel pipettor, mix the liquid in the wells by pipetting it up and down 3 times. Transfer 100 µL to the antibody-coated wells. Incubate **5 minutes** at room temperature 18–30°C (64–86°F), mixing for the first **10–20 seconds** by sliding back and forth on a flat surface. Discard the red-marked mixing wells.
7. The initial reaction is now complete. Shake out the contents of the antibody wells.
8. Fill each antibody well with deionized or distilled water and dump them out. Repeat this step 5 times, then turn the wells upside down and tap out on a paper towel until the remaining water has been removed.
9. Pipette the needed volume of substrate from the green-labeled bottle into the green-labeled reagent boat and, with new tips, pipette 100 µL of substrate into the wells. Incubate **5 minutes** at room temperature 18–30°C (64–86°F), mixing for the first **10–20 seconds** by sliding back and forth on a flat surface. Discard remaining substrate and rinse the reagent boat with water.
10. Pipette Red Stop solution from the red-labeled bottle (same volume as prepared for substrate) into the red-labeled reagent boat. Using the same pipette tips as were used to dispense substrate, add 100 µL Red Stop to each well and mix by sliding back and forth on a flat surface. Discard the tips.
11. Wipe the bottom of the microwells with a dry cloth or towel and read in a microwell reader using a 650 nm filter. Air bubbles should be eliminated, as they could affect analytical results. Results should be read within **20 minutes** after the addition of Red Stop.
12. Read and calculate results using Neogen's Stat Fax microwell reader, or equivalent. If using an EL301 reader or other strip/plate reader, calculate results using Neogen's Veratox software for Windows.

TEST PROCEDURE (FOR SCREENING)

Allow reagents to warm to room temperature (18–30°C, 64–86°F), prior to use.

1. Remove 1 red-marked mixing well for each sample to be tested plus 1 red-marked well for the control, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells that will not be used immediately to the foil pack with desiccant and reseal the pack to protect the antibody. Mark one end of strip with a “1”, and place strip in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Place 100 µL of conjugate from the blue-labeled bottle in each red-marked mixing well.
5. Choose one control level to screen samples against for each test performed. Using a new pipette tip for each, transfer 100 µL of the samples and chosen control to the red-marked mixing wells as shown below. Thoroughly mix the sample and control with the conjugate as you add them to the well by depressing the plunger 5 times. Do not use more than 6 wells at one time.

Control	S1	S2	S3	S4	S5
---------	----	----	----	----	----

6. Using a new tip for each, transfer 100 µL from each mixing well to the corresponding antibody-coated well. Mix by sliding the microwell holder back and forth on a flat surface for **10–20 seconds** without splashing the reagents from the wells. Incubate **5 minutes** at room temperature 18–30°C (64–86°F), mixing for the first **10–20 seconds** by sliding back and forth on a flat surface.
7. The initial reaction is now complete. Shake out the contents of the antibody wells.
8. Fill each antibody well with distilled or deionized water and dump them out. Repeat this step 5 times, then turn the wells upside down and tap out on a paper towel until the remaining water has been removed.
9. Pipette 100 µL of substrate into each well. Incubate **5 minutes** at room temperature 18–30°C (64–86°F), mixing for the first **10–20 seconds** by sliding back and forth on a flat surface.
10. Pipette 100 µL Red Stop to each well and mix by sliding back and forth on a flat surface. Discard tip.
11. Wipe the bottom of the microwells with a dry cloth or towel.
12. Microwells may be read visually or using a 650 nm filter. If the sample well is blue or more blue than the control well, the sample contains less toxin than the control. If the sample well shows less blue color (more red color) than the control, the sample contains more toxin than the control. For optimum observation of color differences, place the wells on a white surface and read looking down through the solution.

RETESTING

If positives occur in commodities not previously tested, confirm with an additional approved method prior to taking action.

PERFORMANCE CHARACTERISTICS

Limit of detection: < 10 ppb of T-2 or HT-2 or a combination of both (Determined by the mean average of 10 T-2/HT-2 toxin-free samples plus 2 standard deviations.)

Limit of quantitation: 25 ppb (Described as the lowest concentration point on the calibration curve that this test can reliably detect T-2/HT-2 toxins.)

Range of quantitation: 25–250 ppb (For quantitating samples above 250 ppb contact Neogen Technical Services for dilution instructions.)

Validated matrix: Barley, corn, corn flour, corn gluten*, corn steep, DDGs wet cake, oats, oat hulls (whole)*, rice (brown), rice flour (white), rice gluten, rice hulls, rye, pea fiber, potato (white), soy, soybean meal, tapioca, wheat, wheat bran*, wheat flour, wheat gluten

NOTE: Neogen continues to validate new commodities. Please contact a representative for the latest validated commodity list.

Cross-reactivity: 100% T-2 toxin and 100% HT-2 toxin. No cross-reactivity with any other trichothecene mycotoxins.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts

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