Read instructions carefully before starting test

Reveal[®] Q+

for Aflatoxin

ENDPOINT

for use with Raptor[®] Solo Readers

THE TOXIN

Aflatoxin is a toxic and carcinogenic substance produced by certain strains of the molds *Aspergillus flavus* and *A. parasiticus*. There are four principle types of aflatoxin: B₁, B₂, G₁ and G₂. Aflatoxin B₁ is the most frequently encountered of the group and the most toxic. The commodities most affected by aflatoxin are corn, peanuts, cottonseed, milo, and the majority of tree nuts.

The effects on animals ingesting excessive amounts of the toxin range from chronic health and performance problems to death. Aflatoxin has been shown to cause liver damage or cancer, decreased milk and egg production, immune suppression, and interference with reproductive efficiency.

Many countries have set maximum allowable levels of aflatoxin in food and feed. Therefore, accurate determination of the toxin's presence is of major importance to those monitoring the quality of food and feed in which aflatoxin may occur. Testing these commodities for the toxin requires careful sampling, chemical extraction, sanitation, and quantitative analysis.

The FDA has issued regulatory levels for aflatoxin as follows:

For	Level	Commodities
Humans	20 ppb	All food except milk
All animal species	20 ppb	All feed (exceptions below)
Exceptions:		
Breeding cattle, breeding swine, mature poultry	100 ppb	Corn
Finishing swine (> 100 lbs.)	200 ppb	Corn
Finishing beef cattle	300 ppb	Corn
Finishing beef cattle, swine, poultry	300 ppb	Cottonseed meal

The European Union has issued regulations for aflatoxin (total B₁, B₂, G₁ and G₂) as follows: (as referenced in the Commission Regulation (EC) No. 1881/2006)

Foodstuffs	Levels
Groundnuts (peanuts) and other oilseeds, subject to sorting and treatment before human consumption or use as an ingredient in foodstuffs Almonds, pistachios, hazelnuts, Brazil nuts and apricot kernels subject to sorting and treatment before human consumption or use as an ingredient in foodstuffs	15 ppb
Tree nuts, other than those listed above, subject to sorting and treatment before human consumption or use as an ingredient in foodstuffs; almonds, pistachios, hazelnuts, Brazil nuts and apricot kernels for direct human consumption or use as an ingredient in foodstuffs Maize, rice and dried fruits subject to sorting and treatment before human consumption or as an ingredient in foodstuffs	10 ppb
Groundnuts (peanuts) and other oilseeds for direct human consumption or use as an ingredient in foodstuffs Dried fruits and processed products thereof, for direct human consumption or use as an ingredient in foodstuffs	4 ppb

INTENDED USE/USER

Reveal Q+ for Aflatoxin ENDPOINT is intended for the quantitative analysis of aflatoxin. The test kit is designed for use by quality control personnel and others familiar with commodities possibly contaminated with aflatoxin.

ASSAY PRINCIPLES

Reveal Q+ for Aflatoxin is a single-step lateral flow immunoassay. This kit can only be used with the Raptor® Solo Integrated Analysis Platform Reader. In the test, the extract is wicked through the reagent zone, which contains antibodies specific for aflatoxin conjugated to colloidal gold particles (gold complex). If aflatoxin is present, it will be captured by the gold complex. The gold complex, along with any free gold-complex, is then wicked onto a membrane, which contains a zone of aflatoxin conjugated to a protein carrier. This zone captures any unbound gold complex, allowing the particles to concentrate and form a visible line. As the level of aflatoxin in a sample increases, free aflatoxin will bind with the gold complex allowing the less gold complex to be captured in the test zones. Therefore, as the concentration of aflatoxin in the sample increases, the test line density decreases. Algorithms programmed into the readers convert the line intensities into a quantitative result displayed in parts per billion (ppb). The membrane also contains

a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of aflatoxin, ensuring the strip is functioning properly.

STORAGE REQUIREMENTS

Store kit components at room temperature (18–30°C, 64–86°F) to ensure full shelf life. Test strips should remain capped in their original tubes until used to ensure optimal performance.

MATERIALS PROVIDED

1. 25 Reveal Q+ for Aflatoxin test strips
2. 25 red sample dilution cups
3. 25 clear sample cups
4. 1 bottle of sample diluent
5. Instructions for use

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Extraction Materials
 - a. 65% ethanol solution (NEOGEN item 8073, 8074)
 - b. Sample collection cups with lids 125 mL (NEOGEN item 9428, 9428B)
 - c. Sample collection tubes with caps (NEOGEN item 9421, 9421B)
 - d. Filter syringes (NEOGEN item 9272, 9860)

OR

 - e. Centrifuge, mini (NEOGEN item 9330)
 - f. Micro-centrifuge tubes (NEOGEN item 9172)
 - g. Dispensing pump or graduated cylinder (NEOGEN item 9448, 9447)
2. Agri-Grind grinder or equivalent (NEOGEN item 9401, 9453)
3. Scale capable of weighing 5–50 g ± 0.1 g (NEOGEN item 9427)
4. Timer (NEOGEN item 9426, 9452)
5. Reveal sample cup rack (NEOGEN item 9475)
6. Raptor Solo Integrated Analysis Platform (NEOGEN item 9696)
7. Pipettor, 100 µL (NEOGEN item 9272, 9860)
8. Pipette tips, 100 µL (NEOGEN item 9407, 9410, 9417)
9. Pipettor, 500 µL (NEOGEN item 9336, 9291)
10. Pipette tips, 200–1000 µL (NEOGEN item 9464, 9487, 9293, 9292)


PRECAUTIONS

1. The test strips must remain inside the stay-dry tube before use.
2. Ethanol is highly flammable. Keep container tightly closed and away from heat, sparks, open flame and those who are smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
3. Store test kit at room temperature (18–30°C, 64–86°F) when not in use. Do not freeze.
4. Do not use kit components beyond expiration date.
5. Treat all used liquids, including sample extract, and labware as if contaminated with aflatoxin. Gloves and other protective apparel should be worn at all times.
6. To avoid cross-contamination, use clean glassware for each sample, and thoroughly wash all glassware between samples.
7. Ensure the device lot number and the curve details match the lot ID number selected on the reader.
8. Commodity extracts should have a pH of 6-8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH, contact a NEOGEN representative or Technical Services.

SOLO READER SETUP

The Solo reader must be in Standard mode to read the Reveal Q+ for Aflatoxin ENDPOINT test strips.

1. From the Main Menu, select the  in the upper left hand corner
2. Select Track Replacement

Testing track removal tool

3. Return to the Home screen by pressing the  in the lower left hand corner
4. Insert the track replacement tool with the arrow facing towards you into the track in the unit
5. Gently pull up to remove the track
6. Insert the Standard ENDPOINT testing track into the unit (track with the blue dot on the right hand side)

Testing track – ENDPOINT

SAMPLE PREPARATION

The sample to be tested should be collected according to accepted sampling techniques (see FGIS sampling protocol or contact your NEOGEN representative). Obtain a representative sample (minimum 100 g). Grind the sample so at least 95% of the ground material passes through a 20 mesh sieve (about the particle size of fine espresso).

SAMPLE EXTRACTION

1. If not using NEOGEN's prepared solution, prepare a 65% ethanol solution by mixing 6.5 parts ethanol with 3.5 parts distilled or deionized water for each sample.
2. Weigh out 10 g \pm 0.1 g of sample into an extraction cup.
3. Add 50 mL of 65% ethanol to the extraction cup.
4. Vigorously shake, using hand or mechanical means, for 3 minutes, or blend for 1 minute.
5. Allow the sample to settle, then filter with a filter syringe or Whatman #4 filter paper to collect a minimum of 3 mL filtrate into a sample collection tube. Or, you may also pipette 1 mL of sample into a 2.0 mL micro-centrifuge tube, and centrifuge for 1 minute using a micro-centrifuge (approx. 2,000 x g).

TEST PROCEDURE

1. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary.
2. Add 100 µL of sample extract to each **red** sample dilution cup.
3. Add 500 µL of sample diluent to the red sample dilution cup containing the sample extract. Mix by pipetting up and down 5 times.
4. Transfer 100 µL of diluted sample extract into a new **clear** sample cup.
5. Place a new Reveal Q+ for Aflatoxin test strip with the sample end down into the sample cup and set timer for **6 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
6. Remove the strip from the sample cup after it has developed for exactly **6 minutes**. Immediately (within 30 seconds) read the test strip in the Raptor Solo reader.

READING A TEST STRIP

1. From the home screen:
 - Select Category: Mycotoxin Q+
 - Select Test: Q+ Afla Endpoint
2. Press Run Test
 - Select Lot
 - (If not listed, add the lot by scanning the QR code located on the device Tube)
3. Enter Sample ID as needed – press accept
4. Fully insert the Reveal Q+ Aflatoxin ENDPOINT test strip into the NEOGEN labeled cartridge adaptor with the sample end first
5. Insert cartridge into the track on the Solo reader – system will read automatically
6. Results are displayed on the unit.
 - Samples greater than 100 ppb must be diluted and re-tested.
7. To read the next test strip press Change at the bottom of the screen – this allows the user to enter the next Sample ID. Press Accept and insert cartridge containing the next test strip

DILUTION PROCEDURE

Samples greater than 100 ppb will need to be diluted and retested.

Example dilution:

1. Add 100 µL sample filtrate to a sample collection tube. Add 100 µL 65% ethanol. Mix well.
2. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary.
3. Add 100 µL of diluted sample extract from step 1 to each **red** sample dilution cup.
4. Add 500 µL of sample diluent to the red sample dilution cup containing the sample extract. Mix by pipetting up and down 5 times.
5. Transfer 100 µL of sample into a new clear sample cup.
6. Place a new Reveal Q+ for Aflatoxin test strip with the sample end down into the sample cup and set timer for 6 minutes. Ensure the test strip comes into contact with liquid and begins to wick.
7. Remove the strip from the sample cup after it has developed for exactly 6 minutes. Immediately (within 30 seconds) read the test strip in the Raptor Solo reader.
8. Final result displayed will need to be multiplied by 2.

READING TEST RESULTS

NOTE: The Reveal Q+ for Aflatoxin ENDPOINT can only be read using the Raptor Solo reader. Test strips should be read within 30 seconds of the completion of the 6-minute incubation.

NOTES:

1. Ensure device is fully inserted into cartridge.
2. Readings should be made between **6 and 6.5 minutes**. Readings after 6.5 minutes may be inaccurate due to overdevelopment of the device.

PERFORMANCE CHARACTERISTICS

1. Limit of detection: 2 ppb
2. Range of detection: 3–100 ppb

NOTE: Samples greater than 100 ppb must be diluted and re-tested. Results below the range of detection should be reported as less than 3 ppb.

VALIDATED/VERIFIED MATRICES

NEOGEN continues to validate new commodities. Please contact a representative for the latest validated commodity list.

CUSTOMER SERVICE

NEOGEN Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all NEOGEN test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of NEOGEN's Food Safety test kits, on NEOGEN's web site at NEOGEN.com, or by calling NEOGEN at 800.234.5333 or 517.372.9200.

TERMS AND CONDITIONS

For NEOGEN's full terms and conditions, please visit neogen.com/terms-and-conditions/

WARRANTY

NEOGEN Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, ergot alkaloids, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, coconut, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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