*Read instructions carefully before starting test*

# Veratox<sup>®</sup> **for Sesame Allergen**

## **Quantitative Test**

*REFRIGERATE AT 2–8°C (35–46°F) • DO NOT FREEZE*

### **SESAME ALLERGEN**

Food allergens are proteins in food that can create an immune response in sensitive individuals. Once ingested, food allergens can cause a number of reactions, ranging in severity from hives and itching to anaphylaxis. Anaphylaxis is a severe allergic reaction, involving vomiting, diarrhea, difficulty breathing, swelling of the mouth and tongue, and a rapid drop in blood pressure.

An estimated 3.5 to 4 percent of adults, and 6 to 8 percent of children, are sensitive in some degree to food allergens. More than 12 million people in the United States alone are known to have a food allergy.

Food manufacturers protect those with food allergies by clearly labeling their products with a list of ingredients. Testing for the presence of sesame components ensures food manufacturers that an unlabeled—and potentially dangerous—ingredient did not make its way into a food product.

### **INTENDED USE**

Veratox for Sesame Allergen is intended for the quantitative analysis of sesame protein residue in food products such as bakery products, spices, cereal products and clean-in-place rinses.

### **INTENDED USER**

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by sesame or sesame products. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed the Neogen training.

## **ASSAY PRINCIPLES**

Veratox for Sesame Allergen is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Sesame residues are extracted from samples with a buffered salt solution (PBS) by shaking in a heated water bath, followed by centrifugation or filtration. Extracted sesame residue is sampled and added to antibody-coated wells (capture antibody) where it binds to the antibody during an incubation. Any unbound sesame residue is washed away and a second antibody (detector antibody), which is enzyme labeled, is added. The detector antibody binds to the already bound sesame residue. After a second wash, the substrate is added. Color develops as a result of the presence of bound detector antibody. Red Stop reagent is added and the color of the resulting solution is observed. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form a standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of sesame residue.

## **STORAGE REQUIREMENTS**

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). Do not freeze.

## **MATERIALS PROVIDED**

1. 48 antibody-coated microwells
2. 48 red-marked transfer wells
3. 5 yellow-labeled bottles of 0, 2.5, 5, 10, 25 ppm sesame controls
4. 2 blue-labeled bottles of enzyme-labeled antibody conjugate
5. 1 green-labeled bottle of K-Blue® Substrate
6. 1 red-labeled bottle of Red Stop solution
7. Foil pouch of 10 mM PBS dry powder extraction solvent. Each pouch is enough to prepare 1 L in distilled or deionized water (pH 7.4)
8. 40 mL of 10 mM PBS-Tween washing reagent in a wide mouth bottle. Each bottle is enough to prepare 1 L in distilled or deionized water (pH 7.4)
9. 50 g of extraction additive in a specimen cup
10. Plastic scoop to measure extraction additive

## **MATERIALS RECOMMENDED BUT NOT PROVIDED**

1. Allergen Extraction Kit (Neogen item 8429)
  - a. 20 disposable plastic extraction bottles
  - b. 20 sample collection tubes (12 x 75 mm) with caps
2. Allergen Environmental Swabbing Kit (Neogen item 8432S)
  - a. 100 sterile swabs
  - b. 100 dropper tips
3. Shaker water bath capable of maintaining  $60 \pm 1^\circ\text{C}$  with clamps to hold 250 mL disposable plastic bottles
4. Whatman #4 filters or equivalent (Neogen item 9429)
5. Centrifuge (optional)
6. Pipettor, adjustable 50–200  $\mu\text{L}$  (Neogen item 9276)
7. Pipettor, 12-channel (Neogen item 9273)
8. Pipette tips (Neogen item 9410, 9407, 9417)
9. Timer (Neogen item 9426)
10. Microwell reader with a 650 nm filter (Neogen item 9303)
11. 1 L bottle to prepare washing solution (Neogen item 9472)
12. 1 L heat safe bottle to prepare extract solution (Neogen item 9472)
13. Paper towels or equivalent absorbent material
14. Microwell holder (Neogen item 9402)
15. Waterproof marker
16. Wash bottle (Neogen item 9400)
17. Distilled or deionized water
18. 3 reagent boats for 12-channel pipettor (Neogen item 9435)
19. Graduated cylinder capable of measuring 125 mL (Neogen item 9368)
20. Scale capable of weighing  $5 \pm 0.1$  g (Neogen item 9427)

## **PRECAUTIONS**

1. Sesame seeds are particulates. To allow extraction of sesame seed protein it is important that samples are ground prior to testing.
2. Components of Veratox for Sesame Allergen, such as controls and extraction reagents, may contain one or more of the following potentially allergic materials: casein, sesame protein and soy protein. If allergic to any of these compounds, please use caution when using this product.
3. Concentrated food additives, colors and flavors may cause interferences on ELISA test methods. Contact Neogen's technical services for validation information.
4. Hydrolyzed and fermented proteins may not be detected using ELISA methods for allergen testing. Due to the breakdown of the proteins to small peptides or amino acids, they may become undetectable by this assay, but still could be allergenic and cause an allergic reaction.
5. Sponges should not be used for sample collection and allergen testing. Sample collection swabs other than Neogen swabs should be validated prior to use. General sponges and swabs may contain solutions or materials that may interfere with the test kit.
6. Store test kit between  $2\text{--}8^\circ\text{C}$  ( $35\text{--}46^\circ\text{F}$ ) when not in use. Do not freeze test kits, and avoid prolonged storage of kits at ambient temperatures.
7. Bring kits to room temperature  $18\text{--}30^\circ\text{C}$  ( $64\text{--}86^\circ\text{F}$ ) prior to use.
8. Cool sample extracts to room temperature prior to testing.
9. Do not use kit components beyond expiration date.
10. Do not mix reagents from one kit serial with reagents from a different kit serial.
11. Do not run more than 24 wells per test when running the quantitative method.
12. When running the screening procedure, do not run more than 6 wells unless using a multichannel pipettor.

13. Veratox for Sesame assay cross reacts with almond, macadamia nut, linseed and sunflower seed when tested at 100%. Please see Veratox for Sesame validation report for additional information.
14. Follow proper pipetting techniques (e.g., prime tips and use clean tips).
15. Use only incubation times specified. Others may give inaccurate results.
16. Use clean pipette tips and glassware for each sample to avoid cross-contamination. Thoroughly wash all glassware between samples.

## PROCEDURAL NOTES

1. **Substrate:** K-Blue Substrate is ready for use. The substrate should be clear to light blue—discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until needed.
  2. **Conjugate.** The conjugate supplied with this kit is ready to use. One bottle is enough for 24 wells. Cover the reagent boat to keep the conjugate protected from direct light and contaminants.
  3. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.
  4. **Extraction solution.** Prepare extraction solution by adding a foil pouch of extraction solvent, 10 mM PBS, to 1 L distilled or deionized water. Swirl to mix thoroughly. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).
  5. **Wash buffer.** Prepare the wash buffer solution by pouring all the wash buffer concentrate into an empty 1 L container. Rinse the wash buffer concentrate bottle with distilled or deionized water and pour into the 1 L container to ensure all the concentrate is used. Fill the 1 L container with additional distilled or deionized water, and swirl to assure thorough mixing. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).
- NOTE:** Discard unused portions of extraction solution and wash buffer when the test kit has been used completely.

## SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques (see Neogen's Food Allergen Handbook). The sample should be ground and thoroughly mixed prior to proceeding with the extraction procedure. For collecting and extracting environmental swabs, reference the Allergen Environmental Swabbing Kit (Neogen item 8432S) instructions.

1. Prepare the extraction solution as described in the procedural notes.
2. Preheat extraction solution to 60°C (140°F) by immersing the bottle containing the solution into the water bath and allowing it to reach 60°C.
3. Using your sampling and collection procedure, obtain a representative sample and grind it to a fine particle size.
4. Transfer 5 g of sample or 5 mL of liquid sample into a disposable extraction bottle.
5. Add one level scoop of the extraction additive to the sample bottle. Do not use the extraction additive from another allergen test kit.
6. Pour 125 mL of the 60°C (140°F) extraction solution to the sample bottle.
7. Cap the sample bottle to prevent contents from splashing during the extraction. Shake to mix before placing in the water bath.
8. Extract by shaking (150 rpm) in a water bath at 60°C (140°F) for 15 minutes. Remove the bottle from the bath.
9. Let material settle for 5 minutes to enable some of the sample to settle before proceeding to the next step.

10. Filter the extract by pouring at least 5 mL through a Whatman #4 filter and collecting the filtrate as a sample.  
**ALTERNATIVE:** Centrifuge at 14,000 rpm for 5 minutes (20 minutes for lower speeds). Use the clear supernatant as a sample.
11. Allow extracts to cool to room temperature before beginning analysis.
12. Discard extracts after completion of analysis.

### TEST PROCEDURE FOR QUANTITATION

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

1. Remove 1 red-marked mixing well for each sample to be tested plus 5 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of the strip with a "1", and place strip in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Using a new pipette tip for each, transfer 150 µL of controls and sample extracts to the red-marked transfer wells as shown in the template below. Only run up to two 12-well strips at a time.

0	2.5	5	10	25	S1	S2	S3	S4	S5	S6	S7	Strip 1
S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19	Strip 2

5. Place tips on the 12-channel pipettor and transfer 100 µL of the controls and sample extracts to the antibody-coated wells. Mix for 20 seconds by sliding the well holder back and forth on a flat surface.
6. Incubate microwells **10 minutes** at room temperature 18–30°C (64–86°F). Discard the red-marked transfer wells.
7. Empty the contents of the wells into a sink. With a wash bottle fill each antibody well with the wash buffer solution and dump out. Repeat the washing 4 times, then turn the wells upside down and tap out on a paper towel until all washing solution is removed.
8. Pour the needed volume of conjugate from the blue-labeled bottle into a clean reagent boat.
9. Using the 12-channel pipettor and new tips, transfer 100 µL of the conjugate into all the wells and mix for 20 seconds by sliding the well holder back and forth on a flat surface.
10. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
11. Wash all wells with the wash buffer solution as described in step 7.
12. Pour the needed volume of substrate solution from the green-labeled bottle into a clean reagent boat.
13. Place new tips on the 12-channel pipettor and transfer 100 µL of substrate into each well and mix for 20 seconds. Do not eject tips.
14. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
15. Pour the needed volume of Red Stop solution from the red-labeled bottle into a clean reagent boat.
16. With the same tips used to dispense the substrate, transfer 100 µL of Red Stop into each well and mix for 20 seconds.
17. Wipe the bottom of the microwells and read in a microwell reader with a 650 nm filter. Air bubbles should be eliminated, as they could affect analytical results. Results should be read within **20 minutes** after the addition of Red Stop solution.
18. Interpret the test's results using Neogen's microwell reader, or an equivalent strip reader. If using a strip reader, calculate the results using Neogen's Veratox for Windows® software.

## TEST PROCEDURE FOR SCREENING

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

1. Remove 1 well for each sample to be tested plus 1 well for the control, and place into the well holder.
2. Choose the 5 ppm yellow-labeled control bottle to serve as the screening level for the test.
3. Mix each reagent by swirling its bottle prior to use.
4. Add 100 µL from the yellow-labeled control bottle to the first well. Add 100 µL from each sample extract to a respective well as indicated in the template below. For environmental swabs, add 3 drops from swab tube with dropper tip. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.

Control	S1	S2	S3	S4	S5
---------	----	----	----	----	----

5. Incubate microwells **10 minutes** at room temperature 18–30°C (64–86°F).
6. Shake out the contents of the wells. Using a wash bottle filled with wash buffer solution, fill each well and shake out. Repeat 4 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
7. Add 100 µL from the blue-labeled conjugate bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.
8. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
9. Shake out the contents of the wells. Using a wash bottle filled with wash buffer solution, fill each well and shake out. Repeat 4 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
10. Add 100 µL from the green-labeled substrate bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.
11. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
12. Add 100 µL from the red-labeled Red Stop bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface. The results are now ready to be interpreted.
13. Visually compare the color of a sample well to the color of the control well. If the sample well has more blue color than the control well, the sample tests positive for sesame contamination of more than the control used. If the sample well has less blue color, or more red color, than the control well, the sample contains less than the control used of sesame contamination.

**Alternative:** Read wells (wipe the bottom of wells with a dry cloth or towel first) in a microwell reader with a 650 nm filter. If the sample well has an optical density (OD) higher than the control well, the sample is positive for sesame contamination of more than the control used. If the sample well has an OD lower than the control well, the sample contains less than the control used of sesame contamination.

## PERFORMANCE CHARACTERISTICS

Limit of quantitation: 2.5 ppm (Described as the lowest concentration point on the calibration curve that this test can reliably detect sesame.)

Range of quantitation: 2.5–25 ppm (For quantitating samples above 25 ppm, contact a Neogen representative for dilution instructions.)

Allergen detection: This test detects sesame protein, and the results are expressed as ppm of total sesame.

Cross-reactivity: Cross reacts with almond, macadamia nut, linseed and sunflower seed when tested at 100%. Please see Veratox for Sesame validation report for additional information.

**CUSTOMER SERVICE**

Neogen Customer Assistance and Technical Service can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

**SDS INFORMATION AVAILABLE**

Safety data sheets (SDS) are available for this test kit, and all of Neogen's Food Safety test kits, at [foodsafety.neogen.com](http://foodsafety.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

**TERMS AND CONDITIONS**

For Neogen's full terms and conditions, please visit [www.neogen.com/Corporate/termsconditions.html](http://www.neogen.com/Corporate/termsconditions.html).

**WARRANTY**

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## TESTING KITS AVAILABLE FROM NEOGEN

### Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

### Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

### Genetic modification

- CP4 (Roundup Ready®)

### Ruminant by-products

- Meat and bone meal, feed

### Species Identification

- Raw and cooked meat samples

#### North America

##### Neogen Headquarters

800/234-5333 (USA/Canada)

foodsafety@neogen.com

foodsafety.neogen.com

#### Europe, Middle East and Africa

##### Neogen Europe

+ 44 (0) 1292 525 600

info\_uk@neogeneurope.com

www.neogeneurope.com

#### Mexico

##### Neogen Latinoamerica

+52 (55) 5254-8235

informacion@neogenlac.com

www.neogenlac.com

#### Brazil

##### Neogen do Brasil

+55 19 3935.3727

info@neogendobrasil.com.br

www.neogendobrasil.com.br

#### China

##### Neogen Bio-Scientific Technology

+86 21 6271 7013

info@neogenchina.com.cn

www.neogenchina.com.cn

#### India

##### Neogen Food and Animal Security

+91 484 2306598, 2301582

info@neogenindia.com

www.neogenindia.com

©Neogen Corporation, 2017. Neogen, K-Blue and Veratox are registered trademarks of Neogen Corporation. All other trademarks are the property of their respective companies. Patent: <http://www.neogen.com/Corporate/patents.html>

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