U.S. English

Product Number:
8024

Reveal[®] Q+

for d9-THC in Hemp

Endpoint test for use with Raptor[®] Solo

Reveal[®] Q+ for d9-THC in Hemp

Product Number: 8024

Introduction — Cannabis and Delta-9-Tetrahydrocannabinol (Δ 9-THC)

Delta-9-tetrahydrocannabinol (Δ 9-THC) is the principal psychoactive constituent of cannabis. Cannabis is characterized as either hemp or marijuana based on the content of total Δ 9-THC in the plant. Total THC is defined as the sum of both Δ 9-THC and the acid form of Δ 9-THC (tetrahydrocannabinolic acid; THCA). THCA is converted to Δ 9-THC with heat.

The Reveal[®] Q+ for d9-THC in Hemp rapid endpoint test was designed to quantify and monitor Δ 9-THC concentrations in hemp throughout the growing season. The test was developed for use on the farm. It removes the need to send samples out to an analytical lab, reducing expenses and time to result. However, this test does not replace the final Δ 9-THC quantification by a certified analytical lab prior to harvest.

Intended Use and User

Reveal Q+ for d9-THC in Hemp is intended for the quantitative analysis of total Δ 9-THC in hemp. The test kit is designed for use by hemp farmers and others familiar with commodities containing Δ 9-THC.

Assay Principles

Reveal Q+ for d9-THC in Hemp is a single-step lateral flow immunoassay. It is used to test dried hemp flower samples that are decarboxylated to convert the THCA to total Δ 9-THC. The decarboxylated hemp samples are ground and extracted with isopropanol to solubilize the Δ 9-THC.

The test strip is added to a small sample cup containing the diluted hemp extract. The extract is wicked through the reagent zone, which contains antibodies specific for Δ 9-THC bound to a gold nanoparticle. If Δ 9-THC is present, it will be captured by the antibody. The extract is then wicked onto a membrane, which contains a zone of Δ 9-THC conjugated to a protein carrier. This zone captures any unbound gold nanoparticles, allowing the particles to concentrate and form a visible line. As the level of Δ 9-THC in a sample increases, free Δ 9-THC will bind with the gold nanoparticle allowing less gold complex to be captured in the test zones. Therefore, as the concentration of Δ 9-THC in the sample increases, the test line density decreases. The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of Δ 9-THC, ensuring the strip is functioning properly.

This assay is analyzed exclusively with the Raptor[®] Solo Integrated Analysis Platform. The reader quantifies the amount of total Δ 9-THC in a sample based on the ratio of the test line intensity and the control line intensity. Algorithms programmed into the Raptor Solo reader convert the line intensities into a quantitative result displayed in percent (%) THC on the reader's screen.

Storage Requirements

Store kit components at ambient temperature of 18–30°C (64–86°F) to ensure full shelf life. Test strips should remain capped in their original tubes until used to ensure optimal performance.

Materials Provided

1. 25 Reveal Q+ test strips
2. 25 prefilled sample dilution bottles containing 5.0 mL of sample diluent
3. 25 clear sample cups

NEOGEN.com

Materials Required, But Not Provided

1. Raptor Solo Integrated Analysis Platform Reader (NEOGEN item 9696)
2. Decarboxylator
3. THC extraction materials
 - a. Cups with lids for sample collection (125 mL; NEOGEN item 9428)
 - b. 70% isopropyl rubbing alcohol
4. Agri-grind grinder or equivalent (NEOGEN item 9401, 9453)
5. Scale capable of weighing $1 \text{ g} \pm 0.1 \text{ g}$ (NEOGEN item 9427)
6. Fixed volume 100 μL pipet (NEOGEN item 9860)
7. 200 μL pipet tips (NEOGEN item 9407, 9410, 9417)
8. 25 mL graduated cylinder (NEOGEN item 90447)
9. Food dehydrator for drying down sample prior to decarboxylation

Materials Recommended, But Not Provided

1. Timer (NEOGEN item 9426, 9452)
2. Reveal sample cup rack (NEOGEN 9475)

Precautions

1. The test strips must remain inside the stay-dry tube before use.
2. Store test kit at ambient temperature of 18–30°C (64–86°F) when not in use.
3. Do not use kit components beyond expiration date.
4. To avoid cross-contamination, use clean clear sample cups for each test.
5. Ensure the test strip lot number and the curve details match the lot ID number selected on the reader.

Raptor Solo Test Setup

The Raptor Solo reader must be in standard mode to read the Reveal Q+ for d9-THC in Hemp test strips.

1. From the main menu, select the  in the upper left-hand corner.
2. Select track replacement.
3. Return to the home screen by pressing the  in the lower left-hand corner.
4. Insert the track replacement tool with the arrow facing towards you into the track in the unit.
5. Gently pull up to remove the Raptor cartridge track.
6. Insert the standard endpoint testing track into the unit (track with the blue dot on the right-hand side).

Sample Preparation

The flower sample to be tested should be collected according to accepted local sampling techniques (contact your NEOGEN representative). Obtain a representative flower sample. To address total $\Delta 9$ -THC by dry weight, hemp samples will need to be dried and decarboxylated, converting THCA to $\Delta 9$ -THC, prior to testing.

Drying

Weigh out a 5.0–7.0 g sample and dry in a food dehydrator on high setting (not to exceed 150°F/66°C) until sample is very dry and brittle.

Decarboxylation With Ardent FX Decarboxylator

1. Weigh out a 1.5–2.0 g dry sample and place into decarboxylator.
2. Put silicone cover on decarboxylation chamber and place chamber into main unit. Place cap on main unit.
3. Plug in decarboxylator, button should have a green light.
4. Press M button to switch to A1 mode.
5. Press power button to begin 1.5 to 2-hour cycle. Red light will show while running.
6. When cycle is complete light will turn green.
7. Empty sample into grinder. Grind for 10 seconds. Ideally, sample should be the particle size of fine espresso.
8. Clean and dry decarboxylator between uses.

NEOGEN.com

Sample Extraction

1. Weigh out 1.0 g of decarboxylated hemp sample into a sample collection cup.
2. Using graduated cylinder add 25 mL of 70% isopropyl alcohol to collection cup containing sample.
3. Put lid on collection cup and shake vigorously side to side for at least 1 minute to extract THC from sample. Shaking in an up and down motion may result in leaking out of the cup.

Test Procedure

1. Place the appropriate number of clear sample cups into a sample cup rack. Label cups if necessary.
2. Transfer 100 µl of the extracted THC solution to a pre-filled sample dilution tube.
3. Cap dilution tube and mix by inverting at least 10 times.
4. Transfer 100 µl of diluted sample extract (from step 3) into a new clear sample cup.
5. Place a new Reveal Q+ test strip into the sample cup with the sample end down and set timer for 5 minutes. Ensure the test strip comes into contact with liquid and begins to wick.
6. Remove the strip from the sample cup after it has developed for exactly 5 minutes. Immediately (within 30 seconds) read the test strip on the Raptor Solo reader.

Reading a Test Strip

1. From the home screen.
Select category: Cannabinoid Q+.
Select test: THC in Hemp.
2. Press run test.
Select lot.
(If not listed, add the lot by scanning the QR code located on the test strip tube).
3. Enter sample ID as needed — press the ✓ in the lower right corner of the reader.
4. Fully insert the Reveal Q+ for d9-THC in Hemp test strip into the NEOGEN labeled cartridge adaptor with the sample end first.
5. Insert cartridge into the track on the Raptor Solo reader — the system will read automatically.
6. Results are displayed on the unit.
7. To read the next test strip press change at the bottom of the screen — this allows the user to enter the next sample ID. Press ✓ and insert cartridge containing the next test strip.

Notes

1. Ensure test strip is fully inserted into cartridge.
2. The Reveal Q+ for d9-THC in Hemp test strips can only be read using the Raptor Solo reader. Test strips should be read within 30 seconds of completing the 5-minute incubation to ensure the most accurate results.

Performance Characteristics

1. Range of detection: 0.1 to 1.6% Δ9-THC.

Note: Results below the range of detection should be considered as less than 0.1% Δ9-THC.

Validated and Verified Matrices

NEOGEN continues to validate new commodities. Please contact a representative for the latest validated commodity list.

Customer Service

NEOGEN Customer and Technical Services can be contacted through [NEOGEN.com](https://www.neogen.com) and product training is available by request.

SDS Information Available

Safety data sheets are available for all test kits at [NEOGEN.com](https://www.neogen.com) or by calling 800.234.5333 or 517.372.9200.

Terms and conditions

NEOGEN's full terms and conditions are available [online](#).

Warranty

NEOGEN makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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