- 
-  **(EN)** **Molecular Detection Assay 2 - *Salmonella* Enteritidis / *Salmonella* Typhimurium screening kit**
  -  **(ES)** **Ensayo de Detección Molecular 2 - Kit de detección de *Salmonella* Enteritidis/*Salmonella* Typhimurium**
  -  **(PT)** **Ensaio de Detecção Molecular 2 – Kit de detecção de *Salmonella* Enteritidis/*Salmonella* Typhimurium**

2

## Product Instructions

# Neogen® Molecular Detection Assay 2 - *Salmonella* Enteritidis / *Salmonella* Typhimurium screening kit

### Product Description and Intended Use

The Neogen® Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium kit is used with the Neogen® Molecular Detection System for rapid and specific detection of two *Salmonella* subtypes, *Salmonella* enterica subsp. *enterica* serotype Enteritidis (SE) and *Salmonella* enterica subsp. *enterica* serotype Typhimurium (including monophasic variants (ST)), in enriched poultry and poultry-related samples. The screening kit contains two separate reagent tubes, one to detect SE serotype strains and the other to detect ST serotype strains. The Neogen® Molecular Detection System Software reports results separately for each of the assays.

The assay uses loop-mediated isothermal amplification to rapidly amplify nucleic acid sequences with high specificity and sensitivity, combined with bioluminescence to detect the amplification. Presumptive positive results are reported in real-time while negative results are displayed after the assay is completed. Presumptive positive results should be confirmed using your preferred method or as specified by local regulations<sup>(1, 2, 3)</sup>.

The Neogen® Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium is intended for use in a laboratory environment by professionals trained in laboratory techniques. Neogen has not documented the use of this product in industries other than food or beverage. For example, Neogen has not documented this product for testing pharmaceutical, cosmetics, clinical or veterinary samples. The Neogen® Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium has not been evaluated with all possible food products, processes, testing protocols or with all possible strains of bacteria. **As with all test methods, the source, formulation, and quality of enrichment medium can influence the results.** Factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may also influence results. Neogen recommends evaluation of the method including enrichment medium, in the user's environment using a sufficient number of samples with particular foods and microbial challenges to ensure that the method meets the user's criteria.

Neogen has evaluated the Neogen Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium with Buffered Peptone Water<sup>(1,2)</sup> and Buffered Peptone Water ISO<sup>(3)</sup>.

The Neogen® Molecular Detection Instrument is intended for use with samples that have undergone heat treatment during the assay lysis step, which is designed to destroy organisms present in the sample. Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should NOT be inserted into the Neogen Molecular Detection Instrument.

Neogen Food Safety is certified to ISO (International Organization for Standardization) 9001 for design and manufacturing.

The Neogen® Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium test kit contains 48 tests of each of SE and ST reagents, described in Table 1.

**Table 1. Kit Components**

Item	Identification	Quantity	Contents	Comments
Lysis Solution (LS) tubes	Pink solution in clear tubes	96 (12 strips of 8 tubes)	580 µL of LS per tube	Racked and ready to use
<i>Salmonella</i> Enteritidis (SE) Reagent tubes	Beige tubes	48 (2 pouches; each containing 3 strips of 8 tubes)	Lyophilized specific amplification and detection mix	Ready to use
<i>Salmonella</i> Typhimurium (ST) Reagent tubes	Grey tubes	48 (2 pouches; each containing 3 strips of 8 tubes)	Lyophilized specific amplification and detection mix	Ready to use
Extra caps	Beige caps	96 (12 strips of 8 caps)		Ready to use
Extra caps	Grey caps	96 (12 strips of 8 caps)		Ready to use
Reagent Control (RC)	Clear flip-top tubes	16 (2 pouches of 8 individual tubes)	Lyophilized control DNA, amplification, and detection mix	Ready to use

The Negative Control, not provided in the kit, is a sterile enrichment medium, e.g., BPW or BPW ISO. Do not use water as a Negative Control.

**Safety**

The user should read, understand, and follow all safety information in the instructions for the Neogen® Molecular Detection System and the Neogen® Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium. Retain the safety instructions for future reference.

**⚠ WARNING:** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

**⚠ CAUTION:** Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury and/or property damage.

**NOTICE:** Indicates a potentially hazardous situation which, if not avoided, could result in property damage.

**⚠ WARNING**

**Do not use the Neogen® Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium in the diagnosis of conditions in humans or animals.**

**The user must train its personnel in current proper testing techniques: for example, Good Laboratory Practices, ISO/IEC 17025<sup>(4)</sup>, or ISO 7218<sup>(5)</sup>.**

**To reduce the risks associated with a false-negative result leading to the release of contaminated product:**

- Follow the protocol and perform the tests exactly as stated in the product instructions.
- Store the Neogen® Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium as indicated on the package and in the product instructions.
- Always use the Neogen® Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium by the expiration date.
- Use the Neogen Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium with foods that have been validated internally or by a third party.

**To reduce the risks associated with exposure to chemicals and biohazards:**

- Perform pathogen testing in a properly equipped laboratory under the control of trained personnel. Incubated enrichment media and equipment or surfaces that have come into contact with incubated enrichment media may contain pathogens at levels sufficient to cause risk to human health.
- This procedure uses/detects pathogenic microorganisms and/or their metabolic products. Care should be taken to avoid ingestion or inhalation of potentially infectious aerosols or contact with the skin. Always follow standard laboratory safety practices, including wearing appropriate protective apparel and eye protection while handling reagents and contaminated samples.
- Avoid contact with the contents of the enrichment media and reagent tubes after amplification.
- Dispose of enriched samples according to current industry standards.
- Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should NOT be inserted into the Neogen Molecular Detection Instrument.

**To reduce the risks associated with cross-contamination while preparing the assay:**

- Always wear gloves (to protect the user and prevent introduction of nucleases).

**To reduce the risks associated with environmental contamination:**

- Follow current industry standards for disposal of contaminated waste.

**⚠ CAUTION**

**To reduce the risks associated with exposure to hot liquids:**

- Do not exceed the recommended temperature setting on heater.
- Do not exceed the recommended heating time.
- Use an appropriate, calibrated thermometer to verify the Neogen® Molecular Detection Heat Block Insert temperature (e.g., a partial immersion thermometer or digital thermocouple thermometer, not a total immersion thermometer). The thermometer must be placed in the designated location in the Neogen Molecular Detection Heat Block Insert.

**NOTICE**

**To reduce the risks associated with cross-contamination while preparing the assay:**

- Change gloves prior to reagent pellet hydration.
- Use of sterile, aerosol barrier (filtered), molecular biology grade pipette tips is recommended.
- Use a new pipette tip for each sample transfer.
- Use Good Laboratory Practices to transfer the sample from the enrichment to the lysis tube. To avoid pipettor contamination, the user may choose to add an intermediate transfer step. For example, the user can transfer each enriched sample into a sterile tube.
- Use a molecular biology workstation containing germicidal lamp where available. Periodically decontaminate laboratory benches and equipment (pipettes, cap/decap tools, etc.) with a 1-5% (v:v in water) household bleach solution or DNA removal solution.

**To reduce the risks associated with a false-positive result:**

- Never open tubes post amplification.
- Always dispose of the contaminated tubes by soaking in a 1-5% (v:v in water) household bleach solution for 1 hour and away from the assay preparation area.
- Never autoclave reagent tubes post amplification.

Consult the Safety Data Sheet for additional information and local regulations for disposal.

If you have questions about specific applications or procedures, please visit our website at [www.neogen.com/foodsafety](http://www.neogen.com/foodsafety) or contact your local Neogen representative or distributor.

**User Responsibility**

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at [www.neogen.com](http://www.neogen.com), or contact your local Neogen representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, laboratory technique and the sample itself may influence results.

It is the user’s responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user’s criteria.

It is also the user’s responsibility to determine that any test methods and results meet its customers’ and suppliers’ requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

To help customers evaluate the method for various food matrices, Neogen has developed the Neogen® Molecular Detection Matrix Control kit. When needed, use the Matrix Control (MC) to determine if the matrix has the ability to impact the Neogen Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium results. Test several samples, representative of the matrix, i.e. samples obtained from different origin, during any validation period when adopting the Neogen method or when testing new or unknown matrices or matrices that have undergone raw material or process changes.

A matrix can be defined as a type of product with intrinsic properties such as composition and process. Differences between matrices may be as simple as the effects caused by differences in their processing or presentation, for example, raw vs. pasteurized; fresh vs. dried, etc.

**Limitation of Warranties / Limited Remedy**

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen Food Safety representative or authorized Neogen distributor for any further questions.

**Limitation of Neogen Liability**

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen’s liability under any legal theory exceed the purchase price of the product alleged to be defective.

**Storage and Disposal**

Store the Neogen® Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium at 2-8°C. Do not freeze. Keep the kit away from light during storage. After opening the kit, check that the foil pouch is undamaged. If the pouch is damaged, do not use it. After opening, unused reagent tubes should always be stored in the re-sealable pouch with the desiccant inside to maintain stability of the lyophilized reagents. Store resealed pouches at 2-8°C for no longer than 60 days. Do not use Neogen® Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium past the expiration date. Expiration date and lot number are noted on the outside label of the box. After use, the enrichment medium and the Neogen® Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium tubes can potentially contain pathogenic materials. When testing is complete, follow current industry standards for the disposal of contaminated waste. Consult the Safety Data Sheet for additional information and local regulations for disposal.

## Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Periodically decontaminate laboratory benches and equipment (pipettes, cap/decap tools, etc.) with a 1- 5% (v:v in water) household bleach solution or DNA removal solution.

The user should complete the Neogen® Molecular Detection System operator qualification (OQ) training, as described in the “Installation Qualification (IQ) / Operational Qualification (OQ) Protocols and Instructions for Neogen Molecular Detection System” document<sup>(6)</sup>.

See Section “Specific Instructions for validated methods” for specific requirements:

Table 3 for enrichment protocols according to AOAC® Performance Tested<sup>SM</sup> Certificate #122302.

### Sample Enrichment

Tables 2 or 3 present guidance for enrichment protocols for poultry and poultry-related samples. It is the user’s responsibility to validate alternate sampling protocols or dilution ratios to ensure this test method meets the user’s criteria.

### Foods

1. Allow BPW or BPW ISO enrichment medium to equilibrate to ambient laboratory temperature. See Table 2 or 3.
2. Aseptically combine the enrichment medium and sample. For all meat and highly particulate samples, the use of filter bags is recommended.
3. Homogenize thoroughly by blending, stomaching, or hand mixing for 2 ±0.2 minutes until the enrichment suspension is homogeneous<sup>(7,8)</sup>. For all meat and highly particulate samples, the use of filter bags is recommended. Incubate as outlined in the appropriate protocol table (see Table 2 or 3).

**Table 2.** General enrichment protocols.

Sample Matrix	Sample Size	Enrichment Broth Volume	Enrichment Temperature	Enrichment Time (hr)	Sample Analysis Volume <sup>(a)</sup>
Raw ground chicken	325 g	1625 mL BPW	35 ± 2 °C	20 - 24	20 µL
Chicken rinse (400 mL nBPW)	30 mL	30 mL BPW	35 ± 2 °C	20 - 24	20 µL
Raw ground chicken	25 g	225 mL BPW ISO	41.5 ± 1 °C	18 - 24	20 µL
Cooked breaded chicken (chicken nuggets)	25 g	225 mL BPW ISO	37 ± 1 °C	18 - 24	20 µL

(a) Volume of sample transferred to Lysis Solution tubes. Refer to step 5.6 of Lysis section.

### Specific Instructions for Validated Methods

AOAC® Performance Tested<sup>SM</sup> Certificate #122302

In the AOAC Research Institute PTM<sup>SM</sup> program, the Neogen Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium was found to be an effective method for the detection of two *Salmonella* subtypes, *Salmonella* enterica subsp. *enterica* serotype Enteritidis (SE) and *Salmonella* enterica subsp. *enterica* serotype Typhimurium (including monophasic variant 1,4,[5],12;i:-) (ST). The matrices tested in the study are shown in Table 3.

**Table 3.** Enrichment protocols according AOAC PTM<sup>SM</sup> Certificate #122302.

Sample Matrix	Sample Size	Enrichment Broth Volume	Enrichment Temperature	Enrichment Time (hr)	Sample Analysis Volume <sup>(a)</sup>
Raw ground chicken	325 g	1625 mL BPW	35 ± 2 °C	20 - 24	20 µL
Chicken rinse (400 mL nBPW)	30 mL	30 mL BPW	35 ± 2 °C	20 - 24	20 µL
Raw ground chicken	25 g	225 mL BPW ISO	41.5 ± 1 °C	18 - 24	20 µL
Cooked breaded chicken (chicken nuggets)	25 g	225 mL BPW ISO	37 ± 1 °C	18 - 24	20 µL

(a) Volume of sample transferred to Lysis Solution tubes. Refer to step 5.6 of Lysis section.

**Preparation of the Neogen® Molecular Detection Speed Loader Tray**

1. Wet a cloth or disposable towel with a 1-5% (v:v in water) household bleach solution and wipe the Neogen Molecular Detection Speed Loader Tray.
2. Rinse the Neogen Molecular Detection Speed Loader Tray with water.
3. Use a disposable towel to wipe the Neogen Molecular Detection Speed Loader Tray dry.
4. Ensure the Neogen Molecular Detection Speed Loader Tray is dry before use.

**Preparation of the Neogen® Molecular Detection Chill Block Insert**

Place the Neogen Molecular Detection Chill Block Insert directly on the laboratory bench: The Neogen Molecular Detection Chill Block Tray is not used. Use the block at ambient laboratory temperature (20-25°C).

**Preparation of the Neogen® Molecular Detection Heat Block Insert**

Place the Neogen Molecular Detection Heat Block Insert in a dry double block heater unit. Turn on the dry block heater unit and set the temperature to allow the Neogen Molecular Detection Heat Block Insert to reach and maintain a temperature of 100 ±1°C.

**NOTE:** Depending on the heater unit, allow approximately 30 minutes for the Neogen Molecular Detection Heat Block Insert to reach temperature. Using an appropriate, calibrated thermometer (e.g., a partial immersion thermometer, digital thermocouple thermometer, not a total immersion thermometer) placed in the designated location, verify that the Neogen Molecular Detection Heat Block Insert is at 100 ±1°C.


**Preparation of the Neogen® Molecular Detection Instrument**

1. Launch the Neogen® Molecular Detection Software and log in. Contact your Neogen Food Safety representative to ensure you have the most updated version of the software.
2. Turn on the Neogen Molecular Detection Instrument.
3. Create or edit a run with data for each sample. Refer to the Neogen Molecular Detection System User Manual for details.
  - 3.1. For the run of the SE and ST assays, select the icon for each assay in the software to run the assays for each sample.
  - 3.2. The SE and ST assays can also be run separately for each sample.
  - 3.3. NC is set up for each of the reagent tubes and one RC is set up for the kit.

**NOTE:** The Neogen Molecular Detection Instrument must reach and maintain a temperature of 60°C before inserting the Neogen Molecular Detection Speed Loader Tray with reaction tubes. This heating step takes approximately 20 minutes and is indicated by an ORANGE light on the instrument’s status bar. When the instrument is ready to start a run, the status bar will turn GREEN.

**Lysis**


1. Remove the bottom of Neogen Lysis Solution Rack with a screwdriver or spatula before placing it in the Neogen Molecular Detection Heat Block Insert.


2. Allow the LS tubes to warm up by setting the rack at room temperature (20-25 °C) overnight (16-18 hours). Alternatives to equilibrate the LS tubes to room temperature are to set the LS tubes on the laboratory bench for at least 2 hours, incubate the LS tubes in a 37 ±1°C incubator for 1 hour or place them in a dry double block heater for 30 seconds at 100°C.
3. Invert the capped tubes to mix. Proceed to the next step within 4 hours.
4. Remove the enrichment broth from the incubator.
5. One LS tube is required for each sample and the Negative Control (NC) sample (sterile enrichment medium).
  - 5.1 LS tube strips can be cut to desired LS tube number. Select the number of individual LS tubes or 8-tube strips needed. Place the LS tubes in an empty rack.
  - 5.2 To avoid cross-contamination, decap one LS tubes strip at a time and use a new pipette tip for each transfer step.
  - 5.3 Transfer enriched sample to LS tubes as described below:

Transfer each enriched sample into an individual LS tube **first**. Transfer the NC **last**.

- 5.4 Use the Neogen® Molecular Detection Cap/Decap Tool-Lysis to decap one LS tube strip – one strip at a time.
- 5.5 Discard the LS tube cap – If lysate will be retained for retest, place the caps into a clean container for re-application after lysis.
- 5.5.1 For processing of retained lysate, see Appendix A.
- 5.6 Transfer 20 µL of sample into a LS tube unless otherwise indicated in the protocol table.
6. Repeat step 5.2 until each individual sample has been added to a corresponding LS tube in the strip.

7. Repeat steps 5.1 to 5.6 as needed, for the number of samples to be tested.
8. When all samples have been transferred, transfer 20 µL of NC (sterile enrichment medium e.g., BPW) into a LS tube. Do not use water as a NC.
9. Verify that the temperature of the Neogen Molecular Detection Heat Block Insert is at 100 ±1°C.
10. Place the uncovered rack of LS tubes in the Neogen Molecular Detection Heat Block Insert and heat for 15 ±1 minutes. During heating, the LS solution will change from pink (cool) to yellow (hot).  
Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should NOT be inserted into the Neogen Molecular Detection Instrument.
11. Remove the uncovered rack of LS tubes from the heating block and allow to cool in the Neogen Molecular Detection Chill Block Insert at least 5 minutes and a maximum of 10 minutes. The Neogen Molecular Chill Block Insert, used at ambient temperature without the Molecular Detection Chill Block Tray, should sit directly on the laboratory bench. When cool, the lysis solution will revert to a pink color.
12. Remove the rack of LS tubes from the Neogen Molecular Detection Chill Block Insert.


### Amplification


1. One Neogen Molecular Detection Assay 2 – *Salmonella* Enteritidis (SE) and one Neogen Molecular Detection Assay 2 – *Salmonella* Typhimurium (ST) are required for each sample and the NC.
  - 1.1 The SE and ST assays can be run in the same rack for each sample, or they can be run in a separate rack.
  - 1.2 Reagent tube strips can be cut to desired tube number. Select the number of individual Reagent tubes or 8-tube strips needed.
  - 1.3 Place SE tubes in an empty rack in one column
  - 1.4 Place ST tubes in the adjacent right column
  - 1.5 Avoid disturbing the reagent pellets from the bottom of the tubes.
2. Select one Reagent Control (RC) tube and place it in the rack.
3. For NC lysate, select one SE Reagent Tube and one ST Reagent Tube and place in the same rack, if both serotypes are running together, or in separate racks, if running them separate.
4. To avoid cross-contamination, decap one Reagent Tube strip at a time and use a new pipette tip for each transfer step.
5. Transfer each of the lysates to a SE and ST Reagent Tube as described below.
  - 5.1 First, transfer each sample lysate to a SE Reagent Tube as described in 5.5 and 5.6.
  - 5.2 Second, transfer each of the same sample lysate to a ST Reagent Tube in the adjacent right column as described in 5.5 and 5.6.

**NOTE: Use new pipette tip for each transfer. Do not use same pipette tip for transfer to SE and ST reagent tubes from the same lysate sample.**

- 5.3 After all sample lysate transfer, add NC lysate to each of SE and ST Reagent Tube.
- 5.4 Transfer NC lysate last to Reagent Control Tube.
- 5.5 Use the Neogen® Molecular Detection Cap/Decap Tool-Reagent to decap the SE and ST Reagent tubes –one Reagent tube strip at a time. Discard cap.
- 5.6 **Transfer 20 µL of sample lysate from the upper ½ of the liquid (avoid precipitate) in the Neogen Lysis Solution Tube into corresponding SE and ST Reagent Tube. Dispense at an angle to avoid disturbing the pellets. Mix by gently pipetting up and down 5 times.**

- 5.7 Repeat step 5.6 until all individual Sample lysate has been added to a corresponding SE and ST Reagent tubes in the strip.
  - 5.8 Cover the Reagent tubes with the provided extra caps and use the rounded side of the Neogen Molecular Detection Cap/Decap Tool-Reagent to apply pressure in a back and forth motion ensuring that the cap is tightly applied.
  - 5.9 Repeat steps 5.6 to 5.8 as needed, for the number of samples to be tested for both SE and ST Reagent Tubes.
  - 5.10 When all sample lysates have been transferred, repeat 5.6 to 5.8 to transfer 20 µL of NC lysate into each SE and ST Reagent tube.
  - 5.11 Transfer **20 µL of NC lysate into a RC tube**. Dispense at an angle to avoid disturbing the pellets. Mix by gently pipetting up and down 5 times.
6. Load capped tubes into a clean and decontaminated Neogen Molecular Detection Speed Loader Tray. Close and latch the Neogen Molecular Detection Speed Loader Tray lid.

- 7. Review and confirm the configured run in the Neogen Molecular Detection Software.
- 8. Click the Start button in the software and select instrument for use. The selected instrument’s lid automatically opens.
- 9. Place the Neogen Molecular Detection Speed Loader Tray into the Neogen Molecular Detection Instrument and close the lid to start the assay. Results are provided within 60 minutes, although positives may be detected sooner.
- 10. After the assay is complete, remove the Neogen Molecular Detection Speed Loader Tray from the Neogen Molecular Detection Instrument and dispose of the tubes by soaking in a 1-5% (v:v in water) household bleach solution for 1 hour and away from the assay preparation area.

**NOTICE:** To minimize the risk of false positives due to cross-contamination, never open reagent tubes containing amplified DNA. This includes Neogen Molecular Detection Assay 2 – *Salmonella* Enteritidis (SE) and Neogen Molecular Detection Assay 2 –*Salmonella* Typhimurium (ST) Reagent, Neogen Reagent Control, and Neogen Matrix Control Tubes. Always dispose of sealed reagent tubes by soaking in a 1-5% (v:v in water) household bleach solution for 1 hour and away from the assay preparation area.

**Results and Interpretation**

An algorithm interprets the light output curve resulting from the detection of the nucleic acid amplification. Results are analysed automatically by the software and are color-coded based on the result. A Positive or Negative result is determined by analysis of a number of unique curve parameters. Presumptive positive results are reported in real-time while Negative and Inspect results will be displayed after the run is completed.

The Neogen Molecular Detection System Software reports results separately for each of the assays (SE and ST).

Presumptive positive samples should be confirmed as per the laboratory standard operating procedures or by following the appropriate reference method confirmation, beginning with transfer from the primary enrichment to secondary enrichment broth(s), followed by subsequent plating and confirmation of isolates using appropriate biochemical and serological methods.

**NOTE:** Even a negative sample will not give a zero reading as the system and Neogen Molecular Detection Assay 2 – *Salmonella* Enteritidis / *Salmonella* Typhimurium amplification reagents have a “background” relative light unit (RLU) reading.

In the rare event of any unusual light output, the algorithm labels this as “Inspect.” Neogen recommends the user to repeat the assay for any Inspect samples. If the result continues to be Inspect, proceed to confirmation test using your preferred method or as specified by local regulations.

### References:

1. US Food and Drug Administration Bacteriological Analytical Manual. Chapter 5: *Salmonella*, September 2023 Version.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 4.14. Isolation and Identification of *Salmonella* from Meat, Poultry, Pasteurized Egg, Carcass, and Environmental Sponges. Effective Date: 05 June 2023.
3. ISO 6579. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of *Salmonella* spp.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
6. Neogen Installation Qualification (IQ) / Operational Qualification (OQ) Protocols and Instructions for Neogen Molecular Detection System. Contact your Neogen Food Safety representative to obtain a copy of this document.
7. ISO 6887-1:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – part 1: General rules for the preparation of the initial suspension and decimal dilutions.
8. ISO 6887-2:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – part 2: Specific rules for the preparation of meat and meat products.

### Explanation of Symbols

[info.neogen.com/symbols](http://info.neogen.com/symbols)

# Neogen Food Safety

## Neogen Corporation

620 Leshar Place  
Lansing, MI 48912 USA  
Neogen.com

## Neogen Europe Ltd.

The Dairy School  
Auchincruive  
Ayr, KA6 5HU  
Scotland, UK

## Neogen Ireland, Ltd.

Bray Business Park, Bray  
Co. Wicklow  
A98YV29, Ireland

## Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA  
www.neogen.com

© Neogen Corporation 2024. All rights reserved.

Neogen is a trademark of Neogen Corporation.

All other trademarks are the property of their respective companies.

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### Referencias:

1. US Food and Drug Administration Bacteriological Analytical Manual. Capítulo 5: *Salmonella*, versión de septiembre de 2023.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 4.14. Isolation and Identification of *Salmonella* from Meat, Poultry, Pasteurized Egg, Carcass, and Environmental Sponges. Fecha de vigencia: 5 de junio de 2023.
3. ISO 6579. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of *Salmonella* spp.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
6. Neogen Installation Qualification (IQ) / Operational Qualification (OQ) Protocols and Instructions for Neogen Molecular Detection System. Póngase en contacto con su representante de Seguridad Alimentaria de Neogen para obtener una copia de este documento.
7. ISO 6887-1:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – part 1: General rules for the preparation of the initial suspension and decimal dilutions.
8. ISO 6887-2:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – part 2: Specific rules for the preparation of meat and meat products.

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### Referências:

1. US Food and Drug Administration Bacteriological Analytical Manual. Chapter 5: *Salmonella*, versão de setembro de 2023.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 4.14. Isolation and Identification of *Salmonella* from Meat, Poultry, Pasteurized Egg, Carcass, and Environmental Sponges. Data de vigência: 05 de junho de 2023.
3. ISO 6579. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of *Salmonella* spp.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
6. Neogen Installation Qualification (IQ) / Operational Qualification (OQ) Protocols and Instructions for Neogen Molecular Detection System. Entre em contato com seu representante da Neogen Food Safety para obter uma cópia deste documento.
7. ISO 6887-1:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – part 1: General rules for the preparation of the initial suspension and decimal dilutions.
8. ISO 6887-2:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – part 2: Specific rules for the preparation of meat and meat products.

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[info.neogen.com/symbols](http://info.neogen.com/symbols)

# Neogen Food Safety

## Neogen Corporation

620 Leshar Place  
Lansing, MI 48912 EUA  
Neogen.com

## Neogen Europe Ltd.

The Dairy School  
Auchincruive  
Ayr, KA6 5HU  
Scotland, Reino Unido

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© Neogen Corporation 2024. Todos os direitos reservados.
Neogen é marca comercial da Neogen Corporation.