U.S. English

*Product Numbers:
9850, 9851*

colitagTM

Water Test System

Colitag™

Product Numbers: 9850, 9851

Intended Use

Colitag is a selective and differential medium for the simultaneous determination of the presence or absence or enumeration of *E. coli* and other total coliforms in water. Colitag is for use in the Environmental Protection Agency's survey and monitoring programs under the Safe Drinking Water Act*. The test system is designed for use by personnel familiar with appropriate aseptic techniques for the identification of total coliforms and *E. coli*. Training, which is available through NEOGEN®, is recommended for those without a basic knowledge of microbiology.

Storage Requirements

Store kit components at 2–30°C (preferably 4–7°C), away from light and moisture to assure full shelf life. The product's expiration date is printed on the box and each individual packet.

Materials Provided

1. Colitag test kits, 100 mL format (NEOGEN items 9850 or 9851)

Materials Required but not Provided

1. Incubator set to 35 ± 0.5°C
2. Long wavelength (365 nanometer) UV lamp
3. Sterile sample containers large enough to contain 100 mL water samples (NEOGEN items 9855, 9855N, 9856, 9859, 9859P)
4. UV protective safety glasses

Optional

1. Colitag Comparator (NEOGEN item 9857)
2. Circulating, covered water bath capable of maintaining 35°C or 44.5°C for prewarming samples
3. Control bacteria strains (see recommended quality control strains below)

Precautions

1. Do not use kit components beyond expiration date.
2. Use of incubation times or temperatures other than those specified may cause erroneous results.
3. Good microbiological laboratory practices should be used where appropriate.

Presence/Absence (P/A) Test Methodology for Colitag Water Test System

For a video demonstration of the instructions, please refer to the training section of the product [page](#).

P/A Methodology

1. Aseptically add the contents of one Colitag packet to a 100 mL water sample in a sterile vessel.
2. Cap vessel and shake until dissolved. Some particles may not fully dissolve but this will not affect test performance.
NOTE: If reading at 16–24 hours: If results are needed in less than 24 hours and the sample is not already at 33–38°C, prewarm the sample by placing vessel in a 35°C water bath for 20 minutes or alternatively, 44.5°C for 7–10 minutes, then go to step 3.
NOTE: If reading at 24–48 hours: Go directly to step 3.
3. Incubate sample at 35 ± 0.5°C or according to local regulation temperature for the remainder of the desired incubation time; between 16–48 hours.
4. Read results according to the result interpretation listed below.

Interpretation of P/A Results

For determination of results, post-incubation samples should be compared to one or more of the following standards: a blank, negative and positive bacterial controls, and/or a Colitag Comparator.

1. Visually check each sample vessel for yellow color. If yellow color is observed in the sample, coliform bacteria are present.
2. If for any reason the chromogenic response is ambiguous (notable turbidity or color cannot be discerned), the sample may be inverted two to four times to mix and incubated for another 2–4 hours. Under these conditions, a lightly positive result will continue to intensify in color.
3. If samples are yellow, examine for blue fluorescence using a 6-watt, 365 nm UV lamp in a darkened environment within five inches of the sample. If bright blue fluorescence is observed, *E. coli* are present.

Enumeration Methodology

1. Aseptically add the contents of one Colitag packet to a 100 mL water sample in a sterile vessel.
2. Cap vessel top and shake until dissolved. Some particles may not fully dissolve — which will not affect test performance.
3. For enumeration, Colitag may be used in different formats:
 - If using multi-well tray, follow multi-well tray product insert directions.
 - If using the multiple tube technique, use one of the tube formats described in section 9221 of the Standard Methods for the Examination of Water and Wastewater.
4. Incubate the sample at 35 ± 0.5°C or according to local regulation temperature for the desired incubation time.
5. Read results according to the sample interpretation method appropriate for enumerative MPN (Most Probable Number) procedure.

Enumeration Results

Colitag may be used in Most Probable Number (MPN) testing formats to determine total coliforms and *E. coli*. This includes MPN formats described in Standard Methods 9221.

1. If using multi-well tray, find results from the MPN table provided in the multi-well tray product insert.
2. If using multiple tube format, find results from the Standard Methods 9221C, Table 9221:II or Table 9221:III.

Result Interpretation

Using Blanks, Controls or Colitag Comparator	Result
Equal or lighter yellow than blank, non-coliform control or lighter yellow than Comparator	Negative for total coliforms and <i>E. coli</i>
Darker yellow than blank or negative control, or equal or darker yellow than Comparator	Positive for total coliforms
Fluorescence greater than blank, non- <i>E. coli</i> coliform control, or equal or greater than Comparator	Positive for <i>E. coli</i>

Recommended Quality Control Strains

Category	Control Strains			Expected Results	
	Organism	ATCC#	NCTC#	ONPG (Yellow)	MUG (Fluorescent)
Non-Coliforms	<i>Salmonella enterica</i> or	700720	-	-	-
	<i>Pseudomonas aeruginosa</i>	10145	10332	-	-
Total Coliforms	<i>Klebsiella pneumoniae</i> or	13883	9633	+	-
	<i>Citrobacter freundii</i>	8090	9750	+	-
<i>Escherichia coli</i>	<i>Escherichia coli</i>	11775	9001	+	+
		25922	12241	+	+

Disposal

Dispose (autoclave, bleach, etc.) Colitag test samples in accordance with all applicable local, state, and federal regulations.

Customer Service

NEOGEN Customer and Technical Services can be contacted through NEOGEN.com and product training is available by request.

SDS Information Available

Safety data sheets are available for all test kits at NEOGEN.com or by calling 800.234.5333 or 517.372.9200.

Terms and conditions

NEOGEN's full terms and conditions are available [online](#).

Warranty

NEOGEN makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

* Reference U.S. EPA: Federal Register/Vol. 74, No. 216/ Nov. 2009 regarding Modified Colitag™ and Federal Register/Vol. 86, No. 100/May 2021 regarding Modified Colitag™ Version 2.0.

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