**Read instructions carefully before starting test**

# Reveal<sup>®</sup> 2.0

## for *Salmonella*

### INTENDED USE

The Reveal<sup>®</sup> 2.0 for *Salmonella* test system provides for rapid recovery of *Salmonella* in food, animal feed, and environmental samples allowing detection and presumptive identification of the test organism, generally within 24 hours.

In an AOAC Research Institute Performance Tested Method<sup>SM</sup> study, the Reveal 2.0 for *Salmonella* test system was found to be an effective procedure for detection of *Salmonella enterica* in the following sample types: chicken carcass rinse, raw ground turkey, raw ground beef, hot dogs, raw shrimp, ready-to-eat meat products, dry pet food, ice cream, fresh spinach, cantaloupe, peanut butter, stainless steel surfaces, and sprout irrigation water.

Antibodies used in the Reveal 2.0 for *Salmonella* test are reactive with *Salmonella enterica* serovars of somatic groups A–E. This includes the most common serovars from both food and non-food sources.<sup>1,2</sup> Serovars that are not within somatic groups A–E will not be detected.

### ASSAY PRINCIPLES

The system utilizes Revive medium, which provides *Salmonella* with readily available nutrients and other factors required for its recovery from a stressed or injured condition. After a brief enrichment in Revive, selective enrichment in Rappaport-Vassiliadis (RV) then favors *Salmonella* growth to levels detectable by the Reveal test device.

For certain matrices, the initial nonselective enrichment step may be bypassed and the sample introduced directly into the RV medium.

Another enrichment option enables the Reveal test device to detect *Salmonella* spp. in matrices that are more difficult due to viscosity, pigment, or containing an extremely high level of nonviable cells. M-broth was selected for post enrichment as it is established in several standard procedures recommended in enzyme assays used in the food industry.

A portion (200 µL) of the enrichment culture is placed into the sample cup. The test device is placed into the sample and allowed to develop at ambient temperature for 15 minutes. The sample is wicked through a reagent zone that contains specific anti-*Salmonella* antibodies conjugated to colloidal gold particles. If antigens are present in the sample, they will bind to the gold-conjugated antibodies. This antigen-antibody complex then leaves the reagent zone and travels through the nitrocellulose membrane which contains a zone of anti-*Salmonella* antibodies. The immune complex with gold conjugate is captured and aggregates in this zone, thus displaying a visible line. The remainder of the sample continues to migrate to the end of the membrane where it is deposited into a waste reservoir.

The reagent zone also contains gold conjugate of a proprietary antigen, which is eluted by the sample regardless of the presence of *Salmonella* antigen. The gold-conjugated control indicator migrates through the membrane to the negative control capture zone (antibody to the proprietary antigen), where it is captured and aggregated to form a visible line. Regardless of the presence or absence of the *Salmonella* antigen, the control line will form in the control zone, ensuring the test is working properly.

## INTENDED USER

The Reveal 2.0 for *Salmonella* test system is designed for use by personnel familiar with the appropriate aseptic techniques for the isolation and identification of *Salmonella*. Training, which is available through Neogen, is recommended for those without a basic knowledge of microbiology.

## MATERIALS PROVIDED

1. 20 Reveal 2.0 for *Salmonella* test devices
2. 20 disposable transfer pipettes
3. 20 reaction cups

## TEST SYSTEMS REQUIRED BUT NOT PROVIDED

**NOTE:** The Reveal 2.0 for *Salmonella* device requires the use of certain enrichment media and supplies. These can be purchased as part of a test system or individually.

1. **Reveal 2.0 one-step system** (contains devices, 1x RV selective medium and sampling bags for 20 tests) (Neogen item 9802)
2. **Reveal 2.0 complete system with RV** (contains devices, Revive, 2x RV selective medium and sampling bags for 20 tests) (Neogen item 9803)
3. **Reveal 2.0 environmental system** (contains devices, Revive, 2x RV, sampling bags and environmental sample kits for 20 tests) (Neogen item 9804)
4. **Reveal 2.0 complete system with RV and M-broth** (contains devices, Revive, 2x RV selective medium, M-broth and sampling bags for 20 tests) (Neogen item 9805)

## MEDIA REQUIRED BUT NOT PROVIDED (DEPENDING ON SAMPLE)

**NOTE:** Media can be purchased individually or as part of a test system.

1. 20 foil pouches of Revive medium (Neogen item 9705)
2. 20 foil pouches 2x concentrated Rappaport-Vassiliadis (2x RV) (Neogen item 9715)
3. 20 foil pouches 1x Rappaport-Vassiliadis (1x RV) (Neogen item 9729)
4. 20 bottles M-Broth (Neogen item 9722)
5. Lactose Broth (Neogen item 7141)

## MATERIALS REQUIRED BUT NOT PROVIDED

1. Scale capable of weighing a minimum of 25 g (Neogen item 9427)
2. Timer (Neogen items 9452, 9426)
3. Sample cup rack (Neogen item 9475)
4. Sterile deionized or distilled water
5. Incubator capable of maintaining  $42 \pm 1^\circ\text{C}$  (Neogen item 9735)
6. Incubator capable of maintaining  $36 \pm 1^\circ\text{C}$  (Neogen item 9735)
7. Graduated cup
8. Stomacher-type bags

## OPTIONAL MATERIALS

1. AccuScan® reader
2. Stomacher 400 machine or equivalent

## STORAGE

Store Reveal devices and dry media at room temperature (15–30°C) when not in use. Store M-broth between 2–8°C.

## TESTING OF DIFFERENT COMMODITIES

Please contact a Neogen representative for an updated list of the commodities for which the test has been validated, or for validation studies for specific commodities.

## PRECAUTIONS

1. Do not use expired media or test devices.
2. Use rehydrated media within 6 hours of preparation.
3. Use of incubation times and temperatures other than those specified may lead to erroneous results.
4. Sample bags must be closed loosely to allow air exchange during incubation, which is vital for sufficient and detectable growth.
5. Sterile water should be warmed to 42°C when rehydrating Revive and  $36 \pm 1^\circ\text{C}$  when rehydrating 2x RV and 1x RV.
6. Test enriched cultures within 6 hours post-incubation.
7. Good microbiological laboratory practices should be used.

## SAMPLE PREPARATION AND ENRICHMENT

**Sample enrichment 1 (Revive/RV) – for most foods and environmental samples. Validated for<sup>3</sup>: cantaloupe, ice cream, ready-to-eat meat products, stainless steel surfaces. For use with Neogen items 9803 and 9804.**

1. Transfer contents of 1 foil pouch of unitized Revive (Neogen item 9705) or 7.2 g of bulk Revive (Neogen item 9708) into a Stomacher-type bag. Using the graduated cup, add 200 mL of sterile-purified water prewarmed to 42°C. Grasp the bag tightly 2–3 inches from the top and mix vigorously until dissolved.
2. Place 25 g of food sample (sample must be at room temperature) or environmental sponge sample into the Stomacher-type bag containing the Revive medium. Grasp bag tightly at top and knead sample until mixed. Shake bag vigorously using a side-to-side motion to ensure complete mixing.  
**ALTERNATIVE:** Place bag in Stomacher apparatus and mix for **30 seconds** at normal speed.
3. Loosely close the bag and place in suitable rack or support. Incubate  $36 \pm 1^\circ\text{C}$  for **4 hours**.
4. Reconstitute 2x RV in a Stomacher-type bag by adding 1 foil pouch of unitized 2x RV concentrated (Neogen item 9715) or 10.6 g of bulk RV (Neogen item 9716). Using the cup provided, add 200 mL of sterile purified water prewarmed to  $36 \pm 1^\circ\text{C}$  to the bag. Mix vigorously until dissolved. Hold prepared 2x RV at 42°C until use.
5. Remove the Revive sample bag from the  $36 \pm 1^\circ\text{C}$  incubator and place in a suitable rack or support.
6. Add the 200 mL of selective 2x RV enrichment prewarmed to 42°C to the entire Revive culture (200 mL) in the sample bag. Grasp tightly 2–3 inches from top and mix gently using a side-to-side motion.
7. Loosely close bag and place in a suitable rack or support. Incubate at  $42 \pm 1^\circ\text{C}$  for **16–24 hours**.

**Sample enrichment 2A (Direct RV) – for select foods including raw ground beef, raw ground poultry, raw seafood, and poultry rinses. Validated for<sup>3</sup>: chicken carcass rinse, raw ground turkey, raw ground beef, raw shrimp. For use with Neogen item 9802.**

1. Transfer contents of 1 foil pouch of unitized 1x RV (Neogen item 9729) or 5.3 g of bulk RV (Neogen item 9716) into a Stomacher-type bag. Using the graduated cup, add 200 mL of sterile purified water prewarmed to 42°C. Grasp the bag tightly 2–3 inches from the top and mix vigorously until dissolved.
2. Place 25 g of sample (sample must be at room temperature) or 30 mL poultry carcass rinse into the Stomacher-type bag containing the 1x RV medium. Grasp bag tightly at top and knead sample until mixed. Shake bag vigorously using a side-to-side motion to ensure complete mixing.  
**ALTERNATIVE:** Place bag in Stomacher apparatus and mix for **30 seconds** at normal speed.
3. Loosely close the bag and place in suitable rack or support. Incubate  $42 \pm 1^\circ\text{C}$  for **20–24 hours**.

**Sample enrichment 2B (Direct RV) – for deli meats and other ready-to-eat meat products. Validated for<sup>3</sup>: hot dogs. For use with Neogen item 9802.**

1. Add 68.9 g of bulk RV (Neogen item 9716) into a large Stomacher-type bag.
2. Add 2.6 L of sterile-purified water prewarmed to 42°C. Grasp the bag tightly 2–3 inches from the top and mix vigorously until dissolved.
3. Place 325 g of sample (sample must be at room temperature) into the Stomacher-type bag containing the 1x RV medium. Grasp bag tightly at the top and knead sample until mixed. Shake bag vigorously using a side-to-side motion to ensure complete mixing.  
**ALTERNATIVE:** Place bag in Stomacher apparatus and mix for **30 seconds** at normal speed.
4. Loosely close the bag and place in a suitable rack or support. Incubate at  $42 \pm 1^\circ\text{C}$  for **20–24 hours**.

**Sample enrichment 2C (Direct RV) – Validated for sprout irrigation water<sup>3</sup>. For use with Neogen item 9706.**

1. Add 19.9 g bulk RV (Neogen item 9716) to a large Stomacher-type bag. Add 375 mL sterile-purified water prewarmed to 42°C. Grasp the bag tightly 2–3 inches from the top and mix vigorously until dissolved.
2. Add 375 mL sprout irrigation water to the Stomacher-type bag containing the 2x RV medium. Grasp bag tightly at the top and knead sample until mixed. Shake bag vigorously using a side-to-side motion to ensure complete mixing.  
**ALTERNATIVE:** Place bag in Stomacher apparatus and mix for **30 seconds** at normal speed.
3. Loosely close the bag and place in a suitable rack or support. Incubate at  $42 \pm 1^\circ\text{C}$  for **20–24 hours**.

**Sample enrichment 3 (Revive/RV – M-Broth) – for matrices that may contain residual antigen, highly viscous or highly pigmented material. Validated for<sup>3</sup>: dry pet food. For use with Neogen item 9805.**

1. Transfer contents of 1 foil pouch of unitized Revive (Neogen item 9705) or 7.2 g of bulk Revive (Neogen item 9708) into a Stomacher-type bag. Using the graduated cup, add 200 mL of sterile-purified water prewarmed to 42°C. Grasp the bag tightly 2–3 inches from the top and mix vigorously until dissolved.
2. Place 25 g of sample (sample must be at room temperature) into the Stomacher-type bag containing the Revive medium. Grasp bag tightly at top and knead sample until mixed. Shake bag vigorously using a side-to-side motion to ensure complete mixing. **ALTERNATIVE:** Place bag in Stomacher apparatus and mix for **30 seconds** at normal speed.
3. Loosely close the bag and place in suitable rack or support. Incubate at 36 ± 1°C for **4 hours**.
4. Reconstitute 2x RV in a Stomacher bag by adding 1 foil pouch of unitized 2x RV concentrated (Neogen item 9715) or 10.6 g of bulk RV (Neogen item 9716). Using the cup provided, add 200 mL of sterile-purified water prewarmed to 36 ± 1°C to the bag. Mix vigorously until dissolved. Hold prepared 2x RV at 42°C until use.
5. Remove the Revive sample bag from the 36 ± 1°C incubator and place in a suitable rack or support.
6. Add the 200 mL of selective 2x RV enrichment prewarmed to 42°C to the entire Revive culture (200 mL) in the sample bag. Grasp tightly 2–3 inches from top and mix gently using a side-to-side motion.
7. Loosely close bag and place in a suitable rack or support. Incubate at 42 ± 1°C for **16–24 hours**.
8. Rehydrate 1 bottle of M-broth (Neogen item 9722) by adding 10 mL of sterile-purified water prewarmed to 36 ± 1°C. Cap tightly and shake to dissolve medium.
9. Remove 1 mL aliquot from the Revive/RV sample and transfer into the bottle of rehydrated M-broth.
10. Incubate at 36 ± 1°C for **6 hours**.

**Sample enrichment 4 (RV – M-broth) – for leafy greens. Validated for<sup>3</sup>: fresh spinach. For use with Neogen item 9802 and M Broth.**

1. Transfer contents of 1 foil pouch of unitized 1x RV (Neogen item 9729) or 5.3 g of bulk RV (Neogen item 9716) into a Stomacher-type bag. Using the graduated cup, add 200 mL of sterile-purified water prewarmed to 42°C. Grasp the bag tightly 2–3 inches from the top and mix vigorously until dissolved.
2. Place 25 g of sample (sample must be at room temperature) into the Stomacher-type bag containing the 1x RV medium. Grasp bag tightly at top and knead sample until mixed. Shake bag vigorously using a side-to-side motion to ensure complete mixing. **ALTERNATIVE:** Place bag in Stomacher apparatus and mix for **30 seconds** at normal speed.
3. Loosely close the bag and place in suitable rack or support. Incubate at 42 ± 1°C for **20–24 hours**.
4. Rehydrate 1 bottle of M-broth (Neogen item 9722) by adding 10 mL of sterile-purified water prewarmed to 36 ± 1°C. Cap tightly and shake to dissolve medium.
5. Remove 1 mL aliquot from the RV sample and transfer into the bottle of rehydrated M-broth.
6. Incubate at 36 ± 1°C for **6 hours**.

**Sample enrichment 5 (Lactose - RV) – Validated for peanut butter<sup>3</sup>. For use with Neogen item 9802 and Lactose Broth**

1. Homogenize 25 g sample in 225 mL Lactose Broth. Incubate at 36 ± 1°C for **22–24 hours**.
2. Remove 0.1 mL aliquot of Lactose Broth culture and add 10 mL RV broth. Incubate at 42 ± 1°C for **22–24 hours**.

**REVEAL TEST PROCEDURE**

1. Remove enriched sample from designated incubator. Mix sample well and transfer 200 µL or 8 drops to the Reveal sample cup.
2. Remove the required number of Reveal 2.0 for *Salmonella* test devices from container.
3. Place Reveal device with sample arrows facing down into sample cup containing sample and incubate at ambient temperature for **15 minutes**.
4. Observe and record the test results immediately (within 60 seconds) after **15 minutes**.

## INTERPRETATION OF RESULTS

### Visual interpretation

1. Line in both control and test zone in 15 minutes is considered positive.
2. Line in control zone only at 15 minutes is considered negative.
3. If no line appears in the control zone, the test is considered invalid and the sample should be retested with another device.

### Optional electronic interpretation

Place the Reveal test device into the AccuScan® reader according to the instrument instructions and follow the on-screen prompts to interpret and record the device's test result.

**NOTE:** The Reveal device forms a distinct line in the test region when *Salmonella* are present; the intensity of the line may be variable based on serovar and/or concentration. If a distinct visible line forms, regardless of intensity, the sample should be considered positive. Migration of any blue pigment from the RV enrichment to the device does not impact results.

### CONFIRMATION

Neogen recommends that presumptively positive Reveal enrichment cultures be verified by plating a sample of the culture following methods described in BAM<sup>4</sup> or MLG<sup>5</sup>.

**NOTE:** Testing a replicate subsample using BAM or MLG procedures may not yield the same result. Bacteria are not evenly distributed with a lot of material hence the subsample may not contain the organism at the onset.

### DISPOSAL

Dispose (autoclave, bleach, etc.) of Reveal test devices, pipettes, and media in accordance with all applicable local, state and federal regulations.

### CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

### SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at [foodsafety.neogen.com](http://foodsafety.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

### TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit [www.neogen.com/en/terms-and-conditions](http://www.neogen.com/en/terms-and-conditions).

### WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### PRODUCTS AVAILABLE

#### REFERENCES

<sup>1</sup> Galanis, E. et al. (2006) *Emerg. Infect. Dis.* **12**, 381-388

<sup>2</sup> [www.who.int/salmsur](http://www.who.int/salmsur)

<sup>3</sup> In an AOAC Research Institute Performance Tested Method<sup>SM</sup> study

<sup>4</sup> <http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/JCM070149>

<sup>5</sup> [http://www.fsis.usda.gov/Science/Microbiological\\_Lab\\_Guidebook/index.asp](http://www.fsis.usda.gov/Science/Microbiological_Lab_Guidebook/index.asp)

- 9706 **Reveal 2.0 for *Salmonella*** – a one-step assay for detection of *Salmonella*. Kit contains 20 devices
- 9705 **Revive for *Salmonella*** – patented resuscitation medium to speed up enrichment cultivation for *Salmonella*. Contains 20 foil pouches and 20 Stomacher-type bags for 20 tests
- 9708 **Bulk Revive for *Salmonella*** – 500 g
- 9729 **1x RV selective media** – Rappaport-Vassiliadis broth 1x media. Kit contains prefilled foil pouches for 20 tests and 20 Stomacher-type bags
- 9715 **2x RV selective media** – Rappaport-Vassiliadis broth 2x concentrate media. Kit contains prefilled foil pouches for 20 tests and 20 Stomacher-type bags. To be used in combination with Revive
- 9716 **Bulk RV media** – 500 g
- 9724 ***Salmonella* confirmation kit** – contains *Salmonella* magnetic beads, buffer packets, and Rambach chromogenic agar and supplement for 20 tests
- 9724B **Bulk *Salmonella* confirmation kit** – contains *Salmonella* magnetic beads for 50 tests
- 7649 **Rambach agar for *Salmonella*** – chromogenic agar for *Salmonella* (4 x 250 mL)
- 9427 **Electronic digital scale**
- 9735 **Incubator**
- 9426 **Timer** (3 channel)
- 9415 **Pipettes** – sterile 10 mL serological (500)
- 9722 **M-broth** – media kit contains prefilled bottles for 20 tests
- 7141 **Lactose Broth**
- 9475 **Rack for Reveal sample cups**
- 9802 **Reveal 2.0 one step system** – contains devices, 1x RV selective medium and sampling bags for 20 tests
- 9803 **Reveal 2.0 complete system with RV** – contains devices, Revive, 2x RV selective medium and sampling bags for 20 tests
- 9804 **Reveal 2.0 environmental system** – contains devices, Revive, 2x RV, sampling bags and environmental sample kits for 20 tests
- 9805 **Reveal 2.0 complete system with RV and M-broth** – contains devices, Revive, 2x RV selective medium, M-broth and sampling bags for 20 kits

#### North America

##### Neogen Headquarters

800/234-5333 (USA/Canada)  
 foodsafety@neogen.com  
 foodsafety.neogen.com/en

#### Europe, Middle East and Africa

##### Neogen Europe

+ 44 (0) 1292 525 600  
 info\_uk@neogeneurope.com  
 foodsafety.neogen.com/uk

#### Mexico

##### Neogen Latinoamerica

+52 (55) 5254-8235  
 informacion@neogenlac.com  
 foodsafety.neogen.com/sp

#### Brazil

##### Neogen do Brasil

+55 19 3935.3727  
 info@neogendobrasil.com.br  
 foodsafety.neogen.com/pt

#### China

##### Neogen Bio-Scientific Technology

+86 21 6271 7013  
 info@neogenchina.com.cn  
 www.neogenchina.com.cn

#### India

##### Neogen Food and Animal Security

+91 484 2306598, 2301582  
 info@neogenindia.com  
 www.neogenindia.com

©Neogen Corporation, 2019. Neogen, AccuScan and Reveal are registered trademarks of Neogen Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.

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