Read instructions carefully before starting test

Reveal[®] 2.0

for
Group D1 *Salmonella*
(including *Salmonella* Enteritidis)

INTENDED USE

Reveal[®] 2.0 for Group D1 *Salmonella* provides for the rapid recovery and detection of Group D1 *Salmonella* spp. serovars (including *S. Enteritidis*) in layer or broiler house environmental samples, raw shell eggs, chicken carcass rinses, and poultry feed. Following plating for confirmation of positive Reveal results, identification of presumptive isolates as *S. Enteritidis* is accomplished using conventional biochemical and serological procedures. This test is approved by the National Poultry Improvement Plan (NPIP) for testing of poultry house environmental samples.

In an AOAC Performance Tested MethodsSM validation study, the Reveal 2.0 for Group D1 *Salmonella* test was found to be an effective method for detection of group D1 *Salmonella*, including *S. Enteritidis*, in raw shell eggs, chicken carcass rinse, and poultry feed.

ASSAY PRINCIPLES

The sample is wicked through a reagent zone in the test device, which contains specific antibodies (anti-SE) conjugated to colloidal gold particles. If antigens are present in the sample, they will bind to the gold-conjugated antibodies. This antigen-antibody complex then leaves the reagent zone and travels through the membrane, which contains a zone of anti-SE antibody. The immune complex with gold conjugate is captured and aggregates in this zone, thus displaying a visible line. The remainder of the sample continues to migrate to the end of the membrane where it will eventually be deposited into a waste reservoir.

The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of Group D1 *Salmonella* spp., ensuring the strip is working properly.

INTENDED USER

The Reveal 2.0 test is designed for use by personnel familiar with the appropriate aseptic techniques for the isolation and identification of *Salmonella* spp. Training, which is available through Neogen, is recommended for those without a basic knowledge of microbiology.

MATERIALS PROVIDED

Reveal 2.0 for Group D1 *Salmonella* (including *S. Enteritidis*)

1. 20 Reveal 2.0 for Group D1 *Salmonella* (including *S. Enteritidis*) test devices

MATERIALS REQUIRED BUT NOT PROVIDED

1. Tetrathionate Broth Base, Hajna (TT Hajna) (Neogen item 7740A) or equivalent
2. Rappaport-Vassiliadis Medium Semisolid Modified (MSRV) (Neogen item 7511A) or equivalent
3. Tryptic Soy Broth (TSB) (Neogen Item 7164A) or equivalent
4. Buffered Peptone Water (BPW) (Neogen item 7418A) or equivalent
5. Reveal 2.0 for *S. Enteritidis* Media (Neogen Item 9865A)
6. Novobiocin (Neogen item 7985) or equivalent
7. Pipette, 20–200 μL (Neogen item 9276) or equivalent
8. Filtered pipette tips, 100 μL (Neogen item 9389) or equivalent
9. Tubes, 13 x 100 mm glass test tubes (box of 250) (Neogen item 9438A) or equivalent
10. Caps for 13 x 100 mm glass tubes (bag of 250) (Neogen item 9509A) or equivalent
11. Potassium iodide (Sigma no. 221945) or equivalent
12. Iodine (Sigma no. 207772) or equivalent
13. Incubator capable of maintaining $37 \pm 1^\circ\text{C}$ (Neogen item 9735) or equivalent
14. Incubator capable of maintaining $42 \pm 0.2^\circ\text{C}$ (Neogen item 9735) or equivalent
15. Scale capable of weighing a minimum of 25 g (Neogen item 9427) or equivalent
16. Stomacher-type bags or equivalent
17. Stomacher or equivalent
18. Pipette pump for serological pipettes, 10 mL (Neogen item 9277) or equivalent
19. Pipettes, sterile serological, 10 mL, pack of 100 (Neogen item 9415) or equivalent
20. Timer (3-channel) (Neogen item 9426) or equivalent
21. Test tube rack, 13 mm (Neogen item 9440) or equivalent

OPTIONAL MATERIALS (FOR ADDITIONAL CONFIRMATION)

1. Reveal *Salmonella* Enteritidis confirmation bead kit (Neogen item 9723)
2. Phosphate Buffered Saline (PBS) (Neogen item 8425) or equivalent

STORAGE

Store Reveal 2.0 devices and dry media at room temperature (15–30°C).

PRECAUTIONS

Good microbiological laboratory practices should be used where appropriate

MEDIA PREPARATION

I. TT Hajna (for environmental samples)

1. Follow manufacturer's instructions for media preparation.

II. TSB (for pooled raw shell eggs)

1. Follow manufacturer's instructions for media preparation.

III. BPW (for chicken carcass rinse)

1. Follow manufacturer's instructions for media preparation.

IV. Reveal 2.0 for *S. Enteritidis* Media (for chicken carcass rinse and poultry feed)

1. Rehydrate 22.5 g of Reveal 2.0 for *S. Enteritidis* Media with 1 L sterile water prewarmed to 37°C. Do not autoclave.

V. MSRV media with novobiocin supplement preparation

1. Suspend 31.6 g of the MSRV medium in 1 L of purified water. Heat with frequent agitation and boil for **1 minute** to completely dissolve the medium. **NOTE:** Do not autoclave.
2. Cool medium to 45–50°C and aseptically add 10 mL (20 mg) of novobiocin supplement.

VI. MSRV tube preparation

1. Add 2.5 mL of MSRV medium with novobiocin supplement (cooled to 45–50°C) to a 13 x 100 test tube.
2. Cap and allow to cool undisturbed for **1 hour**.

SAMPLE PREPARATION AND ENRICHMENT

A. Environmental Samples (drag swabs, boot swabs, chick papers, fluff, rinse, etc.)

NOTE: Not part of the AOAC PTM study or approved matrices.

1. Place environmental sample into a Stomacher-type bag.
2. Aseptically add 100 mL of TT Hajna (prepared earlier) to a Stomacher-type bag containing the environmental sample. Shake the bag vigorously in a side-to-side motion at least 10 times to mix.
3. Incubate the samples at $36 \pm 1^\circ\text{C}$ for **24 ± 2 hours**.
4. Pipette 100 μL of the enriched sample into the bottom of the center of the MSR/V tube. Slowly dispense the sample as the pipette tip is drawn up through the media, making sure to dispense completely before the tip is removed from the top of the media. Cap the tube.
5. Incubate the samples in the capped MSR/V tube(s) at $42 \pm 0.2^\circ\text{C}$ for **24 ± 2 hours**.

B. Pooled Raw Shell Eggs

1. Prepare a pool of 20 shell eggs as described in Bacteriological Analytical Manual, chapter 5 [ref. 2].
2. Add 2 L TSB and mix thoroughly.
3. Incubate the samples at $36 \pm 1^\circ\text{C}$ for **24 ± 2 hours**.
4. Pipette 100 μL of the enriched sample into the bottom of the center of the MSR/V tube. Slowly dispense the sample as the pipette tip is drawn up through the media, making sure to dispense completely before the tip is removed from the top of the media. Cap the tube.
5. Incubate the samples in the capped MSR/V tube(s) at $42 \pm 0.2^\circ\text{C}$ for **24 ± 2 hours**.

C. Chicken Carcass Rinse

1. Prepare the carcass rinse as described in Microbiology Laboratory Guidebook, chapter 4.08 [ref. 3].
2. To 30 mL rinse, add 30 mL Reveal 2.0 for SE Media containing 5 mg/L novobiocin.
3. Incubate the samples at $36 \pm 1^\circ\text{C}$ for **24 ± 2 hours**.
4. Pipette 100 μL of the enriched sample into the bottom of the center of the MSR/V tube. Slowly dispense the sample as the pipette tip is drawn up through the media, making sure to dispense completely before the tip is removed from the top of the media. Cap the tube.
5. Incubate the samples in the capped MSR/V tube(s) at $42 \pm 0.2^\circ\text{C}$ for **24 ± 2 hours**.

D. Poultry Feed

1. Combine 25 g poultry feed with 225 mL Reveal 2.0 for *S. Enteritidis* Media containing 5 mg/L novobiocin.
2. Incubate the samples at $36 \pm 1^\circ\text{C}$ for **24 ± 2 hours**.
3. Pipette 100 μL of the enriched sample into the bottom of the center of the MSR/V tube. Slowly dispense the sample as the pipette tip is drawn up through the media, making sure to dispense completely before the tip is removed from the top of the media. Cap the tube.
4. Incubate the samples in the capped MSR/V tube(s) at $42 \pm 0.2^\circ\text{C}$ for **24 ± 2 hours**.

TESTING WITH REVEAL

1. Remove test tubes containing MSR/V from the $42 \pm 0.2^\circ\text{C}$ incubator.
2. Holding a Reveal 2.0 device at the top with the arrows facing down, place the device into the middle of the tube and incubate at room temperature for **10 minutes**.
NOTE: It is important to ensure the strip is inserted into the inoculated portion of the tube.
3. Observe and record test results immediately (within 30 seconds) after **10 minutes**.

INTERPRETATION OF RESULTS

Line in C (top line) but no line in T zone **Negative**

Line in C zone and T zone in 10 minutes..... **Positive**

If no line appears in the C zone regardless if there is a line in the T zone..... **Invalid**

NOTE: The Reveal 2.0 device forms a distinct line in the test zone in the presence of Group D1 *Salmonella* spp. If a distinct visual line forms in the test window, regardless of intensity, the sample should be considered positive.

CONFIRMATION

For environmental samples, Neogen recommends that presumptively positive enrichment cultures be verified by plating a sample of the enriched culture onto media described in the NPIP procedure [1]. Cultural confirmation can be continued from the MSR/V tubes using a loop for plating onto selective media plates, XLT-4 and BGN. Continue with identification using standard biochemical and serological procedures (MLG 4.06). For raw shell eggs and poultry feed, follow BAM chapter 5 [ref. 2]. For chicken carcass rinse, follow MLG chapter 4.08 [ref. 3]. Immunomagnetic separation using serotype D1-specific magnetic particles (Neogen item 9723) may be used prior to plating to enrich the target organism population.

OPTIONAL: SE IMS BEAD ISOLATION PROTOCOL

To use this protocol, the Reveal *Salmonella* Enteritidis confirmation beads kit (Neogen item 9723) is needed.

Preparation of wash buffer

1. Empty powder from buffer packet into a sterile 1 L bottle.
2. Add 1 L of sterile water to the powder. Cap bottle and shake to dissolve powder.
3. Store buffer at 2–8°C for no longer than 6 months.
4. To avoid buffer contamination, do not work directly from the stock buffer bottle. Place the approximate amount of buffer needed into a smaller sterile container and dispose of any unused portion.

Immunomagnetic separation

1. Remove the magnetic plate from the magnetic particle concentrator (MPC) and load the necessary number (one per sample) of 1.5 mL Eppendorf tubes into the rack.
2. Resuspend the anti-SE beads by vortexing until the pellet in the bottom of the tube disappears. Pipette 20 µL of bead solution into each Eppendorf tube.
3. Add 1 mL PBS to the 20 µL bead suspension.
4. Using a sterile disposable inoculation loop, transfer 10 µL of media from the MSR/V tube to the PBS containing SE specific beads.
5. Follow the instructions in the Reveal *Salmonella* Enteritidis confirmation bead kit insert beginning at step 4 in the **Immunomagnetic separation** section.

References

[REF 1] http://www.aphis.usda.gov/animal_health/animal_dis_spec/poultry/downloads/salmonella_procedures.pdf

[REF 2] US FDA (2014) Bacteriological Analytical Manual, chapter 5
<http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm070149.htm>

[REF 3] USDA-FSIS (2014) Microbiology Laboratory Guidebook, chapter 4.08.
<http://www.fsis.usda.gov/wps/wcm/connect/700c05fe-06a2-492a-a6e1-3357f701f52/MLG-4.pdf?MOD=AJPERES>

DISPOSAL

Dispose (autoclave, bleach, etc.) of Reveal test devices, pipettes, and media in accordance with all applicable local, state and federal regulations.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/en/terms-and-conditions.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

PRODUCTS AVAILABLE

- 9809 **Reveal 2.0 for Group D1 *Salmonella* (including *Salmonella* Enteritidis) – contains 20 devices for 20 tests**
- 7164A **Tryptic Soy Broth, 500 g**
- 9865A **Reveal 2.0 for SE Media, 500 g**
- 7740A **Tetrathionate Broth Base, 500 g**
- 7511A **Rappaport-Vassiliadis Medium Semisolid Modified, 500 g**
- 7418A **Buffered Peptone Water, 500g**
- 7985 **Novobiocin Supplement, 10 mL vial – 10 vials/pkg (lyophilized)**
- 9706 **Reveal 2.0 for *Salmonella* – contains 20 devices for 20 tests**
- 9802 **Reveal 2.0 One-Step *Salmonella* Complete System – contains devices, RV selective medium and sampling bags for 20 tests**
- 9803 **Reveal 2.0 for *Salmonella* Complete System with RV – contains devices, Revive, RV selective medium and sampling bags for 20 tests**
- 9805 **Reveal 2.0 for *Salmonella* Complete System with M-Broth and RV – contains devices, Revive, RV selective medium, M-Broth and sampling bags for 20 tests**
- 9804 **Reveal 2.0 for *Salmonella* Environmental Sampling System – contains devices, Revive, RV selective medium and environmental sampling kits for 20 tests**
- 9724 ***Salmonella* confirmation kit – contains *Salmonella* magnetic beads, buffer packets and Rambach chromogenic agar and supplement for 20 tests**
- 9723 **Bulk *Salmonella* Enteritidis (SE) confirmation beads – contains SE magnetic beads for 50 tests**

North America

Neogen Headquarters
800/234-5333 (USA/Canada)
foodsafety@neogen.com
foodsafety.neogen.com/en

Europe, Middle East and Africa

Neogen Europe
+ 44 (0) 1292 525 600
info_uk@neogeneurope.com
foodsafety.neogen.com/uk

Mexico

Neogen Latinoamerica
+52 (55) 5254-8235
informacion@neogenlac.com
foodsafety.neogen.com/sp

Brazil

Neogen do Brasil
+55 19 3935.3727
info@neogendobrasil.com.br
foodsafety.neogen.com/pt

China

Neogen Bio-Scientific Technology
+86 21 6271 7013
info@neogenchina.com.cn
www.neogenchina.com.cn

India

Neogen Food and Animal Security
+91 484 2306598, 2301582
info@neogenindia.com
www.neogenindia.com

©Neogen Corporation, 2019. Neogen, AccuScan, and Reveal are registered trademarks of Neogen Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.

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