Read instructions carefully before starting test

INTENDED USE

The ANSR[®] *Listeria* Right Now™ assay is intended to be used as a fast and easy screening procedure for the detection of *Listeria* spp. on visibly clean surfaces post-cleaning/pre-production without the need for enrichment. This environmental test can be used on stainless steel, sealed concrete, ceramic tile, rubber and plastic. It has been validated by the AOAC for stainless steel, sealed concrete, ceramic tile, rubber and plastic surfaces. It is not designed for soiled samples or food matrices testing.

ASSAY PRINCIPLES

ANSR *Listeria* Right Now is an isothermal, amplified nucleic acid assay. The method is based on nicking enzyme amplification reaction (NEAR™) technology preceded by the reverse transcription of 23S ribosomal RNA. Target complementary DNA is amplified through a mechanism of polymerization from the ends of nicks created in double-stranded DNA by the action of a specific endonuclease. Amplified target sequences are detected in real time using fluorescent molecular beacon probes.

Listeria Right Now is a novel test for the detection of *Listeria* spp. from environmental surfaces without enrichment. The entire collected contents of the swab are subjected to sample processing. After expression of the swab in the lysis buffer, a 2-stage lysis reaction is performed, first at 37 ± 2°C for 10 minutes, then at 80 ± 2°C for 20 minutes. Next, a portion of the lysed sample is transferred to a strip tube containing lyophilized ANSR reagents. The tubes are sealed and incubated at 56 ± 2°C on the ANSR reader. Results are generated by the reader and displayed in the ANSR software within 18 minutes. Each tube of ANSR reagents contains an internal positive control, ensuring that the reagents are functioning properly.

INTENDED USER

The ANSR *Listeria* Right Now test is designed for use by personnel with appropriate training in microbiology. Training in the use of the ANSR test system is available through NEOGEN.

MATERIALS PROVIDED

1. 96 microbiological swabs in 1 mL Lethen Broth
2. 12 strips of 8 cluster tubes, 1.2 mL
3. 12 strips of 8 reaction tubes, 200 μ L, containing lyophilized ANSR *Listeria* Right Now reagents in 2 sealed foil pouches with desiccant pack
4. 12 strips of 8 permanent caps for reaction tubes
5. 2 bottles lysis reagent suspension buffer, 60 mL
6. 6 vials containing lyophilized lysis reagents
7. 1 kit insert

EQUIPMENT REQUIRED

1. ANSR reader (NEOGEN item 9828)
2. Computer and software for connection to ANSR reader (NEOGEN item 9832)
3. 1 dual heater block with aluminum block inserts for 1.2 mL cluster tubes, $80 \pm 2^\circ\text{C}$ and $37 \pm 2^\circ\text{C}$ (Neogen items 9836-DUAL48, 9829-48) **OR** two single heater blocks with aluminum block inserts for 1.2 mL cluster tubes, $80 \pm 2^\circ\text{C}$ and $37 \pm 2^\circ\text{C}$ (NEOGEN items 9836-48D, 9829-48)
4. Pipettor, 100–1000 μ L (NEOGEN item 9463)
5. Pipette tip rack, 100–1000 μ L, sterile, 96 tips (NEOGEN item 9487)
6. Pipettor, 10–100 μ L, 8-channel (NEOGEN item 9388)
7. Pipette tips, 100 μ L, sterile, filtered, 96 tips (NEOGEN item 9389)
8. Vortex mixer, adjustable speed (NEOGEN item 9494)
9. 3 thermometers, traceable (NEOGEN item 9518)
10. Timer, 3-channel (NEOGEN item 9426)
11. Optional-for-use heater block with 0.2 mL reaction tube aluminum block insert, $56 \pm 2^\circ\text{C}$ (NEOGEN items 9836-48D, 9829-48)
12. Webcam (NEOGEN item WEBCAM)
13. ANSR Ethernet cable (NEOGEN item 9835)
14. 10 mL pipette pump (NEOGEN item 9277)
15. 4.7 inch pipette tip, 100–1,000 μ L, sterile, 96 tips (NEOGEN item 9487-LONG)
16. Pipettes, sterile serological 10 mL (NEOGEN item 8686)
17. 40-slot, 20 mm test tube rack, autoclavable (NEOGEN item 9553)
18. Pipettor, 20–200 mL (NEOGEN item 9276)

STORAGE

Store ANSR reagents at 2–8°C. After removing reaction tubes from the foil pouch, promptly re-seal the pouch. Leave the desiccant pack in the pouch at all times.

PRECAUTIONS

1. Use good microbiology laboratory practices.
2. **Dispose of used pipette tips in a covered container containing a fresh solution of 10% bleach. The 10% bleach solution should be made fresh each day. Use 9 parts water with 1 part commercially available, household-strength bleach to make the 10% bleach solution. Undiluted stock solutions of bleach should be used within 30 days after opening.**

3. **Discard bleach solution and tips as regular waste at the end of each day.**
4. *Listeria monocytogenes* is a known hazard to pregnant women and immunocompromised individuals. Consult with your facility safety director for specific instructions.
5. Do not use reagents beyond expiration date.
6. Remove reaction tubes from the foil pouch just before use and keep covered until heating process begins. Reseal the pouch containing the remaining reaction tubes to **avoid prolonged exposure to light**. More than 15 minutes of total exposure time may lead to erroneous results.
7. Do **NOT**, under any circumstance, remove caps from reaction tubes **AFTER** the assay has been started. This is essential in order to prevent accidental contamination of the environment with amplification products.
8. Exercise care in all pipetting steps to avoid cross-contamination of samples.
9. Complete all assay steps in sequence, avoiding delays between steps.
10. Tap reaction tubes on bench top to make sure lyophilized reagents are at bottom of the tube prior to adding lysed sample.
11. Excess sanitizer must be rinsed from sample site prior to sample collection.

ENVIRONMENTAL SAMPLE COLLECTION

1. Collect environmental sample using the provided microbiological swabs that are pre-moistened with Lethen Broth after labeling swab tube with sample ID.
2. Sample each surface by swabbing horizontally, vertically and diagonally enough times in each direction to provide full coverage of the test area (up to an area of 4" x 4"). 1" x 1" surface areas were validated in an AOAC study.
3. Pour off the remaining Lethen Broth from the sample collection tubes.
4. After swabbing, place the swab back into the original tube without Lethen Broth.
5. Keep the swab sample collection tubes at 2–8°C.
6. The sample should be tested the same day after collection, not to exceed 8 hours after collection. If the collected sample is stored at 2–8°C within 4 hours of collection, testing can take place within 24 hours from sample collection.
7. If a negative control is desired, use a swab that has not been used in sample collection and process as described by removing the Lethen Broth and adding the lysis buffer.

LYSIS REAGENT SOLUTION PREPARATION

Reconstitute 1 vial of lyophilized lysis reagents with 18 mL of lysis reagent suspension buffer by adding the buffer to the reagent vial. Swirl gently to mix.

NOTE: 1 vial of lysis reagents is enough for approximately 18 samples. Prepared lysis reagent solution can be stored at 2–8°C for up to 30 days.

ANSR TEST PROCEDURE

Prior to starting the assay:

1. Preheat one lysis heater block to $80 \pm 2^\circ\text{C}$. Preheat the second lysis heater block to $37 \pm 2^\circ\text{C}$. If using the optional single heater, preheat to $56 \pm 2^\circ\text{C}$. Use the thermometer for the temperature reading.
2. Remove the foil pouch containing the reaction tubes from the refrigerator and **allow them to warm at room temperature for 15 minutes**. To avoid excess light exposure, leave reaction tubes in foil pouch until they are needed. **NOTE:** The lysis reagent solution should be at room temperature for 15 minutes before use.
3. Connect the ANSR reader to the computer via USB or Ethernet and turn the computer on.
4. Turn on the ANSR reader. The reader will preheat to $56 \pm 2^\circ\text{C}$.
5. Start the ANSR software and click the Connect button. Input sample IDs, lot number and user information.
NOTE: For instructions on using the reader and software, see the user guide that came with the ANSR reader.

SAMPLE LYSIS

1. Label the appropriate number of cluster tubes to match the sample ID on the tubes with the swabs.
2. Pipette 1 mL of prepared lysis reagent solution into the tube containing the sampled swab **without** the Lethen Broth. Avoid touching the sides of the swab tube with the exception of the pipette tips (Neogen item 9487).
NOTE: Return the lysis reagent solution to the refrigerator after use (within 1 hour).
3. Secure tube cap to the tube. Vortex vigorously for 5–7 seconds, at medium/high speed (2500 RPM), to express the sample from the swab into the lysis solution.

ASSAY PROCEDURE

OPTIONAL: Pre-pierced plate sealing film (BIO RAD Item 3600040) can be placed over the 1.2mL cluster tubes if desired.

1. Using the long pipette tips (Neogen item 9487-LONG), pipette 500 μL of the lysis reagent solution from the labeled swab tube into a labeled cluster tube. Quickly place the used swab back into the tube after pipetting the lysis reagent solution, and secure the cap.
NOTE: Save remaining 500 μL for additional analysis as needed.
2. After dispensing the pipette tip, take a towel that has been sprayed with 10% bleach (sodium hypochlorite) along the shaft of the pipettor to ensure there is no cross-contamination.
3. Incubate the cluster tube(s) at $37 \pm 2^\circ\text{C}$ for **10 minutes**.
4. Immediately transfer the cluster tube(s) to the $80 \pm 2^\circ\text{C}$ heater block and incubate for **20 minutes**.
NOTE: The $80 \pm 2^\circ\text{C}$ incubation time may be extended to a total of **60 minutes** for the purpose of managing staggered assay start times.
5. **For 3–5 minutes** before the end of the lysis step, preheat the ANSR reagents to $56 \pm 2^\circ\text{C}$ by placing the reaction tubes in the ANSR reader. Optional: A separate heat block can be used. It should be heated to $56 \pm 2^\circ\text{C}$.
NOTE: The strip of reaction tubes may be cut to provide the number of tubes needed. Keep all unused tubes in the sealed foil pack. Ensure the pellet in the reaction tube(s) is at the bottom by tapping the tubes gently on the bench top.

6. After the completion of the 20–60 minute lysis incubation, remove and discard the caps from the reaction tube(s) in the ANSR reader.
IMPORTANT: Proceed with steps 7–9 without delay. **The transfer of the sample from the lysis tubes to the reaction tubes should be completed within 1 minute.**
7. Using an 8-channel pipette and 100 μ L filter tips, carefully transfer 50 μ L from the top third of the lysed sample(s) in the cluster tube(s) to the reaction tube(s). Debris may accumulate at the bottom of the lysis tube(s) that will interfere with assay performance. Avoid transfer of debris by aspirating from the top third of the lysis tube(s). Do not prime the pipette tips and do not mix before aspirating. **Place the provided permanent cap(s) on the reaction tube(s).**
NOTE: Lysed sample may be transferred from the same cluster tube a maximum of 3 times.
8. Remove the strips(s) of tubes from the reader (or $56 \pm 2^\circ\text{C}$ heat block if one was used) and vortex briefly (about 2 seconds) at medium/high speed (2500 RPM). Quickly visually check each reaction tube to assure that no bubbles are on the bottom or in the middle of the sample. A quick tap of the tubes should release any bubbles from the bottom or middle to the top. Then place into the reader without delay. Close the reader's lid.
NOTE: The reader will not provide accurate results if the lid is open. Keep the lid closed at all times while the assay is running. Contamination may occur if the permanent caps are not placed on the reaction tubes and/or if the permanent caps are removed.
9. Click **START** in the ANSR software to begin the **18 minute** assay.
10. Results will be displayed as positive, negative or invalid once the assay is finished. If the result displays an invalid, the test must be repeated. A single repeat test can be run starting at step 6. Alternatively, start over from step 1 if lysis from step 6 has been at $80 \pm 2^\circ\text{C}$ for more than 60 minutes.

INTERPRETATION OF RESULTS

The ANSR software will indicate the test results as positive or negative for the presence of *Listeria* spp. in the environmental swab sample. In addition, the real-time fluorescence curve generated from the assay can be viewed.

FOLLOW-UP TESTING

If desired, follow-up testing of areas producing positive *Listeria* Right Now results can be conducted by re-sampling the areas, prior to any additional cleansing, and by performing one or more of the following analyses: additional *Listeria* Right Now testing, enrichment-based rapid *Listeria* spp. testing, and/or conventional microbiological testing.

DISPOSAL

DO **NOT** AUTOCLAVE.

Do **NOT** remove the permanent caps, for **ANY** reason, from the ANSR reaction tubes once the assay has started, even when disposing of them. Reaction tubes can be disposed of as non-biohazardous waste. It is recommended that they be placed in a sealable container or plastic bag with 10% bleach and immediately disposed of to protect against accidental opening.

NEOGEN ANSR Molecular Diagnostic for Foodborne Pathogen Detection – limited use label license

SYTO® 82

SYTO® 82 contained within this product is provided under an intellectual property license from Life Technologies Corporation, Eugene, OR, and may be used for *in vitro* detection and analysis of (i) food, feeds and beverages, including nutraceuticals, (ii) ingredients for food, feeds and beverages, (iii) process samples from food, feed and beverage preparation, distribution and delivery, and (iv) water from any source for human consumption, all for the purpose of safety and quality assurance. The buyer must not sell or otherwise transfer this product or its components for any other use, including but not limited to: human *in vitro*, veterinary, identity or paternity testing, forensics, or *in vivo* detection of nucleic acid sequences in living beings, or cells. For information on purchasing a license for SYTO 82 for purposes other than food, beverage and water safety, and quality assurance, contact Life Technologies Corporation at outlicensing@lifetech.com.

Molecular beacon probes

One or more molecular beacon probes contained within this product is sold under license from PHRI Properties and may be used under PHRI Properties patent rights only for tests on food products, feeds, beverages and water.

NEAR™ Technology

This product utilizes the patented NEAR isothermal technology and is sold under license from Ionian Technologies, San Diego, CA, and may be used under Ionian Technologies patent rights only for tests on food, beverage and water safety.

CUSTOMER SERVICE

NEOGEN Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all NEOGEN test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of NEOGEN's test kits, on NEOGEN's website at neogen.com, or by calling NEOGEN at 800.234.5333 or 517.372.9200.

TERMS AND CONDITIONS

For NEOGEN's full terms and conditions, please visit: neogen.com/terms-and-conditions/.

WARRANTY

NEOGEN Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

References

[1] USDA-FSIS (2015) Microbiology Laboratory Guidebook, chapter 41.04
<http://www.fsis.usda.gov/wps/wcm/connect/0273bc3d-2363-45b3-befb-1190c25f3c8b/MLG-41.pdf?MOD=AJPERES>

PRODUCTS AVAILABLE

9465	Environmental sampling sponges (with Dey Engley (DE) neutralizing broth) – pack of 100
36002	Sampling sponges, pre-moistened, HiCap – pack of 100
9822	ANSR for <i>E. coli</i> O157:H7 – up to 96 tests
9824	ANSR for <i>Listeria monocytogenes</i> – up to 96 tests
9870	ANSR for <i>Salmonella</i> – up to 96 tests
9871	ANSR for <i>Listeria</i> – up to 96 tests
9872	ANSR for <i>Campylobacter</i> – up to 96 tests
9873	ANSR for <i>Listeria Right Now</i> [™] – up to 96 tests
9837	Complete ANSR Equipment Kit – Reader, computer, pipette set, dry bath heaters and inserts, thermometers, vortex, timer, rack

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