Read instructions carefully before starting test

GeneQuence[®] **for *Salmonella***

**DNA hybridization test
for detection of *Salmonella* spp.**

AOAC Official Method
No. 2007.02 (24-hour enrichment)

INTENDED USE

This test is intended to qualitatively detect *Salmonella* spp. in food samples, and should be used by personnel with appropriate laboratory training in microbiology.

TESTING DIFFERENT COMMODITIES

GeneQuence for *Salmonella* (48-hour enrichment) is an effective tool for a variety of commodities. Through evaluation by the AOAC Research Institute, this test has been shown to detect *Salmonella* in the following food products: raw whole egg, dried whole egg, nonfat dry milk, chocolate, refrigerated dough, soy flour, egg noodles, cheese powder, cake mix, walnuts, food dye, black pepper, dried fruit, frozen fruit, mushrooms, surimi, raw ground beef, raw pork, raw fish, raw turkey, raw chicken, beef franks, dry pet food, and peanut butter.

GeneQuence for *Salmonella* (24-hour enrichment) has been shown to detect *Salmonella* spp. in the following food products: raw poultry, pasteurized egg products (dried, liquid and frozen), milk chocolate, and dry pet food.

BIOLOGICAL PRINCIPLES OF THE TEST

This DNA hybridization test employs *Salmonella*-specific DNA probes directly labeled with horseradish peroxidase, and a colorimetric endpoint to detect as little as one *Salmonella* CFU in a 25 g food sample when the specified enrichment protocols are used. A sample is considered negative for the presence of *Salmonella* spp. if the absorbance value (A_{450}) obtained is less than 0.10. A sample is considered positive for the presence of *Salmonella* spp. if the absorbance value obtained is greater than or equal to 0.10. Positive samples should be confirmed by standard culture procedures, followed by biochemical and serological identification.

MATERIALS PROVIDED

1. 2 bottles of lysis reagent concentrate (labeled **1a**)
2. 1 bottle of 12 mL lysis reagent buffer (labeled **1b**)
3. 1 bottle of 18 mL hybridization solution (labeled **2**)
4. 1 bottle of 6 mL *Salmonella* probe solution (labeled **3**)
5. 1 bottle of 50 mL wash solution 20X concentrate (labeled **4**)
6. 1 bottle of 15 mL substrate chromogen solution (labeled **5**)
7. 1 bottle of 5 mL stop solution (labeled **6**)
8. 1 bottle of 5 mL positive control (labeled **+**)
9. 1 bottle of 5 mL negative control (labeled **-**)
10. 1 96-well microwell plate in divisible strips
11. Hybridization/probe mixture chart

MATERIALS REQUIRED BUT NOT PROVIDED

24-hour and 48-hour Enrichments

Equipment

1. Angled thermometer for digital heat block (Neogen item 9382)
2. Apparatus for plate washing with vacuum source (Neogen item 9348 (manual), Neogen items 6718 or 6718A (automated))
3. Automated plate wash assembly: vacuum tubing assembly (Neogen item 9387), side arm flask assembly (Neogen item 6822), vacuum pump (Neogen item 6824)
4. Carboy with spigot (2) (Neogen item 9349)
5. Digital heat block (Neogen item 9386)
6. Heated block cover (Neogen item 9383)
7. Incubators at $35 \pm 1^\circ\text{C}$ and $42 \pm 0.2^\circ\text{C}$ (Neogen item 9735)
8. Microwell reader to read at 450 nm with discrimination of 0.01 absorbance unit (Neogen items 6704, 6704A or 9303)
9. Modular heating block titer plate (Neogen item 9384)

Optional equipment

1. Dry heater base incubator (Neogen item 9412)
2. 8-channel pipettor and tips to dispense 50–300 μL (Neogen items 9385, 9407)
3. Heater block at $37 \pm 1^\circ\text{C}$ (Neogen item 9411)

Lab supplies and consumables

1. Micropipette and tips to dispense 20–200 μL volumes (Neogen items 9276, 9407)
2. Micropipette and tips to dispense 200–1000 μL volumes (Neogen items 9342, 9343)
3. Minute timer (Neogen item 9426)
4. Test tubes, 12 x 75 mm borosilicate glass, and rack (Neogen items 9439, 9440)
5. Stomacher bags (100) (Neogen item 9736)
6. 50 mL culture tubes for sample pre-enrichment (Neogen item 9381)
7. 1000 mL graduated cylinder (Neogen item 9345)
8. 500 mL wash bottle (Neogen item 9366)
9. 10 mL pipette pump (Neogen item 9277)
10. 10 mL sterile serological pipettes, 0.1 mL graduations (Neogen item 9415)
11. Eppendorf reservoir (Neogen item 4018)
12. Centrifuge tube racks (Neogen item 9479)

Media

1. Lactose Broth and/or other pre-enrichment media (Neogen item 7141)
2. Rappaport-Vassiliadis (RV) Broth (ref. 2, Neogen item 7512)
3. Gram Negative (GN) Broth (ref. 3, Neogen item 7218)

24-hour enrichment only

1. Buffered Peptone Water (ref. 1, Neogen item 7418)
2. Brain Heart Infusion (BHI) Broth (ref. 1, Neogen item 7116)

48-hour enrichment only

1. Tetrathionate (TT) Broth (ref. 2, Neogen item 7241)

Starter kits

1. GeneQuence Starter Kit, includes all items (Neogen item 9268)
2. GeneQuence Starter Kit, without water bath (Neogen item 9268a)
3. GeneQuence Consumables Starter Kit (Neogen item 9266)
4. GeneQuence Small Wash Starter Kit (Neogen item 9267)

Other equipment needed

24-hour and 48-hour enrichments

1. Homogenizer (e.g., Stomacher)
2. Small orbital platform shaker capable of 150 rpm
3. Sterile cotton swabs and absorbent paper
4. 2 mL disposable graduated pipettes
5. Diagnostic reagents as necessary for culture confirmation of positive DNA hybridization assay (ref 2.)

STORAGE REQUIREMENTS

1. Store microwell plate, positive control, negative control, and reagents 1a, 2, 3, and 5 at 2–8°C.
2. Other reagents may be stored at 2–25°C.
3. The entire kit may be stored at 2–8°C.
4. Once reagent 1a has been reconstituted, it must be stored at -20°C and is stable for 60 days in this form.
5. Refrigerated kit reagents should be equilibrated to room temperature before use, but should not be left unrefrigerated for long periods of time.

PRECAUTIONS

1. Reagents are for laboratory use only.
2. Stop solution contains 4.0 N sulfuric acid. Hybridization solution contains formamide. Avoid contact with skin and mucous membranes. Refer to the Material Safety Data Sheet available from Neogen for more information.
3. Reagents from different kit lots should not be interchanged. GeneQuence for *Salmonella* reagents are not interchangeable with other GeneQuence assay reagents.
4. Reagents should not be used beyond their expiration dates.
5. Enriched cultures should be handled/disposed of as potentially infectious material.
6. The preferred method for disposal of contaminated materials, including cultures, pipettes, etc., is autoclaving.
7. Items that cannot be autoclaved should be decontaminated by treatment with a disinfectant solution and rinsing with water.

- This test should be performed in a normal laboratory environment with respect to humidity, lighting, etc. Steps requiring room temperature incubation should be performed at 18–30°C.

SAMPLE PREPARATION AND ENRICHMENT

NOTE: Food samples should be obtained and handled according to standard practices appropriate to analysis for *Salmonella* (ref. 4). **All enrichment media should be prewarmed to room temperature before inoculation.** Samples must be tested when the final GN or BHI incubation period is complete. Do not hold cultures for testing with GeneQuence at a later time.

A. 48-HOUR enrichment procedures for select foods

• For raw foods with a high microbial load

- Homogenize food sample in 9 volumes Lactose Broth. Incubate **22 ± 2 hours** at 35 ± 1°C.
- Remove pre-enriched culture from incubation and mix well. Transfer 1 mL of pre-enrichment culture to 10 mL Tetrathionate (TT) Broth. Transfer another 0.1 mL of pre-enrichment culture to 10 mL Rappaport-Vassiliadis (RV) Broth. Incubate TT and RV cultures **16–18 hours** at 42 ± 0.2°C.
- Remove TT and RV cultures from incubation and mix well. Transfer 1 mL of TT culture to 10 mL Gram Negative (GN) Broth. Transfer 1 mL of RV culture to 10 mL GN. Incubate both GN cultures **6 hours** at 35 ± 1°C. Save TT and RV cultures for possible confirmation.
- Perform GeneQuence for *Salmonella* using a pooled sample containing 0.2 mL of each GN culture. Save GN cultures for possible confirmation.

• For all foods other than raw foods or foods with a high microbial load

- Prepare pre-enrichment culture in accordance with BAM/AOAC guidelines (ref. 4). Incubate **22 ± 2 hours** at 35 ± 1°C.
- Remove pre-enrichment culture from incubation and mix well. Transfer 1 mL of pre-enrichment culture to 10 mL Tetrathionate (TT) Broth. Transfer another 0.1 mL of pre-enrichment culture to 10 mL Rappaport-Vassiliadis (RV) broth. Incubate TT culture **6 hours** at 35 ± 1°C and RV culture **6 hours** at 42 ± 0.2°C.
- Remove TT and RV cultures from incubation and mix well. Transfer 1 mL TT culture to 10 mL Gram Negative (GN) Broth. Transfer 1 mL of RV culture to a second 10 mL of GN. Incubate GN cultures **16–18 hours** at 35 ± 1°C. Save TT and RV cultures for possible confirmation.
- Perform GeneQuence for *Salmonella* using a pooled sample containing 0.2 mL of each GN culture. Save GN cultures for possible confirmation.

B. 24-HOUR enrichment procedures for select foods

• For dried, liquid and liquid frozen pasteurized eggs

- Homogenize sample in 9 volumes Buffered Peptone Water. **NOTE:** The standard sample size as specified by USDA-FSIS is 100 g (ref. 3). For dried egg samples, add BPW to the sample gradually and mix to produce a homogeneous suspension free of lumps. Incubate **18–20 hours** at 35 ± 1°C.
- Remove BPW culture from incubation and mix well. Transfer 10 mL BPW culture to 10 mL **pre-warmed** to 35–37°C **double-strength** Gram Negative Broth (GN) in a 50 mL tube. Incubate **6–7 hours** at 35 ± 1°C.
- Remove GN culture from incubation and mix. Perform GeneQuence for *Salmonella* using 0.4 mL of GN culture. Save GN culture at 4°C for possible confirmation.

• For dry pet food

1. Prepare pre-enrichment culture in accordance with BAM/AOAC guidelines (ref. 2). Incubate **18–20 hours** at $35 \pm 1^\circ\text{C}$.
2. Remove pre-enrichment culture from incubation and mix well. Transfer 10 mL pre-enrichment culture to 10 mL **pre-warmed** to $35\text{--}37^\circ\text{C}$ **double-strength** Gram Negative broth (GN) in a 50 mL tube. Incubate **6–7 hours** at $35 \pm 1^\circ\text{C}$.
3. Remove GN culture from incubation and mix. Perform GeneQuence for *Salmonella* using 0.4 mL of GN culture. Save GN culture at 4°C for possible confirmation.

• For milk chocolate

1. Prepare pre-enrichment culture in accordance with BAM/AOAC guidelines (ref. 2). Incubate **22–24 hours** at $35 \pm 1^\circ\text{C}$.
2. Remove pre-enrichment culture from incubation and mix well. Transfer 10 mL pre-enrichment culture to 10 mL **pre-warmed** to $35\text{--}37^\circ\text{C}$ **double-strength** Brain Heart Infusion Broth (BHI) in a 50 mL tube. Incubate **6–7 hours** at $35 \pm 1^\circ\text{C}$.
3. Remove BHI culture from incubation and mix. Perform GeneQuence for *Salmonella* using 0.4 mL of BHI culture. Save BHI culture at 4°C for possible confirmation.

• For raw meats and poultry

1. Homogenize sample in 9 volumes Rappaport-Vassiliadis Broth (RV). Incubate **18–20 hours** at $42 \pm 0.2^\circ\text{C}$.
2. Remove RV culture from incubation and mix well. Transfer 1 mL RV culture to 10 mL **pre-warmed** to $35\text{--}37^\circ\text{C}$ Gram Negative Broth (GN) in a 50 mL tube. Incubate **6–7 hours** at $35 \pm 1^\circ\text{C}$.
3. Remove GN culture from incubation and mix. Perform GeneQuence *Salmonella* using 0.4 mL of GN culture. Save GN culture at 4°C for possible confirmation.

PRIOR TO STARTING THE TEST

1. Turn on the water bath or heater block and adjust to $65 \pm 1^\circ\text{C}$. Fill a water bath to a level of approximately 1.5 inches, or fill the heater block wells about 1/3 with deionized water.
2. Prepare pretreatment reagent by adding 6 mL of lysis reagent buffer (bottle **1b**) directly to one bottle of lysis reagent concentrate (bottle **1a**). Dissolve contents by gently swirling, and place the bottle on ice. **NOTE:** Lysis reagent is stable in the reconstituted form for 60 days when stored at -20°C . To thaw, place bottle at room temperature. When thawed, place on ice (freeze-thawing does not adversely affect the reconstituted reagent).
3. For each sample to be tested, label a 12 x 75 mm glass test tube with the appropriate sample designation and place in a rack. Include tubes for 1 positive control and 1 negative control per experimental run.
4. Prepare the wash solution by mixing entire contents of wash solution concentrate (bottle 4) with 950 mL of distilled or deionized water (if washing manually with a 500 mL wash bottle, use 25 mL of concentrate with 475 mL of water). Fill the buffer reservoir of the plate-washing device (see manufacturer's instructions for set-up and use). **NOTE:** Wash solution can be stored in a closed bottle at room temperature for up to 60 days.
5. Without touching the bottoms of the wells, place the appropriate number of coated microwells in the plate frame, filling the frame left to right and front to back in rows of 8. Include wells for the reagent blank, negative control and positive control. **NOTE:** If the last row has fewer than 8 wells and a plate-washing device is used, fill the last row with colored wells (colored wells are available free of charge from Neogen).

TEST PROCEDURE

NOTE: This test is adaptable to automation. Contact Neogen for further information.

- 1a. **If using the 48-hour enrichment**, combine the two GN cultures for each sample by adding 0.2 mL of each to the appropriate tube (i.e., 0.2 mL of each of the two GN cultures for each sample pooled into one tube).
- 1b. **If using the 24-hour enrichment**, add 0.4 mL of each sample to the appropriate tubes.
2. Mix the positive and negative control solutions by inverting the bottles several times. Add 0.4 mL of each control to the appropriate tube.
3. Add 0.1 mL of reconstituted lysis reagent (bottle **1a**) to each tube. Mix by gently shaking the rack of tubes by hand for **5 seconds**. The resulting solution should be blue. If any tubes are not blue, check for proper reagent addition. Incubate the rack of tubes in the 65°C water bath or heater block for **5 minutes**.
4. Prepare a 4:1 hybridization/probe mixture by mixing hybridization solution (bottle **2**) and probe solution (bottle **3**) in a plastic or glass container. For mixture guidelines, refer to the mixing chart provided with this test kit or use the formula below. (N = number of test samples + controls)
Volume hybridization solution (bottle 2) = $[(N \times 0.1) + 1.6]$ mL
Volume probe solution (bottle 3) = $[(N \times 0.025) + 0.4]$ mL
5. Remove the tubes from the heat source and shake briefly. Transfer 0.150 mL of each lysed sample, including the controls, to a designated microwell. The first well should be reserved for the reagent blank and receives no sample. The second well should be used for the negative control, and the third for the positive control.
6. Vigorously mix the hybridization/probe solution prepared earlier. Add 0.125 mL to each microwell, except the reagent blank microwell.
7. Shake the plate on an orbital shaker at 150 rpm for **2 minutes** at room temperature.
NOTE: If using a multichannel pipettor, the contents of the wells may be mixed by drawing up and dispensing back into the wells at least 5 times.
8. Incubate the plate at 45 ± 1°C for **60 minutes** on the covered heater block.
9. Wash the wells using one of the following procedures:
 - a) **Plate-washing device:** Wash 5 times at room temperature. For each wash, process one 8-well strip at a time by aspirating the liquid, filling the wells, and then proceeding to the next strip. After the last wash, aspirate the liquid from the wells, then remove any residual liquid by inverting the plate and tapping it onto absorbent paper. Hold the plate by gently squeezing on the sides of the frame to keep the strips in place. Be sure **not** to touch the bottom or sides of wells with platemasher.
 - b) **Manual:** Thoroughly wash the microwells at least 5 times using a plastic wash bottle filled with the wash solution prepared earlier. Wash by emptying the wells into a suitable container, tapping the inverted plate on absorbent paper, completely filling the wells with wash solution, and vigorously shaking out the contents. **NOTES:** All air bubbles must be removed before proceeding to the next step. If air bubbles remain, repeat vigorous striking of the wells onto an absorbent paper towel on a flat surface until eliminated.
10. Add 0.150 mL of substrate chromogen solution (bottle **5**) to each microwell, including the reagent blank microwell. Incubate the plate at room temperature for **20 minutes**.

11. Add 0.05 mL of stop solution (bottle **6**) to each microwell, including the blank microwell.
12. Gently tap the side of the frame a few times to ensure mixing.
13. Read absorbance at 450 nm using a plate or strip reader according to the manufacturer's instructions. Blank using the first microwell that contains the mixture of substrate chromogen and stop solution (do not blank with air). If the plate or strip reader doesn't automatically blank against the first microwell, calculate the values for the negative control, positive control, and all other sample wells, by subtracting the absorbance value of the reagent blank (well **A1**) from the absorbance values registered in all other wells.

INTERPRETATION OF RESULTS

Control values: The absorbance value for the negative control must be ≤ 0.15 , and absorbance value for the positive control must be ≥ 1.00 . If either control falls out of the acceptable range, the test is invalid and should be repeated.

Negative criterion: Tests producing absorbance values < 0.10 are negative for the presence of *Salmonella* spp. in the test samples.

Positive criterion: Tests producing absorbance values ≥ 0.10 are positive for the presence of *Salmonella* spp. in the test samples. A positive test result should be confirmed by standard culture procedures.

RECOMMENDED CONFIRMATION PROCEDURE

Neogen recommends that positive results be verified by plating the sample's GN or BHI culture and/or selective enrichment cultures onto media as described in BAM or FSIS procedures (refs. 2, 3), depending upon the sample type.

LIMITATIONS

The DNA probes used in this kit are not reactive with serovars of the species *Salmonella bongori*.

REFERENCES

1. Atlas, R.M. 1997. Handbook of Microbiological Media, 2nd ed., CRC Press, Boca Raton, FL, pg. 606.
2. US FDA. 2003. Bacteriological Analytical Manual online, chapter 5. <http://www.cfsan.fda.gov/~ebam/bam-5.html>
3. USDA-FSIS. 2004. Microbiology Laboratory Guidebook, chapter 4. http://www.fsis.usda.gov/PDF/MLG_4_03.pdf
4. US FDA. 2003. Bacteriological Analytical Manual online, chapter 1. <http://www.cfsan.fda.gov/~ebam/bam-5.html>

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's Food Safety test kits, at foodsafety.neogen.com, or by calling Neogen at 800.234.5333 or 517.372.9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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