*Read instructions carefully before starting test*

# NeoSEEK™ STEC

## **INTENDED USE**

NeoSeek™ STEC is a multiplex molecular method for detection and identification of select Shiga toxin-producing *Escherichia coli* (STEC), including enterohemorrhagic *E. coli* (EHEC), in enriched samples. The method does not require culture confirmation. NeoSeek STEC can also be used for identification of STEC in pure culture. NeoSeek STEC is a single-source, service-based method performed by Neogen Genomics in Lincoln, Nebraska, United States.

## **ASSAY PRINCIPLES**

The NeoSeek STEC method utilizes the Agena BioScience MassARRAY® platform. Using over 80 targets associated with STEC virulence and serogroup, the test generates a molecular profile of select STEC, including serogroups O26, O45, O103, O111, O121, O145 and O157. Due to the high discriminating power of the multiplex assay and the targeting of multiple alleles of virulence determinants — including *stx*, *eae*, *nle*, and *espA*, among others — as well as STEC-associated single nucleotide polymorphisms, the test can detect the presence of potentially virulent STEC, including EHEC, in mixed populations. Targeting genes for somatic (O group) and flagellar (H) antigens provides for identification of the “Big 6” STEC serogroups as well as *E. coli* O157:H7.

Enriched samples, or a crude DNA extract prepared therefrom, are sent to Neogen via courier. The NeoSeek procedure begins with further sample processing, followed by PCR amplification. Primer extension reactions using the PCR products as templated generate allele-specific DNA products of different mass. Chip-based mass spectrometry analysis of the extension products provides allele-specific determination of target genes present in the sample. The resulting molecular profile allows identification of STEC of specific serogroups, without the need for culture confirmation. Results are available as soon as the next day.

## EQUIPMENT AND SUPPLIES

1. Stomacher, or equivalent
2. Balance, 2,000 g capacity, 0.1 g sensitivity
3. Incubator,  $42 \pm 1^\circ\text{C}$
4. Heat block with insert for 1.5 mL microcentrifuge tubes,  $100 \pm 4^\circ\text{C}$
5. Microcentrifuge, 12,000  $\times$  g
6. Microcentrifuge tubes, 1.5 mL
7. Microcentrifuge tube cap locks
8. Stomacher-type bags for sample enrichment, with filters
9. Serological pipettes, 1–10 mL
10. Micropipettor and tips, 1–1,000  $\mu\text{L}$
11. Parafilm® sealing film
12. STEC sample shipping kit, for packaging and shipping samples (Neogen item 9287, or equivalent)

## MEDIA AND REAGENTS

1. Modified tryptone soya broth with casamino acids (mTSB, ref. 1, 2; Neogen item 7467)
2. Phosphate buffered saline (PBS)
3. Water, sterile

## SAMPLE ENRICHMENT

**NOTE:** Sample enrichment is performed in the customer's laboratory.

1. Place  $325 \pm 32.5$  g raw beef trim sample into a sterile bag with mesh filter. Add  $975 \pm 19.5$  mL mTSB. Pummel, blend or hand massage until well mixed.
2. Incubate enrichment culture for **18–20 hours** at  $42 \pm 1^\circ\text{C}$ .

## SAMPLE PREPARATION

**NOTE:** Sample preparation is performed by the customer or Neogen.

### Enrichment aliquot

1. Turn on heat block to  $100 \pm 4^\circ\text{C}$ .
2. Transfer 1 mL of enrichment aliquot to a 1.5 mL microcentrifuge tube.
3. Centrifuge the enrichment aliquot at 12,000  $\times$  g for **30–60 seconds**.
4. Remove supernatant. Discard supernatant into a suitable waste container to be bleached prior to final disposal. Add 1 mL 1x PBS (or use same volume as original aliquot) to the tube and re-suspend pellet.
5. Centrifuge again at 12,000  $\times$  g for **30–60 seconds** to pellet cells. Remove supernatant into suitable waste container. Add 1 mL sterile ultra-pure water (or same volume as original aliquot) and completely re-suspend pellet.
6. Slip tube cap lock on the tube lid. The tube cap lock prevents the tube from opening during heating. Incubate tubes at  $100 \pm 4^\circ\text{C}$  for **10 minutes** in heat block.
7. Carefully remove tubes from the heat block and let them cool for **15–20 minutes** at room temperature. One cooled, press on lids to ensure tubes are fully shut, then remove the lid locks.
8. Centrifuge the lysate at 12,000  $\times$  g for **1 minute** to pellet cell debris. The supernatant is used as the sample.

### Isolated colony

1. Turn on heat block to  $100 \pm 4^\circ\text{C}$ .
2. From a plate or slant, pick one isolated colony with a sterile inoculating loop. Transfer colony material to a 1.5 mL microcentrifuge tube containing 100  $\mu\text{L}$  of sterile water. Re-suspend thoroughly.
3. Heat at  $95\text{--}100^\circ\text{C}$  for **10 minutes**. Secure tube using a tube cap lock. The tube cap lock

- prevents the tube from opening during heating.
- Remove tube from the heat block and allow to cool for **15–20 minutes** at room temperature. Once cooled, press on the lid to ensure that the tube is fully shut, then remove the lid lock.
  - Centrifuge in a microcentrifuge for **1 minute** at approximately 12,000 x g to pellet cell debris. The supernatant is used as the sample

### **PACKAGING AND SHIPPING**

**NOTE:** Enrichment aliquots or isolates, or crude DNA extracts prepared therefrom, are shipped to Neogen for analysis following the directions below.

#### **Using Neogen's STEC Shipping Kit (Neogen item 9287)**

- Remove the foam tube rack from the plastic bag.
- Place the sealed tube(s) into the foam tube rack. Ensure that they are secure.
- Place the foam rack back into the plastic bag with the absorbent pad, then seal the bag.
- Place the sealed bag into the small sample box and close the box.
- Place the closed box into the insulated foam shipping box along with the two frozen ice bricks.
- Replace the foam shipping box lid.
- Place a copy of the NeoSeek submission form on the top of the foam lid.
- Peel the back off of the UN 3373 hazard label and place on the outside of the cardboard box.  
**NOTE:** UN 3373 is required for live cultures (e.g., if sending sample enrichment or following isolated colony process from above). UN 3373 is recommended, but not required, for heat-killed organisms (e.g., enrichment aliquot process from above).
- Seal the box with tape and ship to:

**NeoSeek**  
**4131 N. 48<sup>th</sup> Street**  
**Lincoln, NE 68504 USA**

#### **Using an alternative shipping kit**

- Place the sealed tube(s) into a foam tube rack or other rack that will hold the tubes securely.
- Place the rack into a plastic bag with an absorbent pad, then seal the bag.
- Place the sealed bag into a small box or canister and close.
- Place the closed box or canister into an insulated shipping box along with two frozen ice bricks.
- Replace the shipping box lid.
- Place a copy of the NeoSeek submission form on the top of the foam lid.
- Peel the back off of the UN 3373 hazard label and place on the outside of the cardboard box.  
**NOTE:** UN 3373 is required for live cultures (e.g., if sending sample enrichment or following isolated colony process from above). UN 3373 is recommended, but not required, for heat-killed organisms (e.g., enrichment aliquot process from above).
- Seal the box with tape and ship to:

**NeoSeek**  
**4131 N. 48<sup>th</sup> Street**  
**Lincoln, NE 68504 USA**

**NOTE:** Samples must be shipped by expedited service for next-day morning delivery. International shipments require U.S. Department of Agriculture documents on company letterhead. Department of Transportation Division 6.2 Category B Biological Substances do not require hazardous materials shipping declaration papers.

## INTERPRETATION OF RESULTS

A result of **“STEC”** indicates that one of the seven target O groups was detected and that the molecular profile obtained is consistent with the organism detected being a potential EHEC (minimally stx+ and eae+).

A result of **“Non”** indicates that one of the seven target O groups was detected but that the molecular profile of the organism detected does not meet the minimal definition of an EHEC.

A result of **“Not Detected”** indicates that a particular O group was not detected.

## CONFIRMATION

The NeoSeek STEC method does not require independent confirmation of results. If desired, results may be confirmed by culture from the original enrichment broths following the most current version of the USDA-FSIS reference procedures for *E. coli* O157:H7<sup>[1]</sup> and non-O157 STEC<sup>[2]</sup>.

---

## References

<sup>[1]</sup> USDA-FSIS (2015) Microbiology Laboratory Guidebook, chapter 5.09

<https://www.fsis.usda.gov/wps/wcm/connect/51507fdb-dded-47f7-862d-ad80c3ee1738/MLG-5.pdf?MOD=AJPERES>

<sup>[2]</sup> USDA-FSIS (2014) Microbiology Laboratory Guidebook, chapter 5B.05

<https://www.fsis.usda.gov/wps/wcm/connect/7ffc02b5-3d33-4a79-b50c-81f208893204/MLG-5B.pdf?MOD=AJPERES>

**CUSTOMER SERVICE**

Neogen Customer Assistance and Technical Service can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

**SDS INFORMATION AVAILABLE**

Safety data sheets (SDS) are available for this test kit, and all of Neogen's Food Safety test kits, at [foodsafety.neogen.com](http://foodsafety.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

**TERMS AND CONDITIONS**

For Neogen's full terms and conditions, please visit [www.neogen.com/en/terms-and-conditions](http://www.neogen.com/en/terms-and-conditions).

**WARRANTY**

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## PRODUCTS AVAILABLE

NSSTEC **NeoSeek STEC method** (service performed by Neogen)

9287 **STEC Shipping Kit**

7167 **Modified tryptic soy broth with acid digest of casein**

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